

User Guide

Standard Quantitative Cartridge Kit (C105201-Q)

A. Specifications

Specifications	Description
Sample Volume Requirement	≥ 2 μ L
DNA Sizing Range	20-5000 bp
DNA Sizing Precision*	2% CV
Fragment DNA Concentration Range**	0.2 ng/ μ L to 50 ng/ μ L
Smear DNA Concentration Range**	1 ng/ μ L to 50 ng/ μ L
DNA Quantitative Precision*	4% CV

* Determined by the DNA ladder (15-622 DNA Size Marker, C109200) as sample.

** Best Quantitative Range: 0.5-5 ng/ μ L (Fragment DNA); 2-20 ng/ μ L (Smear DNA)

B. Kit components and storage conditions

Contents	Storage Conditions
Standard Cartridges (S2)	4°C~27°C (Do Not Freeze)
Quantitative Marker (20 bp & 1500 bp, C109109-500Q)*	-20 °C
Size Marker (15-622bp, C109200-100A)*	-20 °C
Separation Buffer (C104406)	4°C~27°C
Dilution Buffer (C104405)	4°C~27°C
Mineral Oil (C104407)	4°C~27°C

*Short-Term Storage (\leq 3 months): 4°C ~ 27°C ; Long-Term Storage: -20°C

C. Additional reagents (optional)

Contents	Storage Conditions
Quantitative Marker (20 bp & 5000 bp, C109102-500Q)*	-20 °C
Size Marker (50 bp plus, C109300-500A)	-20 °C

*Short-Term Storage (\leq 3 months): 4°C ~ 27°C ; Long-Term Storage: -20°C

D. Cartridge preparation

New cartridge: please follow cartridge unpacking guide to unpack cartridges and do HV check.

§ Note: Please mix 5 μ L Quantitative Marker (20 bp-1500 bp, C109109-500Q) and 15 μ L dilution buffer in 0.2 mL PCR tube for cartridge calibration. Please place Quantitative marker mixture at corresponding position.

Instrument	Position
<i>Qsep₄₀₀</i>	AM01, 04, 07, 10
<i>Qsep₁₀₀</i>	MD1
<i>Qsep₁</i>	8-wells: M 12-wells: M

E. Sample preparation

- When Sample size is between 20 bp to 1500 bp, please follow the instruction of Section E-1.
- When Sample size is between 1500 bp to 5000 bp, please follow the instruction of Section E-2.

E-1: Sample size is between 20 bp to 1500 bp:

Material requirement:

- Quantitative Marker (20 bp & 1500 bp, C109109-500Q)
- 15-622 bp Size Marker (C109200-100A)
- Dilution Buffer (C104405)

Size Marker preparation:

Mix Quantitative Marker with Size Marker according to below table:

Materials	Volume (μ L)
Quantitative Marker (20 & 1500 bp, C109109-500Q)	5
Size Marker (15-622 bp) (C109200-100A)	10
Dilution Buffer (C104405)	5
Total	20

Place Size Marker mixture at corresponding position.

Instrument	Position
<i>Qsep₄₀₀</i>	SM01, 04, 07, 10
<i>Qsep₁₀₀</i>	MA2
<i>Qsep₁</i>	8-wells: S4 12-wells: S8


Qsep₄₀₀



Qsep₁₀₀



Qsep₁



8-wells sample tray

12-wells sample tray

§ Note:

Please set up alignment marker as 20-1.5K and S4/S8 position as Size Marker and C109200 on *Qsep₁*.

Sample preparation:

Mix Quantitative Marker with samples according to below table:

Materials	Volume (μL)
Quantitative Marker (20& 1500 bp, C109109-500Q)	5
Sample	X (2~15)
H ₂ O	15-X (0~13)
Total	20

§ **Note:** Size Marker should be placed at regular 0.2 mL PCR tube. For Qsep₁ 12-wells sample tray, Size Marker and samples should be loaded at 0.1 mL strip tube (C104252).

E-2: Sample size is between 1500 bp to 5000 bp:

Material requirement:

- Quantitative Marker (20 bp & 5000 bp, C109102-500Q)
- 50 bp plus Size Marker (C109300-500A)
- Dilution Buffer (C104405)

Size Marker preparation:

Mix Quantitative Marker with size marker according to below table:

Materials	Volume (μL)
Quantitative Marker (20& 5000 bp, C109102-500Q)	5
Size Marker (50 bp plus) (C109300)	10
Dilution Buffer (C104405)	5
Total	20

Place size marker mixture at corresponding position.

Instrument	Position
Qsep ₄₀₀	SM01, 04, 07, 10
Qsep ₁₀₀	MB2
Qsep ₁	8-wells: S4 12-wells: S8


Qsep₄₀₀




Qsep₁₀₀



Qsep₁



8-wells sample tray



12-wells sample tray

§ **Note:**

Please set up alignment marker as 20-5K and S4/ S8 position as Size Marker and C109300 on Qsep₁.

Sample preparation:


Mix Quantitative Marker with samples according to below table:

Materials	Volume (μL)
Quantitative Marker (20& 5000 bp, C1091102-500Q)	5
Sample	X (2~15)
H ₂ O	15-X (0~13)
Total	20

§ **Note:** Size Marker mixture should be placed at regular 0.2 mL PCR tube. For Qsep₁ 12-wells sample tray, Size Marker and samples should be loaded at 0.1 mL strip tube (C104252).

F. Software Operation

1. Place sample and select corresponding position.
2. Enter Sample information (optional).
3. Set up analytic method.
 - 3-1. Select analytic method "Quantitative".
 - 3-2. Enter sample volume (X).
 - 3-3. Select Quantitative Marker (20-1.5K or 20-5K) which you mixed with samples and Size Markers.
 - 3-4. Select method
 - If samples are DNA fragments, please use method "M4-10-08-160".
 - If samples are smear DNA, please use method "gDNA (NGS)".
 - 3.5. Click OK.



4. Click Run to start analysis.

