

User Guide

High Sensitivity Quantitative Cartridge Kit (C105105-Q/ C105205-Q)

A. Specifications

Specifications	Description
Sample Volume Requirement	≥ 2 μL
DNA Sizing Range	20-1500 bp
DNA Sizing Precision*	2% CV
Fragment DNA Concentration Range**	10 pg/ μL to 500 pg/ μL
Smear DNA Concentration Range**	50 pg/ μL to 1000 pg/ μL
DNA Quantitative Precision*	4% CV

* Determined by the DNA ladder (15-622 DNA Size Marker, C109200) as sample.

** Best Quantitative Range: 10-200 pg/μL (Fragment DNA); 50-500 pg/ μL (Smear DNA)

B. Kit components and storage conditions

Contents	Storage Conditions
High Sensitivity Cartridge (N1)	4°C~27°C (Do Not Freeze)
Quantitative Marker (20 bp & 1500 bp, C109109-60Q)*	-20 °C
Size Marker (15-622 bp, C109200-50)*	-20 °C
Separation Buffer (C104406)	4°C~27°C
Dilution Buffer (C104405)	4°C~27°C
Mineral Oil (C104407)	4°C~27°C

*Short-Term Storage (≤ 3 months): 4°C ~ 27°C ; Long-Term Storage: -20°C

C. Cartridge preparation

New cartridge: please follow cartridge unpacking guide to unpack cartridges and do HV check.

§ **Note:** Please mix 5 μL Quantitative Marker (20 bp-1500 bp, C109109-60Q) and 15 μL dilution buffer in 0.2 mL regular PCR tube for cartridge calibration. Place Quantitative marker mixture at corresponding position.

Instrument	Position
<i>Qsep</i> ₄₀₀	AM01, 04, 07, 10
<i>Qsep</i> ₁₀₀	MD1
<i>Qsep</i> ₁	8-wells: M 12-wells: M

D. Sample preparation

- Sample size must be between 20 bp to 1500 bp.

Material requirement:

- Quantitative Marker (20 bp & 1500 bp, C109109-60Q)
- 15-622 bp Size Marker (C109200-50)
- Dilution Buffer (C104405)

Size Marker preparation:



- Dilute Quantitative Marker (C109109-60Q) and Size Marker (C109200-50) 10x by dilution buffer (C104405).

- Mix 10x diluted Quantitative Marker and 10x diluted Size Marker according to below table:

Materials	Volume (μL)
10x diluted Quantitative Marker (20 & 1500 bp, C109109-60Q)	5
10x diluted Size Marker (15-622 bp) (C109200-50)	10
Dilution Buffer (C104405)	5
Total	20

Place size marker mixture at corresponding position.

Instrument	Position
<i>Qsep</i> ₄₀₀	SM01, 04, 07, 10
<i>Qsep</i> ₁₀₀	MA2
<i>Qsep</i> ₁	8-wells: S4 12-wells: S8



8-wells sample tray

12-wells sample tray

§ Note:

Please set up alignment marker as 20-1.5K and S4/ S8 position as Size Marker and C109200 on *Qsep*₁.

Sample preparation:


Mix 10x diluted Quantitative Marker with samples according to below table:

Materials	Volume (μL)
10x diluted Quantitative Marker (20& 1500 bp, C109109-60Q)	5
Sample	X (2~15)
H ₂ O	15-X (0~13)
Total	20

§ **Note:** Size Marker should be placed at regular 0.2 mL PCR tube. For *Qsep*₁ 12-wells sample tray of *Qsep*₁, size marker and samples should be loaded at 0.1 mL strip tube (C104252).

E. Software Operation

1. Place sample and select corresponding position.
2. Enter Sample information (optional).
3. Set up analytic method.
 - 3-1. Select analytic method "**Quantitative**".
 - 3-2. Enter sample volume (X).
 - 3-3. Select Quantitative Marker (20-1.5K)
 - 3-4. Select method "cfDNA".
 - 3.5. Click OK.



4. Click Run to start analysis.

