


## Human biology-on-a-chip • Simplified

AIM Biotech offers an easy-to-use, modular platform for incorporating 3D cell culture into your neurovascular unit research.



### Connect with AIM








#aimbiotech


www.aimbiotech.com info@aimbiotech.com



# 3D BBB Model

AIM 3D Cell Culture Chips offer a new perspective on the study of the BBB.






Schematic of 3D BBB present in each site of idenTX40 plate; Immunofluorescence staining of the 3D BBB


- Self-organization of the cells into 3D matrix
- Positive signals identified for specific brain cell types
- Live-cell as well as fixed cell staining
- Disease-modelling of neurological dysfunctions

## — Permeability Studies ——



- 2×10-6
  0 1×10-6
  0 10 25
  TNF-α (ng/ml)
  24 h treatment
  24 h treatment + 24 h recovery
- Low permeability values = tight 3D BBB model
- Perturbed BBB recovered post treatment with inflammatory cues
- Presence of interstitial flow which is useful for studying (drug-induced) toxicity, inflammation, quantification of drug transport into the CNS.

# — qPCR of BBB biomarkers ——



- Off-plate analysis of RNA, DNA etc. can be conducted with qPCR analysis
- BBB specific biomarkers are identified
- EC transporters e.g., Pgp, ABCG2, MRP4 identified

qPCR data of some markers in the 3D BBB: (a) AC & PC markers in 4Day vs 7Day BBB; (b) Transporters identified in BBB

#### References

- [1] Winkelman MA., Dai G.; Sci Adv, 2023, 9, eaaz9499
- [2] Wan Z., Zhong AX., Zhang S., Pavlou G., Coughlin MF., Shelton SE., Nguyen HT., Lorch JH., Barbie DA., Kamm RD., Small Methods, 2022, 6, 2200143
- [3] Hajal C., Offenddu GS., Shin Y., Zhang S., Morozova O., Hickman D., Knutson CG., Kamm RD.; Nat Prot, 2022, 17, 125-128
- [4] Campisi M, Shin YJ, Osaki T, Hajal C, Chiono V, Kamm RD.; Biomat. 2018, 180, 117-119

