


PRODUCT INFORMATION

GENERAL INFORMATION

- Product Name :** pET Expression Vector pETIA
Code No. : DV215
Size : 15 µg (lyophilized plasmid contains salt of TE buffer)
Storage : This product is shipped at ambient temperature. Upon receipt, store at - 20 °C
Reconstitution : Resuspend the lyophilized pETIA with 15 µl of sterile water to make 1 µg/µl plasmid in 1 × TE buffer. After reconstitution, store at - 20 °C

Product Description : pETIA is a medium copy number, ampicillin resistant, stringent controllable T7 bacterial expression vector. The T7 expression system is one of the strongest expression systems and has been widely used with a coupling of BL21 (DE3) *E. coli* cell. T7 RNA polymerase gene is integrated in a genome of BL21(DE3) under control of lacUV5 promoter. Upon addition of isopropyl-1-thio-β-D-galactopyranoside (IPTG), T7 RNA polymerase is expressed in the BL21(DE3) cells harboring pETIA vector, and it induces a high-level protein expression from T7 promoter of pETIA. The pETIA has a lacI gene, which represses T7 RNA polymerase gene in the absence of IPTG. The regulation with lac repressor is beneficial to repress a basal level protein expression and to maintain a recombinant plasmid in BL21 (DE3) cell.

Plasmid Map:



T7 promoter :	213-229
T7 transcription start :	230
His Tag :	305-322
T7 terminator :	498-545
lacI :	2891-3973
pMB1 ori :	1695-2314
Amp (bla) :	680-1540

PRODUCT INFORMATION

Features of T7 expression vectors

BioDynamics Laboratory Inc. provides 6 kinds of T7 expression vectors, pETUA, pETBA, pETIA, pETUK, pETBK, and pETIK. These vectors have the same multicloning site and specific feature of each vector is below:

	Plasmid copy number	Replicon	Antibiotic resistance	Feature and recommendation
pETUA	high copy	pUC	ampicillin	for non-toxic protein expression
pETBA	medium copy	pMB1	ampicillin	general expression
pETIA	medium copy	pMB1	ampicillin	stringent regulation with lac repressor
pETUK	high copy	pUC	kanamycin	for non-toxic protein expression
pETBK	medium copy	pMB1	kanamycin	general expression
pETIK	medium copy	pMB1	kanamycin	stringent regulation with lac repressor

pETIA Sequence

```

GTTTGACAGC TTATCATCGA CTGCACGGTG CACCAATGCT TCTGGCGTCA GGCAGCCATC 60
GGAAGCTGTG GTATGGCTGT GCAGGTCGTA AATCACTGCA TAATTCGTGT CGCTCAAGGC 120
GCACTCCCGT TCTGGATAAT GTTTTTTGGC CCGACATCAT AACGGTCTCTG GCAAATATTC 180
TGAAATGAGC TGAGATCTCG ATCCCAGCAA ATTAATACGA CTCACTATAG GGAGACCACA 240
ACGGTTTCCC TCTAGAAATA ATTTTGTTTA ACTTTAAGAA GGAGATATAC ATATGCGGGG 300
TTCTCATCAT CATCATCATC ATGGTATGGC TAGCATGACT GGTGGACAGC AAATGGGTCG 360
GGACGATGAC GATAAGGATC CCCGGGTACC GAGCTCGAAT TCGATTTCTG CGACAAGCTT 420
AGCGGCCGCC GTTTAATCCG GCTGCTAACA AAGCCCAGAA GGAAGCTGAG TTGGCTGCTG 480
CCACCGCTGA GCAATAACTA GCATAACCCC TTGGGGCCTC TAAACGGGTC TTGAGGGGTT 540
TTTTGCTGAA AGGAGGAACT ATATCCGGAT GCGTTTCTAC AAACCTTTTT GTTTATTTTT 600
CTAAATACAT TCAAATATGT ATCCGCTCAT GAGACAATAA CCCTGATAAA TGCTTCAATA 660
ATATTGAAAA AGGAAGAGTA TGAGTATTC AATTTCCGT GTCGCCCTTA TTCCCTTTTT 720
TGCGGCATTT TGCCTTCCTG TTTTGTCTCA CCCAGAAACG CTGGTGAAAG TAAAAGATGC 780
TGAAGATCAG TTGGGTGCAC GAGTGGGTTA CATCGAACTG GATCTCAACA GCGGTAAGAT 840
CCTTGAGAGT TTTCGCCCCG AAGAACGTTT TCCAATGATG AGCACTTTTA AAGTTCTGCT 900
ATGTGGCGCG GTATTATCCC GTGTTGACGC CGGGCAAGAG CAACTCGGTC GCCGCATACA 960
CTATTCTCAG AATGACTTGG TTGAGTACTC ACCAGTCACA GAAAAGCATC TTACGGATGG 1020
CATGACAGTA AGAGAATTAT GCAGTGCTGC CATAACCATG AGTGATAACA CTGCGGCCAA 1080
CTTACTTCTG ACAACGATCG GAGGACCGAA GGAGCTAACC GCTTTTTTGC ACAACATGGG 1140
GGATCATGTA ACTCGCCTTG ATCGTTGGGA ACCGGAGCTG AATGAAGCCA TACCAAACGA 1200
CGAGCGTGAC ACCACGATGC CTACAGCAAT GGCAACAACG TTGCGCAAAC TATTAAGTGG 1260
CGAACTACTT ACTCTAGCTT CCCGGCAACA ATTAATAGAC TGGATGGAGG CGGATAAAGT 1320
TGCAGGACCA CTTCTGCGCT CGGCCCTTCC GGCTGGCTGG TTTATTGCTG ATAAATCTGG 1380
AGCCGGTGAG CGTGGGTCTC GCGGTATCAT TGCAGCACTG GGGCCAGATG GTAAGCCCTC 1440
CCGTATCGTA GTTATCTACA CGACGGGGAG TCAGGCAACT ATGGATGAAC GAAATAGACA 1500

```

PRODUCT INFORMATION

GATCGCTGAG	ATAGGTGCCT	CACTGATTAA	GCATTGGTAA	CTGTCAGACC	AAGTTTACTC	1560
ATATATACTT	TAGATTGATT	TAAAACTTCA	TTTTTAAATTT	AAAAGGATCT	AGGTGAAGAT	1620
CCTTTTTGAT	AATCTCATGA	CCAAAAATCCC	TTAACGTGAG	TTTTTCGTTCC	ACTGAGCGTC	1680
AGACCCCGTA	GAAAAAGATCA	AAGGATCTTC	TTGAGATCCT	TTTTTTCTGC	GCGTAATCTG	1740
CTGCTTGCAA	ACAAAAAAAC	CACCGCTACC	AGCGGTGGTT	TGTTTGCCGG	ATCAAGAGCT	1800
ACCAACTCTT	TTTCCGAAGG	TAACTGGCTT	CAGCAGAGCG	CAGATACCAA	ATACTGTCCT	1860
TCTAGTGTAG	CCGTAGTTAG	GCCACCACTT	CAAGAACTCT	GTAGCACCGC	CTACATACCT	1920
CGCTCTGCTA	ATCCTGTTAC	CAGTGGCTGC	TGCCAGTGGC	GATAAGTCGT	GTCTTACCGG	1980
GTTGGACTCA	AGACGATAGT	TACCGGATAA	GGCGCAGCGG	TCGGGCTGAA	CGGGGGGTTC	2040
GTGCACACAG	CCCAGCTTGG	AGCGAACGAC	CTACACCGAA	CTGAGATACC	TACAGCGTGA	2100
GCTATGAGAA	AGCGCCACGC	TTCCCGAAGG	GAGAAAGGCG	GACAGGTATC	CGGTAAGCGG	2160
CAGGGTCGGA	ACAGGAGAGC	GCACGAGGGA	GCTTCCAGGG	GGAAACGCCT	GGTATCTTTA	2220
TAGTCCTGTC	GGGTTTCGCC	ACCTCTGACT	TGAGCGTCTGA	TTTTTGTGAT	GCTCGTCAGG	2280
GGGGCGGAGC	CTATGGAAAA	ACGCCAGCAA	CGCGGCCTTT	TTACGGTTCC	TGGCCTTTTG	2340
CTGGCCTTTT	GCTCACATGT	TCTTTCCTGC	GTTATCCCCT	GATTCTGTGG	ATAACCGTAT	2400
TACCGCCTTT	GAGTGAGCTG	ATACCGCTCG	CCGCAGCCGA	ACGACCGAGC	GCAGCGAGTC	2460
AGTGAGCGAG	GAAGCGGAAG	AGCGCCTGAT	GCGGTATTTT	CTCCTTACGC	ATCTGTGCGG	2520
TATTTACACAC	CGCATATATG	GTGCACTCTC	AGTACAATCT	GCTCTGATGC	CGCATAGTTA	2580
AGCCAGTATA	CACTCCGCTA	TCGCTACGTG	ACTGGGTTCAT	GGCTGCGCCC	CGACACCCGC	2640
CAACACCCGC	TGACGCGCCC	TGACGGGCTT	GTCTGCTCCC	GGCATCCGCT	TACAGACAAG	2700
CTGTGACCGT	CTCCGGGAGC	TGCATGTGTC	AGAGGTTTTT	ACCGTCATCA	CCGAAACGCG	2760
CGAGGCAGCA	GATCAATTCG	CGCGCGAAGG	CGAAGCGGCA	TGCATTTACG	TTGACACCAT	2820
CGAATGGTGC	AAAACCTTTC	GCGGTATGGC	ATGATAGCGC	CCGGAAGAGA	GTCAATTCAG	2880
GGTGGTGAAT	GTGAAACCAG	TAACGTTATA	CGATGTCGCA	GAGTATGCCG	GTGTCTCTTA	2940
TCAGACCGTT	TCCCGCGTGG	TGAACCAGGC	CAGCCACGTT	TCTGCGAAAA	CGCGGGAAAA	3000
AGTGGAAGCG	GCGATGGCGG	AGCTGAATTA	CATTCCCAAC	CGCGTGGCAC	AACAACGGC	3060
GGGCAAACAG	TCGTTGCTGA	TTGGCGTTGC	CACCTCCAGT	CTGGCCCTGC	ACGCGCCGTC	3120
GCAAATTGTC	GCGGCGATTA	AATCTCGCGC	CGATCAACTG	GGTGCCAGCG	TGGTGGTGTG	3180
GATGGTAGAA	CGAAGCGGCG	TCGAAGCCTG	TAAAGCGGCG	GTGCACAATC	TTCTCGCGCA	3240
ACGCGTCAGT	GGGCTGATCA	TTAACTATCC	GCTGGATGAC	CAGGATGCCA	TTGCTGTGGA	3300
AGCTGCCTGC	ACTAATGTTT	CGGCGTTATT	TCTTGATGTC	TCTGACCAGA	CACCCATCAA	3360
CAGTATTATT	TTCTCCCATG	AAGACGGTAC	GCGACTGGGC	GTGGAGCATC	TGGTCGCATT	3420
GGGTACCCAG	CAAATCGCGC	TGTTAGCGGG	CCCATTAAGT	TCTGTCTCGG	CGCGTCTGCG	3480
TCTGGCTGGC	TGGCATAAAT	ATCTCACTCG	CAATCAAATT	CAGCCGATAG	CGGAACGGGA	3540
AGGCGACTGG	AGTGCCATGT	CCGGTTTTCA	ACAAACCATG	CAAATGCTGA	ATGAGGGCAT	3600
CGTTCCCACT	GCGATGCTGG	TTGCCAACGA	TCAGATGGCG	CTGGGCGCAA	TGCGCGCCAT	3660
TACCGAGTCC	GGGCTGCGCG	TTGGTGCGGA	TATCTCGGTA	GTGGGATACG	ACGATACCGA	3720
AGACAGCTCA	TGTTATATCC	CGCCGTCAAC	CACCATCAAA	CAGGATTTTC	GCCTGCTGGG	3780
GCAAACCAGC	GTGGACCGCT	TGCTGCAACT	CTCTCAGGGC	CAGGCGGTGA	AGGGCAATCA	3840
GCTGTTGCCC	GTCTCACTGG	TGAAAAGAAA	AACCACCCTG	GCGCCCAATA	CGCAAACCGC	3900
CTCTCCCCGC	GCGTTGGCCG	ATTCATTAAT	GCAGCTGGCA	CGACAGGTTT	CCCGACTGGA	3960
AAGCGGGCAG	TGAGCGCAAC	GCAATTAATG	TGAGTTAGCG	CGAATTGATC	TG	4012

PRODUCT INFORMATION

PRODUCT USAGE

Cloning of a gene to pETIA:

Below is the multiple cloning site of pETIA. To express a recombinant protein correctly, it is necessary to clone the gene of interest in frame with an N-terminal peptide of pETIA. The start codon of pETIA is boxed ATG in the below figure. Digest pETIA completely with appropriate restriction enzyme(s) to form DNA ends which can be ligated to the gene of interest. If only one restriction enzyme is used, dephosphorylation of a vector is often performed. Ligation of processed pETIA and the gene of interest can be performed by the standard procedure. The following transformation procedure should be done with non-expression hosts such as DH5 α or JM109. Recombinant plasmids derived from pETIA are selected by colony-PCR, enzyme digestion of prepared plasmids, or other methods. Sequencing of cloning portion and an insert region on the obtained plasmid is recommended to determine the correct recombinant plasmids for expression experiments.

```

                T7 promoter                                XbaI
GATCCCCGCA AATTAATACG ACTCACTATA GGGAGACCAC AACGGTTTCC CTCTAGAAAT 259
AspProAlaL ysLeuIleAr gLeuThrIle GlyArgProG lnArgPhePr oSerArgAsn
                NdeI                                6×His
AATTTTTGTTT AACTTTAAGA AGGAGATATA CATATGCGGG GTTCTCATCA TCATCATCAT 319
AsnPheVal* **Leu***Gl uGlyAspIle HisMetArgG lySerHisHi sHisHisHis
                EK                                BamHI
CATGGTATGG CTAGCATGAC TGGTGGACAG CAAATGGGTC GGGACGATGA CGATAAGGAT 379
HisGlyMetA laSerMetTh rGlyGlyGln GlnMetGlyA rgAspAspAs pAspLysAsp
                KpnI                                EcoRI                                SalI                                NotI                                ↑
CCCCGGGTAC CGAGCTCGAA TTCGATTTTCG TCGACAAGCT TAGCGGCCGC CGTTTAATCC 439
                SmaI                                SacI                                HindIII
ProArgValP roSerSerAs nSerIleSer SerThrSerL euAlaAlaAl aVal***Ser

```

EK: Enterokinase recognition sequence (AspAspAspAspLys↓)

ATG: start codon

TAA: stop codon

Protein Expression Procedure :

The following protocol is a general guide for the protein expression by using T7 expression vectors, pETUA, pETBA, pETIA, pETUK, pETBK, and pETIK, coupling with an expression host *E. coli* cell, BL21(DE3) cells or BL21(DE3)pLysS cells.

PRODUCT INFORMATION

- Before starting:

Transform BL21(DE3) or BL21(DE3)pLysS cells with the prepared expression plasmid by the standard procedure.

- ‡ Notes for transformation

1. Sometimes, expression may vary among transformants. If large and small colonies are observed in the same plate, the expressed protein may affect the growth of the *E. coli* cells.

2. If the expressed protein is toxic to *E. coli* cells, transformants may not be obtained.

In this case, repression of a basal level expression by T7 promoter may work, see “Notes for expression.”

- Expression:

1. Following transformation, pick a colony and inoculate it into 3 ml of LB medium containing the appropriate antibiotic with shaking at 37°C, overnight. For the BL21(DE3)pLysS strain, it is preferable to add chloramphenicol at a final concentration of 34 µg/ml in the overnight culture to maintain pLysS.

2. The next morning, transfer 0.5 ml of the overnight culture to a new 10 ml of LB medium containing the appropriate antibiotic to select the expression plasmid. Grow the culture with shaking at 37°C until the OD₆₀₀ reaches 0.5 (approximately 2 hrs but this depended on the expression plasmids).

When using BL21(DE3)pLys, chloramphenicol is not usually required in the short-period culture.

3. When the OD₆₀₀ reaches 0.5, transfer an aliquot (e.g., 1 ml) of the culture to a new centrifuge tube and centrifuge it to harvest cells. Store the cells at -80°C until analysis.

Add IPTG to a final concentration of 1 mM to the rest of the culture and grow the culture with shaking at 37°C for 3 hours.

The IPTG concentration and induction time are general values. It is recommended to determine the optimal condition for the target gene expression.

4. After the induction, harvest the cells. To analyze the expression, before harvesting the cells, transfer an aliquot of the culture (e.g., 1 ml) and centrifuge it to precipitate the cells.

- Analysis

1. Suspend the precipitated cells (from the 1 ml culture) in 200 µl of 1× PBS buffer.


2. Mix an aliquot of the suspension (e.g., 100 µl) with an equal volume of 2 × SDS sample buffer.

3. Heat the mixture at 85°C for 5 min, then centrifuge at 10,000 g for 10 min. Subject the supernatant (e.g., 5-25 µl) to SDS-PAGE. Western blot will help analyzing the expression of the target protein.

- 2 × SDS sample buffer : 2 % sodium dodecyl sulfate, 5 % 2-mercaptoethanol, 20 % glycerol, 0.02 % BPB, 62.5 mM Tris-HCl, pH6.8

- 1× PBS buffer.: 20 mM sodium phosphate, 150 mM sodium chloride, pH7.4

PRODUCT INFORMATION



An arrow shows the expressed 70 KDa proteins. Only induced cells expressed 70 KDa proteins.

Figure of protein expression from pETIA

A gene of 70 KDa protein was cloned into pETIA (pETIA/70K). BL21(DE3) cell was transformed with the pETIA/70K, one of colonies were cultured overnight and transferred to two tubes (#1, #2) containing culture medium. IPTG was added to only tube #2 when the OD₆₀₀ reaches 0.5. At each stage, OD₆₀₀ of the culture was determined and the same amount of cells were lysed and subjected to 10 % polyacrylamide gel SDS electrophoresis.

Lane 1: DynaMarker Protein Eco (#DM610)

Lane 2, 3 : Cells from tubes #1 and 2 before induction.

Lane 4 : Cells (tubes #1), two hours after OD0.5.

Lane 5 : Cells (tubes #2), two hours after induction

Lane 6 : Cells (tubes #1), 4 hours after OD0.5.

Lane 7 : Cells (tubes #2), 4 hours after induction

‡ Notes for expression:

1. As the T7 expression method is a high-level protein expression system, some basal level expression of the target protein will occur in uninduced cells. This is likely problematic in cases in which the target protein is toxic to *E. coli* cells. In this case, it may be necessary to decrease the basal level expression as follows:

- Use a lower-copy number T7 expression vector, pETBA, pETBK, but not pETUA, pETUK
- Use a stringent regulated expression vector, pETIA, pETIK.
- Use liquid medium and agar plates supplemented with glucose (0.5 -1 %).
Glucose is known to decrease a basal expression from *lacUV5* promoter²⁾.
- Use BL21(DE3)pLysS strain but not BL21(DE3) strain.

The T7 Lysozyme encoded in a pLysS plasmid reduces the basal level of T7 RNA polymerase Expression³⁾. This leads to suppression of the basal level expression of the target protein.

2. When expressing proteins in BL21(DE3) cells, if it takes longer time (5 hrs or more) to reach 0.5 at OD₆₀₀ after inoculating the overnight culture (0.5 ml) to a new LB medium (10 ml), the expressed protein is likely toxic to *E. coli* cells.

3. When BL21(DE3) cells lyse after induction with IPTG, the expressed protein is likely toxic to *E. coli* cells.

PRODUCT INFORMATION

Reference:

- 1) Studier, F.W. and Moffatt, B.A., *J. Mol. Biol.* 189 (1986) 113–130.
- 2) Pan, S. and Malcom, B.A., *BioTechniques* 29 (2000), 1234–1238
- 3) Moffatt, B.A. and Studier, F.W., *Cell* 49 (1987) 221-227

General reference in this Product Information

Sambrook, J. and Russell, D.W. (2001) *Molecular Cloning: A Laboratory Manual*, 3rd ed., Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

Related Products:

DV210	pET Expression Vector pETBA	DV215	pET Expression Vector pETIA
DV220	pET Expression Vector pETUK	DV230	pET Expression Vector pETBK
DV235	pET Expression Vector pETIK	DS110	DNA Ligation Kit ver. 2
DS210	Competent Cell JM109	DS220	Competent Cell DH5 α
DS225	Jet Competent Cell (DH5 α)	DS240	Competent Cell BL21
DS255	Zip Competent Cell BL21(DE3)	DS260	Competent Cell BL21(DE3)pLysS
DS500	QuickBlue Protein Staining Solution		

● Purchaser Notification

This product is manufactured based on the T7 expression system which is the subject of US patent applications assigned to Brookhaven Science Associates, LLC (BSA). The product must be used only outside the United States and its territories. Neither this product nor materials prepared used by the T7 expression system are allowed to be distributed in the US and its territories without license of BSA. Information about license regarding the T7 expression system may be obtained from the Office of Intellectual Property and Sponsored Research, Brookhaven National Laboratory, Building 185, P.O. Box 500, Upton, New York 11973-5000, USA.