

DATA SHEET for BSEP Vesicle Uptake Assay Reagents

BSEP ProVesicle Uptake Assay Reagents

Cat. #PB21001; Lot #PN-B-1001

Contents: Kit contains reagents for 200 assays

- A. BSEP Uptake Buffer: 10 mM HEPES, pH 7.5; 50 mM sucrose, 100 mM KNO₃; 12.5 mM Mg(NO₃)₂
- B. 10 x Wash Buffer: 100 mM Tris-HCl, pH 7.4; 1 M KNO₃; 500 mM sucrose
- C. 200 mM ATP
- D. 200 mM AMP

Store at 4C upon receipt.

Representative data using PB21001 showing time dependence, concentration dependence and competitive inhibition of uptake at hBSEP when using 2',7'-dichlorodihydrofluorescein diacetate (H_2DCFDA), a fluorescent substrate:



ATP-dependent H₂DCFDA (10 μ M) uptake: 1000 pmol/mg/min IC₅₀ for reference inhibitor rifampicin: 30 μ M

Vesicle Uptake Assay Protocol:

- Incubate a 95 μL reaction containing 50 μg vesicles and H₂DCFDA in BSEP Uptake Buffer (10 mM HEPES, pH 7.5; 50 mM sucrose, 100 mM KNO₃; 12.5 mM Mg(NO₃)₂) for 5 min at 37C.
- 2. Uptake was started by adding MgATP (or AMP) at a final of 5 mM and incubating for 20 min at 37C.
- 3. Reaction was terminated by transferring vesicle samples to a filter plate and washing the filter plate 6 x with ice cold washing buffer (10 mM Tris-HCl, pH 7.4; 100 mM KNO₃; 50 mM sucrose).
- 4. Filter plate was dried and fluorescent counts measured with a spectrophotometer.

Calculating ATP-Dependent Transport:

ATP-dependent transport = uptake in the presence of ATP – uptake in the presence of AMP.

This product is strictly for laboratory research use only.