

DATA SHEET for MRP/BCRP Vesicle Uptake Assay Reagents

MRP/BCRP ProVesicle Uptake Assay Reagents


Cat. #PB21002; Lot #PN-B-1001

Contents: ***Kit contains reagents for 200 assays***

- A. MRP/BCRP Uptake Buffer: 10 mM Tris-HCl, pH 7.4; 250 mM sucrose, 10 mM MgCl₂
- B. 10 x Wash Buffer: 400 mM MOPS-Tris, pH 7.0; 700 mM KCl
- C. 300 mM GSH
- D. 200 mM ATP
- E. 200 mM AMP

Store at 4°C upon receipt.

Representative data using PB21002 showing time dependence, concentration dependence and competitive inhibition of uptake at hMRP2 when using CDCF, a fluorescent substrate:



ATP-dependent CDCF (10 μ M) uptake: 40 pmol/mg/min
IC₅₀ for reference inhibitor benzobromarone: 65 μ M

Vesicle Uptake Assay Protocol:

1. Incubate a 95 μ L reaction containing 50 μ g vesicles and CDCF in MRP/BCRP Uptake Buffer (10 mM Tris-HCl, pH 7.4; 250 mM sucrose, 10 mM MgCl₂ and 3 mM GSH) for 5 min at 37°C.
2. Uptake was started by adding MgATP (or AMP) at a final of 5 mM and incubating for 20 min at 37°C.
3. Reaction was terminated by transferring vesicle samples to a filter plate and washing the filter plate 6 x with ice cold washing buffer (40 mM MOPS-Tris, pH 7.0; 70 mM KCl).
4. Filter plate was dried and fluorescent counts measured with a spectrophotometer.

Calculating ATP-Dependent Transport:

ATP-dependent transport = uptake in the presence of ATP – uptake in the presence of AMP.

This product is strictly for laboratory research use only.

