



# Mouse TGF-β1 ELISA

**660 050 096**      **1 x 96 tests**  
**660 050 192**      **2 x 96 tests**

## INTENDED USE

The mTGF-β1 ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of mouse transforming growth factor beta-1 levels in cell culture supernatants, serum, plasma, or other body fluids. **The mTGF-β1 ELISA is for research use only. Not for use in diagnostic or therapeutic procedures.**

## SUMMARY

Transforming growth factor-β (TGF-β) belongs to a family of dimeric 25 kDa polypeptides that are ubiquitously distributed in tissues and synthesized by many different cells (2). Three isoforms of transforming Growth Factor-β (TGF-beta1, beta-2 and beta-3) exist in mammals. They play critical roles in growth regulation and development. Each isoform is encoded by a unique gene on different chromosomes. All three of these growth factors are secreted by most cell types, generally in a latent form, requiring activation before they can exert biological activity. The TGF-betas possess three major activities: they inhibit proliferation of most cells, but can stimulate the growth of some mesenchymal cells; they exert immunosuppressive effects and they enhance the formation of extracellular matrix. Two types of membrane receptors possessing kinase activity are involved in signal transduction. The TGF-betas are involved in wound repair processes and in starting inflammatory reaction and then in the resolution through chemotactic attraction of inflammatory cells and fibroblast (3).

TGF-β1 is the first recognized transforming growth factor (1), its subunits of each 12.5 kDa are bound via disulphide bridges. TGF-β1 is inhibitive to T- and B cell proliferation as well as to maturation and activation of macrophages. It furthermore inhibits activity of natural killer cells and lymphokine activated killer cells and blocks production of cytokines.


## PRINCIPLES OF THE TEST

An anti-mTGF-β1 coating antibody is adsorbed onto microwells.

mTGF-β1 present in the sample or standard binds to antibodies adsorbed to the microwells; a HRP-conjugated monoclonal anti-mTGF-β1 antibody is added and binds to mTGF-β1 captured by the first antibody.

Following incubation unbound enzyme conjugated anti-mTGF-β1 is removed during a wash step and substrate solution reactive with HRP is added to the wells.

A coloured product is formed in proportion to the amount of mTGF-β1 present in the sample. The reaction is terminated by addition of acid and absorbance is measured at 450 nm. A standard curve is prepared from seven mTGF-β1 standard dilutions and mTGF-β1 sample concentration determined.



**REAGENTS PROVIDED**

REAGENTS (store at 2-8°C)	QTY	QTY	RECONSTITUTION
	1 plate	1 plates	
96-wells precoated microtiter plate	1	2	Ready-to-use
Plate covers	2	4	
mTGF-β1 Standard	2 vials	4 vials	See label on the vial
HRP conjugate monoclonal antibody to mTGF-β1	1 vial	2 vials	(0.1 ml)
Assay Buffer Concentrate*	1 vial	2 vials	(5 ml) 20X concentrate. Dilute in distilled water
Wash Buffer Concentrate*	1 bottle	2 bottles	(50 ml) 20X concentrate. Dilute in distilled water
Substrate Solution	1 vial	2 vials	(15 ml) Ready-to-use
Stop Solution (1 M Phosphoric acid)	1 vial	2 vials	(12 ml) Ready-to-use
1N HCl (pretreatment of samples)	1 vial	2 vials	(3 ml) Ready to use
1N NaOH (pretreatment of samples)	1 vial	2 vials	(3 ml) Ready to use
Blue Dye	1 vial	2 vials	(0.4 ml) Make a 1/250 dilution in the appropriate diluent
Green Dye	1 vial	2 vials	(0.4 ml) Make a 1/100 dilution in the appropriate diluent

\* contains preservative

**MATERIALS REQUIRED BUT NOT PROVIDED**

- 5 ml and 10 ml graduated pipettes
- 5 µl to 1000 µl adjustable single channel micropipettes with disposable tips
- 50 µl to 300 µl adjustable multichannel micropipette with disposable tips
- Multichannel micropipette reservoir
- Beakers, flasks, cylinders necessary for preparation of reagents
- Device for delivery of wash solution (multichannel wash bottle or automatic wash system)
- Microwell strip reader capable of reading at 450 nm (620 nm as optional reference wave length)
- Glass-distilled or deionized water
- Statistical calculator with program to perform linear regression analysis

**SAFETY**

- Reagents are intended for research use only and are not for use in diagnostic or therapeutic procedures.
- Do not mix or substitute reagents with those from other lots or other sources.
- Do not use kit reagents beyond expiration date on label.
- Do not expose kit reagents to strong light during storage or incubation.
- Do not pipette by mouth.
- Do not eat or smoke in areas where kit reagents or samples are handled.
- Avoid contact of skin or mucous membranes with kit reagents or specimens.
- Rubber or disposable latex gloves should be worn while handling kit reagents or specimens.
- Reagents containing thimerosal as preservative may be toxic if ingested.
- Avoid contact of substrate solutions with oxidizing agents and metal.
- Avoid splashing or generation of aerosols.
- In order to avoid microbial contamination or cross-contamination of reagents or specimens which may invalidate the test use disposable pipette tips and/or pipettes.
- Use clean, dedicated reagent trays for dispensing the conjugate and substrate reagents.
- Exposure to acids will inactivate the conjugate.
- Glass-distilled water or deionized water must be used for reagent preparation.
- Substrate solutions must be at room temperature prior to use.
- Decontaminate and dispose specimens and all potentially contaminated materials as if they could contain infectious agents. The preferred method of decontamination is autoclaving for a minimum of 1 hour at 121.5°C.
- Liquid wastes not containing acid and neutralized waste may be mixed with sodium hypochlorite in volumes such that the final mixture contains 1.0 % sodium hypochlorite. Allow 30 minutes for effective decontamination. Liquid waste containing acid must be neutralized prior to the addition of sodium hypochlorite.

## STORAGE INSTRUCTIONS

Store kit reagents between 2° and 8°C. Immediately after use remaining reagents should be returned to cold storage (2° to 8°C). Expiry of the kit and reagents is stated on labels. The expiry of the kit components can only be guaranteed if the components are stored properly, and if, in case of repeated use of one component, the reagent is not contaminated by the first handling.

## SPECIMEN COLLECTION, PROCESSING AND STORAGE

Cell culture supernatants, murine serum and plasma, or other biological samples will be suitable for use in the assay. Remove serum from the clot or red cells, respectively, as soon as possible after clotting and separation.

Samples containing a visible precipitate must be clarified prior to use in the assay. Do not use grossly haemolysed or lipemic specimens.

Clinical samples should be kept at 2° to 8°C and separated rapidly (within half an hour) before storing at -20°C to avoid loss of bioactivity. Avoid repeated freeze-thaw cycles.

## PREPARATION OF REAGENTS

### 1. Wash Buffer

If crystals have formed in the Wash Buffer Concentrate, warm it gently until they have completely dissolved.

Pour entire contents (50 ml) of the **Wash Buffer Concentrate** into a clean 1,000 ml graduated cylinder. Bring final volume to 1,000 ml with glass-distilled or deionized water. Mix gently to avoid foaming. The pH of the final solution should adjust to 7.4.

Transfer to a clean wash bottle and store at 2° to 25°C. Please note that Wash Buffer is stable for 30 days. Wash Buffer may be prepared as needed according to the following table:

Number of Strips	Wash Buffer Concentrate (ml)	Distilled Water (ml)
1 - 6	25	475
1 - 12	50	950

### 2. Assay Buffer

Mix the contents of the bottle well. Add contents of **Assay Buffer Concentrate** (5.0 ml) to 95 ml distilled or deionized water and mix gently to avoid foaming. Store at 2° to 8°C. Please note that the Assay Buffer is stable for 30 days. Assay Buffer may be prepared as needed according to the following table:

Number of Strips	Assay Buffer Concentrate (ml)	Distilled Water (ml)
1 - 6	2.5	47.5
1 - 12	5.0	95.0

### 3. Preparation of HRP-Conjugate

The **HRP-Conjugate** must be diluted 1:100 with **Assay Buffer** (reagent B) just prior to use in a clean plastic test tube.

Please note that the HRP-Conjugate should be used within 30 minutes after dilution. HRP-Conjugate may be prepared as needed according to the following table:

Number of Strips	HRP-Conjugate(ml)	Assay Buffer (ml)
1 - 6	0.03	2.97
1 - 12	0.06	5.94

**4. Preparation of Standard**

Dilute mTGFβ1 **Standard** by addition of **Assay Buffer**. Dilution volume is stated on the label of the standard vial. Mix gently to insure complete homogeneity.

**5. Addition of Colour-giving Dyes**

In order to help our customers to avoid any mistakes in pipetting, DIACLONE now offers a new tool that helps to monitor the addition of even very small volumes of a solution to the reaction well by giving distinctive colours to each step of the ELISA procedure.

This procedure is optional, does not in any way interfere with the test results, and is designed to help the customer with the performance of the test, but can also be omitted, just following the instruction booklet.

Alternatively, the dye solutions from the stocks provided (**Blue-Dye, Green-Dye**) can be added to the reagents according to the following guidelines:

**A. Assay Buffer:**

Before sample dilution add the **Blue-Dye** at a dilution of 1:250 (see table below) to the appropriate diluent (1x) according to the test protocol. After addition of **Blue-Dye**, proceed according to the instruction booklet.

5 ml Assay Buffer	20 µl <b>Blue-Dye</b>
12 ml Assay Buffer	48 µl <b>Blue-Dye</b>
50 ml Assay Buffer	200 µl <b>Blue-Dye</b>

**B. HRP-Conjugate:** Before dilution of the concentrated conjugate, add the **Green-Dye** at a dilution of 1:100 (see table below) to the Assay Buffer (1x) used for the final conjugate dilution. Proceed after addition of **Green-Dye** according to the instruction booklet, preparation of HRP-conjugate.

3 ml Assay Buffer	30 µl <b>Green-Dye</b>
6 ml Assay Buffer	60 µl <b>Green-Dye</b>

**TEST PROTOCOL**

- a. Prepare your samples before starting with the test procedure. Dilute serum, plasma and cell culture samples 1:10 with Assay Buffer (180 µl Assay Buffer + 20 µl sample). Add 20 µl 1N HCl to 200 µl of prediluted sample, mix and incubate for 1 hour at room temperature. Neutralize by addition of 20 µl 1N NaOH.
- b. Mix all reagents thoroughly without foaming before use.
- c. Determine the number of microwell strips required to test the desired number of samples plus appropriate number of wells needed for running blanks and standards. Each sample, standard, blank, and optional control sample should be assayed in duplicate. Remove sufficient **Microwell Strips coated with Antibody** to mouse TGF-β1 from the aluminium pouch immediately prior to use. Return any wells not required for this assay with desiccant to the pouch. Seal tightly.
- d. Wash the microwell strips twice with approximately 400 µl **Wash Buffer** per well with thorough aspiration of microwell contents between washes. Allow the Wash Buffer to sit in the wells for about 10-15 seconds before aspiration. Take care not to scratch the surface of the microwells.  
  
After the last wash, tap microwell strips on absorbent pad or paper towel to remove excess Wash Buffer. Use the microwell strips immediately after washing or place upside down on a wet absorbent paper for not longer than 15 minutes. Do not allow wells to dry.
- e. Add 100 µl of **Assay Buffer**, in duplicate, to all standard wells. Prepare standard dilutions by pipetting 100 µl of **mTGF-β1 Standard**, in duplicate, into well A1 and A2 (see Figure 1 and 2). Mix the contents by repeated aspiration and ejection, and transfer 100 µl to well B1 and B2, respectively. Take care not to scratch the inner surface of the microwells. Mix the contents by repeated aspiration and ejection, and transfer 100 µl to well C1 and C2, respectively. Continue this procedure four times, creating two rows of mTGF-β1 standard dilutions ranging from 30 to 0.5 ng/ml. Discard 100 µl of the contents from the last microwells (G1, G2) used.

Figure 1. Preparation of mTGF-β1 standard dilutions:




Figure 2. Diagram depicting an example of the arrangement of blanks, standards and samples in the microwell strips:

	Standard Concentrations ng/mL		Sample wells									
	1	2	3	4	5	6	7	8	9	10	11	12
A	30	30										
B	15	15										
C	7.5	7.5										
D	3.75	3.75										
E	1.88	1.88										
F	0.94	0.94										
G	0.47	0.47										
H	Blank	Blank										

- f. Add 100 µl of **Assay Buffer**, in duplicate, to the blank wells.
  - g. Add 100 µl of each pretreated **Sample**, in duplicate, to the designated wells.
  - h. Prepare **HRP-Conjugate**. (Refer to preparation of reagents)
  - i. Add 50 µl of diluted **HRP-Conjugate** to all wells.
  - j. Cover with a **Plate Cover** and incubate at room temperature (18°C to 25°C) for 4 hours on a **rotator set at 100 rpm**
  - k. Remove **Plate Cover** and empty wells. Wash microwell strips 3 times according to point d. of the test protocol. Proceed immediately to the next step.
  - l. Pipette 100 µl of **TMB Substrate Solution** to all wells, including the blank wells.
  - m. Incubate the microwell strips at room temperature (18° to 25°C) for about 10 minutes, if available on a rotator set at 100 rpm. Avoid direct exposure to intense light.  
**The colour development on the plate should be monitored and the substrate reaction stopped (see point o. of this protocol) before positive wells are no longer properly recordable.**
- It is recommended to add the stop solution when the highest standard has developed a dark blue colour. Alternatively the colour development can be monitored by the ELISA reader at 620 nm. The substrate reaction should be stopped as soon as an OD of 0.6 – 0.65 is reached.
- n. Stop the enzyme reaction by quickly pipetting 100 µl of **Stop Solution** into each well, including the blank wells. It is important that the Stop Solution is spread quickly and uniformly throughout the microwells to completely inactivate the enzyme. Results must be read immediately after the Stop Solution is added or within one hour if the microwell strips are stored at 2 - 8°C in the dark.
  - o. Read absorbance of each microwell on a spectro-photometer using 450 nm as the primary wave length (optionally 620 nm as the reference wave length; 610 nm to 650 nm is acceptable). Blank the plate reader according to the manufacturer's instructions by using the blank wells. Determine the absorbance of both, the samples and the mTGF-β1 standards.


### CALCULATION OF RESULTS

- Calculate the average absorbance values for each set of duplicate standards and samples. Duplicates should be within 20 per cent of the mean.
- Create a standard curve by plotting the mean absorbance for each standard concentration on the ordinate against the mTGF-β1 concentration on the abscissa. Draw a best fit curve through the points of the graph.
- To determine the concentration of circulating mTGF-β1 for each sample, first find the mean absorbance value on the ordinate and extend a horizontal line to the standard curve. At the point of intersection, extend a vertical line to the abscissa and read the corresponding mTGF-β1 concentration.
- **For samples which have been diluted according to the instructions given in this manual (e.g. 1:12) the concentration read from the standard curve must be multiplied by the dilution factor (e.g. x12).**

**Note:** Calculation of samples with an O.D. exceeding the range of the standard curve may result in incorrect low mTGF-β1 levels. Such samples require further dilution with Assay Buffer in order to precisely quantitate the actual mTGF-β1 level.

- It is suggested that each testing facility establishes a control sample of known mTGF-β1 concentration and runs this additional control with each assay. If the values obtained are not within the expected range of the control, the assay results may be invalid.
- A representative standard curve is shown in Figure 3. This curve cannot be used to derive test results. Every laboratory must prepare a standard curve for each group of microwell strips assayed.

Figure 3. Representative standard curve for mTGF-β1 ELISA. Recombinant soluble mTGF-β1 was diluted in Assay Buffer; symbols represent the mean of three parallel titrations. Do not use this standard curve to derive test results. A standard curve must be run for each group of microwell strips assayed.



Typical data using the mTGF-β1 ELISA

Measuring wavelength: 450 nm  
 Reference wavelength: 620 nm

Standard	mTGF-β1 Concentration (ng/ml)	O.D. (450 nm)	O.D. Mean	C.V. (%)
1	30	2.165	2.1515	0.6
	30	2.138		
2	15	1.364	1.379	1.1
	15	1.394		
3	7.5	0.823	0.822	0.1
	7.5	0.821		
4	3.75	0.462	0.4695	1.6
	3.75	0.477		
5	1.88	0.244	0.245	0.4
	1.88	0.246		
6	0.94	0.123	0.126	2.4
	0.94	0.129		
7	0.47	0.070	0.073	4.1
	0.47	0.076		
Blank	0	0.020	0.0185	8.1
	0	0.017		

**LIMITATIONS OF THE PROCEDURE**

- Since exact conditions may vary from assay to assay, a standard curve must be established for every run.
- Bacterial or fungal contamination of either screen samples or reagents or cross-contamination between reagents may cause erroneous results.
- Disposable pipette tips, flasks or glassware are preferred, reusable glassware must be washed and thoroughly rinsed of all detergent before use.
- Improper or insufficient washing at any stage of the procedure will result in either false positive or false negative results. Completely empty wells before dispensing fresh Wash Buffer, fill with Wash Buffer as indicated for each wash cycle and do not allow wells to sit uncovered or dry for extended periods.

**PERFORMANCE CHARACTERISTICS**

**Sensitivity**

The limit of detection for recombinant mTGF-β1, spiked into normal human serum, defined as the analyte concentration resulting in an absorption significantly higher than the absorption of the dilution medium (mean plus two standard deviations) was determined to be 48 pg/ml (mean of three independent assays).

**Reproducibility**

**a. Intra-assay**

Intra-assay variability was determined by 6 replicates of 7 serum samples. The average coefficient of variation was 5.4%.

**b. Inter-assay**

Inter-assay variability was determined by 18 replicates of 7 serum samples. The average coefficient of variation was 10.6%.

**Specificity**

The ELISA was shown to specifically detect mTGF-β1. No cross reactivity was apparent with serum factors spiked into mouse serum, notably there was no cross reactivity with murine TGF-β2.

**Spike Recovery**

Spike recovery was determined to be 88% in average of recombinant mTGF-β1 spiked into various samples.

**Dilution Linearity**

Linearity of dilution was measured in various samples. Recovery was 117% in average.

**BIBLIOGRAPHY**

1. Chambaz EM, Souchelnitskiy S, Pellerin S, Defaye G, Cochet C, Feige JJ. Transforming growth factors-beta s: a multifunctional cytokine family. Implication in the regulation of adrenocortical cell endocrine functions. Horm Res 1996;45(3-5):22-226.
2. Kropf J, JO Schurek, A Wollner, and AM Gressner. Immunological measurement of transforming growth factor-beta I (TGF-β1) in blood; assay development and comparison. Clinical Chemistry 1997;43(10):1965-1974.
3. Lawrence DA. Transforming growth factor-beta: a general review Eur Cytokine Netw 1996 Sep;7(3):363-374.

**REAGENT PREPARATION SUMMARY**

**A. Wash Buffer** Add **Wash Buffer Concentrate** 20 x (50 ml) in to 950 ml distilled water.

<b>B. Assay Buffer</b>	Number of Strips	Assay Buffer Concentr. (ml)	Distilled Water (ml)
	1 - 6	2.5	47.5
	1 - 12	5.0	95.0

<b>C. HRP-Conjugate</b>	Number of Strips	HRP- Conjugate (ml)	Assay Buffer (ml)
	1 - 6	0.03	2.97
	1 - 12	0.06	5.94

**D. Standard** Dilute Standard by addition of **Assay Buffer** according to vial label.

**TEST PROTOCOL SUMMARY**

- Pretreat Samples
- Wash Microwell Strips twice with **Wash Buffer**
- Add 100 µl **Assay Buffer**, in duplicate, to all standard wells
- Pipette 100 µl **mTGF-β1 Standard** into the first standard wells and create standard dilutions ranging from 30 to 0.5 ng/ml by transferring 100 µl from well to well; discard 100µl from the last wells.
- Add 100 µl **Assay Buffer**, in duplicate, to the blank wells
- Add 100 µl pretreated **Sample** to designated wells
- Prepare **HRP-Conjugate**



- Add 50 µl diluted **HRP-Conjugate** to all wells
- Cover microwell strips and incubate 4 hours at room temperature (18° to 25°C) **on a rotator**
- Empty and wash microwell strips 3 times with **Wash Buffer**
- Add 100 µl of **TMB Substrate Solution** to all wells including blank wells
- Incubate the microwell strips for about 10 minutes at room temperature (18° to 25°C)
- Add 100 µl **Stop Solution** to all wells including blank wells
  
- Blank microwell reader and measure colour intensity at 450 nm

Note: For samples which have been diluted according to the instructions given in this manual 1:12 the concentration read from the standard curve must be multiplied by the dilution factor (x12x). Calculation of samples with an O.D. exceeding 2.0 may result in incorrect low mTGF-β1 levels. Such samples require further dilution with Assay Buffer in order to precisely quantitate the actual mTGF-β1 level.