Murine VEGF-A ELISA Kit

	4.5	•	
Inetri	uctions	? t∩r	пед
เมเงนเ	มนเบบ เจ	וטו כ	usc

Catalogue numbers: 1x96 tests: 660.140.096

2x96 tests: 660.140.192

For research use only

Fast Track Your Research.....

Table of Contents

1.	Intended use	2
2.	Introduction	2
2.1.	Summary	2
2.2.	Principle of the method	.3
3.	Reagents provided and reconstitution	4
4.	Materials required but not provided	4
5.	Storage Instructions	5
6.	Specimen collection, processing & storage	5
7.	Safety & precautions for use	6
8.	Assay Preparation	7
8.1.	Assay Design	7
8.2.	Preparation of Wash Buffer	7
8.3.	Preparation of Conjugate Diluent Buffer	8
8.4.	Preparation of Standard	8
8.5.	Preparation of Biotin-Conjugate	8
8.6.	Preparation of Streptavidin-HRP	9
8.7.	Addition of Colour-giving Reagents: Blue-Dye, Green-Dye, Red-Dye	9
9.	Method	.10
10.	Data Analysis	.11
11.	Assay limitations	.12
12.	Performance Characteristics	.13
12.1.	Sensitivity	.13
12.2.	Specificity	.13
12.3.	Precision	.13
12.4.	Dilution Parallelism	14
12.5.	Spike Recovery	.14
12.6.	Stability	.15
13.	Assav Summary	16

Murine VEGF-A ELISA KIT

1. Intended use

The mouse VEGF-A ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of mouse VEGF-A. The mouse VEGF-A ELISA is for research use only. Not for diagnostic or therapeutic procedures.

2. Introduction

2.1. Summary

Normal tissue function depends on a regular supply of oxygen through the blood vessels. Understanding the formation of blood vessels has become the focus of a major research effort throughout the last decade. Vasculogenesis in the embryo is the process by which new blood vessels are generated *de novo* from primitive precursor cells. Angiogenesis is the process of new blood vessel formation from pre-existing vasculatures. It plays an essential role in development, normal tissue growth, wound healing, the female reproductive cycle (placental development, ovulation, corpus luteum) and also plays a major role in various diseases. Special interest is focused on tumor growth, since tumors cannot grow more than a few millimeters in size without developing a new blood supply. This process is described as tumor angiogenesis which is also essential for the spread and growth of tumor cell metastasis.

One of the key molecules for angiogenesis and for the survival of the endothelium is vascular endothelial growth factor (VEGF-A). It is a specific endothelial cell mitogen and a strong vascular permeability factor (VPF). VEGF-A is a heparin-binding glycoprotein, secreted as a homodimer of 45 kDa by many different cell types. VEGF-A also causes vasodilation through the nitiric oxide synthase pathway in endothelial cells and can activate migration in monocytes. Many different splice variants of VEGF-A have been described, but VEGF₁₆₅ is the most predominant protein and anchors with its heparin binding domain to extracellular matrix and to heparin sulfate. During the past few years, several other members of the VEGF family have been cloned, including VEGF-B, -C- and -D. In terms of vascular angiogenesis, which mainly is regulated by VEGF-A, lymphangiogenesis is mainly regulated by VEGF-C and -D.

VEGF-A transcription is highly activated by hypoxia and by oncogenes like H-ras and several transmembrane tyrosine kinases, such as epidermal growth factor receptor and ErbB2. Together these pathways account for a marked upregulation of VEGF-A in tumors compared to normal tissues and are often of prognostic importance and relevance.

VEGF-A can be detected in both plasma and serum samples of patients, with much higher levels in serum. Extremely high levels can be detected in the cystic brain fluid of brain tumor patients or in ascites fluid of patients. Platelets release VEGF-A upon aggregation and may be another major source of VEGF-A delivery to tumors. Several other studies have shown that association of high serum levels of VEGF-A with poor prognosis in cancer patients may be correlated with an elevated platelet count. Tumors can release cytokines and growth factors that stimulate the production of megakaryocytes in the marrow and elevate the platelet count. This can result in another, indirect increase of VEGF-A delivery to tumors.

Furthermore, VEGF-A is implicated in several other pathological conditions associated with enhanced angiogenesis or enhanced vascular permeability. Examples where VEGF-A plays an important role are psoriasis and rheumatoid arthritis, as well as the ovarian hyperstimulation syndrome. Diabetic retinopathy is associated with high intraocular levels of VEGF-A, and inhibition of VEGF-A function may result in infertility by blockage of corpus luteum function. Direct demonstration of the importance of VEGF-A in tumor growth has been achieved using dominant negative VEGF receptors to block in vivo proliferation, as well as blocking antibodies to VEGF or to one of the VEGF receptors. Interference with VEGF-A function has therefore become of major interest for drug development to block angiogenesis and metastasis. More than 110 pharmaceutical companies world-wide are involved in the development of such antagonists. Their

approaches include antagonists of VEGF-A or its receptors, selective tyrosine kinase inhibitors, targeting of drugs and toxins to VEGF receptors and gene therapy regulated by the same hypoxia pathway that controls VEGF-A production. Targeting the VEGF signalling may be of major therapeutic importance for many diseases and serves as a basis for the design of future (anti)-angiogenic treatments.

2.2. Principle of the method

An anti-mouse VEGF-A coating antibody is adsorbed onto microwells.

Figure 1




Figure 2

Mouse VEGF-A present in the sample or standard binds to antibodies adsorbed to the microwells.

First Incubation




Figure 3

Following incubation unbound biological components are removed during a wash step and a biotin-conjugated antimouse VEGF-A antibody is added and binds to mouse VEGF-A captured by the first antibody.







Figure 4






Following incubation unbound biotin-conjugated anti-mouse VEGF-A antibody is removed during a wash step. Streptavidin-HRP is added and binds to the biotin-conjugated anti-mouse VEGF-A antibody.

Following incubation unbound Streptavidin-HRP is removed during a wash step, and substrate solution reactive with HRP is added to the wells.


Figure 5

Fourth Incubation



A coloured product is formed in proportion to the amount of mouse VEGF-A present in the sample or standard. The reaction is terminated by addition of acid and absorbance is measured at 450 nm. A standard curve is prepared from 7 mouse VEGF-A standard dilutions and mouse VEGF-A sample concentration determined.

Figure 6



3. Reagents provided and reconstitution

(Note: Quantity shown is for the 1x96 format, for the 2x96 format all reagents will be supplied in duplicate where applicable)

- 1 aluminium pouch with a Microwell Plate coated with polyclonal antibody to mouse VEGF-A
- 1 vial (120 µl) **Biotin-Conjugate** anti-mouse VEGF-A polyclonal antibody
- 1 vial (150 µl) Streptavidin-HRP
- 2 vials mouse VEGF-A **Standard** lyophilized, 4 ng/ml upon reconstitution
- 1 vial (10 ml) Sample Diluent
- 1 vial (12 ml) Assay Diluent
- 1 vial (5 ml) Conjugate Diluent Concentrate 20x
- 1 bottle (50 ml) Wash Buffer Concentrate 20x (PBS with 1% Tween 20)
- 1 vial (15 ml) **Substrate Solution** (tetramethyl-benzidine)
- 1 vial (15 ml) **Stop Solution** (1M Phosphoric acid)
- 1 vial (0.4 ml) Blue-Dye
- 1 vial (0.4 ml) Green-Dye
- 1 vial (0.4 ml) **Red-Dye**
- 6 Adhesive Films

4. Materials required but not provided

- Microtitre plate reader fitted with appropriate filters (450nm required with optional 620nm reference filter)
- Microplate washer or wash bottle
- 10, 50, 100, 200 and 1,000µl adjustable single channel micropipettes with disposable tips
- 50-300ul multi-channel micropipette with disposable tips
- Multichannel micropipette reagent reservoirs
- Distilled water
- Vortex mixer
- Miscellaneous laboratory plastic and/or glass, if possible sterile

5. Storage Instructions

Store kit reagents between 2° and 8°C. Immediately after use remaining reagents should be returned to cold storage (2° to 8°C). Expiry of the kit and reagents is stated on labels.

Expiry of the kit components can only be guaranteed if the components are stored properly, and if, in case of repeated use of one component, this reagent is not contaminated by the first handling.

6. Specimen collection, processing & storage

Cell culture supernatant and serum* were tested with this assay. Other body fluids might be suitable for use in the assay. Remove serum from the clot as soon as possible after clotting.

Samples containing a visible precipitate must be clarified prior to use in the assay. Do not use grossly hemolyzed or lipemic specimens.

Samples should be aliquoted and must be stored frozen at -20°C to avoid loss of bioactive mouse VEGF-A. If samples are to be run within 24 hours, they may be stored at 2° to 8°C. Avoid repeated freeze-thaw cycles. Prior to assay, the frozen sample should be brought to room temperature slowly and mixed gently.

* Pay attention to a possibly elevated serum level of mouse VEGF-A due to VEGF-A release by platelets during platelet activation (sampling process).

Do not thaw samples in a 37°C water bath. Do not vortex or sharply agitate samples.

7. Safety & precautions for use

- Handling of reagents, serum or plasma specimens should be in accordance with local safety procedures , e.g.CDC/NIH Health manual : " Biosafety in Microbiological and Biomedical Laboratories" 1984
- · Laboratory gloves should be worn at all times
- Avoid any skin contact with H₂SO₄ and TMB. In case of contact, wash thoroughly with water
- Do not eat, drink, smoke or apply cosmetics where kit reagents are used
- Do not pipette by mouth
- When not in use, kit components should be stored refrigerated or frozen as indicated on vials or bottles labels
- All reagents should be warmed to room temperature before use. Lyophilized standards should be discarded after use
- Once the desired number of strips has been removed, immediately reseal the bag to protect the remaining strips from deterioration
- Cover or cap all reagents when not in use
- Do not mix or interchange reagents between different lots
- Do not use reagents beyond the expiration date of the kit
- Use a clean disposable plastic pipette tip for each reagent, standard, or specimen addition in order to avoid cross contamination, for the dispensing of H₂SO₄ and substrate solution, avoid pipettes with metal parts
- Use a clean plastic container to prepare the washing solution
- Thoroughly mix the reagents and samples before use by agitation or swirling
- All residual washing liquid must be drained from the wells by efficient aspiration or by decantation followed by tapping the plate forcefully on absorbent paper. Never insert absorbent paper directly into the wells
- The TMB solution is light sensitive. Avoid prolonged exposure to light. Also, avoid contact of the TMB solution with metal to prevent colour development. Warning TMB is toxic avoid direct contact with hands. Dispose off properly
- If a dark blue colour develops within a few minutes after preparation, this indicates that the TMB solution has been contaminated and must be discarded. Read absorbance's within 1 hour after completion of the assay
- When pipetting reagents, maintain a consistent order of addition from well-to-well. This will ensure equal incubation times for all wells
- Follow incubation times described in the assay procedure
- Dispense the TMB solution within 15 min of the washing of the microtitre plate

8. Assay Preparation

Bring all reagents to room temperature before use

8.1. Assay Design

Determine the number of microwell strips required to test the desired number of samples plus appropriate number of wells needed for running zeros and standards. Each sample, standard and zero should be tested in duplicate. Remove sufficient Microwell Strips for testing from the aluminium pouch immediately prior to use. Return any wells not required for this assay with desiccant to the pouch. Seal tightly and return to 2-8°C storage.

Example plate layout (example shown for a 7 point standard curve)

	Stand	dards		Sample Wells								
	1	2	3	4	5	6	7	8	9	10	11	12
Α	4000	4000										
В	2000	2000										
С	1000	1000										
D	500	500										
E	250	250										
F	125	125										
G	62.5	62.5										
Н	Zero	Zero										

All remaining empty wells can be used to test samples in duplicate

8.2. Preparation of Wash Buffer

Pour entire contents (50 ml) of the Wash Buffer Concentrate (20x) into a clean 1000 ml graduated cylinder. Bring to final volume of 1000 ml with glass-distilled or deionized water. Mix gently to avoid foaming.

Transfer to a clean wash bottle and store at 2° to 25°C. Please note that Wash Buffer (1x) is stable for 30 days.

Wash Buffer (1x) may also be prepared as needed according to the following table:

Number of Strips	Wash Buffer Concentrate (20x)	Distilled Water		
	(ml)	(ml)		
1 - 6	25	475		
1 - 12	50	950		

8.3. Preparation of Conjugate Diluent Buffer

Pour the entire contents (5 ml) of the **Conjugate Diluent Concentrate** (20x) into a clean 100 ml graduated cylinder. Bring to final volume of 100 ml with distilled water. Mix gently to avoid foaming.

Store at 2° to 8°C. Please note that the Conjugate Diluent (1x) is stable for 30 days.

Conjugate Diluent (1x) may also be prepared as needed according to the following table:

Number of Strips	Conjugate Diluent Concentrate (20x)	Distilled Water	
	(ml)	(ml)	
1 - 6	2.5	47.5	
1 - 12	5.0	95.0	

8.4. Preparation of Standard

Reconstitute **mouse VEGF-A standard** by addition of Sample Diluent (for subsequent measurement of **serum** or **plasma samples**) or distilled water (for subsequent measurement of **cell culture supernatant samples**).

Reconstitution volume is stated on the label of the standard vial. Allow the reconstituted standard to sit for 10-30 minutes. Swirl or mix gently to insure complete and homogeneous solubilization (concentration of reconstituted standard = 4000 pg/ml).

Label 6 tubes, one for each standard point. S2, S3, S4, S5, S6, S7. Then prepare 1:2 serial dilutions for the standard curve as follows: Pipette 150 μ l of Sample Diluent (for subsequent measurement of serum or plasma samples) or Assay Diluent (for subsequent measurement of cell culture supernatant samples) into each tube. Pipette 150 μ l of reconstituted standard (serves as highest standard S1 = 4 ng/ml) into the first tube, labelled S2 and mix (concentration of standard 2 = 2 ng/ml). Pipette 150 μ l of this dilution into the second tube and mix thoroughly before the next transfer. Repeat serial dilutions 4 more times thus creating the points of the standard curve ranging from 4000 to 62.5 pg/ml.

Sample Diluent (serum or plasma samples) or Assay Diluent (cell culture supernatant samples) serve as blank.

Alternativley this can be conducted directly on the microtitre plate.

8.5. Preparation of Biotin-Conjugate

Please note that the Biotin-Conjugate should be used within 30 minutes after dilution.

Make a 1:100 dilution of the concentrated **Biotin-Conjugate** solution with Conjugate Diluent (1x) in a clean plastic tube as needed according to the following table:

Number of Strips	Biotin-Conjugate	Conjugate Diluent (1x)
	(ml)	(ml)
1 - 6	0.06	5.94
1 - 12	0.12	11.88

8.6. Preparation of Streptavidin-HRP

Please note that the Streptavidin-HRP should be used within 30 minutes after dilution.

Make a 1:100 dilution of the concentrated **Streptavidin-HRP** solution with Conjugate Diluent (1x) in a clean plastic tube as needed according to the following table:

Number of Strips	Streptavidin-HRP	Conjugate Diluent (1x)
•	(ml)	(ml)
1 - 6	0.06	5.94
1 - 12	0.12	11.88

8.7. Addition of Colour-giving Reagents: Blue-Dye, Green-Dye, Red-Dye

In order to help our customers to avoid any mistakes in pipetting

Platinum ELISAs, eBioscience offers a tool that helps to monitor the addition of even very small volumes of a solution to the reaction well by giving distinctive colours to each step of the ELISA procedure.

This procedure is optional, does not in any way interfere with the test results, and is designed to help the customer with the performance of the test, but can also be omitted, just following the instruction booklet.

Alternatively, the dye solutions from the stocks provided (*Blue-Dye, Green-Dye, Red-Dye*) can be added to the reagents according to the following guidelines:

1. Diluent: Before standard dilution add the **Blue-Dye** at a dilution of 1:250 (see table below) to the appropriate diluent (1x) according to the test protocol. After addition of **Blue-Dye**, proceed according to the instruction booklet.

5 ml	20 µl Blue-Dye
Sample Diluent (for serum or plasma samples)	
Assay Diluent (for cell culture supernatant samples)	
12 ml	48 μΙ <i>Blue-Dye</i>
Sample Diluent (for serum or plasma samples)	
Assay Diluent (for cell culture supernatant samples)	
50 ml	200 μΙ <i>Blue-Dye</i>
Sample Diluent (for serum or plasma samples)	
Assay Diluent (for cell culture supernatant samples)	

2. Biotin-Conjugate: Before dilution of the concentrated Biotin-Conjugate, add the *Green-Dye* at a dilution of 1:100 (see table below) to the Conjugate Diluent used for the final conjugate dilution. Proceed after addition of *Green-Dye* according to the instruction booklet: Preparation of Biotin-Conjugate.

3 ml Conjugate Diluent	30 µl Green-Dye
6 ml Conjugate Diluent	60 µl Green-Dye
12 ml Conjugate Diluent	120 µl Green-Dye

3. Streptavidin-HRP: Before dilution of the concentrated Streptavidin-HRP, add the *Red-Dye* at a dilution of 1:250 (see table below) to the Conjugate Diluent used for the final Streptavidin-HRP dilution. Proceed after addition of *Red-Dye* according to the instruction booklet: Preparation of Streptavidin-HRP.

6 ml Conjugate Diluent	24 µl Red-Dye
12 ml Conjugate Diluent	48 μΙ Red-Dye

9. Method

We strongly recommend that every vial is mixed thoroughly without foaming prior to use.

Prepare all reagents as shown in section 8.

Note: Final preparation of Biotin conjugate (section 8.5) and Streptavidin-HRP (section 8.6) should occur immediately before use.

ssay Step	Details
Wash	a) Dispense 0.4 ml of 1x washing solution into each well b) Aspirate the contents of each well c) Repeat steps a and b Do not allow wells to dry before use
Addition	Add 50μl of Assay diluent to all wells
Preparation	Prepare Standard curve as shown in section 8.4
Addition	Add 50µl of each Standard, Sample and appropriate blank (sample or assay diluent) in duplicate to appropriate number of wells
Incubation	Cover with a plastic plate cover and incubate at room temperature (18 to 25°C) for 2 hour(s) on a microplate shaker set at 100 rpm
Wash	Remove the cover and wash the plate as follows: a) Aspirate the liquid from each well b) Dispense 0.4 ml of 1x washing solution into each well c) Aspirate the contents of each well d) Repeat step b and c another five times
Addition	Add 100μl of diluted biotin-conjugate to all wells
Incubation	Cover with a plastic plate cover and incubate at room temperature (18 to 25°C) for 1 hour(s) on a microplate shaker set at 100 rpm
Wash	Repeat wash step 6.
Addition	Add 100µl of Streptavidin-HRP solution into all wells
Incubation	Cover with a plastic plate cover and incubate at room temperature (18 to 25°C) for 1 hour on a microplate shaker set at 100 rpm
Wash	Repeat wash step 6.
Addition	Add 100µl of ready-to-use TMB Substrate Solution into all wells
Incubation	Incubate in the dark for 30 minutes * at room temperature. Avoid direct exposure to light by wrapping the plate in aluminium foil.
Addition	Add 100µl of Stop Reagent into all wells
	Addition Preparation Addition Incubation Wash Addition Uncubation Wash Addition Incubation Wash Addition Incubation Uncubation Uncubation Uncubation Uncubation

nm as the primary wavelength and optionally 620 nm as the reference wave length (610 nm to 650 nm is acceptable).

Note: In case of incubation without shaking the obtained O.D. values may be lower than indicated below. Nevertheless the results are still valid.

^{*}Incubation time of the substrate solution is usually determined by the ELISA reader performance. Many ELISA readers only record absorbance up to 2.0 O.D. Therefore the colour development within individual microwells must be observed by the analyst, and the substrate reaction stopped before positive wells are no longer within recordable range

10. Data Analysis

Calculate the average absorbance values for each set of duplicate standards and samples. Duplicates should be within 20 per cent of the mean value.

Create a standard curve by plotting the mean absorbance for each standard concentration on the ordinate against the mouse VEGF-A concentration on the abscissa. Draw a best fit curve through the points of the graph (a 5-parameter curve fit is recommended).

To determine the concentration of circulating mouse VEGF-A for each sample, first find the mean absorbance value on the ordinate and extend a horizontal line to the standard curve. At the point of intersection, extend a vertical line to the abscissa and read the corresponding mouse VEGF-A concentration.

Calculation of samples with a concentration exceeding standard 1 may result in incorrect, low mouse VEGF-A levels. Such samples require further external predilution according to expected mouse VEGF-A values with Sample Diluent (serum or plasma samples) or Assay Diluent (cell culture supernatant samples) in order to precisely quantitate the actual mouse VEGF-A level.

It is suggested that each testing facility establishes a control sample of known mouse VEGF-A concentration and runs this additional control with each assay. If the values obtained are not within the expected range of the control, the assay results may be invalid.

Figure 7
Representative standard curve for mouse VEGF-A ELISA. Mouse VEGF-A was diluted in serial 2-fold steps in Sample Diluent (for serum or plasma samples) or Assay Diluent (for cell culture supernatant samples). Do not use this standard curve to derive test results. A standard curve must be run for each group of microwell strips assayed.




Table 1

Typical data using the mouse VEGF-A ELISA

Measuring wavelength: 450 nm Reference wavelength: 620 nm

	Mouse VEGF-A		Mean	
	Concentration	O.D. at	O.D. at	C.V.
Standard	(pg/ml)	450 nm	450 nm	(%)
1	4000.0	2.401	2.394	0.3
		2.388		
2	2000.0	1.266	1.290	1.9
		1.314		
3	1000.0	0.706	0.714	1.1
		0.722		
4	500.0	0.364	0.371	1.9
		0.379		
5	250.0	0.236	0.235	0.6
		0.233		
6	125.0	0.143	0.142	0.9
		0.141		
7	62.5	0.106	0.105	0.5
		0.104		
Blank	0	0.066	0.065	0.7
		0.065		

The OD values of the standard curve may vary according to the conditions of assay performance (e.g. operator, pipetting technique, washing technique or temperature effects). Furthermore shelf life of the kit may affect enzymatic activity and thus colour intensity. Values measured are still valid.

11. Assay limitations

Do not extrapolate the standard curve beyond the maximum standard curve point. The dose-response is non-linear in this region and good accuracy is difficult to obtain. Concentrated samples above the maximum standard concentration must be diluted with Standard diluent or with your own sample buffer to produce an OD value within the range of the standard curve. Following analysis of such samples always multiply results by the appropriate dilution factor to produce actual final concentration.

The influence of various drugs on end results has not been investigated. Bacterial or fungal contamination and laboratory cross-contamination may also cause irregular results.

Improper or insufficient washing at any stage of the procedure will result in either false positive or false negative results. Completely empty wells before dispensing fresh Washing Buffer, fill with Washing Buffer as indicated for each wash cycle and do not allow wells to sit uncovered or dry for extended periods.

Disposable pipette tips, flasks or glassware are preferred, reusable glassware must be washed and thoroughly rinsed of all detergents before use.

As with most biological assays conditions may vary from assay to assay therefore **a fresh standard curve must be prepared and run for every assay.**

12. Performance Characteristics

12.1. Sensitivity

The limit of detection of mouse VEGF-A defined as the analyte concentration resulting in an absorbance significantly higher than that of the dilution medium (mean plus 2 standard deviations) was determined to be 20.0 pg/ml (mean of 6 independent assays).

12.2. Specificity

The assay detects both natural and recombinant mouse VEGF-A. The cross reactivity and interference of circulating factors of the immune system was evaluated by spiking these proteins at physiologically relevant concentrations into a mouse VEGF-A positive serum. There was no cross reactivity or interference detected, notably not with mouse VEGF-B, VEGF-C, VEGF-D and PfGF.

Interference was detected for VEGF-R1 at concentrations > 200 pg/ml, and not for VEGF-R2.

12.3. Precision

Intra Assay

Reproducibility within the assay was evaluated in 3 independent experiments. Each assay was carried out with 6 replicates of 6 serum samples containing different concentrations of mouse VEGF-A. 2 standard curves were run on each plate. Data below show the mean mouse VEGF-A concentration and the coefficient of variation for each sample (see Table). **The calculated overall intra-assay coefficient of variation was 3.8%**.

Table 2The mean mouse VEGF-A concentration and the coefficient of variation for each sample

		Mean Mouse VEGF-A	Coefficient of Variation
Sample	Experiment	Concentration (pg/ml)	(%)
1	1	2932	3.6
	2	3158	1.7
	3	3148	1.3
2	1	1747	1.4
	2	1920	2.1
	3	1971	4.2
3	1	1214	3.7
	2	1306	4.2
	3	1349	3.9
4	1	614	5.1
	2	690	6.1
	3	710	2.8
5	1	494	3.3
	2	543	2.3
	3	547	1.0
6	1	121	5.7
	2	143	8.9
	3	123	6.6

Inter-assay

Assay to assay reproducibility within one laboratory was evaluated in 3 independent experiments. Each assay was carried out with 6 replicates of 6 serum samples containing different concentrations of mouse VEGF-A. 2 standard curves were run on each plate. Data below show the mean mouse VEGF-A concentration and the coefficient of variation calculated on 18 determinations of each sample (see Table). The calculated overall inter-assay coefficient of variation was 6.4%.

Table 3The mean mouse VEGF-A concentration and the coefficient of variation of each sample

Sample	Mean Mouse VEGF-A Concentration (pg/ml)	Coefficient of Variation (%)	
1	3079	4.2	
2	1879	6.3	
3	1289	5.4	
4	672	7.5	
5	528	5.7	
6	129	9.4	

12.4. Dilution Parallelism

Serum, plasma and cell culture supernatant samples with different levels of mouse VEGF-A were analysed at serial 2 fold dilutions with 4 replicates each. For recovery data see Table .

Table 4

Sample matrix	Recovery of Exp. Val.		
	Range (%)	Mean (%)	
Serum	84 - 112	102	
Plasma (EDTA)	98 - 119	106	
Cell culture supernatant	78 - 105	92	

12.5. Spike Recovery

The spike recovery was evaluated by spiking 4 levels of mouse VEGF-A into serum, plasma and cell culture supernatant. Recoveries were determined in 3 independent experiments with 4 replicates each. The amount of endogenous mouse VEGF-A in unspiked serum was subtracted from the spike values. For recovery data see Table 5.

Table 1

Sample matrix *	Spike high (%)	Spike medium (%)	Spike low (%)
Serum	95	87	95
Plasma (EDTA)	88	104	nd
Cell culture supernatant	88	97	75

^{*} Due to high endogen mouse VEGF-A levels data for low spikes are not indicated

12.6. Stability

Freeze-Thaw Stability

Aliquots of serum samples (spiked or unspiked) were stored at -20°C and thawed 5 times, and the mouse VEGF-A levels determined.

There was a significant decrease of mouse VEGF-A immunoreactivity detected. Therefore samples should be stored in aliquots at -20°C and thawed only once.

Storage Stability

Aliquots of serum samples (spiked or unspiked) were stored at -20°C, 2-8°C, room temperature (RT) and at 37°C, and the mouse VEGF-A level determined after 24 h. There was no significant loss of mouse VEGF-A immunoreactivity detected during storage at -20°C and 2-8°C.

A significant loss of mouse VEGF-A immunoreactivity was detected during storage at RT and 37°C after 24 h.

13. Assay Summary

Total procedure length: 4h30mn

Wash two times

↓

Add sample and diluted standard
↓

Incubate 2 hours at room temperature
↓

Wash six times
↓

Add 100µl of Biotin-Conjugate
↓

Incubate 1hour at room temperature
↓

Wash six times

Add 100µl of Streptavidin-HRP
↓

Incubate 1hour at room temperature
↓

Wash six times

Add 100µl of ready-to-use TMB

Protect from light. Let the color develop for 30 mn.
↓

Add 100 µl Stop Reagent

Read Absorbance at 450 nm

Supplier:

Diaclone SAS 1 Boulevard A.Fleming 25020 Besancon Cedex France Tel +33 381413838 Fax +33 381413636

Email: info@diaclone.com