

# Clover Direct TM

tRNA Reagents for Site-Directed Protein Functionalization








**Questions about Related Products** 

Product Description	
About Clover <i>Direct</i> TM	2
Principle	4
Product Contents	6
Brief Protocol	7
Points to Note	8
Product List	9
Site-Directed Fluorescence Labeling	
Clover <i>Direct</i> ™ CR110-X-AF-tRNA	
Clover <i>Direct</i> ™ HiLyte Fluor™ 488-AF-tRNA	
Clover <i>Direct</i> ™ TAMRA-X-AF-tRNA	
Clover <i>Direct</i> ™ ATTO 633-AF-tRNA	
Clover <i>Direct</i> ™ ATTO 655-X-AF-tRNA	
Site Directed Pictin Labeling	
Site-Directed Biotin Labeling Clover Direct ™ Biotin-AF-tRNA	
Clover Direct TM Biotin-X-AF-tRNA	
Clover Direct TM Biotin-XX-AF-tRNA	
Clovel Direct ···· Blottil-AA-AF-tRIVA	
Site-Directed Post-translational Modification	
Clover <i>Direct</i> ™ Tyr(PO <sub>3</sub> H <sub>2</sub> )-tRNA	
Clover <i>Direct</i> ™ Lys(Me)-tRNA	
Clover <i>Direct</i> TM Lys(Me <sub>2</sub> )-tRNA	
Clover <i>Direct</i> ™ Lys(Ac)-tRNA	
Site-Directed Unnatural Mutagenesis	
PEGylated amino acids	
Clover Direct TM PEG4-AF-tRNA	
Clover Direct ™ PEG8-AF-tRNA	
Clover Direct TM PEG12-AF-tRNA	
CIOVEIDINOC I EGIZ / II II II II II	
Cross-linking amino acids	
Clover <i>Direct</i> ™ BPA-tRNA	
Clover <i>Direct</i> <sup>™</sup> AcPhe-tRNA	
Photo-isomerizable amino acid	
Clover <i>Direct</i> ™ azoAla-tRNA	
Custom Service	11
Custom Services for Unnatural Aminoacyl-tRNAs	
Custom Service for Proteins with Unnatural Amino Acids	
Applications	12
Expression of Labeled Proteins	12
Applications of Site-Directly Fluorescent Labeled Proteins	
Applications of Site-Directly Fluorescent Labeled Ploteins	
References	14




# [About Clover Direct TM]

Clover *Direct* TM tRNA Reagents for Site-Directed Protein Functionalization allow the incorporation of unnatural amino acids at defined positions of proteins using *in vitro* translation. Unnatural amino acids containing fluorescent groups, biotin, PEG, photo-crosslink, etc are available. Proteins with unnatural amino acids will be obtained within a few hours just by adding Clover *Direct* TM reagents and DNA template having an amber stop codon (UAG) or a four-base codon (CGGG) to an *in vitro* translation system. Clover *Direct* TM covers the following four applications. In addition, we provide custom services for the expression of proteins with unnatural amino acids (see page 11 for details).



#### **Site-Directed Fluorescence Labeling**

it is not easy to incorporate fluorescent groups into proteins in a site-direct and quantitative fashion by chemical modification. Clover *Direct* ™ tRNA Reagents for Site-Directed Fluorescence Labeling allow the incorporation of fluorescent unnatural amino acids into proteins in a site-direct and quantitative fashion. Various fluorescent dyes are available including those for 488 nm, 543 nm and 633 nm excitation.









#### Site-Directed Biotin Labeling

Clover *Direct* TM tRNA Reagents for Site-Directed Biotin Labeling allow the incorporation of biotinylated unnatural amino acids into proteins in a site-direct and quantitative fashion. Labeling proteins are available for the oriented immobilization onto avidin-coated plates and beads. The biotinylated amino acids have one or two aminohexyl liners between amino acid and biotin.




#### **Site-Directed Post-Translational Modification**

It is not easy to prepare post-translationally modified proteins (phosphorylation, methylation, etc.). Clover *Direct* TM tRNA Reagents for Site-Directed Post-Translational Modification allow the incorporation of modified amino acids into proteins to obtain post-translationally modified proteins in a site-direct and quantitative fashion.



#### Site-Directed Unnatural Mutagenesis

By incorporation of unnatural amino acids containing functional groups, novel functional proteins can be designed and synthesized. Clover *Direct* <sup>TM</sup> tRNA Reagents for Site-Directed Unnatural Mutagenesis allow the incorporation of unnatural amino acids with PEG, photo-crosslinking, photo-isomerizable groups, etc.






# [Principle of incorporation of unnatural amino acids]

Incorporation position of unnatural amino acids is defined by a UAG amber codon or CGGG four-base codon. An unnatural aminoacyl-tRNA recognizes the UAG amber codon or the CGGG codon during translation. Consequently, the unnatural amino acid is incorporated at the directed site of the protein. By using two tRNAs for amber and four-base codons, dual-labeled proteins can be obtained which are available for fluorescence resonance energy transfer (FRET).

#### [UAG amber codon]


If the UAG codon is recognized by the amber suppressor tRNA, full-length protein containing the unnatural amino acid is successfully synthesized. On the contrary, if the UAG codon is recognized by release factor 1 (RF1) which is one of the termination factors, the protein synthesis is terminated. Therefore, the translation product obtained as a full-length protein contains the unnatural amino acid at 100% efficiency.





#### [CGGG four-base codon]

If the CGGG codon is recognized by the four-base anticodon tRNA, full-length protein containing the unnatural amino acid is successfully synthesized. On the contrary, if the CGG is recognized as a triplet codon by Arg-tRNA, the reading frame shifts to +1 frame and a downstream stop codon terminates the protein synthesis. Therefore, the translation product obtained as a full-length protein contains the unnatural amino acid at 100% efficiency.



Incorporation at N-terminal regions (within 20 amino acid residues from the N-terminus) in response to CGGG codon sometimes results in the production of full-length proteins without unnatural amino acids, possibly because of spontaneous +1 frameshifting. In such case, ProteinExpress recommend the use of ProX<sup>TM</sup> tag, which is original peptide tag developed for the CGGG codon-mediated incorporation of unnatural amino acids.

5'-	AUG	UCU	AAA	CAA	AUC	GAA	GUA	AAC	CGGG	UCU	AAU	GAG	-3'
	Met	Ser	Lys	Gln	lle	Glu	Val	Asn	Xaa	Ser	Asn	Glu	

Sequence of ProXTM tag





# [Product Contents]

Unnatural aminoacyl-tRNA (See note1)
 tRNA dissolving buffer
 X 1

Note 1 : One tube contains unnatural aminoacyl-tRNA sufficient for 300  $\mu$ L of *in vitro* translation reaction. Once thawed, unnatural aminoacyl-tRNA can be stored at -70 °C for 2 months.

#### [Equipment and reagents to be supplied by others]

#### **Protein Expression**

Cell-free translation system(*E.coli*)
 (e.g. RTS100 *E.coli* HY kit; Roche Applied Science, #3186156)

• Expression gene containing UAG or CGGG codon (plasmid DNA, linear DNA, or messenger RNA). (Suitable expression vector should be used for the cell-tree translation system.)

### Purification, buffer exchange, and concentration

Affinity column and buffer for purification
 (e.g. His SpinTrap<sup>™</sup> kit; GE Healthcare Science, #28-9321-71)

· Ultrafiltration membrane and buffer (e.g. ULTRAFREE®-0.5 Centrifugal Filter Devices 10k; Millipore, UFV5BC00)



# [Brief protocol]

#### Step 1 Protein Expression


Dissolve unnatural aminoacyl-tRNA with tRNA buffer, and mix with template DNA and cell-free translation system. The mixture is incubated for 1 hour to synthesize protein containing unnatural amino acid.

#### Step 2 Purification

Reaction mixture includes several proteins derived from cell-free translation system and unnatural amino acid which is not incorporated into protein. Full-length protein containing unnatural amino acid can be isolated by purification for C-terminal tag such as His tag.

#### Step 3 Buffer exchange and concentration (optional)

Buffer exchange and concentration can be done by ultrafiltration membrane. Purified protein containing unnatural amino acid can be directly used for the downstream experiment.



Technical Overview of Clover Direct ™



## [Points to note]

#### 1. Protein expression

Confirm your protein can be expressed in *E.coli* cell-free translation system.

In case of very low expression of a wild-type gene that does not contain UAG codon or CGGG codon, optimization of nucleotide sequence (codon usage, addition of N-terminal tags, etc.) is required to improve the expression. ProteinExpress provides custom service for protein expression using Clover*Direct* TM including gene construction for efficient expression.

#### 2. Incorporation site-dependency

Some unnatural amino acids are allowed to be incorporated only at the N-terminal region (within 20 amino acid residues from the N-terminus). Please check the product list for details.

	Sit	Site-dependency (*1)			Q'tv	Product
Pin-point Fluorescence Labeling	N-te	erminal	Internal, C-terminal	codon	(translation)	No
CR110-X-AF	NH <sub>2</sub> *			Amber	300 µL	CLD1001
[5-CR110-X : Abs/Em = 498/521nm]	<i>\bigg\</i>	0	×	Amber	5 X 300 μL	CLD2001
NH2	NH <sub>2</sub>	•	•	CGGG	5 X 300 μL	CLD2002
	△ : Ava	: Available : Available in some sites : Unavailable				

N-terminal region = within 20 amino acid residues from the N-terminus

Incorporation at N-terminal regions (within 20 amino acid residues from the N-terminus) in response to CGGG codon sometimes results in the production of full-length proteins without unnatural amino acids, possibly because of spontaneous +1 frameshifting. In such case, ProteinExpress recommend the use of ProX<sup>TM</sup> tag, which is original peptide tag developed for the CGGG codon-mediated incorporation of unnatural amino acids.

5'-	AUG	UCU	AAA	CAA	AUC	GAA	GUA	AAC	CGGG	UCU	AAU	GAG	-3'
	Met	Ser	Lys	Gln	lle	Glu	Val	Asn	Xaa	Ser	Asn	Glu	

## Sequence of ProX<sup>TM</sup> tag

#### 3. Custom Services

ProteinExpress provides custom services for the synthesis of unnatural aminoacyl-tRNAs (fluorescent-labeled, functional amino acids, etc.), which allows the expression of proteins with your requested unnatural amino acids.

● Custom Services for Unnatural Aminoacyl-tRNA......p.11

ProteinExpress also provides custom services for the expression of proteins with unnatural amino acids including construction of recombinant DNA, cell-free translation, and purification of proteins.

● Custom Service for Protein Expression using Clover *Direct* TM......p.11





# **Product List**

	Site-dependency (*1)			Q'ty	Product
Site-Directed Fluorescence Labeling	N-terminal	Internal, C-terminal	codon	(translation)	No
CR110-X-AF			Amban	300 µL	CLD1001
[5-CR110-X : Abs/Em = 498/521nm]	0	×	Amber	5 X 300 μL	CLD2001
NH <sub>2</sub> H NH <sub>2</sub> NH <sub>3</sub>			CGGG	5 X 300 μL	CLD2002
HiLyte Fluor™ 488-AF				300 µL	CLD01
[HiLyte Fluor <sup>™</sup> 488 : Abs/Em = 497/525nm]	0	×	Amber	5 X 300 μL	CLD05
Not Available			CGGG	5 X 300 μL	CLD2004
TAMRA-X-AF OH 5				300 µL	CLD02
NH <sub>2</sub>	0	×	Amber	5 X 300 μL	CLD06
[5(6)-TAMRA-X : Abs/Em = 546/575nm]			CGGG	Contact to O	LYMPUS
ATTO633-AF			Amber	300 µL	CLD03
[ATTO633 : Abs/Em = 629/657nm]	0	×		5 X 300 μL	CLD07
Not Available			CGGG	5 X 300 μL	CLD2008
ATTO655-X-AF				300 µL	CLD1009
[ATTO655-X : Abs/Em = 633/684nm]			Amber	5 X 300 μL	CLD2009
	0	×	cccc	300 µL	CLD1010
Not Available			CGGG	5 X 300 μL	CLD2010
Site-Directed Biotin Labeling					
Biotin-AF [Biotin]	0	0	Amber	5 X 300 μL	CLD2101
[DIOUIT]		0	CGGG	5 X 300 μL	CLD2102
Biotin-X-AF [Biotin-X]	0	0	Amber	5 X 300 μL	CLD2103
[DIOTIT-X]			CGGG	5 X 300 μL	CLD2104
Biotin-XX-AF			Anaham	300 µL	CLD04
[Biotin-XX]	0		Amber	5 X 300 μL	CLD08
NH <sub>2</sub> HN H			CGGG	5 X 300 μL	CLD2106

<sup>\*1 :</sup> Incorporation site-dependency

[Please inquire about the price]

Some unnatural amino acids are incorporated in a incorporation site-dependent manner.

(N-terminal region = within 20 amino acid residues from the N-terminus.)

O: Available

 $\triangle$  : Available in some sites

× : Unavailable



# Product List



		Incorporable site			O'to a	Drodust	
Site-Directed Post-Translation	ite-Directed Post-Translational Modification		N-terminal Internal, C-terminal		Q'ty (translation)	Product No	
Tyr(PO <sub>3</sub> H <sub>2</sub> ) [O-phospho-Tyr]	ОН	0	Δ	Amber	5 X 300 μL	CLD2201	
[O-pnospno-Tyr]	NH <sub>2</sub> OH			CGGG	5 X 300 μL	CLD2202	
Lys(Me)	ОН	0	0	Amber	5 X 300 μL	CLD2203	
$[\varepsilon$ -methyl-Lys]	O NH <sub>2</sub>	0	0	CGGG	5 X 300 μL	CLD2204	
Lys(Me <sub>2</sub> )	OH			Amber	5 X 300 μL	CLD2205	
$[\varepsilon$ -dimethyl-Lys]	O NH <sub>2</sub>	0	0	CGGG	5 X 300 μL	CLD2206	
Lys(Ac)	ОН		0	Amber	5 X 300 μL	CLD2207	
$[\varepsilon$ -acetyl-Lys]	NH <sub>2</sub>	0		CGGG	5 X 300 μL	CLD2208	
Site-Directed Unnatural Mutag	enesis						
PEGylated amino acids							
PEG4-AF [Methyl-PEG4]	OH PI	0	0	Amber	5 X 300 μL	CLD2301	
[modify 1 20 t]	NH <sub>2</sub> NH <sub>3</sub>			CGGG	5 X 300 μL	CLD2302	
PEG8-AF [Methyl-PEG8]	OH O(	0	0	Amber	5 X 300 μL	CLD2303	
[Metrly Pr 200]	NH <sub>2</sub> CH <sub>3</sub>			CGGG	5 X 300 μL	CLD2304	
PEG12-AF [Methyl-PEG12]	OH O	0	Δ	Amber	5 X 300 μL	CLD2305	
[Metryl-PEG12]	NH <sub>2</sub> CH <sub>3</sub>			CGGG	5 X 300 μL	CLD2306	
Cross-linking amino acids							
BPA	OH OH			Amber	5 X 300 μL	CLD2321	
[ $ ho$ -benzoyl-phenylalanine]	ŇH <sub>2</sub>	0	0	CGGG	5 X 300 μL	CLD2322	
AcPhe	OH 1 ~ ~			Amber	5 X 300 μL	CLD2323	
[ $ ho$ -acetyl-phenylalanine]	ONH <sub>2</sub>	0	0	CGGG	5 X 300 μL	CLD2324	
Photo-isomerizable amino ac	id	•				1	
azoAla	OH OH	_		Amber	5 X 300 μL	CLD2331	
[ $ ho$ -phenylazophenyl-alanine]	O' NH <sub>2</sub>	0	0	CCCC	5 V 300 ul	CI D2332	

[Please inquire about the price]

5 X 300 μL

CLD2332


CGGG





# [Custom Services for Unnatural Aminoacyl-tRNA]


ProteinExpress provides custom services for the synthesis of unnatural aminoacyl-tRNAs (fluorescent-labeled, functional amino acids, etc.), which allows the expression of proteins with your requested unnatural amino acids.



**Custom Synthesis of Unnatural Aminoacyl-tRNA** 

## [Custom Service for Proteins with Unnatural Amino Acids]

ProteinExpress provides custom services for the expression of proteins with your requested unnatural amino acids at requested positions, including artificial gene synthesis, cell-free translation, and protein purification.



**Custom Service for Protein with Unnatural Amino Acids** 


# **Applications**



# [Expression of site-directly labeled proteins]

Fluorescent- and biotin-labeled unnatural amino acids are incorporated into eight prokaryote and eukaryote proteins. The site-directly fluorescent-labeled proteins can be visualized on SDS-PAGE using a laser-based fluorescence scanner. The proteins are also detectable by an antibody against tag peptide or biotin. A 0.25 ~ 1 µL of translational reaction mixture is sufficient for the detection.

· Labeled unnatural amino acids that are not incorporated into proteins are detected at the bottom of SDS-PAGE gel.



ProX tag : ProX tag is fused to the N-terminus. Applied volume: 0.25 µL of translational reaction mix Fluorescence image (Top) are visualized with Ex and Em

wavelengths listed below: HiLyteFluor488 Ex: 488nm / Em: 520 nm

**TAMRA** Ex: 532nm / Em: 580 nm ATTO633 Ex: 635nm / Em: 670 nm

Western blotting (Bottom) are visualized by anti-His tag antibody (for fluorescent amino acids) and anti-biotin antibody (for biotin).






# [Applications for site-directly fluorescent labeled proteins]

The site-directly fluorescent-labeled proteins are available to the following measurements.

- Interaction analysis using single molecule fluorescence analysis
   Single Molecule fluorescence detection system (MF20 / Fluor Point-Light MF20; OLYMPUS)
- -Conformation analysis of protein by inter- or intra-molecular fluorescence resonance energy transfer (FRET)
- ·Functional analysis in cell imaging
- Interaction analysis in protein array
- ·Expression analysis by in-gel fluorescent detection of SDS-PAGE



# References / Questions about Products



## [References]

- 1) FRET analysis of protein conformational change through position-specific incorporation of fluorescent amino acids Daisuke Kajihara, Ryoji Abe, Issei lijima, Chie Komiyama, Masahiko Sisido, Takahiro Hohsaka *Nature Methods.*, 3, 923-929 (2006).
- 2) Position-specific incorporation of biotinylated non-natural amino acids into a protein in a cell-free translation system Takayoshi Watanabe, Norihito Muranaka, Issei lijima, Takahiro Hohsaka *Biochem. Biophys. Res. Commun.*, 361, 794-799 (2007)
- 3) Comprehensive screening of amber suppressor tRNAs suitable for incorporation of non-natural amino acids in a cell-free translation system

Hikaru Taira, Yosuke Matsushita, Kenji Kojima, Kaori Shiraga, Takahiro Hohsaka *Biochem. Biophys. Res. Commun.*, 374, 304-308 (2008).

4) Efficient Incorporation of Nonnatural Amino Acids with Large Aromatic Groups into Streptavidin in In Vitro Protein Synthesizing Systems

Takahiro Hohsaka, Daisuke Kajihara, Yuki Ashizuka, Hiroshi Murakami, Masahiko Sisido *J. Am. Chem. Soc.*, 121, 34-40 (1999).

# [Questions about Products]

ProteinExpress Co.,Ltd URL: http://www.proteinexpress.co.jp Chiba University Inohana Innovation Plaza 1-8-15, Inohana, Chuo-ku,

Chiba-shi, Chiba 260-0856, Japan E-mail: tech@proteinexpress.co.jp



tRNA Reagents for Site-Directed Protein Functionalization

ProteinExpress Co.,Ltd
URL: http://www.proteinexpress.co.jp
Chiba University Inohana Innovation Plaza
1-8-15, Inohana, Chuo-ku,
Chiba-shi, Chiba 260-0856, Japan
E-mail: tech@proteinexpress.co.jp

