



mutaFISH™ IFNGwt RNA Probes

Catalog Number FP0031

10 Assays

Version: 1.1

Intended for research use only

Table of Contents

Introduction	3
Intended Use.....	3
General Information	3
Materials Supplied.....	3
Storage Instruction.....	3
Note	4
Precautions	4
Protocol.....	5
Cell <i>In Situ</i> Hybridization	5
CTC <i>In Situ</i> Hybridization.....	8

Introduction

Intended Use

mutaFISH™ IFNGwt RNA Probes is designed to detect IFNG gene on single strand RNA in cells using *in situ* rolling-circle amplification technology.

General Information

Materials Supplied

Component	Qty.	Storage temp.
RT IFNG Primer	10 µL	-20°C
mutaFISH™ IFNGwt RNA Probe	10 µL	-20°C
Detection Probe-FITC	10 µL	-20°C

Reagent Required but Not Provided

- ✓ 100% ethanol.
- ✓ 85% ethanol: a solution of 85% ethanol
- ✓ 70% ethanol: a solution of 70% ethanol
- ✓ 3% paraformaldehyde

*mutaFISH™ RNA Accessory Kit (Catalog#: KA4915) contains reagents for cell *in situ* hybridization and CTC *in situ* hybridization. When frozen tissue and FFPE tissue is used with mutaFISH™ probe, additional DEPC-H₂O, DEPC-PBS and PBS-T for immersing the slide should be prepared by the user.

Storage Instruction

Store at -20°C. Aliquot to avoid repeated freezing and thawing.

Note

We recommend mutaFISH™ RNA Accessory Kit (Catalog #: KA4915) which provides necessary reagents and enzymes for *in situ* reverse transcription, RNA digestion, mutaFISH™ hybridization, ligation and amplification prior to mutaFISH™.

Precautions

- General Precautions
 - ✓ For research use only.
 - ✓ Wear gloves and laboratory coats when handling immunodiagnostic materials.
 - ✓ Do not drink, eat or smoke in the areas where immunodiagnostic materials are being handled.
 - ✓ Wear gloves and eye and clothing protection when handling these reagents.
 - ✓ The materials must not be pipetted by mouth.

- General Technical Hints
 - ✓ Reagents with different lot numbers should not be mixed.
 - ✓ Use thoroughly cleaned glassware.
 - ✓ Use deionized (distilled) water, stored in clean containers.
 - ✓ Avoid any contamination among samples and reagents. For this purpose, disposable tips should be used for each sample and reagent.
 - ✓ Dispose consumable materials and unused contents in accordance with applicable national regulatory requirements.

- Technical hints for Cell *In Situ* Hybridization
 - ✓ Attached cells could be prepared by culturing on poly-L-lysine treated cover slide. For convenience, the cover slide with growth cells could be placed on a microscope slide before starting mutaFISH™ procedures.
 - ✓ Suspension or attached cells could be adhered on poly-L-lysine treated microscope slide by Cytospin.
 - ✓ Cell *In Situ* Hybridization protocol is suitable for 240 mm² cell area for each assay, user should make sure reagents applied could cover all of cell area.
 - ✓ Use RNase free and DNase free tips for handling all reagents
 - ✓ Perform immunostaining prior to mutaFISH™ procedures. After immunostaining is completed, follow Cell *In Situ* Hybridization protocol below.

Protocol

Cell *In Situ* Hybridization

1. Fix cells with 100 μ L of 2% paraformaldehyde for 15 min at RT.
2. Wash cells twice with 100 μ L of DEPC-PBS.
3. **Fixation:** Pipette 100 μ L of 3% paraformaldehyde into the slide and incubate the slide at RT for 30 minutes in humidity box.
4. **Wash step:** Pipette 100 μ L of DEPC-PBS to the slide for 1 minute to remove the residual 3% paraformaldehyde.
5. Repeat the wash step.
6. **Dehydration:** Pipette 100 μ L of 70%, 85% and 100% ethanol to the slide for 1 minute each.
7. **Permeabilization:** Wash the slide with 100 μ L of PBS-T followed by 100 μ L of permeabilization buffer at RT for 10 minutes.
8. **Wash step:** Wash the slide with 100 μ L of PBS-T twice.
9. ***In Situ* RT:** The reaction contents are indicated in the **Table 1** below. Pipette 100 μ L of total mixture to the slide and incubate the slide at 37°C for 3 hours in humidity box. .

Table 1.

Components of the reaction mixture	Amount (μL)
DEPC-H ₂ O	60.5
5X RT Buffer	20
BSA	1
dNTP Mix	5
RT IFNG Primer	1
RNase Inhibitor	2.5
RT Enzyme	10
Total reaction volume (μL)	100

10. **Wash step:** Wash the slide twice with 100 μ L of PBS-T.
11. **Postfixation:** Pipette 100 μ L of 3% paraformaldehyde to the slide and incubate the slide at RT for 30 minutes in humidity box.S
12. **Wash step:** Wash the slide twice with 100 μ L of PBS-T.
13. **mutaFISH™ probe hybridization:** The reaction contents are indicated in the **Table 2** below. Pipette 100 μ L of total mixture to the slide and incubate the slide at 37°C for 30 minutes, followed by at 45 °C for 45 min in humidity box.

Table 2.

Components of the reaction mixture	Amount (μL)
DEPC-H ₂ O	33.5
Formamide	20
10X Hybrid Enzyme Buffer	10
1M KCL	5
mutaFISH™ IFNGwt RNA Probe	1
RNase Inhibitor	2.5
RNase H	8
Hybrid Enzyme	20
Total reaction volume (μL)	100

14. **Wash step:** Wash the slide twice with 100 μL of PBS-T.
15. **Amplification:** The reaction contents are indicated in the **Table 3** below. Pipette 100 μL of total mixture to the slide and incubate the slide for 2 hours at 37°C in humidity box.

Table 3.

Components of the reaction mixture	Amount (μL)
DEPC-H ₂ O	61.5
50% Glycerol	10
10X DNA Polymerase Buffer	10
BSA	1
dNTP Mix	5
RNase Inhibitor	2.5
DNA Polymerase	10
Total reaction volume (μL)	100

16. **Wash step:** Wash the slide twice with 100 μL of PBS-T.
17. **Detection:** The reaction contents are indicated in the **Table 4** below. Pipette 100 μL of total mixture to the slide and incubate the slide for 15 minutes at 37°C in humidity box.

Table 4.

Components of the reaction mixture	Amount (μL)
Detection Buffer	99
Detection Probe-FITC	1
Total reaction volume (μL)	100

18. **Wash step:** Wash the slide twice with 100 μL of PBS-T.
19. **Dehydration:** Pipette 100 μL of 70%, 85% and 100% ethanol to the slide for 30 seconds each.
20. **Counter Stain:** Apply 100 μL of 300nM DAPI (Add 1 ul of 100 μM DAPI and DEPC-PBS to 333 μL) to the slide at RT for 3 minutes.
21. Wash the slide with 100 μL of DEPC-PBS twice. The slide is now ready for examine under fluorescence microscope.

CTC In Situ Hybridization



Bring Components and mixtures prepared below to room temperature before applying to the chip

1. Wash cells twice with 100 μ L of DEPC-PBS.
2. **Fixation:** Pipette 100 μ L of 3% paraformaldehyde into the chip and incubate the chip at RT for 30 minutes in humidity box.
3. **Wash step:** Pipette 100 μ L of DEPC-PBS into the chip for 1 minute to remove the residual 3% paraformaldehyde.
4. Repeat the wash step.
5. **Dehydration:** Pipette 100 μ L of 70%, 85% and 100% ethanol into the chip for 1 minute each.
6. **Permeabilization:** Wash the chip with 100 μ L of PBS-T followed by 100 μ L of permeabilization buffer at RT for 10 minutes.
7. **Wash step:** Wash the chip with 100 μ L of PBS-T twice.
8. **In Situ RT:** The reaction contents are indicated in the **Table 5** below. Pipette 100 μ L of total mixture into the chip and incubate the chip at 37°C for 3 hours in humidity box. Make sure the inlet and outlet of the chip are consistently covered with DEPC-PBS throughout incubation.

Table 5.

Components of the reaction mixture	Amount (μ L)
DEPC-H ₂ O	60.5
5X RT Buffer	20
BSA	1
dNTP Mix	5
RT IFNG Primer	1
RNase Inhibitor	2.5
RT Enzyme	10
Total reaction volume (μL)	100

9. **Wash step:** Wash the chip twice with 100 μ L of PBS-T.
10. **Postfixation:** Pipette 100 μ L of 3% paraformaldehyde into the chip and incubate the chip at RT for 30 minutes in humidity box.
11. **Wash step:** Wash the chip twice with 100 μ L of PBS-T.
12. **mutaFISH™ probe hybridization:** The reaction contents are indicated in the **Table 6** below. Pipette 100 μ L of total mixture into the chip and incubate the chip at 37°C for 30 minutes, followed by at 45 °C for 45 min in humidity box. Make sure the inlet and outlet of the chip are consistently covered with DEPC-PBS throughout incubation.

Table 6.

Components of the reaction mixture	Amount (μL)
DEPC-H ₂ O	33.5
Formamide	20
10X Hybrid Enzyme Buffer	10
1M KCL	5
mutaFISH™ IFNGwt RNA Probe	1
RNase Inhibitor	2.5
RNase H	8
Hybrid Enzyme	20
Total reaction volume (μL)	100

13. **Wash step:** Wash the chip twice with 100 μL of PBS-T.
14. **Amplification:** The reaction contents are indicated in the **Table 7** below. Pipette 100 μL of total mixture into the chip and incubate the chip for 2 hours at 37°C in humidity box. Make sure the inlet and outlet of the chip are consistently covered with DEPC-PBS throughout incubation.

Table 7.

Components of the reaction mixture	Amount (μL)
DEPC-H ₂ O	61.5
50% Glycerol	10
10X DNA Polymerase Buffer	10
BSA	1
dNTP Mix	5
RNase Inhibitor	2.5
DNA Polymerase	10
Total reaction volume (μL)	100

15. **Wash step:** Wash the chip twice with 100 μL of PBS-T.
16. **Detection:** The reaction contents are indicated in the **Table 8** below. Pipette 100 μL of total mixture into the chip and incubate the chip for 15 minutes at 37°C in humidity box. Make sure the inlet and outlet of the chip are consistently covered with DEPC-PBS throughout incubation.

Table 8.

Components of the reaction mixture	Amount (μL)
Detection Buffer	99
Detection Probe-FITC	1
Total reaction volume (μL)	100

17. **Wash step:** Wash the chip twice with 100 μL of PBS-T.
18. **Dehydration:** Pipette 100 μL of 70%, 85% and 100% ethanol into the chip for 30 seconds each.
19. **Counter Stain:** Apply 100 μL of 300nM DAPI (Add 1 μl of 100 μM DAPI and DEPC-PBS to 333 μL) at RT for 3 minutes.
20. Wash the chip with 100 μL of DEPC-PBS twice. The chip is now ready for examine under fluorescence microscope.



Abnova (Taiwan) Corporation

Customer Service: 02-87511888 8:30~17:30

(Except for Sat. ~ Sun. and holidays)

Email: sales@abnova.com