

Product Specification

AIE™ pH



Product Description

- This product can indicate pH.
- As compared with the widely-used pH dyes on the market, the product has an especially broad response range of pH.

• Demonstrations




Figure 1: Aggregation-induced emission properties of AIE™ pH




Figure 2: pH responses in aqueous phase




Figure 3: Cell staining.

Recommended storage condition

Store away from sunlight at 2-8 °C

Product parameters

Purpose	pH sensor
Color:	Red
Imaging platform:	Fluorescence microscope Confocal microscope Flow cytometry and et, al
Pack size and quantity:	10 µmol
Detection method:	Fluorescence
Excitation/Emission (nm):	Channel 1: 405/440 - 525 Channel 2: 488/561 - 659
Recommended transport condition:	Room temperature
Product declaration:	Only used for research. Do not apply to any detection procedure.

AI Egen Probe for pH Sensing

Introduction

- This product can indicate pH value from 1 to 13 with red to blue fluorescent emission shifting.
- This product can also be used for living cell pH mapping
- After incubation with this product in cell, the pH of the stained cell can be observed under UV illumination

Channel 1: Excitation/Emission = 405 / 440 - 525 nm

Channel 2: Excitation/Emission = 488 / 561 - 659 nm

Material Preparation and Visualization Recommendation

- **Stock solution prepare:** AIE™ pH (10 mM) stock solution is prepared with the 10 μ mol of AIE™ pH in 1 mL DMSO. Further dilute with 1000 times to get the working concentration at 10 μ M.

Before Your Experiment, You might NEED

- | | |
|-------------------|---------------------------|
| 1. DMSO | 3. Millipore water |
| 2. PBS 7.4 buffer | 4. Suitable culture media |

Protocol (Recommended)

Cell Culture

HeLa cells were cultured in minimum essential medium containing 10 % fetal bovine serum and antibiotics (100 units/mL penicillin and 100 µg/mL streptomycin) in a 5 % CO₂ humidity incubator at 37 °C.

Cell Staining and Imaging

1. HeLa cells are seeded in a 35-mm petri dish with a glass cover slide.
2. After overnight culture, cells were stained with the working concentration of 10 µM AIETM pH for 2 hours.
3. Wash the cells with phosphate buffered saline (PBS) (pH 7.4) solution for three times. Ready to do imaging using confocal microscope.
4. Confocal imaging condition is recommended as follows:

Channel 1: Excitation/Emission = 405 / 440 - 525 nm

Channel 2: Excitation/Emission = 488 / 561 - 659 nm

Reference

1. S. Chen, and B. Z. Tang "An AIE-active Hemicyanine Fluorogen with Stimuli-responsive Red/Blue Emission: Extending the pH Sensing Range by 'Switch + Knob' Effect" *Chem. Sci.*, **2012**, 3, 1804-1809
2. S. Chen and B. Z. Tang "Full-Range Intracellular pH Sensing by an Aggregation-Induced Emission-Active Two-Channel Ratiometric Fluorogen" *J. Am. Chem. Soc.*, **2013**, 135 (13), 4926-4929
3. Optical information and suggested storage conditions:

Item	Ex/Em	Qty	Storage Condition*
AIE TM pH	405/440-525 nm 488/561-659 nm	10 µmol	<ul style="list-style-type: none">• 2-8 °C (Upon receive this product)• Avoid Light• Keep Dry

* Remember to warm up to room temperature upon opening the vial