



SensoLyte® anti-Human MOG (1-125) Human IgG Specific Quantitative ELISA Kit

Colorimetric

Revision number: 1.2

Last updated: 05/23/2014

Catalog #	55153-H
Kit Size	One 96-well strip plate

This kit is optimized to detect anti-human MOG (1-125) IgG in human samples. Wells are pre-coated with recombinant human MOG (1-125) protein and pre-blocked with proprietary solution. The amount of anti-human MOG IgG in serum or plasma is quantified using ELISA. Human anti-Human MOG (1-125) standard is included. Ample materials and reagents are provided to perform 96 assays.

- **Convenient Format**
 - Pre-coated and pre-blocked 96-well strip plate
 - Ready-to-use substrate solution and other assay components
 - 2-3 hours assay time at room temperature
- **Minimal Sample Size**
 - Requires only 5-10 µl of serum or plasma to perform the assay
- **High Sensitivity**
 - Detects as low as 1 ng/ml anti-human MOG (1-125) IgG
- **Broad Dynamic Range**
 - 1.56-100 ng antibody/ml serum

Kit Components and Handling

Component	Description	Quantity
Component A	Human-MOG (1-125) coated and blocked 8-well strips	12 strips
Component B	Human anti-Human MOG (1-125) IgG standard	500 µl (100ng/ml) x 3
Component C	1X Sample Dilution Buffer	30 ml
Component D	10X Wash Buffer	50 ml
Component E	TMB color substrate solution	10 ml
Component F	Stop Solution	10 ml
Component G	Secondary antibody, Goat anti-Human IgG-HRP	30 µl (0.1mg/ml)

Other Materials Required (but not provided)

- Microplate reader: Capable of reading absorbance at 450 nm
- Rocking platform or shaker
- Strip ejector (to eject strips for later assay if not all strips are used in one experiment)
- Computer software: Capable of plotting Four Parameter Logistic Curve Fit (4-PL)
- Plate washer (optional)

Shipment and Storage

- Kit is shipped on blue ice. Upon receipt, store all kit components at 2-8°C for up to 12 months.
- Store **Component B** at -20 °C upon arrival

For Research Use Only. Not for Diagnostics.

©AnaSpec, Inc. 34801 Campus Drive, Fremont, CA 94555
Tel: (800)-452-5530 | service@anaspec.com | www.anaspec.com

Introduction

Myelin oligodendrocyte glycoprotein (MOG) is a member of the immunoglobulin superfamily and is expressed exclusively in the central nervous system.¹⁻³ Human MOG (1-125) is able to induce autoantibody production and relapsing-remitting neurological disease causing extensive plaque-like demyelination.¹⁻³ Autoantibody response to human MOG (1-125) has been observed in induced experimental autoimmune encephalomyelitis (EAE) in DA and Lewis rats, C57/BL6 and SJL mice, and common marmoset.¹⁻³ However, the exact pathological role and action of anti-human MOG (1-125) autoantibody is not known and is currently under vigorous investigation.¹⁻³

The SensoLyte® Anti-Human MOG (1-125) IgG Quantitative ELISA Kit provides a convenient quantitative assay for anti-human MOG (1-125) autoantibody. This kit is useful to researchers for determining the amount of anti-human MOG (1-125) antibody present in human serum or plasma, and can help provide information on the role it plays in the development and treatment of multiple sclerosis.

Safety Precautions!

Human source materials used as a standard in this kit have been tested and found negative for antibodies to HIV1, HIV2, and Hepatitis C, Hepatitis B antigen, and Syphilis. Nevertheless, all human source materials should be handled as potentially infectious. Therefore, wear gloves, protective clothing, and safety goggles while working with human source products.

Protocol

Please Note:

- a) Allow kit components to warm up to room temperature before starting the assay
- b) Spin down all components with volume less than 100 µl before use
- c) Mix thoroughly wash buffer before diluting to dissolve any precipitated salt
- d) More Sample Dilution Buffer can be made by adding 1% BSA into 1 X Wash Buffer

1. ELISA assay:

- 1.1 Establish dilution range of serum samples: we recommend 1:40 dilution for human serum or plasma. Use 1X Sample Dilution Buffer (Component C) to dilute samples. Depending on the amount of anti-Human MOG antibody present, the dilution range can be further adjusted.
- 1.2 Arrange and label strips (Component A) based on the number of wells with standard and samples. An example is shown in [Table 1](#). Although diluted standard and samples can be run as single points, duplicates are recommended.

Table 1. An example of four samples layout in duplicates using 3 strips.

	Standard [ng/ml]	Standard [ng/ml]	3
A	100	100	Sample A 1:40
B	50	50	Sample A 1:40
C	25	25	Sample B 1:40
D	12.5	12.5	Sample B 1:40
E	6.25	6.25	Sample C 1:40
F	3.125	3.125	Sample C 1:40
G	1.56	1.56	Sample D 1:40
H	Blank	Blank	Sample D 1:40

- 1.3 Dilute anti- human MOG (1-125) IgG standard (Component B) with 1X Sample Dilution Buffer (Component C) according to the Table 2.

Table 2. Serial dilution of anti-human MOG (1-125) IgG standard.

Step	Concentration [ng/ml]	Anti-MOG IgG standard	Sample Dilution Buffer (Component C)
1	100.00	Component B	0
2	50.00	250 µl from step 1	250 µl
3	25	250 µl from step 2	250 µl
4	12.5	250 µl from step 3	250 µl
5	6.25	250 µl from step 4	250 µl
6	3.125	250 µl from step 5	250 µl
7	1.56	250 µl from step 6	250 µl

- 1.4 Add 100 µl of the diluted standards and 100 µl of 1X Sample Dilution Buffer (Component C) as a blank into appropriate wells in duplicates.
- 1.5 Add diluted samples into appropriate wells (depends on the number of samples to be tested). After adding the standards and samples to the wells, cover the plate and incubate at room temperature for 60 min with gentle shaking.
- 1.6 Prepare 1X working wash buffer by diluting the 10X Wash Buffer (Component D) with DI H₂O.
- 1.7 Wash wells four times at 250 µl/well with 1X wash buffer. Pat dry.
- 1.8 Dilute goat anti-Human IgG-HRP (Component G) secondary antibody (2nd Ab) with Sample Dilution Buffer (Component C): working solution at 1:2,000 dilution (based on 0.1 mg/ml concentration). Add 100 µl of the diluted 2nd Ab into each well; incubate plate at room temperature for 45-60 min with gentle shaking.
- 1.9 Wash wells five times with 250 µl per well of 1X wash buffer. Pat dry. Clean the outside bottom of the wells with lens paper if necessary before the next step (this ensures accurate absorbance reading).
- 1.10 Add 100 µl of the TMB color substrate solution (Component E) into each well. Tap plate gently and incubate at room temperature until blue gradient is clearly observed across the wells (1-5 min). It may be necessary to adjust color development time so that absorbance values are within the detection range.
- 1.11 Add 50 µl of the Stop Solution (Component F) into each well and tap plate gently (blue color will turn to yellow). Measure absorbance (OD) at 450 nm using a microplate absorbance reader within 20 minutes after adding stop solution.

2. Calculate concentration of the samples.

- 2.1 Determine the average values (if replicates are used) for the standard and sample absorbance readings. Plot calibration curve using Four Parameter Logistic (4-PL) curve-fit. R² should be higher than 0.98. There should be at least 5 standard concentrations in the calculation to ensure statistical significance.
- 2.2 Choose absorbance values for the samples that are within the range used in the standard curve, and calculate the concentration of anti-human MOG (1-125) IgG using 4-PL curve-fit.
- 2.3 Example of calculation of Human anti-human MOG (1-125) IgG concentrations:

Please note, new standard curve must be generated each time the assay is run.


Table 3. An example of the assay with 4 samples.

	1	2	3
A	4.178	4.059	0.932
B	3.565	3.305	0.844
C	2.768	2.502	1.413
D	2.032	1.837	1.382
E	1.123	1.126	1.237
F	0.628	0.604	1.242
G	0.370	0.364	1.192
H	0.024	0.019	1.344

Note: Columns 1 and 2 are duplicate Human anti-Human MOG (1-125) IgG standards 100, 50, 25, 12.5, 6.25, 3.125, 1.56, and 0 ng/ml (Row A ~ H). Human Multiple Sclerosis Patients Serum: Sample-1, 3A-B; Sample-2, 3C-D; Sample-3, 3E-F; Sample-4, 3G-H at 1:40 dilution in duplicates.

2.3.1 Four-parameter logistic curve-fit (4-PL) based on the average absorbance reading values:

anti-Human MOG (1-125) IgG Standard Curve



$$Y=[(A-D)/(1+\{x/C\}^B)]+D, \underline{A}=-0.21 \underline{B}=0.787 \underline{C}=29.09 \underline{D}=5.687 \underline{R}^2=0.999$$

2.3.2 From 4-PL curve-fit data table (not shown) generated by computer software the following concentrations for human samples were obtained (based on the average absorbance readings):

	Absorbance@ 450nm, mean value	Calculated Concentration, [ng/ml]	Dilution Factor	Adjusted Serum Concentration, [ng/ml]
Sample1	0.866	4.34	1:40	173.6
Sample2	1.375	8.21	1:40	328.3
Sample3	1.217	6.86	1:40	274.6
Sample4	1.246	7.24	1:40	289.6

References:

1. Von Büdingen, H-C. et al. (2001) *J. Clin. Immunol.* 21 (3): 155-170.
2. Marta C.B. et al. (2005) *PNAS* 102 (39): 13992-13997.
3. Lyons, J-A. et. al. (1999) *European Journal of Immunology* 29: 3432-3439