



# SensoLyte<sup>®</sup> MFP Protein Phosphatase Assay Kit

## *\*Fluorimetric\**

Revision Number: 1.1		Last Updated: October 2014	
Catalog #		AS-71104	
Kit Size		500 Assays (96-well) or 1250 Assays (384-well)	

- **Convenient Format:** Complete kit includes all the assay components.
- **Optimized Performance:** Optimal conditions for detecting protein phosphatase activity.
- **Enhanced Value:** Less expensive than the sum of individual components.
- **High Speed:** Minimal hands-on time.
- **Assured Reliability:** Detailed protocol and references are provided.

## Kit Components, Storage and Handling

Component	Description	Quantity
Component A	MFP (1 vial), Ex/Em= 470 nm/510 nm upon phosphate group removal	1 mL
Component B	Assay buffer	60 mL
Component C	10X Lysis buffer	50 mL
Component D	Triton X-100	500 µL
Component E	Stop solution	30 mL
Component F	1 M DTT	100 µl

### Other Materials Required (but not provided)

- 96-well or 384-well microplate: Black, flat-bottom microplates with non-binding surface.
- Fluorescence microplate reader: Capable of detecting emission at 510±20 nm with excitation at 470±20 nm.
- Protease Inhibitors: Aprotinin, Leupeptin, PMSF and Pepstatin A

### Storage and Handling


- Store kit components at -20°C
- Components B, C, D, and E can be stored at room temperature for convenience

---

## Introduction

Protein phosphorylation/dephosphorylation, a potent and versatile mechanism for the regulation of protein activity, plays a key role in signal transduction and cellular function modulations. Consequently, protein phosphatases have received great attention as potential drug-screening targets.

MFP is a proprietary fluorogenic substrate for measuring the activity of protein phosphatases, such as protein tyrosine phosphatases, serine/threonine phosphatases,  $\text{Na}^+/\text{K}^+$  ATPase<sup>1</sup>, and plasma membrane  $\text{Ca}^{2+}$ -ATPase. Upon de-phosphorylation by phosphatases, MFP generates MF, which has bright green fluorescence and can be detected at excitation/emission=470 nm/510 nm. The fluorescence of MF is pH-insensitive, its maximal fluorescence intensity is at a wide pH range of between 4 to 10 - an important feature since a variety of phosphatases have different optimal pH requirements. For those phosphatases that prefer buffer with pH <6.0, MFP is the choice for measuring their activity. For phosphatases that prefer buffer with pH ≥6.0, although MFP is still good for them, but FDP-based phosphatase assay (Cat# 71100) will be more sensitive.



**Figure 1.** The assay sensitivity in measuring protein tyrosine phosphatase 1B (PTP 1B).

Recombinant PTP 1B was serially diluted in assay buffer, pH 6.5, and its activity was measured according to the protocol. The assay can detect as low as 0.5 ng of PTP 1B.

---

## Protocol

Note 1: Warm all kit components to room temperature before starting the experiment.

Note 2: Since MFP is a generic phosphatase substrate, in order to measure the activity of the phosphatase of interest, your sample has to be purified by immuno-affinity or other methods before measuring its activity using MFP.

### 1. Prepare protein phosphatase-containing sample.

- 1.1 For protein phosphatase-containing biological sample, please refer to **Appendix I** for the preparation of cell extract or tissue extract.
- 1.2 For purified protein phosphatase, dilute the enzyme in assay buffer (Component B) to the appropriate concentration.

Note: The activity of protein phosphatase can be preserved better if the purified enzyme is diluted with assay buffer containing 1 mg/mL of bovine serum albumin. Keep enzyme on ice before the experiment. Avoid vigorously mixing of the enzyme.

## 2. Prepare MFP reaction solution according to Table 1.

**Table 1. MFP reaction solution for one 96-well plate (100 assays)**

Components	Volume
MFP (100X, Component A)	50 $\mu$ L
Assay Buffer (Component B)	4.935 mL
1 M DTT (Component F)	15 $\mu$ L
Total volume	5 mL

Note: The assay buffer (Component B) works for protein tyrosine phosphatase. Since some protein phosphatases require their unique assay buffer, you may use your own formulated buffer. Please refer to Appendix II for some references.

## 3. Start the protein phosphatase detection.

- 3.1 Add 50  $\mu$ L/well (black 96-well plate) or 20  $\mu$ L/well (black 384-well plate) of the protein phosphatase-containing sample. Include a non-phosphatase-containing sample as a negative control.
- 3.2 Add 50  $\mu$ L/well (96-well plate) or 20  $\mu$ L/well (384-well plate) of the MFP reaction solution. Mix the reagents by gently shaking the plate for 30 sec.
- 3.3 Measure fluorescence signal:
  - For kinetic reading: Immediately start measuring fluorescence intensity at Ex/Em=470  $\pm$ 20 nm/510  $\pm$  20 nm continuously and record data every 5 min for 30 to 60 min.
  - For end-point reading: Incubate the reaction at the desired temperature for 30 to 60 min, and keep away from light. Optional: Add 50  $\mu$ L/well (96-well plate) or 20  $\mu$ L/well (384-well plate) of stop solution (Component E). Measure fluorescence intensity at Ex/Em=470  $\pm$ 20 nm/510  $\pm$  20 nm.

## Appendix I. Preparation of cell and tissue extract

### Prepare cell extract for protein phosphatase assay.

- Prepare 1X lysis buffer by adding 1 mL of 10X lysis buffer (Component C) to 9 mL of deionized water.
- Wash cells with 1X lysis buffer twice gently.
- Add 20  $\mu$ L of Triton X-100 (Component D) to 10 mL of 1X lysis buffer. Add protease inhibitors to a final concentration of 10  $\mu$ g/ml Aprotinin, 10  $\mu$ g/mL Leupeptin, 100  $\mu$ M PMSF and 10  $\mu$ g/ml Pepstatin A.  
Note: Protease inhibitors are not provided.
- Add an appropriate amount of the above lysis buffer to the cells or cell pellet. Scrape off the adherent cells or resuspend the cell pellet. Collect the cell suspension in a microcentrifuge tube.
- Incubate the cell suspension at 4°C for 10 min with agitation.
- Centrifuge the cell suspension at 2500 X g for 10 min at 4°C.
- Collect the supernatant to perform the protein phosphatase assay.

### **Prepare tissue extract for protein phosphatase assay.**

- Prepare 1X lysis buffer by adding 1 mL of 10X lysis buffer (Component C) and 20 µL of Triton X-100 (Component D) to 9 mL of deionized water. Add protease inhibitors to a final concentration of 10 µg/ml Aprotinin, 10 µg/mL Leupeptin, 100 µM PMSF, and 10 µg/ml Pepstatin A.

Note: Protease inhibitors are not provided.

- Add an appropriate amount of 1X lysis buffer to the tissue sample, and homogenize.
- Centrifuge the tissue sample at 10,000 X g for 10 min at 4°C.
- Collect the supernatant for protein phosphatase assay.

### **Appendix II. References for protein phosphatase assay buffer**

Phosphatases	Assay Buffer
CD45, PTP1B	50 mM Bis-tris, pH 6.5, 2 mM EDTA, 5 mM DTT, 0.05% Brij35 <sup>2</sup>
PP1	100 mM Tris-HCl, pH 7.5, 4 mM DTT, 0.2 mM EDTA, 0.5 mM MnCl <sub>2</sub> , 0.4 mg/mL BSA <sup>3</sup>
PP2A	40 mM Tris-HCl, pH 8.4, 34 mM MgCl <sub>2</sub> , 4 mM EDTA, 4 mM DTT <sup>4</sup>
Na <sup>+</sup> /K <sup>+</sup> ATPase	80 mM Tris-HCl, pH 7.2, 4 mM MgCl <sub>2</sub> , 0.5 mM EGTA, 5 mM creatinine phosphate, activated by 10 mM KCl
PTEN	100 mM Tris-HCl, pH 8, 10 mM DTT <sup>5</sup>
PP2Cα	50 mM Tris-HCl, pH 7.5, 0.1 mM EDTA, 0.5 mM DTT, 60 mM MgCl <sub>2</sub> <sup>6</sup>

### **References**

1. Johansson, M. et al. *J. Clin. Endocrinol. Metab.* **88**, 2831 (2003).
2. Huang, Z. et al. *J. Biomol. Screen.* **4**, 327 (1999).
3. Heresztyn, T. et al. *Environ. Toxicol.* **16**, 242 (2001).
4. Takai, A. et al. *Biochem. J.* **287** (Pt 1), 101 (1992).
5. Maehama, T. et al. *Anal. Biochem.* **279**, 248 (2000).
6. Marley, A. E. et al. *Biochem. J.* **320** (Pt 3), 801 (1996).