



SensoLyte[®] 520 Generic Caspase Assay Kit *Fluorimetric*

Revision number: 1.1	Last updated: 10/20/14
Catalog #	AS-72211
Kit Size	100 Assays

- **Optimized Performance:** This kit is optimized for detecting generic caspase activity in cells.
- **Enhanced Value:** It provides ample reagents to perform 100 assays.
- **High Speed:** The cells can be stained and quantified within 2 hours.
- **Assured Reliability:** Detailed protocol and references are provided.

Kit Components, Storage and Handling

Component	Description	Quantity
Component A	FAM-VAD-FMK Reagent (Lyophilized powder)	4 vials
Component B	Dimethyl sulfoxide (DMSO)	1 mL
Component C	10X Phosphate buffered saline (10X PBS), pH 7.4	60 mL
Component D	Hoechst 33342 solution	1 mL, 0.2 mg/mL
Component E	Propidium Iodide (PI) solution	1 mL, 0.25 mg/mL
Component F	Fixative solution	12 mL

Other Materials Required (but not provided)

- Slides, coverslips, and fluorescence microscope
- 96-well black microplate and fluorescence microplate reader (capable of detecting excitation at 490±20 nm with emission at 520 ±20 nm)

Storage and Handling

- Store Component A, D, E and F at 4°C. Protect Component A from light and moisture.
- Components B, C can be stored at room temperature for convenience.

Introduction

Apoptosis is involved in a variety of physiological and pathological events,¹ ranging from normal fetal development to development of diseases, such as cancer,² organ failure and neurodegenerative diseases. Caspases are a large family of cysteinyl aspartate-specific proteases. They are activated in response to diverse cell death stimuli and ultimately dismantle a cell through restricted proteolysis of numerous cellular proteins, leading to apoptotic cell death.³

The SensoLyte[®] 520 Generic Caspase Assay Kit provides a method for in situ detection of caspase activity in cells undergoing apoptosis using a fluorochrome inhibitor of caspase, FAM-VAD-FMK. This reagent is cell permeable and non-cytotoxic. When added to a population of cells, FAM-VAD-FMK enters cells and covalently binds to activated caspases-1, -3, -4, -5, -7, -8, and -9.⁴ Bound FAM-VAD-FMK is retained within the cell, while any unbound FAM-VAD-FMK diffuses out of the cells and is washed away. Signal in the cells containing the bound FAM-VAD-FMK can be evaluated with a fluorescence plate reader, fluorescence microscopy, or flow cytometry.

Preparations

1. Prepare reagents

- 1.1 FAM-VAD-FMK stock solution: Reconstitute one vial of lyophilized powder (Component A) with 50 µL of DMSO (Component B). If not used immediately, stock solution may be frozen in aliquots in -20°C for up to 6 months.
- 1.2 FAM-VAD-FMK working solution: Dilute stock solution 5 fold in phosphate buffered saline (PBS), immediately before use. The working solution must be prepared fresh for each experiment. Protect from light.
- 1.3 PBS: Dilute 10X PBS (Component C) 10 fold in dH₂O.
- 1.4 Hoechst solution (Component D), PI solution (Component E), and fixative solution (Component F) are ready to use.

2. Prepare apoptotic cells.

- 2.1 Cells should be cultured at a density that is optimal for apoptosis induction.

Note: Optimal cell density will vary depending on the cell line used and experimental conditions. In general, cell density should not exceed 10⁶/mL.

- 2.2 Induce apoptosis according to your protocol.
- 2.3 Concurrently incubate a non-induced negative control cell population at the same density as the induced cell population for each labeling conditions.

Protocol

1. Protocol for fluorescence microplate reader.

- 1.1 Concentrate both induced and non-induced cells to ~1 x 10⁷ cells/mL.
- 1.2 Transfer 300 µL of each cell suspension to sterile tubes.

- 1.3 Add 10 μ L FAM-VAD-FMK working solution and mix cells by slightly flicking the tubes.
- Note: Optimal amount of FAM-VAD-FMK can vary depending on cell line used and experimental conditions.
- 1.4 Incubate cells with FAM-VAD-FMK for 1 hour, at 37 °C under 5% CO₂. Protect tubes from light. Swirl tubes twice during this time to gently resuspend settled cells.
- 1.5 Wash cells by adding 1 mL PBS with gently mixing. Centrifuge cells at 1000 rpm, 1 min at room temperature (RT). Carefully remove and discard supernatant and gently vortex cell pellet to disrupt any cell-to-cell clumping. Repeat washing procedure twice.
- 1.6 Resuspend cells in 400 μ L PBS. Take aliquot from each tube to count cells with a hemocytometer or other method. Compare density of induced and non-induced cell populations.
- 1.7 Place 3 replicates of 100 μ L for each cell suspension into 96-well black microtiter plate.
- 1.8 Measure fluorescence intensity at Ex/Em=490/520 nm.
- 1.9 Analyze results: Normalize RFU based on counted cell number in samples (Fig.1).




Figure 1. Active caspase was detected with FAM-VAD-FMK in Jurkat cells (n=3). Apoptosis was induced with 1 μ M staurosporin for 4 h and control cells were treated with DMSO at 37°C under 5% CO₂. After treatment, FAM-VAD-FMK was incubated with control and induced cells at 37°C for 1 h. Fluorescence signal was measured by a microplate reader (FlexStation II384, Molecular Device, CA) at Ex/Em=490/520 nm.

2. Protocol for fluorescence microscopy

- 2.1 Transfer 300 μ L of each cell suspension to sterile tubes.
- 2.2 Add 10 μ L FAM-VAD-FMK working solution and mix cells by slightly flicking the tubes.
- Note: The amount of FAM-VAD-FM can be adjusted based on cell type and experiment conditions.
- 2.3 Incubate cells for 1 hour at 37°C under 5% CO₂, protecting tubes from light. Swirl tubes twice during this time to gently resuspend settled cells.
- 2.4 Hoechst stain can be used to label the nuclei of cells undergoing apoptosis. If desired, add 1.5 μ L Hoechst solution (Component D, 0.5% v/v). Incubate for 15 minutes at 37°C under 5% CO₂.
- 2.5 Add 1 mL PBS to each tube and gently mix.

- 2.6 Centrifuge cells at 1000 rpm for 1 min at RT.
- 2.7 Carefully remove and discard supernatant, and gently vortex cell pellet to disrupt any cell-to-cell clumping.
- 2.8 Wash cells again with 1 mL PBS.
- 2.9 Resuspend cell pellet in 300 μ L PBS and place cells on ice.
- 2.10 PI stain can be used to exclude dead cells from the analysis. If desired, add 1.5 μ L PI solution (Component E, 0.5% v/v). Incubate for 5 min at RT.
- 2.11 Place one drop of cell suspension onto a microscope slide. Cover cells with a cover slip.
- 2.12 To analyze cells later, the following fixative procedure can be used.
 - Add 100 μ L Fixative solution (Component F) to each 300 μ L cell suspension.
 - Incubate cells for 15 min at RT in the dark.
 - Dry cells onto a microscope slide.
 - Cover cells with mounting media and cover slip.
 - Store slides at 2-8°C protected from light for up to 24 hours.
 - Proceed to step 2.13.
- 2.13 Observe cells under a fluorescence microscope.
 - Observe cytoplasm of FAM-VAD-FMK labeled cells using filters at Ex/Em 490/520 nm
 - Observe nuclei of Hoechst labeled cells using filters at Ex/Em 350/461 nm.
 - Observe nuclei of PI labeled cells using filters at Ex/Em 535/617 nm.

References

1. Thornberry, N.A. and Lazebnik Y. *Science* **281**, 1312 (1998).
2. Reed, J.C. *J. Clin. Oncol.* **17**, 2941 (1999).
3. Timmer, JC. et al. *Cell Death Differ.* **14**, 66 (2007).
4. Thornberry, NA. et al. *J. Biol. Chem.* **272**, 17907 (1997).