

Anti-IGF-1R aptamer, Direct Magnetic AP Kit

Introduction

When a protein is expressed at low levels and is difficult to detect with western blot analysis, aptoprecipitation (AP, Aptamer based protein pull down method) may be the method of choice. An aptoprecipitating reagent has to be specific in order to avoid precipitation of unwanted protein. Furthermore, sufficient affinity is required to pull down the protein and it has to withstand stringent washing steps. AptSci IGF-1R aptamer molecule is a specific affinity ligand and has been proven well suited for pull down experiments of IGF-1R proteins. Most commonly encountered problems with IP approach is interference from antibody heavy and light chains that may comigrate with relevant bands, masking important results. However aptamer as an oligonucleotide will not contribute to protein/peptide background that can interfere with subsequent analysis.

AptSci has developed proprietary protein pull down method using target protein-specific aptamers. The aptamer-coupled magnetic bead included in the kit has low nonspecific binding characteristic and enables convenient magnetic isolation of protein targets and reusable magnetic beads. Mild elution condition enables isolation of non-denatured proteins which can be used for further study.




Fig. 1. Flow cytometry histograms showing the binding of representative IGF-1R aptamer against the target HEK293 cells. Approximately 1×10^6 cells were washed and incubated with FITC-conjugated IGF-1R aptamer (Pink histogram). The untreated cell was used as background fluorescence signal (Gray histogram).

Result of Aptoprecipitation (AP)

As shown in Figure 2, Western blot analysis revealed that IGF-1R aptamer precipitates IGF-1R with high specificity, while no IGF-1R was observed when precipitating with blank beads. This result indicates that IGF-1R aptamer is highly specific to IGF-1R protein and IGF-1R aptamer-coupled magnetic bead efficiently precipitates IGF-1R from a protein complex.




Fig. 2. Aptoprecipitation of IGF-1R protein from HEK293 cells using the AptSci direct IGF-1R AP Kit. HEK293 cell lysates (1mg/lane) were incubated with IGF-1R aptamer (150pmol)-coupled magnetic bead. After washing the beads, the bound protein was eluted with boiling SDS loading buffer. The eluate was loaded on a SDS-PAGE (4-15% gradient gel) and blotted onto a PVDF membrane. Western blot was probed with anti IGF-1R Ab.

Product Information

- Product name: Anti-IGF-1R aptamer, Direct Magnetic AP Kit
- Catalog number: IGF-1R-2293DM
- Content: Magnetic agarose conjugated IGF-1R aptamer molecule and all buffers required to perform small scale AP
- Form: As 25% slurry in 20% ethanol containing 0.04% (w/v) sodium azide.
- Protein source for generation of aptamer: Recombinant protein produced in mammalian cells
- Specificity: Anti-IGF-1R aptamer binds to human IGF-1R. Cross reactivity with other species has not been tested.
- MW: ~26 kDa
- Conjugation yield: > 90% as determined by spectrometer analysis.
- Tested applications: FACS and Aptoprecipitation.
- Storage: At +4°C.
- **Shipping**: At ambient temperature.
- Stability: There is no decrease in performance of the kit after storage for 6 months at ambient temperature.



www.aptsci.com. Aptamer Sciences Inc.

Postech Biotech Center, San31, Hyoja-Dong, Pohang,
Gyeongbuk 790-784, South Korea. TEL +82-54-279-8691
FAX +82-54-279-8245. E-mail: aptamer@aptsci.com.

LIMITATIONS

Warranty: AptSci AptoPrepTM products are warranted to meet stated product specifications and to confirm to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to one year from date of sales for products used, handled and stored according to AptSci's instructions. AptSci's sole liability is limited to replacement of the product or refund of the purchase price. AptoPrepTM products are supplied for research use only. They are not intended for medicinal, diagnostic or therapeutic use. AptoPrepTM products may not be resold, modified for research use of the manufacture commercial products without prior written approval from AptSci.