ExpressArt FFPE RNAready

FFPE RNA isolation kit

Catalogue No. 9000-A100

for 100 RNA isolations from archival paraffin samples

This protocol provides the required laboratory procedures

For additional information about ExpressArt technology and products see http://www.amp-tec.com/products » more information

Kit Contents

Box I	Volume per RNA isolation	Kit contents for 100 isolations
Lysis Buffer (LB)	500 μl	55 ml
Demodification Solution (DeS)	130 µl	15 ml
Wash Buffer 1 (WB 1)	1 ml	70 ml concentrate
Add 35.0 ml ethanol (96-100%) to		
70.0 ml concentrate of WB 1		
Wash Buffer (WB 2)	1.5 ml	35 ml concentrate
Add 140.0 ml ethanol (96-100%) to		
35.0 ml concentrate of WB 2		
RNase-free water	30 µl	5 ml
Spin columns	1 piece	100 pieces

Box II	Volume per RNA isolation	Kit contents for 100 isolations
NucleoGuard	5 μl	0.6 ml
Proteinase K, PCR Grade*	30 µl	powder
RNase-Free DNase**	45 μl (working solution)	powder
DiB, Dilution buffer for DNase	39.13 µl	2x 2 ml

* Dissolve proteinase K (powder) in 4 ml water (25 mg/ml); store at -20°C.

**Prepare "DNase working solution"
First step:
Dissolve DNase (powder) in 600 µl water; store at -20℃.
Second step:
DNase working solution (prepare according to your needs)
The example is for 51 isolations:
Mix 0.3 ml DNase + 2.0 ml Dilution Buffer (DiB)
Unused Aliguots are stored at -20℃

<u>Please note, these reagents are not included in the kit</u> Xylene Ethanol, abs.

Required instrument

Thermomixer (Eppendorf # 5350 000.013 or "Thriller" from PEQLAB #91-7010) or equivalent

<u>Tissue samples</u>

Up to 5 FFPE slides (10 µm thick with appr. 1-2 cm² of tissue) per isolation

1. Deparaffinisation (in Thermomixer with 1000 rpm)

Insert paraffin sections in reaction tube Add 1.5 ml xylene and incubate for 20 min@37℃ Then centrifuge for 3 min @ 16'000 g and remove xylene (important!) Repeat: Add 1.5 ml xylene and incubate for 20 min@37℃ Then centrifuge for 3 min @ 16'000 g and remove xylene (important!)

Then add 1 ml EtOH abs "Flick tubes" (to dislodge pellet) Centrifuge for 3 min @ 16'000 g and remove ethanol Dry at 37℃ (can take some time – but pellet has to be dry – not "**sticky**"!)

Dissolve dry pellet in 500 μl **Lysis Buffer (LB)** + 30 μl Proteinase K (25 mg/ml) + 5 μl NucleoGuard

2. Lysis: Incubation in Thermomixer (1000 rpm) for 3 h @ 55°C

<u>3. Demodification</u>: Add 130 µl **Demodification Solution (DeS)** (to 535 µl lysate) Vortex, then 20 min @ 80℃ Chill on ice (3 min)

4. RNA Purification

Before you start: Put RNase-free water in thermoblock for heating to 95°C

Centrifuge lysate: 3 min @ 13'000 rpm Transfer 600 μ l of the supernatant in 2 ml reaction tube Add 750 μ l EtOH abs (1.25x Vol.) Mix by pipetting (Volume = 1350 μ l)

Transfer aliquot of 700 µl of the EtOH-supplemented lysate on column Centrifuge: 30 sec @ 10'000 rpm & discard flow-through Transfer the remaining lysate on column & discard flow-through Add 500 µl **Wash Buffer 1 (WB 1)** Centrifuge: 30 sec @ 10'000 rpm & discard flow-through Plus 500 µl **Wash Buffer 2 (WB 2)** Centrifuge: 30 sec @ 10'000 rpm & discard flow-through Centrifuge: 1 min @ Max <u>4.1. DNA-Digestion:</u> Transfer 45 µl of the "DNase Working solution" on column Incubate: 15 min @ R.T.
Plus 500 µl Wash Buffer 1 (WB 1) Incubate: 1 min
Centrifuge: 30 sec @ 10'000 rpm & discard flow-through
Plus 500 µl Wash Buffer 2 (WB 2)
Centrifuge: 30 sec @ 10'000 rpm & discard flow-through
Plus 500 µl Wash Buffer 2 (WB 2)
Centrifuge: 30 sec @ 10'000 rpm & discard flow-through
Plus 500 µl Wash Buffer 2 (WB 2)
Centrifuge: 30 sec @ 10'000 rpm & discard flow-through
Plus 500 µl Wash Buffer 2 (WB 2)

<u>4.2. Elution</u>: Transfer column in fresh 1.5 ml reaction tube
Add 30 µl RNase-free water (preheated at 95℃)
Incubate: 2 min
Centrifuge: 1 min @ Maximum speed
Reapply the eluate on column
Incubate: 2 min
Centrifuge: 1 min @ Max
Transfer eluate to fresh tube for storage or for further processing (avoid risk: lid can get contaminated during entrifugation).

RNA is recovered in a total volume of ~30 µl

Expected yield ~10-30 µg RNA (~2-6 µg per slide)

AmpTec GmbH Koenigstrasse 4a | 22767 Hamburg | Germany Tel: +49 (0)40 – 636 747 22 | Fax: +49 (0)40 – 636 747 19

Technical Support: "Dr. Guido Krupp" <krupp@amp-tec.com> "Dr. Peter Scheinert" < scheinert@amp-tec.com > www.amp-tec.com