ExpressArt RNA Clean RNAready

**RNA Clean-up kit** 

## Catalogue No. 9003-A100

## for 100 RNA Clean-up procedures

# Purification of in vitro transcription products, Removal of small RNAs (below ~200nt) from RNA samples, Dye removal from labelled RNA

This protocol provides the required laboratory procedures

For additional information about ExpressArt technology and products see http://www.amp-tec.com/products » more information

#### **Kit Contents**

	Volume per RNA isolation	Kit contents for 100 isolations
Lysis Buffer (LB)	340 µl	40 ml
DeS Buffer	87 µl	10 ml
Wash Buffer 1 (WB 1): add 18.0 ml ethanol (100%)	500 µl	36 ml
to 36.0 ml concentrate of WB 1		concentrate
Wash Buffer (WB 2): add 44.0 ml ethanol (100%)	500 µl	11 ml
to 11.0 ml concentrate of WB 2		concentrate
RNase-free water	30 µl	5 ml
Spin columns with collection tubes	1 piece	100 pieces

#### Please note, these reagents are not included in the kit

Ethanol, abs.

#### **Chemical hazards**

The Lysis Buffer (LB) contains guanidine thiocyanate, which is harmful in contact with skin, when inhaled or swallowed. Guanidine thiocyanate also liberates toxic gas, when mixed with strong acids. Always store and use the Binding Buffer away from food. Always wear gloves, and follow standard safety precautions during handling and make sure to comply with the safety rules of your laboratory.

#### PROCEDURE

Purification of in vitro transcription products, removal of small RNAs (below ~200nt) from RNA samples, dye removal from labelled RNA, etc., with Spin Columns

Before starting: Add 18 ml of 100% ethanol to WB1 concentrate (36 ml) and 44 ml of 100% ethanol to WB2 concentrate (11 ml), as indicated on the bottles.

At the start: Put RNase-free water in thermoblock at 95°C.

RNA Purification Mix		
LB Buffer	340 µl	
DeS Buffer	87 µl	
100% Ethanol (at room temperature; not supplied)	560 µl	

- Add 987 µl of **RNA Purification Mix** to 20 µl RNA sample (adjust smaller volumes by addition of water). Mix thoroughly.
- Transfer an aliquot of 700 µl of this mixture on the RNA spin column, centrifuge for 30 sec@10'000 rpm & discard flow-through
- Transfer the **remaining mixture** on column, centrifuge for 30 sec@10'000 rpm & discard flow-through
- Add 500 µl Wash Buffer 1 (WB 1), centrifuge for 30 sec@10'000 rpm & discard flowthrough
- Add 500 µl Wash Buffer 2 (WB 2), centrifuge for 30 sec@10'000 rpm & discard flowthrough
- Re-insert the column in the same collection tube. Add 500 μl 80% EtOH, centrifuge for 30 sec@10'000 rpm & discard flow-through
- Re-insert the column in the same collection tube, and centrifuge for 1 min@maximum speed to get rid of residual salt on the spin column matrix.
- <u>Elution</u>: Transfer column in fresh 1.5 ml reaction tube and add 30 μl RNase-free water (preheated at 95°C)
- Incubate: 2 min and centrifuge for 1 min@maximum speed
- Reapply the eluate on the column
- Incubate: 2 min and centrifuge for 1 min@maximum speed
- RNA is in a total volume of ~ 30 µl.
- Use immediately or store at -80°C.
- Long-term storage or shipping:
  - a) for concentrations  $\geq 1 \mu g/\mu I$ , storage at -80°C and shipping on dry ice is adequate
  - b) for low concentrations, storage/shipment as EtOH suspension is recommended: add 1/10 volume of sodium acetate (pH 5) and 1 to 2 µl of carrier, like "Pellet Paint" (Novagen), mix thoroughly; then add 2.2 volumes of EtOH (100%).
    Store at -20℃ / ship on dry ice.

RNA is recovered by centrifugation (20min@max speed, at 4 $^{\circ}$ C), remove supernatant with a pipette, wash with 200 µl 80% EtOH, spin again, discard supernatant and air dry the RNA pellet, dissolve in RNase-free water.



### AmpTec GmbH

Koenigstrasse 4a 22767 Hamburg Germany Tel: +49 (0)40 41 364 781 Fax: +49 (0)40 41 364 887 Technical Support: "Dr. Guido Krupp" <krupp@amp-tec.com> "Dr. Peter Scheinert" < scheinert@amp-tec.com > www.amp-tec.com