# **ExpressArt LBR RNAready**

**RNA isolation kit** <u>for solid tissues or for Bacteria</u> in combination with mechanical lysis (beads)

Catalogue No. 9005-A100

# for 100 RNA isolations

# including on-column DNase digestion

This protocol provides the required laboratory procedures

For additional information about ExpressArt technology and products see http://www.amp-tec.com/products » more information

#### **Kit Contents**

Box I	Volume per RNA	Kit contents for
	isolation	100 isolations
1 <sup>st</sup> Lysis Buffer LBR-1 (no detergent)	350 µl	40 ml
2 <sup>nd</sup> Lysis Buffer LBR (with detergent)	300 µl	30 ml
Wash Buffer 1 (WB 1): add 35.0 ml ethanol (100%) to	1 ml	70 ml concentrate
70.0 ml concentrate of WB 1		
Wash Buffer (WB 2): add 140.0 ml ethanol (96-100%)	1.5 ml	35 ml concentrate
to 35.0 ml concentrate of WB 2		
RNase-free water	30 µl	5 ml
Spin columns	1 piece	100 pieces

Box II	Volume per RNA Kit contents for	
	isolation	100 isolations
NucleoGuard	6.5 μl	650 µl
RNase-Free DNase	45 μl	powder
	(working	
	solution)	
DiB, Dilution buffer for DNase	39.13 µl	2x 2 ml

#### Please note, these consumables are not included in the kit

<u>For solid tissues:</u> Precellys CKM tubes (<u>www.peqlab.com</u> Cat.-No. 91-PCS-CKM) <u>For bacterial samples:</u> Precellys VK01 tubes (<u>www.peqlab.com</u> Cat.-No. 91-PCS-VK01) Ethanol, abs.

 $\beta$ -Mercapto-ethanol

#### **Required equipment**

Precellys homogeniser (<u>www.peqlab.com</u>) or FastPrep homogeniser (<u>www.qbiogene.com</u>) Thermomixer (<u>www.eppendorf.com</u>) or Thermoshaker (<u>www.peqlab.com</u>)

#### **Before you start:**

- 1) Dissolve lyophilized DNase enzyme in 600 μl RNase-free water. Unused Aliquots are stored at -20°C.
- 2) Prepare "DNase working solution" (51 isolations) Add 0.3 ml DNase to 2.0 ml Dilution Buffer (DiB) Unused Aliquots are stored at -20°C.
- **3)** Prepare supplemented Lysis Buffers Use supplemented buffer within one day.

<u>Preparation of 10.0 ml LBR-1+ and LBR+</u> 10 ml of Lysis Buffer, add 100 μl <u>each</u>: β-Mercapto-Ethanol and NucleoGuard. Dispense in a fume hood and wear appropriate protective clothing.

4) Preheat RNase-free water in thermoblock to 95°C

### PROCEDURE

### A. Lysis step

### A1. Solid tissues

- A1.1 Use Precellys CKM tube (contains a mixture of 1.4 & 2.8mm ceramic Zirkonium beads) and add <u>350 µl 1<sup>st</sup> Lysis Buffer LBR-1+</u>
- A1.2 Transfer fresh/frozen biopsy (up to ~50 mg) to the CKM tube. Homogenise for 2x 30sec@6.5K (15 sec Pause) Maybe, some optimisation is required.
- A1.3 Add 300 µl 2<sup>nd</sup> Lysis Buffer LBR+
- A1.4 Incubate for 15 min@RT in Thermomixer@1000 rpm.
- A1.5 Transfer ~600 µl supernatant to fresh 1.5ml tube, spin for 5 min@max speed.
- A1.6 Transfer ~580 μl supernatant to fresh 2ml tube <u>Per 580 μl of transferred lysate:</u> Add 150 μl water and 910 μl EtOH abs and mix.

Proceed to step B.

### A2. Bacterial samples

- A2.1 Use Precellys VK01 tube (contains 0.1 mm glass beads) and add <u>350 µl 1<sup>st</sup> Lysis Buffer LBR-1+</u>
- A2.2 Transfer fresh/frozen sample (up to 50 mg) to the VK01 tube. Homogenise for 2x 20sec@5.5K (5 sec Pause) Maybe, some optimisation is required.
- A2.3 Add 300 µl 2<sup>nd</sup> Lysis Buffer LBR+
- A2.4 Incubate for 15 min@RT in Thermomixer@1000 rpm.
- A2.5 Transfer ~600 µl supernatant to fresh 1.5ml tube, spin for 5 min@max speed.
- A2.6 Transfer ~580 μl supernatant to fresh 2ml tube <u>Per 580 μl of transferred lysate:</u> Add 150 μl water and 910 μl EtOH abs and mix.

Proceed to step B.

### **B.** Spin column purification

B1	Transfer an aliquot of 700 µl of the mixture on column Centrifuge: 30 sec @ 10,000 rpm & discard flow-through Transfer the remaining mixture on column Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
B2	Add 500 μl Wash Buffer 1 (WB 1) Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
B3	Plus 500 μl Wash Buffer 2 (WB 2) Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
<u>B4</u>	<u>DNA-Digestion:</u> Transfer 45 μl of the "DNase Working solution" on column Incubate: 15 min @ RT.
B5	Plus 500 μl Wash Buffer 1 (WB 1) Incubate: 1 min Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
<b>B6</b>	Plus 500 μl Wash Buffer 2 (WB 2) Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
B7	Plus 500 μl Wash Buffer 2 (WB 2) Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
<b>B8</b>	Add 500 µl 80% EtOH (at room temperature) Centrifuge: 30 sec @ 10,000 rpm & discard flow-through
<b>B9</b>	Re-insert the column in the same collection tube, centrifuge for 1 min@maximum speed to get rid any residual liquid
<b>B10</b>	Elution: Transfer column in fresh 1.5 ml reaction tube and add 30 $\mu$ l RNase-free water (preheated at 95°C)
	Incubate: 2 min and centrifuge for 1 min@maximum speed
B11	Reapply the eluate on the column Incubate: 2 min and centrifuge for 1 min@maximum speed
DNA	is in a total valume of 20 ul

6 41

RNA is in a total volume of  $\sim 30 \mu l$ . Use RNA immediately or store at -80°C.

### Note for subsequent mRNA amplification

If required, the volume of the eluted RNA can be reduced by limited concentration in a SpeedVac (avoid complete drying!).

### **Further note:**

A modified kit protocol can be used for mRNA amplification with the ExpressArt mRNA amplification kits:

A larger volume of the isolated RNA (~10  $\mu$ l) can be used, if water is omitted in Mix 1 and Mix 2.

Technical Support:	Dr. Guido Krupp	<u>krupp@amp-tec.com</u>
	Dr. Peter Scheinert	t scheinert@amp-tec.com

#### **AmpTec GmbH**

Koenigstrasse 4a | 22767 Hamburg | Germany | Tel:+49 40 636 747 22 | Fax: +49 40 636 747 19