## ReliaBLOT® 12% SDS-PAGE Gels 10 X 10 cm, 12-well



### Catalog No. WB101-01212G

Bethyl Laboratories ReliaBLOT® 12% SDS-PAGE (sodium dodecyl sulfate polyacrylamide **Description:** 

gel electrophoresis) gels are ready-to-use polyacrylamide gels cast in plastic 10(W) X 10(H) cm cassettes. The gels are cast with 12 sample wells, and the gel matrix is 1mm in

thickness.

10 gels **Amount:** 

Bethyl Laboratories ReliaBLOT® 12% SDS-PAGE gels are used to resolve proteins under **Application:** 

denaturing conditions. For optimal results, The ReliaBLOT® Running Buffer and The ReliaBLOT<sup>®</sup> Transfer Buffer should be used for electrophoresis and transfer of proteins in ReliaBLOT® SDS-PAGE Gels. Subsequent to electrophoresis, the gels can be transferred

for western blot analysis or the gel may be stained and dried.

**Preparation:** Remove the gel from the package, drain away the buffer, and rinse the gel and wells with

> distilled water. Place the gel on the buffer core and assemble the electrophoresis cell according to the manufacturer's directions. Prepare a 1X solution of ReliaBLOT® Running Buffer for the inner and outer chambers. Run the gel(s) at 150 V until the dye front reaches

the bottom of the gel (approximately 60 minutes).

**Storage:** Store at 2 - 8°C.

**Stability:** Stable for 2 years from date of manufacturing.

**Expiration Date:** See gel cassette.

 $\label{eq:Reliablot} \begin{aligned} & \text{ReliaBLOT}^{\textcircled{\$}} \text{ SDS-PAGE Gels, 4-8\%, 10 x 10 cm} \\ & \text{ReliaBLOT}^{\textcircled{\$}} \text{ SDS-PAGE Gels, 4-12\%, 10 x 10 cm} \end{aligned}$ **Related Products:** WB101-40812G

WB101-41212G ReliaBLOT<sup>®</sup> SDS-PAGE Gels, 4-20%, 10 x 10 cm WB101-42012G ReliaBLOT<sup>®</sup> SDS-PAGE Gels, 4-8%, 10 x 8 cm WB102-40812G ReliaBLOT<sup>®</sup> SDS-PAGE Gels, 4-12%, 10 x 8 cm WB102-41212G ReliaBLOT® SDS-PAGE Gels, 12%, 10 x 8 cm ReliaBLOT® SDS-PAGE Gels, 4-20%, 10 x 8 cm WB102-01212G WB102-42012G ReliaBLOT® Prestained Protein Marker WB103-600 ReliaBLOT<sup>®</sup> LDS Sample Buffer (4X) WB104-10 ReliaBLOT® Running Buffer (20X) WB105-500 ReliaBLOT® Transfer Buffer (10X)
ReliaBLOT® Nitrocellulose Membrane Filter Sandwiches
ReliaBLOT® DTT Reducing Agent 10X WB106-500 WB107-20 **WB108** ReliaBLOT® Chemiluminescent HRP Substrate WB110 ReliaBLOT® IP/Western Blot Reagents WB120 Goat anti-Rabbit IgG-h+l HRP A120-101P Goat anti-Mouse IgG-h+l HRP

A90-116P

# ReliaBLOT® 12% SDS-PAGE Gels 10 X 10 cm, 12-well



Catalog No. WB101-01212G

#### MATERIAL SAFETY DATA SHEET

#### **Contents Description**

This product contains 3-20% Polyacrylamide, 1-5% Triethanolamine-hydrochloride, 0.01 to 0.2% SDS, <0.02% Acrylamide, and <0.002% N,N'-methylene bisacrylamide.

#### **Hazardous Ingredients**

<u>Polyacrylamide</u> may be harmful if inhaled. May Cause respiratory tract irritation. May be harmful if absorbed through the skin and may cause respiratory tract, skin, and eye irritation. May be harmful if swallowed. Potential cancer hazard. Wash thoroughly after handling; avoid contact with eyes, skin, clothing. Avoid ingestion.

<u>Triethanolamine-hydrochloride</u> has target organ effects and may be harmful if inhaled, absorbed through the skin, or swallowed. May cause respiratory tract, skin, and eye irritation.

#### **Disposal Guidelines**

Observe all Federal, State and Local laws concerning health and pollution.

The above information is believed to be correct but does not purport to be all-inclusive and is intended to be used only as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state, or local laws or regulations.