ReliaBLOT[®] 4 – 20% SDS-PAGE Gels 10 X 8 cm, 12-well



Catalog No. WB102-42012G		BETHYL LABORATORIES, INC	
Description:	Bethyl Laboratories ReliaBLOT [®] 4 -20% SDS-PAGE (sodiu polyacrylamide gel electrophoresis) gels are ready-to-use gra X 8(H) cm cassettes. The gels are cast with 12 sample wells, thickness.	dient gels cast in plastic 10(W)	
Amount:	10 gels		
Application:	Bethyl Laboratories ReliaBLOT [®] 4 -20% SDS-PAGE gels are used to resolve proteins under denaturing conditions. For optimal results, The ReliaBLOT [®] Running Buffer and The ReliaBLOT [®] Transfer Buffer should be used for electrophoresis and transfer of proteins in ReliaBLOT [®] SDS-PAGE Gels. Subsequent to electrophoresis, the gels can be transferred for western blot analysis or the gel may be stained and dried.		
Preparation:	Remove the gel from the package, drain away the buffer, and rinse the gel and wells with distilled water. Place the gel on the buffer core and assemble the electrophoresis cell according to the manufacturer's directions. Prepare a 1X solution of ReliaBLOT [®] Running Buffer for the inner and outer chambers. Run the gel(s) at 150 V until the dye front reaches the bottom of the gel (approximately 60 minutes).		
Storage:	Store at $2 - 8^{\circ}$ C.		
Stability:	Stable for 2 years from date of manufacturing.		
Expiration Date:	See gel cassette.		
Related Products:	ReliaBLOT [®] SDS-PAGE Gels, 4-8%, 10 x 10 cm ReliaBLOT [®] SDS-PAGE Gels, 4-12%, 10 x 10 cm ReliaBLOT [®] SDS-PAGE Gels, 12%, 10 x 10 cm ReliaBLOT [®] SDS-PAGE Gels, 4-20%, 10 x 10 cm ReliaBLOT [®] SDS-PAGE Gels, 4-8%, 10 x 8 cm ReliaBLOT [®] SDS-PAGE Gels, 4-12%, 10 x 8 cm ReliaBLOT [®] SDS-PAGE Gels, 12%, 10 x 8 cm ReliaBLOT [®] Prestained Protein Marker ReliaBLOT [®] LDS Sample Buffer (4X) ReliaBLOT [®] Transfer Buffer (10X) ReliaBLOT [®] Nitrocellulose Membrane Filter Sandwiches ReliaBLOT [®] DTT Reducing Agent 10X ReliaBLOT [®] IP/Western Blot Reagents Goat anti-Rabbit IgG-h+1 HRP Goat anti-Mouse IgG-h+1 HRP	WB101-40812G WB101-41212G WB101-01212G WB101-42012G WB102-40812G WB102-01212G WB102-01212G WB103-600 WB104-10 WB105-500 WB106-500 WB106-500 WB107-20 WB108 WB110 WB120 A120-101P A90-116P	
	25043 West FM 1097 • Montgomery TX 77356		

250+5 West 1 W 1077 • Wongomery, 174 77550			
800-338-9579 • 936-597-6111 • Fax 866-597-6105 • 936-597-6205			
www.bethyl.com • bethyl@bethyl.com			

ReliaBLOT[®] 4 – 20% SDS-PAGE Gels 10 X 8 cm, 12-well



Catalog No. WB102-42012G

MATERIAL SAFETY DATA SHEET

Contents Description

This product contains 3-20% Polyacrylamide, 1-5% Triethanolamine-hydrochloride, 0.01 to 0.2% SDS, <0.02% Acrylamide, and <0.002% N,N'-methylene bisacrylamide.

Hazardous Ingredients

<u>Polyacrylamide</u> may be harmful if inhaled. May Cause respiratory tract irritation. May be harmful if absorbed through the skin and may cause respiratory tract, skin, and eye irritation. May be harmful if swallowed. Potential cancer hazard. Wash thoroughly after handling; avoid contact with eyes, skin, clothing. Avoid ingestion.

<u>Triethanolamine-hydrochloride</u> has target organ effects and may be harmful if inhaled, absorbed through the skin, or swallowed. May cause respiratory tract, skin, and eye irritation.

Disposal Guidelines

Observe all Federal, State and Local laws concerning health and pollution.

The above information is believed to be correct but does not purport to be all-inclusive and is intended to be used only as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state, or local laws or regulations.