

12. Short protocol.

A. Cloning and transformation of donor strain.

1. Ligate the DNA fragment with pUTmini-Tn5 vector, prior digestion with *Not* I.
2. Transform the ligation reactions into host strain λ pir cells and select the transformants by plating onto selective LB plates. Incubate overnight at 37°C.

B. Transfer of the DNA into the recipient strain (triparental mating).

Filter system

1. Mix the 10-50 μ l of overnight cultures of the donor, recipient and helper strains (1:1:1 mixtures of cultures overnight), in 5 ml of 10 mM MgSO₄, vortex for a few seconds, transfer to a 5 ml disposable syringe, and filter through a Millipore membrane (or equivalent).
2. Remove carefully the drained membrane on the agar surface of an LB plate (cell side up), avoiding bubbles formation.
3. Incubate the plates at 30-37 °C for 8-18 hours.

Drop system

1. Mix the 10-50 μ l of overnight cultures of the donor, recipient and helper strains (1:1:1 mixtures of cultures overnight), centrifugate the mixture and discard the supernatant. Resuspend the pellet in 20 μ l of LB medium and place on a LB plate. Let to dry the drop on the plate.
2. Incubate the plates at 30-37 °C for 8-18 hours.

C. Selection of exconjugants.

1. To filter system: Resuspend the filter with the mating mix in 5 ml of 10 mM MgSO₄.
To drop system: Recover the mating mix on the plate and resuspend in 5 ml of 10 mM MgSO₄.
2. Plate 100-500 μ l of this suspension on selective LB agar medium containing tetracycline.
3. Incubate the plates at the optimal growth temperature of the recipient strain until colonies become visible.
4. Confirm that the acquisition of the selected phenotype is due to an authentic transposition event and not to integration of the delivery plasmid in a recipient replicon or its illegitimate replication in the new host, by Southern blot, PCR analysis or growing the exconjugants onto medium containing β -lactam antibiotic.

PANC-1 plus mini-Tn5 package

A collection of Delivery System pUT mini-Tn5 plasmids carrying antibiotic selection markers and donor and helper strains.

For research use only

Contents

1. Product Description	2
2. Vectors and strains	3
2.1. Vectors	3
2.1.1. pUT mini-Tn5 vectors	3
2.1.2. Auxiliary vectors	4
2.2. Strains	5
2.2.1. Donor strains	5
2.2.2. Helper strains	6
2.2.3. Strain revival and long term storage	6
3. About mini-Tn5 System	7
4. Vectors Maps	10
5. Restriction analysis	13
5.1. pUTmini-Tn5 Cm Restriction analysis	13
5.2. pUTmini-Tn5 Km Restriction analysis	13
5.3. pUTmini-Tn5 Sm/Sp Restriction analysis	14
5.4. pUTmini-Tn5 Tc Restriction analysis	14
5.5. pUC18 <i>Not</i> Restriction Enzymes	15
5.6. pUC18 <i>Sfi</i> Restriction analysis	16
6. Protocol	16
6.1. Cloning and transformation of donor strain	16
6.2. Transfer of the DNA into the recipient strain	17
6.3. Selection of exconjugants	19
7. Troubleshooting	21
8. Reagents and recipes	22
9. Appendix: Transformation using TSS	24
9.1. TSS description	24
9.2. Protocol	24
10. References	26
11. Related products	27
12. Short protocol	28

1. Product Description.

Mini-Tn5 system is the generic name of a system derived of mini-Tn5 transposon vectors, that can be used in the analysis, construction and manipulation of complex phenotypes in a wide range of gram-negative bacteria. This system includes cloning vectors (pUTmini-Tn5) and strains for mating experiments (donor and helper strains). The main feature of the cloning vectors of this system is their ease for insertion of one or more segments of heterologous DNA in the chromosome of the strain of interest, by means of a mating between the donor strain and the recipient strain. This system can be used for insertional mutagenesis, promoter probing and analysis of transcriptional terminators.

KIT COMPONENTS (CAT. NO. KT- 3261)

PRODUCT	QUANTITY	STORAGE	CAT. NO.
pUTmini-Tn5 Cm	8 µg	4°C or -20°C.	CV-3223
pUTmini-Tn5 Km	8 µg	4°C or -20°C.	CV-3224
pUTmini-Tn5 Sm/Sp	8 µg	4°C or -20°C.	CV-3225
pUTmini-Tn5 Tc	8 µg	4°C or -20°C.	CV-3226
pUC18Not	8 µg	4°C or -20°C.	CV-3282
pUC18Sfi	8 µg	4°C or -20°C.	CV-3283
<i>E. coli</i> DH5α λpir	Stab	RT	BS-3233
<i>E. coli</i> S17-1 λpir	Stab	RT	BS-3234
<i>E. coli</i> DH5α(pRK2013)	Stab	RT	BS-3236
<i>E. coli</i> DH5α(pRK2073)	Stab	RT	BS-3263
TSS solution	1.5 ml	4 °C	RS-3215
Sequences CD	1 CD	-----	

! NOTE:

- All components of the Kit are shipped at room temperature. Upon arrival, store the kit components according to the directions in the table above.
- The plasmid is supplied dried, the strains in stabs and TSS as solution.
- CD-Rom supplied with the kit contains vectors sequences.

11. Related products.

PRODUCT	DESCRIPTION	Cat. No.		
STRAINS	<i>E. coli</i> DH5α λpir	Propagation of pUT plasmids	BS-3233	
	<i>E. coli</i> CC118 λpir	Propagation of pUT plasmids	BS-3235	
	<i>E. coli</i> S17-1 λpir	Propagation of pUT plasmids	BS-3234	
	<i>E. coli</i> DH5α (pRK2013)	Mating helper strain	BS-3236	
	<i>E. coli</i> DH5α (pRK2073)	Mating helper strain	BS-3263	
PRIMERS	pUToriR6K	Amplification by PCR (selection of positive clones) and sequencing of cloned fragments	PR-3280	
VECTORS	pUT mini-Tn5 Cm	Mini-transposons derivatives vectors (repeat insertion events)	CV-3223	
	pUT mini-Tn5 Km	Mini-transposons derivatives vectors (repeat insertion events)	CV-3224	
	pUT mini-Tn5 Sm/Sp	Mini-transposons derivatives vectors (repeat insertion events)	CV-3225	
	pUT mini-Tn5 Tc	Mini-transposons derivatives vectors (repeat insertion events)	CV-3226	
	pUT mini-Tn5 Tel	Mini-transposons derivatives vectors (repeat insertion events)	CV-3227	
	pUT/lacZ1	Mini-transposons derivatives vectors with a reporter gene	CV-3256	
	pUT/lacZ2	Mini-transposons derivatives vectors with a reporter gene	CV-3257	
	pUT/phoA	Mini-transposons derivatives vectors with a reporter gene	CV-3258	
	pUT/luxAB	Mini-transposons derivatives vectors with a reporter gene	CV-3259	
	pUC18Not	Auxiliary plasmid for cloning	CV-3282	
	pUC18NotSfi	Auxiliary plasmid for cloning	CV-3248	
	ANTIBIOTICS	Apicillin	Test culture	RS-3217
		Chloramphenicol	Test culture	RS-3218
Kanamycin		Test culture	RS-3219	
Streptomycin		Test culture	RS-3221	
Tetracycline		Test culture	RS-3220	
REAGENTS and SOLUTIONS	Potassium Tellurite	Selection of transformants	RS-3222	
	TSS	Competent cells preparation	RS-3215/ RS-3216	

10. References.

1. **De Lorenzo V., Herrero M., Jakubzik U., Timmis K.** Mini-Tn5 Transposon Derivatives for Insertion Mutagenesis, Promoter Probing, and Chromosomal Insertion of Cloned DNA in Gram-Negative Eubacteria. *Journal Bacteriology*, 172, p. 6568-6572.
2. **Prentki P., Binda A., Epstein A.** Plasmid vectors for selecting IS1-promoted deletions in cloned DNA: sequence analysis of the omega interposon. *Gene*, 103, p 17-23.
3. **De Lorenzo V., Timmis K.N.** Analysis and Construction of Stables Phenotypes in Gram-Negative Bacteria with Tn5- and Tn10-Derived Minitransposons. *Methods in Enzimology*, 235, p. 386-405.
4. **Herrero M., De Lorenzo V., Timmis K.N.** Transposon vectors Containing Non-Antibiotic Resistanse Selection Markers for Cloning and Stable Chromosomal Insertion of Foreign Genes in Gram-Negative Bacteria. *Journal Bacteriology*, 172, p.6557-6567.
5. **De Lorenzo V., Herrero M., Sánchez J., Timmis K.N.** Mini-transposons in microbial ecology and environmental biotechnology. *FEMS Microbiology Ecology*, 27, p. 211-224.
6. **Boyd D., Weiss D.S., Chen J.C., Beckwith J.** Towards single-Copy gene expression systems making gene cloning physiologically relevant: Lamda Inch, a simple Escherichia coli plasmid-chromosome shuttle system. *Journal Bacteriology*, 182, p. 842-847.

26

2. Vectors and strains.

2.1. Vectors.

2.1.1. pUT mini-Tn5 vectors

All pUTmini-Tn5 vectors included in this kit contain a mini-Tn5 element, carrying an antibiotic resistance cassettes, flanked by I and O termini (terminal repeat sequences) of Tn5. The cassettes conferring resistance to tetracycline, streptomycin-spectinomycin and chloramphenicol are derived from Ω interposon and the cassette conferring resistance to kanamycin contains the kanamycin resistance interposon originally isolated from Tn5 (1). These four mini-Tn5 carry 156-nt inverted repeats consisting of the T4 phage gene 32 transcriptional terminators (strong transcriptional terminator) flanking the antibiotic resistance gene as well as synthetic translational stop codons in all three reading frames (1,2).

Selection for resistance to	Antibiotic selection marker
Replication Origin	R6K
Transfer Origin	RP4oriT
Antibiotic Resistance	Ampicillin, selection marker
Host strain	<i>E. coli</i> λ pir strains
Copy number	Low

pUTmini-Tn5 Cm.

pUTmini-Tn5 Cm is a 8686 basepair plasmid, that contains the minitransposon mini-Tn5 Cm cloned in the pUT backbone. Mini-Tn5 Cm carries as selection marker the *cat* gene, that confers chloramphenicol resistance.

pUTmini-Tn5 Km.

pUTmini-Tn5 Km is a 7575 basepair plasmid, that contains the minitransposon mini-Tn5 Km cloned in the pUT backbone. Mini-Tn5 Km carries as selection marker a kanamycin resistance gene. Mini-Tn5 Km has one more I end at the left extreme of the kanamycin resistance gene, that itself originates from Tn5.

pUTmini-Tn5 Sm/Sp.

pUTmini-Tn5 Sm/Sp is a 7386 basepair plasmid, that contains the minitransposon mini-Tn5 Sm/Sp cloned in the pUT backbone. Mini-Tn5 Sm/Sp carries as selection marker the *add* gene, that confers streptomycin and spectinomycin resistance.

pUTmini-Tn5 Tc.

pUTmini-Tn5 Tc is a 7477 basepair plasmid, that contains the minitransposon mini-Tn5 Tc cloned in the pUT backbone. Mini-Tn5 Tc carries as selection marker the *tet* gene, that confers tetracycline resistance.

2.1.2. Auxiliary vectors.

pUC18Not and pUC18Sfi are plasmids derived from pUC18 vector, so it can be used like cloning vector in any *E. coli* strain and exhibit the blue/white screening property in conjunction with X-Gal containing plates. These vectors are an excellent auxiliary plasmids to clone fragments which will be introduced like fragments flanked by Not I and Sfi I restriction sites, in other vectors with these restriction sites.

pUC18Not

pUC18Not is a 2715 basepair plasmid that carries pUC18 multiple cloning site polylinker (in the same orientation as M13mp18), flanked by Not I cleavage sequence. This vector is an excellent auxiliary plasmid to clone fragments, that will be introduced as Not fragments in other vectors with this restriction site.

pUC18Sfi

pUC18Sfi is a 2724 basepair plasmid that carries pUC18 multiple cloning site polylinker (in the same orientation as M13mp18), flanked by Sfi I cleavage sequence. This vector is an excellent auxiliary plasmid to clone fragments, that will be introduced like Sfi fragments in other vectors with this restriction sites. Note that because of the degeneracy of the Sfi I cleavage sequence, the Sfi I fragments can only be cloned in vector designed to be compatible with pUC18Sfi. For the same reason, pUC18Sfi - derived fragments can be cloned in only one orientation in others Sfi-vectors.

! NOTE:

- Alternatively, you can use exponential growing cells, obtained in the same day of the transformation directly from fresh colonies.

! IMPORTANT! While the culture is growing place on ice or at 4°C LB medium that you will use in step 6.

4. When the culture reaches the desired OD₆₀₀, divide the cultures in microcentrifuge tubes (1 ml/tube) and place on ice for 5 minutes.
5. Centrifuge at 12000 r.p.m. for 30 seconds at room temperature.
6. Discard the medium and resuspend gently the pellets in 75 µl of cold LB medium. Place on ice for 5 minutes.
7. Add 75 µl of 2xTSS solution and mix gently. Replace on ice for 5 minutes.

! IMPORTANT: Invert the tube of TSS solution several times before use.

Now the cells are ready to perform the transformation!

TRANSFORMATION PROCEDURE (STANDARD PROTOCOL)

8. Add 1-50 ng of DNA (supercoiled DNA) to the competent cells obtained in step 7. We recommend adding 5-10 µl of DNA; the DNA added do not must be more than 15 µl of volume (10% of the competent cell volume). Mix gently and place the transformation reaction on ice for 30 minutes.
9. Heat shock: Incubate the tubes by immersing in a 42°C water bath for exactly 45 seconds. Immediately replaces the tubes on ice.
10. Add 1 ml of LB medium to each tube. Incubate at 37°C for the time required to the resistance antibiotic expression (typical 60 min).
11. For each transformation reactions, we recommend diluting the cells 1:10 and plating out 100 µl of the 1:10 dilutions, 100 µl of the transformation (undiluted cells), and the rest of the cells concentrated for centrifugation, onto LB plates containing the appropriate selection antibiotics. Allow the plates to completely absorb any excess media.
12. Incubate the plates overnight at 37°C (12-14 hours).

9. Appendix: Transformation using TSS.

9.1. TSS description.

TSS is a solution used to prepare competent cells ready for plasmid transformation, using an easy and fast procedure in a single step. You only need an early exponential growing cell. The TSS solution is mixed directly with an equal volume of a culture cells in a tube, in which is also added the plasmid to the transformation. This means no extra pipeting steps so you will get consistent high-efficiency transformations results. For cloning experiments, in which DNA ligation is directly used to transformate, the transformation efficiency could be low.

It is a very convenient system to perform transformations of different strains simultaneously. If you use a collection of different strains in your experimentation, using this system you will not have to repeat long protocols to perform competent cells of each strain. Besides, competent cells stocks are not necessary because you can get competent cells at the moment of the transformation.

This system is very fast, because you only need exponential growing cells, which can be obtained in the same day of the transformation directly from fresh colonies. Alternatively, overnight cultures in rich media can be diluted into fresh liquid media and grow up to early-log phase ($OD_{600}=0.2-0.4$). The competent cells performed using TSS solutions are competent for one day but cannot be stored.

9.2. Protocol.

COMPETENT CELLS PERFORMANCE

1. Streak the strain that you want to transformate onto a LB plate containing the appropriate antibiotic. Grow the plate overnight at 37°C (12-16 hours).
2. Pre-culture: Inoculate a single colony from LB plate, in 3 ml of LB medium, including the corresponding antibiotic. Grow the culture overnight (12-16 hours) at 37°C with shaking (200-225 r.p.m.).
3. Culture: Inoculate 20 ml of LB medium with 200 µl of the pre-culture (in a 100 ml flask). Grow the culture at 37°C with shaking (200-225 r.p.m.) until the OD_{600} reaches 0.2-0.4.

! NOTE: pUC18*Not* and pUC18*Sfi* carry a single mutation that eliminates a C and changes the position of the start codon of *lacZ* gene. This change eliminates the first and second aminoacids of β -galactosidase without changing the structure and function of the *lacZ alpha* product.

2.2. Strains.

2.2.1. Donor strains

E. coli S17-1 λ pir

E. coli S17-1 λ pir strain contains the *pir* gene (λ pir) and is designed for cloning and propagation of plasmids with R6K origin of replication sequence (such as pUTmini-Tn5 vectors). The plasmids containing the R6K origin of replication require the R6K-specified replication protein π (*pir* gene) for propagation. *E. coli* S17-1 λ pir strain has chromosomally integrated conjugal transfer functions (RP4 transfer functions), so when it is used as specific host strain into which the transposon vector DNA is transformed, the transfer occurs by means of a biparental mating (a helper strain is not necessary).

Genotype

$Tp^R S^R recA, thi, pro, hsdR^-M^+RP4: 2-Tc:Mu: Km Tn7 \lambda pir$

E. coli DH5 α λ pir

E. coli DH5 α λ pir strain contains the *pir* gene (lysogenized with λ pir phage) and is designed for cloning and propagation of plasmids with R6K origin of replication sequence (such as pUTmini-Tn5 vectors). The plasmids containing the R6K origin of replication require the R6K-specified replication protein π (*pir* gene) for propagation. To transfer this plasmids directly from DH5 α λ pir donor strain to the recipient strain, triparental mating with a helper strain, such as *E. coli* DH5 α (pRK2013) or *E. coli* DH5 α (pRK2073) is necessary.

Genotype

$supE44, \Delta lacU169(F80lacZ\Delta M15), recA1, endA1, hsdR17, thi-1, gyrA96, relA1, \lambda pir$ phage lysogen

2.2.2. Helper strains.

E. coli DH5 α (pRK2013)

E. coli DH5 α (pRK2013) strain contains the helper plasmid pRK2013, that provides *tra* and *mob* genes that encode for mobilisation and chromosomal transfer. Both of genes are required to transfer the DNA cloned in pUTmini-Tn5 plasmids into the recipient strain, from donor strains when lack of the chromosomally integrated conjugal transfer functions.

Genotype

supE44, Δ *lacU169*(F80*lacZ* Δ M15), *recA1*, *endA1*, *hsdR17*, *thi-1*, *gyrA96*, *relA1*, pRK2013 (Km^R oriColE1 RK2-Mob⁺ RK2-Tra⁺)}

E. coli DH5 α (pRK2073)

E. coli DH5 α (pRK2073) strain contains the helper plasmid pRK2073, that provides *tra* and *mob* genes that encode for mobilisation and chromosomal transfer. Both of genes are required to transfer the DNA cloned in pUTmini-Tn5 plasmids into the recipient strain, from donor strains when lack of the chromosomally integrated conjugal transfer functions.

Genotype

supE44, Δ *lacU169*(F80*lacZ* Δ M15), *recA1*, *endA1*, *hsdR17*, *thi-1*, *gyrA96*, *relA1*, pRK2073 (Sm^{RS} oriColE1 RK2-Mob⁺ RK2-Tra⁺)}

2.2.3. Strain revival and long term storage.

6 Bacterial strains in this kit are delivered as stab cultures, which can be stored for up to 15 days at 4°C. We strongly recommend immediate revival of the strains upon arrival, and preparation of frozen stocks for long term storage. To this end, insert a sterile loop or toothpick into the stab agar, smear a small amount of the culture onto the surface of an LB plate without any antibiotic (in the case of strains S17-1 λ pir and DH5 α λ pir), LB supplemented with 25 mg/l kanamycin (strain DH5 α [pRK2013]) or LB supplemented with 50 mg/l streptomycin (strain DH5 α [pRK2073]), and streak for single colonies. Incubate overnight at 37°C.

For long term strain storage, inoculate 5 ml of liquid LB medium (alone or supplemented with the corresponding antibiotic) with a single colony and incubate overnight at 37 °C with shaking to obtain a saturated culture. Mix this culture with an equal volume of a sterile 30% glycerol solution, aliquot the mix into cold-resistant plastic vials and store at -70 °C. Strain revival from frozen glycerol stocks can be done by scraping off splinters of solid ice with a sterile toothpick or pipette and streaking these splinters onto an LB plate (alone or supplemented with the corresponding antibiotic). Do not allow the contents of the vial to thaw.

Chloramphenicol 25 mg/ml stock solution

Dissolve 125 mg chloramphenicol in 5 ml of ethanol.
Store the stock solution at 4°C.

Biomedal Cat. No. RS-3218

Kanamycin 25 mg/ml stock solution

Dissolve 0.25 g of kanamycin in 9 ml of deionized water.
Bring up to a 10 ml final volume with water deionized.
Filter sterilizes the solution with a 0.22 μ filter.
Store the stock solution at -20°C.

Biomedal Cat. No. RS-3219

Streptomycin 10 mg/ml stock solution

Dissolve 0.1 g of streptomycin in 9 ml of deionised water.
Bring up to a 10 ml final volume with water deionised.
Filter sterilizes the solution with a 0.22 μ filter.
Store the stock solution at -20°C.

Biomedal Cat. No. RS-3221

Tetracycline 25 mg/ml stock solution

Dissolve 200 mg tetracycline in 5 ml water deionised/ 5 ml ethanol mix.
Store the stock solution at -20 °C.

Biomedal Cat. No. RS-3220

8. Reagents and recipes.

A. Media and solutions.

1. MEDIA.

LB (per liter) (liquid)

- 10 g tryptone
- 5 g yeast extract
- 10 g NaCl

Dissolve in 950 ml of deionised water.
Adjust pH to 7.0 with 5 N NaOH.
Add deionised water to 1 l final volume.
Autoclave and let cool to below 55°C.

LB agar (plates)

Exactly like liquid LB but adding 15 g agar after adjusting the pH.

22

2. MATING SOLUTION.

Magnesium phosphate ($MgSO_4 \cdot 7H_2O$) 1M solution

Dissolve 24.6 g of $MgSO_4 \cdot 7H_2O$ in 100 ml of deionised water.
Autoclave sterilizes.

3. ANTIBIOTICS AND ADDITIVES.

Ampicillin 100 mg/ml stock solution

Dissolve 1 g ampicillin in 9 ml of deionised water.
Add deionised water to a 10 ml volume final.
Filter sterilizes the solution with a 0.22 μ filter.
Store the stock solution at -20°C.

Biomedal Cat. No. RS-3217

3. About mini-Tn5 System

Mini-transposons derivatives vectors provide a straightforward tool to clone and insert foreign genes stably into the chromosomes of a variety of Gram negative bacteria such *E. coli*, *Klebsiella*, *Salmonella*, *Proteus*, *Vibrio*, *Bordetella*, *Brucella*, *Actinobacillus*, *Rhizobium*, *Acinetobacter*, *Rhodobacter*, *Agrobacterium*, *Alcaligenes* and several *Pseudomonads* (3). These vectors also simplify the generation of insertion mutants (insertion mutagenesis) and in vivo fusions with reporter genes (promoter probing). Mini-transposons cloning vectors have the advantages of natural transposons, but lack their disadvantages.

All pUTmini-Tn5 plasmids carry a backbone common (pUT backbone) and a mobile element variable (mini-Tn5). The mobile element (mini-Tn5 elements) consists of an *Sfi* I cassette containing a selection marker and a single *Not* I site outside of the cassette that can be used for cloning foreign DNA fragments. Flanking these two features are the 19 basepair I and O termini of Tn5 transposon. In this vectors any heterologous DNA segment can be cloned within the boundaries of a mini-Tn5 element and finally inserted into the chromosome of Gram negative bacteria target. The modular nature of the mini-transposons facilitates the construction of mobile elements *à la carte*.

Hybrid transposons bearing the DNA segment of interest can be easily assembled with the help of the specialised cloning vector as pUC18*Not* or pUC18*Sfi* vector.

The delivery system employed for all mini-Tn5 transposons is the pUT plasmid. The backbone pUT includes the R6K origin of replication sequence, the RP4oriT origin of transfer sequence, the *bla* gene sequence (ampicillin resistance) and the *tnp** gene sequence. The *tnp** is a mutant *tnp* gene of IS50_R (Tn5 transposon), that lacks *Not* I sites and which encodes for the transposase needed for transposition of the mini-Tn5 elements. *tnp** carries a single mismatch that changed the GCG codon specifying Ala-168 of the *tnp* gene to the Ala codon GCC. This change eliminate the *Not* I site without changing the structure of the *tnp* product (4).

The *tnp** gene is presented in cis but external to the mobile element. Due to the loss of the *tnp** gene after insertion, minitransposons are stably inherited and do not favour DNA rearrangements or others forms of genetic instability.

Plasmids having the R6K origin of replication require the R6K-specified replication protein π and can be maintained only in host strains producing this protein. It also carry the origin of transfer $oriT$ of plasmid RP4, which results in its efficient conjugal transfer to recipient strain from donor strains expressing RP4-conjugative functions (5). Delivery plasmids are thus maintained stably in λ_{pir} lysogens or in *E. coli* strains with the *pir* gene recombined in their chromosome, and can be mobilised into target strains through RP4 transfer functions. However, once transferred it is unable to replicate in recipients that lack the π protein. The pUT plasmids are “suicide vectors” since the $oriR6K$ is functional only in the presence of the *pir* gene and they cannot replicate in the recipient cells which have not the π protein. Delivery of the donor plasmids into selected host bacteria is accomplished through mating with the target strain.

One important feature of mini-Tn5 elements is that, as a result of the loss of the transposase-cognate inhibitor along the pUT system after transposition (transposase function is not maintained in target cells), a single recipient strain can be used for repeated insertion events with differentially marked minitransposons. Multiple insertions in the same strain are therefore only limited by the availability of distinct selection markers.

7. Troubleshooting.

PROBLEM	WHY?	SOLUTIONS
Loss of pUT plasmids.	Plasmid derivatives cannot be maintained into delivery λ_{pir} strains	Using of alternative donor strains.
Low or zero yield of exconjugants.	Selection conditions are inappropriate. Zero or low frequency of RP4-mediated transfer.	It is advisable to make and assay with the target cells before the matting/plating is carried out, because the minimal inhibitory concentration of the toxic compound used as marker can vary considerably among genera and strains. Properties of target cells such as restriction system, surface exclusion phenomena may lead to low frequency or zero transfer. If the target strain is unable to act as a recipient for RP4, the pUT plasmid cannot be introduced at high frequency. Also, some changes in the mating protocol can be introduced. For example, avoid the pregrowth of the recipient strains on antibiotic medium prior to mating, since traces of the antibiotic present in the suspension or accumulated by the recipient cells can kill donor cells.
All exconjugants are resistant to ampicillin or others β -lactam antibiotics	The exconjugant phenotype is not due to an authentic transposition events	Shorter times. Lower temperatures (30°C). To <i>E. coli</i> recipient strains use LB medium with 0.1 M citrate to avoid the recipient strain infection by λ_{pir} phage.

2. Plate 100-500 µl of this suspension on selective LB agar medium with tetracycline. The tetracycline concentration must be adequate to ensure maintenance of the delivery plasmid, but you must determine the minimal inhibitory concentration of antibiotic.

! NOTE: When the marker carried by the minitransposon is an antibiotic resistance determinant, optimal concentrations of antibiotics to the selection can vary among strains. Must be determined the minimal inhibitory concentration (MIC) in each case. Can be orientatives the concentrations which appear in the next table.

ANTIBIOTIC	ORIENTATIVE MIC
Chloramphenicol	5 - 50 µg/ml
Kanamycin	25 - 75 µg/ml
Streptomycin	50 -100 µg/ml
Spectinomycin	50 -100 µg/ml
Tetracycline	2.5 - 15 µg/ml

20

3. Incubate the plates at the optimal growth temperature of the recipient strain until colonies become visible.

4. Confirm that the acquisition of the selected phenotype is due to an authentic transposition event and not to integration of the delivery plasmid in a recipient replicon or its illegitimate replication in the new host. The insertion can be confirmed by Southern blot or PCR analysis. Other procedure consists to grow the exconjugants onto medium containing β-lactam antibiotic, because authentic transposition results in the loss of the portion of the delivery plasmid containing *bla* gene and the exconjugants will be sensitive to ampicillin and other β-lactams (see Troubleshooting)


ADVANTAGES of mini-Tn5 System vs OTHER DNA INTEGRATION INTO THE CHROMOSOME SYSTEMS

- Ease for insertion of one or more segments of the heterologous DNA in the chromosome of the strain of interest (any heterologous DNA segments can be inserted into the chromosome of target cells after a few simple genetic manipulations).
- Conditional replication to allow selection for integration into the chromosome ("suicide plasmids"). The mini-Tn5 would be able transpose from the plasmid to the chromosome but when successive cell divisions occur, the plasmid vector would be lost to the population.
- Insert foreign DNA stably into the chromosomes.
- Wide range of inserted fragments size. pUTmini-Tn5 plasmids allow to clone from small fragments until large size fragments. Heterologous DNA of 12 Kb cloned in mini-Tn5 derivatives are transposed at frequencies in the range of those observed with insert-lacking transposon (3). In other system the range of sizes is more restricted. For example, lambda Inch vectors system allow to insert fragments a total of about 7 Kb (6).
- Inserted DNA segments remain stably inherited and produce no burden to the carrier strain. Possibility of multiple insertions in the same strains (only limited by the available selection markers), as a result of the loss of the transposase function, which is not maintained in target cells.
- Easy modification of the cloning vector, because of the modular nature of the constructions.
- Use of cleavage sites of the rare cutters *Not* I and *Sfi* I (it's convenient that the vector has an array of unique cloning sites).
- The pUTmini-Tn5 plasmids are transferable to a variety of bacteria.

4. Vectors Maps.

pUTmini-Tn5 Cm


mini-Tn5 Cm: 1-3467
I termini: 1-19
O termini: 3449-3467
T4 gene transcription terminators (inverted repeats): 64-188/ 3259-3383
cat II coding region (Cm resistance): 1452-2090
R6K origin of replication: 3560-3940
mobRP4 origin of transfer: 3944-5670
bla coding region (ampicillin resistance): 6130-6987
*tnp** coding region (transposase): 7139-8581



pUTmini-Tn5 Km

SEQUENCE REFERENCE POINT

miniTn5 Km: 1-2356
I termini: 1-19 / 831-849 (reverse complement)
O termini: 2338-2356
km coding region (kanamycin resistance): 857-1648
T4 gene 32 transcription terminators: 64-188/ 2148-2272 (inverted repeats)
R6K origin of replication: 2448-2829
mobRP4 origin of transfer: 2833-4559
bla coding region (ampicillin resistance): 5018-5876
*tnp** coding region (transposase): 6031-7458



3. After the filtering, remove carefully the drained membrane on the agar surface of an LB plate (cell side up). You can help yourself with sterile tweezers, preferably with curved tips.

! NOTE: Air bubbles should be avoided between the filter and agar surface.

4. Incubate the plates at 30-37 °C for 8-18 hours.

! NOTE: Several alternatives to mating protocol are possible, according to individual needs. Strains ratios (donor/recipient/helper), temperatures, time of mating and culture medium can be changed to suit specific requirements of the recipients.

Biparental mating

If you use *E. coli* SM10 λ pir or *E. coli* SM17-1 λ pir like donor strains, to set up the biparental mating, you must follow the same protocol that when use triparental mating, but without helper strain. You must mix equal volumes (5-10 μ l) of overnight cultures of the donor and recipient strain, in 5 ml of 10 mM MgSO₄, and proceed as described above.

6.3. Selection of exconjugants.

The selection of clones with the minitransposon inserted in the chromosome is a critical step. The procedure is as follows.

1. If you use filter system:

Resuspend the filter with the mating mix in 5 ml of 10 mM MgSO₄.

If you use drop system:

Recover the matting mix on the plate using a sterile inoculating loop and resuspend in 5 ml of 10 mM MgSO₄.

! NOTE: This suspension can be kept at 4°C for several weeks.

Triparental Mating

1. Grow overnight cultures in LB of each of the three strains: donor strain, recipient strain and helper strain:
 - a. Inoculate a single colony of donor strain DH5 α λ pir carrying the desired pUT derivative plasmid in LB medium containing 100 μ g/ml ampicillin and the adequate concentration of the antibiotic of selection to ensure maintenance of the delivery plasmid. Incubate with shaking at 37°C overnight.
 - b. Grow the recipient strain: Inoculate a single colony under the same conditions (at different temperature, if necessary), but preferably without selection (see Troubleshooting).
 - c. Inoculate a single colony of *E. coli* DH5 α (pRK2013) or *E. coli* DH5 α pRK2073) (if you use pUTmini-Tn5 Km) in LB medium, preferably without selection.
2. Mix the 10-50 μ l (~10⁸ cells) of overnight cultures of each of the three strains (1:1:1 mixtures of cultures overnight), in 5 ml of 10 mM MgSO₄, vortex for a few seconds, transfer to a 5 ml disposable syringe, and filter through a Millipore membrane, 13-25 mm diameter type HA (0.45 μ m) or equivalent, placed on a reusable filter case.

18


! NOTE: To remove antibiotics from donor strain culture, an additional washing step of the culture is necessary.

! IMPORTANT! Alternatively, mix the 10-50 μ l of overnight cultures of helper strain, donor strain and recipient strain, centrifugate the mix and discard the supernatant. Resuspend the pellet in 20 μ l of LB medium or 10 mM MgSO₄ using vortex and place it on LB plate, without dispense it. Let it dry. Go to step 4.

! NOTE:

- You can use alternative culture medium instead of LB, to favour recipient strain viability.
- When you use the filter system, if too many donor, recipient or helper cells are used and the membrane clogs, it is better to refilter a more diluted mix.
- If you need only few insertions, 10 μ l drops of cultures of donor, recipient and helper strains can be mixed and spotted on an LB plate, which is then dried and incubated for several hours before the mating mixture is streaked out on selective medium (go to step 4).


pUTmini-Tn5 Sm/Sp



SEQUENCE REFERENCE POINT

mini-Tn5 Sm/Sp: 1-2167
I termini: 1-19
O termini: 2149-2167
T4 gene transcription terminators (inverted repeats): 64-188/ 1959-2083
aad coding region (Sp/Sm resistance): 443-1450
R6K origin of replication: 2260-2640
mobRP4 origin of transfer: 2644-4370
bla coding region (ampicillin resistance): 4830-5687
tnp* coding region (transposase): 5842-7269

pUTmini-Tn5 Tc




SEQUENCE REFERENCE POINT

mini-Tn5 Tc: 1-2258
I termini: 1-19
O termini: 2239-2258
T4 gene transcription terminators (inverted repeats): 64-188/ 2050-2174
tet coding region (Tc resistance): 489-1676
R6K origin of replication: 2351-2731
mobRP4 origin of transfer: 2735-4461
bla coding region (ampicillin resistance): 4921-5768
tnp* coding region (transposase): 5932-7372

11

pUC18Not



Multiple Cloning Site: 406-480

lacZ: 149-405/481-493 (complementary strand)

lacZ start codon: 493 (complementary strand)

lacZ promoter: 499-572 (complementary strand)

ColE1 origin (pUC origin): 837-1510 (complementary strand)

bla coding region (ampicillin resistance): 1658-2515 (complementary strand)

bla promoter: 2516-2614 (complementary strand)

Binding site of pUC/M13 Forward Sequencing Primer (-20): 379-395

Binding site of pUC/M13 Forward Sequencing Primer (-40): 359-375

Binding site of pUC/M13 Forward Sequencing Primer (-47): 352-375

Binding site of pUC/M13 Reverse Sequencing Primer (-48): 506-529 (complementary strand)

1. Excise the DNA fragment from pUC18Not plasmid by means of *Not* I digest.
2. Digest the pUTmini-Tn5 vector (transfer vector) with *Not* I. We recommend treating the digested plasmid with alkaline phosphatase, to prevent re-ligation.
3. Ligate your DNA fragment and pUTmini-Tn5 backbone.
4. Transform the ligation reactions into λ pir cells. Select the transformants by plating onto LB /ampicillin (100 μ g/ml) / antibiotic that the pUTmini-Tn5 confers resistance (the minimal inhibitory concentration in each case) plates. Incubate at 37°C overnight.

! NOTE: You can use the TSS solution, to perform the transformation.

6.2. Transfer of the DNA into the recipient strain.


To transfer the DNA cloned in pUTmini-Tn5 plasmid into the recipient strain, the genes that code for mobilisation and chromosomal transfer (*tra* and *mob* genes) are necessary. These genes can be present in a host strain or can be provided by way of a helper plasmid. When it is used the chromosomally integrated conjugal transfer functions of a specific host strain, such as *E. coli* SM10 λ pir or *E. coli* S17-1 λ pir, into which the transposon vector DNA is transformed, the transfer is accomplished by means of biparental mating. The transfer can also be achieved by triparental mating mediated by a helper plasmid.

In the triparental mating the conjugatable suicide minitransposon donor is introduced into the recipient strain by a helper plasmid, which self mobilises from its own host into the donor strain and provides the genes necessary for the conjugal transfer of the miniTn5 from the donor host strain to the recipient, provided the vector plasmids contain the specific recognition site for mobilisation.

Since some λ pir RP4 strains like *E. coli* SM10 λ pir or *E. coli* S17-1 λ pir, cannot stably maintain the pUT plasmids, we recommend other strains to propagate this vectors without problems, such as DH5 α λ pir. In this case, to mobilize the delivery plasmid directly from the donor strain into target cells, a triparental mating with a helper strain such as *E. coli* DH5 α (pRK2013) or *E. coli* DH5 α (pRK2073) is necessary.

! NOTE: To propagate pUTmini-Tn5 in λ pir strains, you can use the TSS solution to perform the transformation.

12 pUC18Sfi



Multiple Cloning Site: 406-489

lacZ: 149-405/491-502 (complementary strand)

lacZ start codon: 502 (complementary strand)

lacZ promoter: 508-581 (complementary strand)

ColE1 origin (pUC origin): 846-1519 (complementary strand)

bla coding region (ampicillin resistance): 1667-2525 (complementary strand)

bla promoter: 2525-2623 (complementary strand)

Binding site of pUC/M13 Forward Sequencing Primer (-20): 379-395

Binding site of pUC/M13 Forward Sequencing Primer (-40): 359-375

Binding site of pUC/M13 Forward Sequencing Primer (-47): 352-375

Binding site of pUC/M13 Reverse Sequencing Primer (-48): 515-538 (complementary strand)

5.6. pUC18Sfi Restriction analysis.

RESTRICTION ENZYMES THAT DO NOT CUT pUC18Sfi

Aar I	Bcl I	BsiW I	Bsu36 I	Mfe I	Oli I	Psr I	Stu I
Afe I	BfrB I	Bsm I	Btg I	Mlu I	Pac I	Rsr II	Swa I
Afl II	Bgl II	BsmF I	Btr I	Msc I	PflM I	Sac II	Tih111 I
Age I	Blp I	BspE I	Cla I	Nae I	Pme I	SanD I	Xcm I
Alo I	Bpl I	BsrG I	Dra III	Nco I	Pml I	SexA I	Xho I
Apa I	Bpu10 I	BssH II	EcoN I	NgoM1 V	Ppu10 I	Sgf I	
Asc I	BsaA I	BstB I	EcoR V	Nhe I	PpuM I	SgrA I	
Bae I	BsaB I	BstE II	Fse I	Not I	PshA I	SnaB I	
Bbs I	BseR I	BstX I	FspA I	Nru I	Psi I	Spe I	
BbvC I	Bsg I	BstZ17 I	Hpa I	Nsi I	PspOM I	Srf I	

RESTRICTION ENZYMES THAT CUT ONCE pUC18Sfi

Aat II	Apo I	Bcg I	BstAP I	<u>HinD III</u>	<u>Pst I</u>	Sfo I
Acc I	Ava I	Bpm I	Eag I	Kas I	<u>Sac I</u>	<u>Sma I</u>
Acc65 I	Avr I	Bsa I	Ecl136 II	<u>Kpn I</u>	<u>Sal I</u>	<u>Sph I</u>
Afl III	<u>BamH I</u>	BsaX I	EcoO109 I	Nar I	Sap I	Ssp I
Ahd I	Ban II	BspM I	<u>EcoR I</u>	Nde I	Sbf I	<u>Xba I</u>
AlwN I	Bbe I	BsrF I	Hinc II	Pci I	Sca I	Xma I

! **NOTE:** Underlined multicloning restriction sites.

6. Protocol.

6.1. Cloning and transformation of donor strain.

To insert your DNA fragment in the chromosome of the recipient strain, first of all you must clone it into the pUTmini-Tn5 vector. We recommend that you clone your fragment in pUC18Not auxiliary plasmid, prior the cloning in pUT plasmid. The protocol should be as follows.

5. Restriction analysis.

5.1. pUTmini-Tn5 Cm Restriction analysis.

RESTRICTION ENZYMES THAT DO NOT CUT pUTmini-Tn5 Cm

Aat II	BbvC I	EcoLCR I	Nru I	Sbf I
Acc65 I	BfrB I	EcoN I	Nsi I	Spe I
Afl II	Blp I	Fse I	Pme I	Srf I
Ale I	Bpl I	FspA I	PshA I	Swa I
Apa I	BsiW I	Kpn I	PspOM I	Zra I
Asc I	BstB I	Mlu I	Psr I	
AsiS I	Bsu36 I	Nco I	Rsr II	
Bae I	Cla I	Nde I	Sac I	

RESTRICTION ENZYMES THAT CUT ONCE pUTmini-Tn5 Cm

Age I	BssH II	Dra III	Pac I	Sph I
Bcl I	BssS I	Drd I	Pml I	Stu I
Bmt I	BstAP I	EcoR V	Pvu I	Xba I
Bpu10 I	BstE II	Mfe I	Sca I	Xcm I
BsaX I	BstX I	Msc I	SexA I	Xmn I
BsrG I	BstZ17 I	Nhe I	SnaB I	

5.2. pUTmini-Tn5 Km Restriction analysis.

RESTRICTION ENZYMES THAT DO NOT CUT pUTmini-Tn5 Km

Aat II	Bbs I	BsiW I	Cla I	Fsp AI	Pac I	Psp OMI	Spe I
Afl II	Bbv CI	Bsp EI	Ecl 136II	Mlu I	Pfl MI	Psr I	Srf I
Apa I	Bfr BI	Bsr GI	Eco NI	Nde I	Pme I	Sac I	Stu I
Asc I	Blp I	Bst EII	Eco RV	Nsi I	Ppu 10I	Sbf I	Swa I
Bae I	Bsa XI	Bst XI	Fse I	Oli I	Psh AI	Sgf I	

RESTRICTION ENZYMES THAT CUT ONCE pUTmini-Tn5 Km

Acc 65I	Bsg I	Dra III	Not I	Rsr II	Xba I
Age I	Bsm I	Kpn I	Nru I	Sca I	Xho I
Alo I	Bst BI	Mfe I	Pml I	Sex AI	Xmn I
Bpl I	Bst Z17I	Nco I	Psi I	Sna BI	
Bpu 10I	Bsu 36I	Nhe I	Pvu I	Sph I	

5.3. pUTmini-Tn5 Sm/Sp Restriction analysis.

RESTRICTION ENZYMES THAT DO NOT CUT pUTmini-Tn5 Sm/Sp

<i>Aat</i> II	<i>Bae</i> I	<i>Bsi</i> WI	<i>Bsu</i> 36I	<i>Fse</i> I	<i>Nde</i> I	<i>Pfl</i> MI	<i>Psr</i> I	<i>Spe</i> I
<i>Acc</i> 65I	<i>Bbv</i> CI	<i>Bsp</i> EI	<i>Cla</i> I	<i>Fsp</i> AI	<i>Nru</i> I	<i>Pme</i> I	<i>Rsr</i> II	<i>Srf</i> I
<i>Afl</i> II	<i>Bfr</i> BI	<i>Bsr</i> GI	<i>Ecl</i> 136II	<i>Kpn</i> I	<i>Nsi</i> I	<i>Ppu</i> 10I	<i>Sac</i> I	<i>Stu</i> I
<i>Apa</i> I	<i>Bpl</i> I	<i>Bst</i> BI	<i>Eco</i> NI	<i>Mlu</i> I	<i>Oli</i> I	<i>Psh</i> AI	<i>Sbf</i> I	<i>Swa</i> I
<i>Asc</i> I	<i>Bsa</i> XI	<i>Bst</i> XI	<i>Eco</i> RV	<i>Nco</i> I	<i>Pac</i> I	<i>Psp</i> OMI	<i>Sgf</i> I	

RESTRICTION ENZYMES THAT CUT ONCE pUTmini-Tn5 Sm/Sp

<i>Age</i> I	<i>Bsg</i> I	<i>Bst</i> Z17I	<i>Nhe</i> I	<i>Sca</i> I	<i>Xba</i> I
<i>Alo</i> I	<i>Bsm</i> I	<i>Drd</i> I	<i>Not</i> I	<i>Sex</i> AI	<i>Xho</i> I
<i>Bbs</i> I	<i>Bss</i> SI	<i>Hpa</i> I	<i>Pml</i> I	<i>Sna</i> BI	
<i>Blp</i> I	<i>Bst</i> API	<i>Mfe</i> I	<i>Psi</i> I	<i>Sph</i> I	
<i>Bpu</i> 10I	<i>Bst</i> EII	<i>Msc</i> I	<i>Pvu</i> II	<i>Tth</i> 111I	

5.4. pUTmini-Tn5 Tc Restriction analysis.

RESTRICTION ENZYMES THAT DO NOT CUT pUTmini-Tn5 Tc

<i>Aat</i> II	<i>Bae</i> I	<i>Bsa</i> XI	<i>Bst</i> EII	<i>Fse</i> I	<i>Nsi</i> I	<i>Psp</i> OMI	<i>Sgf</i> I
<i>Acc</i> 65I	<i>Bbv</i> CI	<i>Bsi</i> WI	<i>Bst</i> XI	<i>Kpn</i> I	<i>Oli</i> I	<i>Psr</i> I	<i>Spe</i> I
<i>Afl</i> II	<i>Bfr</i> BI	<i>Bsp</i> EI	<i>Bsu</i> 36I	<i>Mlu</i> I	<i>Pac</i> I	<i>Rsr</i> II	<i>Srf</i> I
<i>Apa</i> I	<i>Blp</i> I	<i>Bsr</i> GI	<i>Cla</i> I	<i>Nco</i> I	<i>Pme</i> I	<i>Sac</i> I	<i>Stu</i> I
<i>Asc</i> I	<i>Bpl</i> I	<i>Bst</i> BI	<i>Ecl</i> 136II	<i>Nde</i> I	<i>Ppu</i> 10I	<i>Sbf</i> I	<i>Swa</i> I

RESTRICTION ENZYMES THAT CUT ONCE pUTmini-Tn5 Tc

<i>Age</i> I	<i>Bsg</i> I	<i>Drd</i> I	<i>Mfe</i> I	<i>Psh</i> AI	<i>Sex</i> AI	<i>Xcm</i> I
<i>Alo</i> I	<i>Bss</i> HII	<i>Eco</i> NI	<i>Msc</i> I	<i>Psi</i> I	<i>Sna</i> BI	<i>Xho</i> I
<i>Bbs</i> I	<i>Bss</i> SI	<i>Eco</i> RV	<i>Not</i> I	<i>Pvu</i> I	<i>Sph</i> I	<i>Xmn</i> I
<i>Bcl</i> I	<i>Bst</i> Z17I	<i>Fsp</i> AI	<i>Nru</i> I	<i>Pvu</i> II	<i>Tth</i> 111I	
<i>Bpu</i> 10I	<i>Dra</i> III	<i>Hpa</i> I	<i>Pml</i> I	<i>Sca</i> I	<i>Xba</i> I	

5.5. pUC18Not Restriction Enzymes.

RESTRICTION ENZYMES THAT DO NOT CUT pUC18Not

<i>Aar</i> I	<i>Blp</i> I	<i>BstE</i> II	<i>Nae</i> I	<i>Psr</i> I
<i>Afe</i> I	<i>BmgB</i> I	<i>BstX</i> I	<i>Nco</i> I	<i>Rsr</i> II
<i>Afl</i> II	<i>Bmt</i> I	<i>BstZ17</i> I	<i>NgoMI</i> V	<i>Sac</i> II
<i>Age</i> I	<i>Bpl</i> I	<i>Bsu36</i> I	<i>Nhe</i> I	<i>SanD</i> I
<i>Ale</i> I	<i>Bpu</i> 10 I	<i>Btg</i> I	<i>Nru</i> I	<i>SexA</i> I
<i>Alo</i> I	<i>BsaA</i> I	<i>Cla</i> I	<i>Nsi</i> I	<i>Sfi</i> I
<i>Apa</i> I	<i>BsaB</i> I	<i>Dra</i> III	<i>Oli</i> I	<i>Sgf</i> I
<i>Asc</i> I	<i>BseR</i> I	<i>EcoN</i> I	<i>Pac</i> I	<i>SgrA</i> I
<i>AsiS</i> I	<i>Bsg</i> I	<i>EcoR</i> V	<i>PflM</i> I	<i>SnaB</i> I
<i>Avr</i> II	<i>BsiW</i> I	<i>Fal</i> I	<i>Pme</i> I	<i>Spe</i> I
<i>Bae</i> I	<i>Bsm</i> I	<i>Fse</i> I	<i>Pml</i> I	<i>Srf</i> I
<i>Bbs</i> I	<i>BsmF</i> I	<i>FspA</i> I	<i>Ppu</i> 10 I	<i>Stu</i> I
<i>BbvC</i> I	<i>BspE</i> I	<i>Hpa</i> I	<i>PpuM</i> I	<i>Sty</i> I
<i>Bcl</i> I	<i>BsrG</i> I	<i>Mfe</i> I	<i>PshA</i> I	<i>Swa</i> I
<i>BfrB</i> I	<i>BssH</i> II	<i>Mlu</i> I	<i>Psi</i> I	<i>Tth</i> 111 I
<i>Bgl</i> II	<i>BstB</i> I	<i>Msc</i> I	<i>PspOM</i> I	<i>Xho</i> I

RESTRICTION ENZYMES THAT CUT ONCE pUC18Not

<i>Acc</i> I	<u><i>Bam</i>H I</u>	<i>BspM</i> I	<u><i>HinD</i> III</u>	<u><i>Pst</i> I</u>	<u><i>Sma</i> I</u>	<i>Zra</i> I
<i>Aat</i> II	<i>Ava</i> I	<i>BsaX</i> I	<i>Hinc</i> II	<i>Pfo</i> I	<i>Sfo</i> I	<i>Xmn</i> I
<i>Acc65</i> I	<i>Ban</i> II	<i>BsrF</i> I	<i>Kas</i> I	<u><i>Sac</i> I</u>	<u><i>Sph</i> I</u>	
<i>Afl</i> III	<i>Bbe</i> I	<i>BstAP</i> I	<u><i>Kpn</i> I</u>	<u><i>Sal</i> I</u>	<i>Ssp</i> I	
<i>Ahd</i> I	<i>Bcg</i> I	<i>Ecl</i> 136 II	<i>Nar</i> I	<i>Sap</i> I	<u><i>Xba</i> I</u>	
<i>AlwN</i> I	<i>Bpm</i> I	<i>Eco</i> 0109 I	<i>Nde</i> I	<i>Sbf</i> I	<i>Xcm</i> I	
<i>Apo</i> I	<i>Bsa</i> I	<u><i>EcoR</i> I</u>	<i>Pci</i> I	<i>Sca</i> I	<i>Xma</i> I	

! NOTE: Underlined multicloning restriction sites.