

## Anti-Ovalbumin IgE (mouse) EIA Kit

Item No. 500840



**Customer Service** 800.364.9897 \* **Technical Support** 888.526.5351

[www.caymanchem.com](http://www.caymanchem.com)

## TABLE OF CONTENTS

<b>GENERAL INFORMATION</b>	3	Materials Supplied
	4	Precautions
	4	If You Have Problems
	4	Storage and Stability
	4	Materials Needed but Not Supplied
<b>INTRODUCTION</b>	5	Background
	5	About This Assay
	6	Description of Immunometric EIAs
<b>PRE-ASSAY PREPARATION</b>	8	Buffer Preparation
	8	Sample Preparation
<b>ASSAY PROTOCOL</b>	9	Preparation of Assay-Specific Reagents
	11	Plate Set Up
	12	Performing the Assay
<b>ANALYSIS</b>	14	Calculations
	15	Performance Characteristics
<b>RESOURCES</b>	20	Troubleshooting
	21	References
	21	Related Products
	22	Warranty and Limitation of Remedy
	23	Plate Template
	24	Notes

## GENERAL INFORMATION

### Materials Supplied

Item Number	Item	96 wells Quantity/Size
400840	Ovalbumin-biotin Conjugate	1 vial/100 dtn
400842	Goat Anti-Mouse IgE Precoated 96-Well Strip Plate	1 plate
400844	Anti-Ovalbumin IgE (mouse) EIA Standard	1 vial/200 ng
400054	Immunoassay Buffer B Concentrate (10X)	2 vials/10 ml
400063	Streptavidin-HRP	2 vials/1.5 ml
400062	Wash Buffer Concentrate (400X)	1 vial/5 ml
400035	Polysorbate 20	1 vial/3 ml
400074	TMB Substrate Solution	1 vial/12 ml
10011355	HRP Stop Solution	1 vial/12 ml
400012	96-Well Cover Sheet	3 covers

If any of the items listed above are damaged or missing, please contact our Customer Service department at (800) 364-9897 or (734) 971-3335. We cannot accept any returns without prior authorization.



**WARNING: This product is for laboratory research use only; not for administration to humans. Not for human or veterinary diagnostic or therapeutic use.**

## Precautions

Please read these instructions carefully before beginning this assay.

The reagents in this kit have been tested and formulated to work exclusively with Cayman's Anti-Ovalbumin IgE (mouse) EIA Kit. This kit may not perform as described if any reagent or procedure is replaced or modified. The Stop Solution provided with this kit is an acid solution. Please wear appropriate personal protection equipment (*e.g.*, safety glasses, gloves, and lab-coat) when using this material.

**For research use only. Not for human or diagnostic use.**

## If You Have Problems

### Technical Service Contact Information

**Phone:** 888-526-5351 (USA and Canada only) or 734-975-3888  
**Fax:** 734-971-3641  
**Email:** techserv@caymanchem.com  
**Hours:** M-F 8:00 AM to 5:30 PM EST

In order for our staff to assist you quickly and efficiently, please be ready to supply the lot number of the kit (found on the outside of the box).

## Storage and Stability

This kit will perform as specified if stored as directed at 4°C and used before the expiration date indicated on the outside of the box.

## Materials Needed But Not Supplied

1. A plate reader capable of measuring absorbance at 450 nm.
2. Adjustable pipettes and a repeating pipettor.
3. A source of pure water; glass distilled or deionized water is acceptable
4. Materials used for **Sample Preparation** (see page 8).

## INTRODUCTION

### Background

Immunization of mice with chicken egg albumin (ovalbumin/OVA) in a precipitate complex with aluminum hydroxide (alum) is a highly effective means of inducing a potent TH2-mediated immune response.<sup>1-6</sup> OVA/alum immunized mice produce anti-OVA antibodies predominantly of the IgG<sub>1</sub> and IgE isotypes that mediate tissue-specific effector functions in multiple mouse models of chronic inflammation, including allergic asthma, allergic rhinitis, and cutaneous hypersensitivity.<sup>1-4</sup>

When using one of these models, it is often desirable to measure anti-OVA antibody levels in the plasma or serum to determine the effectiveness of the immunization, the activity of a drug, or the effect of a specific gene deletion. IgG<sub>1</sub> is the predominant anti-OVA immunoglobulin isotype found in the serum or plasma of mice immunized with OVA/alum; the plasma concentration of OVA-specific IgG<sub>1</sub> is typically 1,000-fold greater than that of OVA-specific IgE.<sup>1</sup> Thus, plate-bound ovalbumin can be used to capture and quantify ovalbumin-specific IgG<sub>1</sub> directly from the plasma or serum using Cayman's Anti-Ovalbumin IgG<sub>1</sub> (mouse) EIA Kit (Item No. 500830). However, it is not possible to use plate-bound ovalbumin to capture IgE; the more abundant IgG<sub>1</sub> occupies all of the available binding sites. For the measurement of ovalbumin-specific IgE, total plasma or serum IgE is selectively captured using a plate-bound antibody specific for mouse IgE. A biotin-labeled ovalbumin reagent is then added to quantify the ovalbumin-specific subset of this captured IgE.

### About This Assay

Cayman's Anti-Ovalbumin IgE (mouse) EIA Kit is an immunometric assay which can be used to measure anti-ovalbumin of the IgE isotype in mouse plasma and serum without prior sample purification. A monoclonal anti-ovalbumin IgE produced from mice immunized with OVA/alum is used as the standard. The standard curve spans the range of 1.56-100 ng/ml, with an LLOQ of 3.12 ng/ml.

## Description of Immunometric EIAs

Each well of the microwell plate supplied in the kit has been coated with Goat anti-mouse IgE polyclonal antibody. IgE antibodies, if present in the biological fluid sample, will bind to the immobilized anti-IgE antibody. A biotin-conjugated ovalbumin reagent is then added to the well and is bound by ovalbumin-specific IgE. An HRP-conjugated streptavidin reagent is then added and binds to the biotin, allowing quantitation of the ovalbumin-specific IgE antibody. Addition of the HRP Substrate 3,3',5,5'-tetramethylbenzidine (TMB), followed by Stop Solution produces a yellow colored product which can be measured spectrophotometrically. The intensity of the color is directly proportional to the amount of bound streptavidin-HRP, which is proportional to the concentration of the ovalbumin-biotin, which is proportional to the anti-ovalbumin IgE antibody.

$\text{Absorbance} \propto [\text{streptavidin-HRP}] \propto [\text{ovalbumin-biotin}] \propto [\text{Anti-ovalbumin IgE antibody}]$

A schematic of this process is shown in Figure 1, on page 7.




Figure 1. Schematic of the Immunometric EIA

## Buffer Preparation

Store all diluted buffers at 4°C; they will be stable for about two months.

### 1. Assay Buffer Preparation

Dilute the contents of one vial of Immunoassay Buffer B Concentrate (10X) (Item No. 400054) with 90 ml of water. Be certain to rinse the vial to remove any salts that may have precipitated. *NOTE: It is normal for the concentrated buffer to contain crystalline salts after thawing. These will completely dissolve upon dilution with water.*

### 2. Wash Buffer Preparation

**5 ml vial Wash Buffer Concentrate (400X) (Item No. 400062):** Dilute to a total volume of 2 L with water and add 1 ml of Polysorbate 20 (Item No. 400035).

*NOTE: Polysorbate 20 is a viscous liquid and cannot be measured by a regular pipette. A positive displacement pipette or a syringe should be used to deliver small quantities accurately.*

## Sample Preparation

Prior to use, it is recommended that serum or plasma samples from OVA/alum-immunized mice be diluted in Assay Buffer at least 1:25 in order to fall within the range of the standard curve (see Table 4 on page 16). In general, mouse serum or plasma (prepared using heparin or EDTA as the anticoagulant) can be used directly in the assay following dilution in Assay Buffer.

## Preparation of Assay-Specific Reagents

### Anti-Ovalbumin IgE (mouse) EIA Standard

Reconstitute the lyophilized purified Anti-Ovalbumin IgE (mouse) EIA Standard (Item No. 400844) with 2.0 ml of Assay Buffer. Mix gently. The concentration of this solution (the bulk standard) is 100 ng/ml. The reconstituted standard is stable for two weeks at 4°C. Enough standard is provided to produce four duplicate-well standard curves for use on different days, if necessary.

To prepare the standard for use in the EIA: Obtain eight clean test tubes or plastic microfuge tubes and label them #1 through #8. Aliquot 250 µl of Assay Buffer into tubes #2-8. Transfer 500 µl of freshly prepared stock standard (100 ng/ml) to tube #1. Serially dilute the standard by removing 250 µl from tube #1 and placing into tube #2. Mix gently. Next, remove 250 µl from tube #2 and place into tube #3; mix gently. Repeat this process for tubes #4-7.




Figure 2. Preparation of the Anti-Ovalbumin IgE (mouse) standards

## Ovalbumin-biotin Conjugate

Reconstitute the lyophilized Ovalbumin-Biotin conjugate (Item No. 400840) with 12 ml Assay Buffer. Mix gently. The reconstituted conjugate is stable for two weeks at 4°C. Enough conjugate is provided for 96 wells.

## Streptavidin-HRP

This reagent is supplied as a concentrated (10X) stock solution of Streptavidin conjugated to HRP. On the day of the assay, prepare a Working Solution by adding 1.2 ml of the Streptavidin-HRP (Item No. 400063) to 10.8 ml Assay Buffer (12 ml total). This Working Solution is stable for 24 hours at 4°C, *protected from light*. In the event that two or more experiments are performed with this kit more than 24 hours apart, two vials of stock solution has been provided to produce an additional 12 ml of the Working Solution.

## Plate Set Up

The 96-well plate(s) included with this kit is supplied ready to use. It is not necessary to rinse the plate(s) prior to adding the reagents. *NOTE: If you do not need to use all of the strips at once, place the unused strips back in the plate packet and store according to the plate insert at 4°C. Be sure the packet is sealed with the desiccant inside.*

Each plate or set of strips must contain an eight point standard curve run in duplicate. *NOTE: Each assay must contain this minimum configuration in order to ensure accurate and reproducible results.* Each sample should be assayed at a minimum of two dilutions and each dilution should be assayed in duplicate. For statistical purposes, we recommend assaying samples in triplicate.

A suggested plate format is shown in Figure 3, below. The user may vary the location and type of wells present as necessary for each particular experiment. The plate format provided below has been designed to allow for easy data analysis using a convenient spreadsheet offered by Cayman (see page 14, for more details). We suggest you record the contents of each well on the template sheet provided (see page 23).

	1	2	3	4	5	6	7	8	9	10	11	12
A	S1	S1	1	1	1	9	9	9	17	17	17	25
B	S2	S2	2	2	2	10	10	10	18	18	18	25
C	S3	S3	3	3	3	11	11	11	19	19	19	25
D	S4	S4	4	4	4	12	12	12	20	20	20	26
E	S5	S5	5	5	5	13	13	13	21	21	21	26
F	S6	S6	6	6	6	14	14	14	22	22	22	26
G	S7	S7	7	7	7	15	15	15	23	23	23	27
H	S8	S8	8	8	8	16	16	16	24	24	24	27

S1-S8 - Standards 1-8  
1-27 - Samples

Figure 3. Sample plate format

## Performing the Assay

### Pipetting Hints

- Use different tips to pipette each reagent.
- Before pipetting each reagent, equilibrate the pipette tip in that reagent (*i.e.*, slowly fill the tip and gently expel the contents, repeat several times).
- Do not expose the pipette tip to the reagent(s) already in the well.

### Addition of Standards and Samples and First Incubation

1. Add 100  $\mu\text{l}$  of the standards or diluted sample to the appropriate wells on the plate. Each sample should be assayed in duplicate, triplicate recommended.
2. Cover the plate with a 96-Well Cover Sheet (Item No. 400012). Incubate for two hours at room temperature on an orbital shaker.

### Addition of Ovalbumin-biotin Conjugate and Second Incubation

1. Empty the wells and rinse four times with Wash Buffer. Each well should be completely filled with Wash Buffer during each wash. Invert the plate between wash steps to empty the fluid from the wells. After the last wash, gently tap the inverted plate on absorbent paper to remove the residual Wash Buffer.
2. Add 100  $\mu\text{l}$  of the Ovalbumin-biotin Conjugate working solution to each well of the plate.
3. Cover the plate with a 96-Well Cover Sheet and incubate for one hour at room temperature on an orbital shaker.

### Addition of Streptavidin-HRP and Third Incubation

1. Empty the wells and rinse four times with Wash Buffer. Each well should be completely filled with Wash Buffer during each wash. Invert the plate between wash steps to empty the fluid from the wells. After the last wash, gently tap the inverted plate on absorbent paper to remove the residual Wash Buffer.
2. Add 100  $\mu\text{l}$  of the Streptavidin-HRP Working Solution to each well of the plate.
3. Cover the plate with a 96-Well Cover Sheet and incubate for 30 minutes at room temperature on an orbital shaker.

### Development of the Plate

1. Empty the wells and rinse four times with Wash Buffer.
2. Add 100  $\mu\text{l}$  of TMB Substrate Solution (Item No. 400074) to each well of the plate.
3. Cover the plate with plastic film and incubate for 15 minutes at room temperature on an orbital shaker. Development of the blue color can be monitored at 650 nm. When the maximum standard (S1) optical density (O.D.) value reaches 0.5-0.6, Stop Solution should be added to the entire plate.
4. DO NOT WASH THE PLATE. Add 100  $\mu\text{l}$  of HRP Stop Solution (Item No. 10011355) to each well of the plate. Blue wells will turn yellow and colorless wells will remain colorless. *NOTE: The Stop Solution in this kit contains an acid. Wear appropriate protection and use caution when handling this solution.*

### Reading the Plate

1. Wipe the bottom of the plate with a clean tissue to remove fingerprints, dirt, etc.
2. Read the plate at a wavelength of 450 nm.

## ANALYSIS

Many plate readers come with data reduction software that plots data automatically. Alternatively a spreadsheet program can be used. *NOTE: Cayman has a computer spreadsheet available for data analysis. Please contact Technical Service or visit our website ([www.caymanchem.com/analysis/immuno](http://www.caymanchem.com/analysis/immuno)) to obtain a free copy of this convenient data analysis tool.*

### Calculations

#### Plotting the Standard Curve and Determining the Sample Concentration

Using computer reduction software, plot absorbance (linear y-axis) for standards (S1-S8) versus concentration (linear x-axis) and fit the data with a quadratic equation. Using the equation of the line, calculate the concentration of anti-ovalbumin IgE in each sample.

## Performance Characteristics

### Sample Data

The standard curve presented here is an example of the data typically produced with this kit; however, your results will not be identical to these. You **must** run a new standard curve. Do not use the data below to determine the values of your samples. Your results could differ substantially. Development of the plate for 15 minutes typically results in an absorbance of >1.0 O.D. units for the 100 ng/ml standard.

Anti-OVA IgE (ng/ml)	Absorbance	
100	2.330	2.264
50	1.312	1.266
25	0.686	0.664
12.5	0.359	0.324
6.25	0.181	0.171
3.12	0.100	0.095
1.56	0.041	0.045
0	0.036	0.026

Table 1. Typical results




Figure 4. Typical standard curve

### Precision:

The intra- and inter-assay CVs have been determined at multiple points on the standard curve. These data are summarized in the graph on page 16 and in the table below.

Anti-OVA IgE (ng/ml)	%CV* Intra-assay variation	%CV* Inter-assay variation
100	2.99	2.50
50	3.73	2.80
25	2.02	4.32
12.5	2.55	6.58
6.25	3.08	2.74
3.12	8.62	10.72
1.56	7.98	30.25

**Table 2. Intra- and inter-assay variation**

\*%CV represents the variation in concentration (not absorbance) as determined using a reference standard curve.

†Outside of the recommended usable range of the assay.

Anti-OVA IgE (ng/ml)	Mean of O.D.	Standard Deviation (S.D.)	O.D. - (1.64 x S.D.)
100	2.485	0.054	2.397
50	1.426	0.045	1.353
25	0.745	0.014	0.722
12.5	0.382	0.009	0.367
6.25	0.195	0.006	0.186
3.12	0.104	0.008	0.090
1.56	0.040	0.003	0.035
0	0.020	0.013	0.042*

\*O.D. + (1.64 x S.D.)

**Table 3. Determination of LLOQ**

The lower limit of quantitation (LLOQ) is defined as the lowest standard concentration in which O.D. - (1.64 x S.D.) is higher than the blank value of O.D. + (1.64 x S.D.). The LLOQ is 3.12 ng/ml.

	Calculated Anti-OVA IgE concentration (ng/ml)							
	Unimmunized Control Mice				OVA/Alum Immunized Mice			
Dilution Factor	1	2	3	4	1	2	3	4
5	BLQ*	BLQ	BLQ	BLQ				
25	-	-	-	-	3060	880	3532	2464
50	-	-	-	-	3584	892	4403	3109
100	-	-	-	-	3904	918	4903	4298
200	-	-	-	-	4394	927	5594	5939
Mean	-	-	-	-	3735	905	4608	3952
Std. Deviation					560	22	868	1527

**Table 4. Reproducibility of the assay over a wide dilution range**

Plasma from four unimmunized BALB/c female mice and four BALB/c female mice immunized with OVA/alum was diluted from 1:5 to 1:200 and the concentration of Anti-OVA IgE determined using Cayman's Anti-Ovalbumin IgE (mouse) EIA Kit.

\*BLQ = Below the limit of quantitation

## Troubleshooting

Problem	Possible Causes	Recommended Solutions
Erratic values; dispersion of duplicates	A. Trace organic contaminants in the water source B. Poor pipetting/technique	A. Replace activated carbon filter or change source of UltraPure water
Poor development (low signal) of standard curve	A. Plate required more development time B. Standard was diluted incorrectly C. Standard is degraded	

## References

- Hogan, S.P., Mould, A., Kikutani, H., *et al.* Aeroallergen-induced eosinophilic inflammation, lung damage, and airways hyperreactivity in mice can occur independently of IL-4 and allergen-specific immunoglobulins. *J. Clin. Invest.* **99(6)**, 1329-1339 (1997).
- Kennedy, J.D., Hatfield, C.A., Fidler, S.F., *et al.* Phenotypic characterization of T lymphocytes emigrating into lung tissue and the airway lumen after antigen inhalation in sensitized mice. *Am. J. Respir. Cell Mol. Biol.* **12(6)**, 613-623 (1995).
- Saito, H., Howie, K., Wattie, J., *et al.* Allergen-induced murine upper airway inflammation: Local and systemic changes in murine experimental allergic rhinitis. *Immunology* **104(2)**, 226-234 (2001).
- Sawada, K., Nagai, H., Basaki, Y., *et al.* The expression of murine cutaneous late phase reaction requires both IgE antibodies and CD4 T cells. *Clin. Exp. Allergy* **27(2)**, 225-231 (1996).
- Pichavant, M., Goya, S., Hamelmann, E., *et al.* Animal models of airway sensitization. *Curr. Protoc. Immunol.* **79**, 15.18.1-15.18.19 (2007).
- Kung, T.T., Jones, H., Adams, G.K.I., *et al.* Characterization of a murine model of allergic pulmonary inflammation. *Int. Arch. Allergy Immunol.* **105(1)**, 83-90 (1994).

## Related Products

Anti-Ovalbumin IgG<sub>1</sub> (mouse) EIA Kit - Item No. 500830  
 Interleukin-4 (mouse) EIA Kit - Item No. 500860  
 Interleukin-5 (mouse) EIA Kit - Item No. 500880  
 Interleukin-13 (mouse) EIA Kit - Item No. 500870

## Warranty and Limitation of Remedy

Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman **warrants only** to the original customer that the material will meet our specifications at the time of delivery. Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, indirect, incidental or consequential damages, even if Cayman is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees.

Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a refund of the purchase price, or at Cayman's option, the replacement, at no cost to Buyer, of all material that does not meet our specifications.

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.

12								
11								
10								
9								
8								
7								
6								
5								
4								
3								
2								
1								
	A	B	C	D	E	F	G	H

## NOTES

This document is copyrighted. All rights are reserved. This document may not, in whole or part, be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form without prior consent, in writing, from Cayman Chemical Company.

©09/15/2014, Cayman Chemical Company, Ann Arbor, MI, All rights reserved. Printed in U.S.A.