

## miRNASelect™ pEP-miR Cloning and Expression Vector

**CATALOG NUMBER:** MIR-EXP-C      **STORAGE:** -80°C

**QUANTITY:** 2 vectors; each contains 100 µL of bacterial glycerol stock

### Components

1. miRNASelect™ pEP-miR Cloning and Expression Vector (Part No. MIR-EXP): One tube
2. miRNASelect™ pEP-miR Null Control Vector (Part No. MIR-NULL): One tube


### Background

MicroRNAs (miRNAs) are 18–24 nucleotide RNA molecules that regulate the stability or translational efficiency of target mRNAs. These regulatory RNAs function by acting as sequence-specific guides which recruit a large protein complex known as the RNA-induced silencing complex (RISC) to target mRNAs which are subsequently silenced. Diverse functions have been attributed to miRNAs including the regulation of cellular differentiation, proliferation, and apoptosis. Moreover, significant evidence has accumulated implicating a fundamental role for miRNAs in the development of cancer.

miRNAs are initially transcribed as long precursor transcripts known as primary microRNAs (pri-miRNAs). Within these transcripts, the mature miRNA sequences are found in ~60–80 nucleotide hairpin structures. Mature miRNAs are generated from pri-miRNAs by sequential processing (Figure 1). Pri-miRNAs are initially recognized in the nucleus by the microprocessor complex which includes as core components the RNase-III enzyme Drosha and its obligate partner DGCR8. This complex excises the hairpin structure containing the mature miRNA sequence. The liberated hairpins, referred to as precursor miRNAs (pre-miRNAs), are recognized by the nuclear export factor exportin 5 which transports them to the cytoplasm. There, the RNase-III enzyme Dicer performs a second cleavage to generate a double-stranded 18–24 nucleotide RNA molecule. The RISC then associates with this RNA duplex and unwinds it. Generally, only one strand is stably incorporated into the RISC; the other is discarded and rapidly degraded. miRNAs guide the RISC to target messages that are subsequently cleaved or translationally silenced.

Synthetic miRNA molecules based on predicted mature miRNA sequence are sometimes used. Despite their optimized design criteria, synthetic miRNAs underscore the importance of primary miRNA in its native expressed form. The primary miRNA contains critical biological components involved in mature miRNA expression and cellular processing, and is often processed into several mature miRNA molecules.


Cell Biolabs' miRNASelect™ Cloning and Expression Vector is designed to clone and express an individual miRNA precursor in its native context while preserving putative hairpin structures to ensure biologically relevant interactions with endogenous processing machinery and regulatory partners, leading to properly cleaved microRNAs. Individual miRNA precursor from any species can be cloned between BamHI and Nhe I sites (Figure 2).



**Figure 1.** miRNA Biogenesis and function


The miRNASelect™ pEP-miR cloning and expression vector contains the following features:

- **miRNA Processing** – miRNA stem loop precursor in its native context is cloned between BamHI and Nhe I sites. To preserve the putative hairpin structure and proper endogenous processing, miRNA stem loop sequence is flanked by its native intron sequence.
- **EF-1 $\alpha$  Promoter** - ensures a high level of expression in mammalian cells
- **Puromycin resistance marker** - to monitor cells positive for expression and stable selection
- **SV40 polyadenylation signal** - enables efficient termination of transcription
- **pUC origin** - for high copy replication and maintenance of the plasmid in *E. coli*
- **Ampicillin resistance gene** - for selection in *E. coli*



**Figure 2.** Schematic representation of pEP-miR cloning and expression vector.

Cell Biolabs' miRNASelect™ pEP-miR Null Control Vector is similar to the pEP-miR vector except it does not contain any miRNA precursor or the BamHI and NheI cloning sites (Figure 3).



**Figure 3.** Schematic representation of pEP-miR null control vector.

miRNA precursor cloning sites: BamHI-1.1 kb -NheI

```

1  tcgaggatcc  cttaagtct  tctgaggagc  atggacctgc  tagactgcaa  gctgctgcat
61  tctgaggagc  agggacctgc  taaactgcaa  gctgccgcag  gcggagattt  ttgcctgttt
121 ggtttatatt  gttaccactg  agtccccagc  tctgagaaca  gagcgaggca  tgtgatcgat
181 gcatcgtaac  tgttgagcaa  atgaatgctg  agctgtgatg  gcaccgacgt  gtccggggag
241 aggacgcacc  cggctgtgtg  gacatgtgcc  cagggcccag  gacagcgcca  cggagaggga
301 cacaccggc  tgtgtggaca  tgtgcccagg  gcccgggaca  gcgccacgga  agaggacgca
361 cccggctgtg  tgcacatgtg  cccagggccc  gggacagcgc  cacggaagag  gacgcacccc
421 gctgtgtgca  catgtgccca  gggcccggga  cagcgccacg  gaagaggacg  caccggctgt
481 tgtgcacatg  tgcccagggc  ccgggacagc  gccacggaag  aggacgcacc  cggctgtgtg
541 gacatgtgcc  cagggcccgg  gacagcgcca  cggagaggga  cgcacaggac  agcgccacgg
601 aagaggacgc  acccggtgtg  gtgcacatgt  gccacgggcc  cgggacagcg  ccacggaaga
661 ggacgcaccc  ggctgtgtgc  acatgtgcc  agggcccggg  acagcgccat  ggaagaggac
721 gcaccggct  gtgtgcacat  gtgccaggg  cccgggacag  cgccacggaa  gaggacgcac
781 ccggctgtgt  gcacatgtgc  ccagggccc  ggacagcgcc  acggaagagg  atgcaccggg
841 ctgtgtggac  atgtgcccag  ggcccgggac  agcgccacgg  aagaggacgc  acccggtgtg
901 gtggacatgt  gccacgggcc  cgggacagcg  ccacggaaga  ggacgcaccc  ggctgtgtgg
961 acatgtgcc  agggcccggg  acagcgccac  ggaagaggac  gcacaggaca  gcgccacgga
1021 agaggacgca  cccggctgtg  tggcagggga  ggttcccaag  agcaggctca  gggctctggg
1081 gtctgagtga  gggctctcca  gccaggctag  ctcgca

```

*Note: pEP-miR Cloning and Expression Vector should NOT be used as a transfection control vector, because the vector itself contains 1.1 kb of mir-941 precursor sequence between the BamHI and NheI sites. For transfection control, we recommend pEP-miR Null Control Vector.*

### miRNA Precursor Cloning

All of our premade human miRNA precursor clones are based on the following design, and the resulting overexpression of the mature miRNA is confirmed by Northern blot. Here we use human let-7a-2 miRNA as an example:



Reverse PCR Primer: tcga-gctagc-aaataccataaaataatcgta

- 3) PCR the miRNA precursor from genomic DNA and clone into the BamHI/NheI sites of the expression vector.

PCR Product of let-7a-2 precursor: let-7a-2 stem-loop sequence is underlined.

```
1 tcgaggatcc gcccaaatag gtgacagcac gatgaatcat tataagacta acttgtaatt
61 tccctgctta agaaatggta gttttccagc cattgtgact gcatgctccc aggttgaggt
121 agtaggttgt atagtttaga attacatcaa gggagataac tgtacagcct cctagctttc
181 cttgggtcct gcactaaaca acatggtgag aacgatcatg attcctccag gccttttctc
241 cctatgaaag gtaagattgg gtacgattat tttatgggat ttgctagctc ga
```

- 4) Validate the insert by DNA sequencing.

Forward Sequencing Primer: tcctcagccg tcgcttcatg

Reverse Sequencing Primer: gtgtggggaa actccatcgc

## Methods

- 1) Bacterial culture: the miRNA cloning and expression vector and the null control vector are provided as bacterial glycerol stocks. Individual colonies can be obtained by culturing in an LB-ampicillin plate.
- 2) Plasmid isolation: we recommend EndoFree Plasmid Kits (QIAGEN).
- 3) Transfection into target cells: we recommend Lipofectamine 2000 (Invitrogen).
- 4) Stable selection: 48 hrs post-transfection, select stable clones in 1-10 µg/mL Puromycin-containing medium.

## References

1. microRNA sequences listed in Sanger's miRBase (<http://microrna.sanger.ac.uk/sequences/>).
2. John, B., C. Sander and D. S. Marks (2006) *Methods Mol Biol* **342**: 101-13.
3. Johnson, S. M., H. Grosshans, J. Shingara, M. Byrom, R. Jarvis, A. Cheng, E. Labourier, K. L. Reinert, D. Brown and F. J. Slack (2005) *Cell* **120**: 635-47.
4. Kim, V. N. (2005) *Mol Cells* **19**: 1-15.
5. Lee, R. C., R. L. Feinbaum and V. Ambros (1993) *Cell* **75**: 843-54.
6. Lee, Y., K. Jeon, J. T. Lee, S. Kim and V. N. Kim (2002) *Embo J* **21**: 4663-70.
7. Yi, R., Y. Qin, I. G. Macara and B. R. Cullen (2003) *Genes Dev* **17**: 3011-6..

## Recent Product Citation

Mazzacurati, L. et al. (2015). Use of miRNA response sequences to block off-target replication and increase the safety of an unattenuated, glioblastoma-targeted oncolytic HSV. *Mol Ther.* **23**:99-107.

## Warranty

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of

this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

*This product is for RESEARCH USE ONLY; not for use in diagnostic procedures.*

**Contact Information**

Cell Biolabs, Inc.  
7758 Arjons Drive  
San Diego, CA 92126  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2008-2015: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.