IngoFlowEx[®] Kit

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Cat.no. / Kat. č.: ED7040

ENGLISH

1. Intended use

IngoFlowEx[®] Kit is intended for quantification of **phagocytic activity of human granulocytes and monocytes** by measuring the ingestion of fluorescently labeled *E. coli* bacteria in human heparinized whole blood using flow cytometry.

2. Introduction

The ability of phagocytosis is a characteristic feature of the so-called professional phagocytic cells (neutrophil and eosinophil granulocytes, monocytes and macrophages). The process includes the binding of foreign particles, e.g. bacteria, to cell surface receptors and the particle engulfment by the phagocytic cell. Differentiation of bound and internalized material is achieved by using trypan blue dye that suppresses the fluorescence of free and surface-bound bacteria. Opsonisation of bacteria, i.e. their coating with plasma proteins, opsonins (C3b, antibodies) facilitates the binding and engulfment. The bacteria used in the test were preopsonised by human AB plasma thus the results are not primarily dependent on the presence of opsonins in the tested blood sample.

3. Test principle

The test principle is based on measuring the intensity of the fluorescence of the phagocytes that engulfed FITC-labeled *E. coli*. A sample of heparinized blood is mixed with fluorescent *E. coli* and the mixture is incubated at 37 °C. For each reaction with *E. coli* a simultaneous control reaction without *E. coli* is performed that is used to assess the discrimination boundary between the phagocytosing and the non-phagocytosing cell population. Quenching solution is added at the end of the incubation period to suppress the fluorescence of non-engulfed bacteria. After washing the trypan blue away, the cells are fixed and RBCs lysed by adding a fixative-lysing reagent. After washing of the residues of erythrocytes the cellular DNA is labeled with propidium iodide. The labeling of cellular DNA helps to filter leukocytes out of debris and clumps of bacteria.

4. Precautions

- Intended for research use only.
- **Blood must be collected into a tube containing heparin.** Anticoagulants EDTA and citrate negatively affects results of the analysis.
- **Blood samples should be measured within 8 hours after collection.**

- The flow cytometer should be calibrated on a routine basis using fluorescent microbeads to ensure the stable sensitivity of detectors.
- Do not use reagents after the expiration date stated on vial labels.
- Avoid contamination of reagents.
- Use protective gloves and follow procedures for handling potentially infectious materials.
- Avoid contact of samples with skin, eyes and mucous membranes.

Warning: 10xLysing Solution (ED7040-3) contains methanol, formaldehyde and 2,2' -oxybisethanol (diethylene glycol).

Harmful

R-phrases

- R 20/21/22 Harmful if inhaled, in contact with skin and if swallowed.
- R 36/37/38 Irritating to eyes, respiratory system and skin.
- R 40 Limited evidence of a carcinogenic effect.
- R 43 May cause sensitisation by skin contact.
- R 68/20/21/22 Harmful: possible risk of irreversible effects if inhaled, in contact with skin and if swallowed.

S-phrases

- S 24 Avoid contact with skin.
- S 36/37 Wear suitable protective clothing and gloves.
- S 53 Avoid exposure – obtain special instructions before use.

5. Reagents provided

- ED7040-1 **E. coli-FITC** (ready to use) - amber glass vial containing 1 ml of a suspension of fluorescein-labeled opsonized *E. coli* strain K-12, vial is intended for 100 reactions.
- ED7040-2 **Quenching Solution** (ready to use) – vial containing 20 ml of buffered solution of trypan blue, dark blue solution, the vial is intended for 200 reactions.
- ED7040-3 **10xLysing Solution** (10x concentrated solution) - vial containing 60 ml of 10x concentrated transparent solution of fixative-lysing agent. Each vial is intended to prepare 600 ml of 1xLysing solution, i.e. 200 reactions.
- ED7040-4 **25xWash buffer** (25x concentrated solution) - vial containing 80 ml of concentrated wash buffer. Each vial is intended to prepare 2 liters of 1xWash buffer, i.e. 222 reactions.

- **ED7040-5 DNA Staining Solution** (ready to use) - amber vial containing 60 ml of pinkish solution of propidium iodide. Each vial is intended for 200 reactions.

6. Necessary material not supplied

- Deionized water (dH₂O), approximately 0.6 liter for one kit
- Phosphate buffered saline (PBS), approximately 2 liters
- Suitable 5ml test tubes for blood staining (e.g. 12 × 75 mm)
- Racks for the test tubes
- Centrifuge with suitable rotor for 5 ml tubes
- Automatic pipettes with disposable tips
- Vortex
- Thermal incubator or water bath (37 °C)
- Container with ice (ice cubes or crushed ice) suitable for placing a rack with tubes
- Flow cytometer - blue laser excitation 488 nm, light emission at 525 nm (FITC channel) and 617 nm (PE channel).
- Cylinders and glass bottles for the dilution of Wash buffer and Lysing solution

7. Storage

Store the IngoFlowEx® Kit at 2-8 °C. Expiration date is stated on reagent labels and on the IngoFlowEx® Kit box.

8. Assay procedure

8.1 Preparation of reagents before measurement

Wash buffer

Dilute 25xWash buffer with PBS to working concentration in a clean glass bottle. Allow the buffer to cool to 2-8 °C. Assign the 1xWash buffer with the date of preparation and store at 2-8 °C for up to 4 weeks. Use the information given in the table below as a quick guidance for dilution volumes.

No. of blood samples	Volume of 25xWash buffer in ml	Volume of PBS in ml	Total volume in ml
10	8	192	200
100	80	1920	2000

Lysing solution

Dilute 10xLysing solution with deionized water to working concentration in a clean glass bottle. Assign the 1xLysing solution with the date of preparation and store at 2-8 °C for up to 4 weeks. Use the information given in the table below as a quick guidance for dilution volumes.

No. of blood samples	Volume of 10xLysing solution in ml	Volume of dH ₂ O in ml	Total volume in ml
10	4	36	40
100	40	360	400

E. coli-FITC

Before use, thoroughly mix the contents of the vial using vortex.

8.2 Procedure

Cool down **1xWash buffer** to 2-8 °C.

Let the **1xLysing solution** reach the room temperature.

Place a rack for tubes into a container with ice.

Place the vials with **E. coli-FITC**, **Quenching solution** and **DNA Staining Solution** on ice.

Prepare **two tubes for each tested blood sample**:

- tube for the reaction with *E. coli* (E)
- tube for the control reaction without *E. coli* (C)

1. From each blood sample pipette 50 µl into the tube for reaction with *E. coli* (E) and 50 µl into the control tube reaction without *E. coli* (C).
2. Place the tubes on the ice and let them cool for 10 minutes.
3. Add 10 µl of the *E. coli* suspension into the tubes intended for reaction with *E. coli*, vortex and return the tubes to the ice.
4. Do not add anything into the control tubes.
5. Transfer all tubes (including the control tubes) into the thermal incubator set to 37 °C (or a water bath set to 37 °C).
6. Incubate for 30 minutes (shorten the incubation to 10 minutes if using the water bath)
7. Place the tubes back into the ice and let them cool down for 1-2 minutes.
8. Add 100 µl of precooled Quenching solution to all the tubes and mix well.
9. Add 3 ml of precooled 1xWash buffer and centrifuge the tubes for 5 min at 200×g at room temperature.
10. Remove the supernatant by aspiration into a container with an appropriate disinfectant; leave about 50 µl of residual volume.
11. Repeat wash (steps 9 and 10).
12. Resuspend the pellet in the residual volume of 1x Wash buffer.
13. Add 2 ml of 1xLysing solution (room temperature) and incubate for 20 minutes at room temperature in the dark.
14. Centrifuge for 5 min at 200×g at room temperature.
15. Remove the supernatant by aspiration into a container with appropriate disinfectant; leave about 50 µl of residual volume.
16. Wash the cells with 3 ml precooled 1xWash buffer and centrifuge for 5 min at 200×g at laboratory temperature.
17. Remove the supernatant by aspiration into a container with an appropriate disinfectant; leave about 50 µl of residual volume.
18. Resuspend in 300 µl DNA Staining Solution and incubate for 10 minutes on ice in the dark.
19. Measure with flow cytometer **within 2 hours** after adding the DNA Staining Solution.

9. Flow cytometric analysis

Setting up the cytometer before the first run: (Optional)

The color compensation for spectral overlap may be adjusted by using the assay procedure for a blood sample and fluorescent *E. coli* only (no quenching solution no DNA staining), a blood sample with Quenching solution only (no *E. coli* and no DNA staining) and a blood sample with DNA Staining Solution (no *E. coli* and no Quenching solution). Use the above mentioned measurements to adjust the cytometer settings.

Analysis of samples

Set up a gate in PE channel to acquire a sufficient number of cellular PE-bright events (10,000 events), the gate excludes debris and clumps of bacteria from further analysis (Figure 1). Visualize the events acquired in the side-scatter (SSC) versus forward-scatter (FSC) dot-plot and set gates around the granulocyte and monocyte subpopulations (Figure 2). Then bring the gated subpopulations to histograms in FITC channel. Use the control reaction histogram to place the discrimination line. The setting of the discrimination line is usually sufficient for all tested samples within the run. Calculate the percentage of positive events (phagocytic activity) (Figure 3A and 3B) and the median fluorescence intensity of positive events (phagocytic index).

10. Interpretation of results and the expected values

There are two basic evaluations: 1) the percentage of cells with engulfed *E. coli* (**phagocytic activity**) and 2) the median fluorescence intensity of cells that phagocytosed *E. coli*, which corresponds to the number of engulfed bacteria per cell (**phagocytic index**). The following table shows the average values of phagocytic activity that were calculated from data obtained by measurement of a set of blood samples from blood donors.

Subpopulation of cells	% phagocytosing cells (average +/- 2SD)
Granulocytes	91 (83-100)
Monocytes	82 (72-94)

11. Performance characteristics

Repeatability was determined by measuring four parallel reactions with a set of blood samples from blood donors and the phagocytic activity was assessed. The variability is expressed as the average coefficient of variation from these measurements.

	Range of the acquired values of phagocytic activity of tested samples (MIN-MAX)	Average CV (%)
% phagocytosing granulocytes	83.2-97.1	1.8
% phagocytosing monocytes	73.2-91.0	4.8

12. Limitations of the assay

- Reproducibility of the assay strongly depends on precise working (incubation times, temperature, pipetting).
- Flow cytometer may produce false results if the device has not been regularly calibrated and maintained appropriately.

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