CellCycleFlowEx[®] Kit

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Cat. No. / Kat. č.: ED7069

ENGLISH

1. Intended use

CellCycleFlowEx[®] Kit is intended for cell cycle analysis using flow cytometry. The kit is suitable for testing suspensions of isolated cells, such as leukocytes isolated from peripheral blood (PBMC) or cells from tissue culture.

2. Introduction

Proliferating cells duplicate their DNA during the cell cycle. The DNA is subsequently segregated equally into the daughter cells. Measurement of cellular DNA content within a proliferating cellular population recognizes the cells in S phase (actively duplicating their DNA), G2/M phase (having double the amount of DNA than resting cells) and G0/G1 phase (the same amount of DNA as resting cells) of the cell cycle. Analysis of the percentage of cells at different stages of the cycle according to their DNA content is used to assess the proliferative response to stimulation with mitogens (lymphocyte transformation test) or to monitor the effect of cell cycle inhibitors.

3. Test principle

The tested cells are washed and resuspended in a buffer to achieve a homogeneous suspension. Then the cells are fixed and permeabilized with 70% ethanol. Subsequently the cellular DNA is stained with propidium iodide, which passes through the permeabilized membrane and intercalates between the bases of nucleic acid. Propidium iodide binds to the DNA stoichiometrically, which means that the amount of bound propidium iodide is proportional to the amount of DNA within the cell. Labeling of DNA takes place in buffer with RNase which digests the contaminating cellular RNA. After laser excitation the cells emit fluorescence of which intensity is proportional to the DNA content. Flow cytometric analysis of the cell suspension distinguishes the cells in G0/G1, S and G2/M phase of the cell cycle.

4. Precautions

- For laboratory research only, not for drug, diagnostic or other use.
- Wear protective gloves.
- Do not use reagents after the expiration date.
- Avoid contamination of reagents.
- Protect the reagent Propidium Iodide from light.
- The flow cytometer should be calibrated on a routine basis using fluorescent microbeads to ensure stable sensitivity of detectors.
- Any deviation from the recommended procedure can affect the test results.

5. Reagents provided

- **ED7069-1 RNase A** – 2 vials, each contains 1 ml of RNase A solution. The amount of the reagent in each vial is sufficient for 100 tests. The solution is provided ready to use.
- **ED7069-2 Propidium Iodide** – 1 vial, contains 2 ml of propidium iodide solution. The amount of the reagent is sufficient for 200 tests. The solution is provided ready to use.
- **ED7069-3 10x Wash Buffer** – 2 bottles, each contains 30 ml of 10x concentrated solution. The amount of the reagent in each bottle is sufficient for 100 tests.

6. Necessary material not supplied

- Suitable 5 ml test tubes for cell staining (e.g. 12 × 75 mm)
- 70% ethanol at least of p.a. purity
- Ice bath / container with crushed ice
- Automatic pipettes with disposable tips
- Vortex
- Centrifuge with rotor suitable for test tubes
- Flow cytometer - blue laser excitation at 488 nm, light emission at 617 nm (PE)

7. Storage

CellCycleFlowEx[®] Kit consists of two parts:

- a) Box with 10x Wash Buffer and Propidium iodide. Store the box at 2-8 °C
- b) Resealable bag with RNase A vials. Store the bag with vials at -20 °C.

Expiration date is printed on each reagent label, on the box and the bag label.

8. Assay procedure

Preparation of reagents before the assay

1. Dilute **10x Wash Buffer** with deionized water (1 part of 10x Wash Buffer and 9 parts of deionized water) and keep it cold at 2-8 °C. Store the diluted solution (**1x Wash Buffer**) at 2-8 °C for up to 6 months.
10x Wash Buffer may contain precipitated salts. If present, place the bottle to room temperature or into a water bath set to 37°C and wait till the salts dissolve. Mix to ensure homogeneity before dilution.

2. Prepare **PI staining solution** (0,5 ml per test) according to the following formulation:

1x Wash Buffer	1 ml
Propidium iodide	0.02 ml
RNase A	0.02 ml

The solution may be stored up to 2 weeks at 2-8 °C.

3. Prepare **70% ethanol** (not provided) in deionized water and keep it cool at -15 to -30° C.

Procedure

Washing the cells

1. Centrifuge the cells at 300 g for 5 minutes at laboratory temperature.
Remove the supernatant.
2. Resuspend the cells in 1 ml of cold 1x Wash Buffer (2-8 °C).
3. Centrifuge the cells at 300 g for 5 minutes at laboratory temperature.
Remove the supernatant.

Resuspend the cells in cold 1x Wash Buffer (2-8 °C) to concentration 2.5-10 x 10⁶ cells/ml.

For one test you need 200 µl of cell suspension.

Ethanol fixation

4. Pipet 200 µl of the cell suspension to a 5 ml test tube and place it on ice.
5. Add dropwise 2 ml of cold 70% ethanol (-15 to -30 °C) while continuously mixing (vortexing) the tube.
Then immediately place the tube back to ice.

The procedure may be stopped at this step. The fixed cells can be transported or stored at 2-8 °C for up to 2 weeks.

Staining of DNA

6. Centrifuge the cells at 300 g for 5 minutes at laboratory temperature.
Remove the supernatant.
7. Resuspend the cells in 1 ml of cold 1x Wash Buffer (2-8 °C).
8. Centrifuge the cells at 300 g for 5 minutes at laboratory temperature.
Remove the supernatant.
9. Resuspend the cells in 0.5 ml of PI staining solution and incubate for 30 minutes at room temperature in the dark.
10. Thoroughly vortex the samples and measure with flow cytometer within 4 hours.

The cells in PI staining solution form aggregates in time. They need to be resuspended by vortexing prior the flow cytometric analysis.

The level of DNA labeling increases in time thus the samples measured in 30 minutes time and those

measured in 4 hour time differ in their G1 peak mean fluorescence of about 10% and the latter have higher G2/G1 fluorescence ratio.

9. Flow cytometric analysis

Analyze stained samples using flow cytometer. Acquire at least 20,000 events per sample. For precise measurements we recommend setting the flow at the lowest rate. High speed of acquisition negatively affects the resolution of the populations.

Visualize the measured data in the side-scatter (**SSC**) versus forward scatter (**FSC**) dot-plot. Set a gate around the target cell population (gate in Figure 1). Display the target population in the PE-Width (**PE-W**) versus PE-Area (**PE-A**) and draw a gate containing singlets to eliminate cell aggregates from further analysis. Display singlets in the histogram, where the X axis represents the fluorescence intensity in PE channel in linear scale (Figure 3). Set the PMT in PE channel to ensure the G1 peak is positioned in at least one third of the scale, see Figure 3.

Set borderlines between the G0/G1, S and G2/M phases of cell cycle, as shown in Figure 3. Calculate the proportion of cells in each phase. Alternatively, use suitable software (e.g., FlowJo) that can provide cell cycle analysis automatically.

10. Limitations of the assay

1. Flow cytometer may produce false results if the device has not been regularly calibrated and maintained appropriately.
2. There is a risk of misinterpretation of the data if the singlet gate and the G0/G1, S and G2/M phase borderlines were not placed correctly.

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