BrdUFlowEx[®] FITC Kit

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REF ED7073

ENGLISH

1. Manufacturer

EXBIO Praha, a.s.
Nad Safinou II 341
252 42 Vestec, Czech Republic
Tel: +420 261 090 666
Fax: +420 261 090 660
E-mail: orders@exbio.cz
www.exbio.cz

2. Intended use

The BrdUFlowEx[®] Kit is designed for labeling and subsequent immunodetection of proliferating cells by flow cytometry. The kit can be used to analyze proliferation of human or mouse leukocytes and variety of human and murine cell lines.

3. Introduction

There are in principle four techniques how to assess cellular proliferation: (I) measurement the level of DNA synthesis; (II) fluorescence dilution assays; (III) detection of unique proliferation specific cellular markers; and (IV) measurement of cell biochemical activity.

Among the mentioned methods the most reliable and accurate assay available today is to quantify the newly synthesized DNA by incorporation of labeled nucleosides or nucleoside analogues.

BrdUFlowEx[®] FITC Kit labels cellular DNA in proliferating cells with the use of thymidine analogue 5-bromo-2-deoxyuridine (BrdU). The labeling can be performed *in vitro* by adding BrdU to culture medium or *in vivo* by exposing experimental animals to BrdU via intraperitoneal injections or drinking water. Long time exposures (several hours to overnight) are used for experiments with identification of cycling versus non-cycling populations, whereas short time exposures allow analysis of cell cycle kinetics (so called pulse-chase experiments). Once BrdU is incorporated within the double stranded DNA helix, it is masked by the pairing bases and inaccessible to anti-BrdU antibodies. There are several ways to unmask incorporated BrdU. Either the cellular DNA is exposed to denaturing agents or DNA cleaving enzymes are used to remove the masking bases. Stretches of single stranded DNA with exposed BrdU are incubated with fluorescently labeled antibody. Often the test is performed as bivariate analysis of BrdU detection and DNA content. Some protocols allow simultaneous detection of surface and intracellular molecules together with BrdU detection.

4. Test principle

BrdUFlowEx[®] FITC Kit unmasking method relies on oxidative attack at the deoxyribose moiety of DNA strand by copper (I) in the presence of oxygen. The unmasking principle so called "atomic scissors" was successfully used to detect BrdU by immunocytochemistry (Ligasova et al., 2012). The atomic scissors reaction causes mild DNA damage and allows simultaneous detection of other cellular nuclear proteins that are usually lost during other BrdU unmasking methods (e.g. PCNA-1). The reaction requires considerable oxygen supply which is achieved by shaking the reaction mixture. Following the DNA scission the cells are stained with fluorescently labeled BrdU specific antibody (Liboska et al., 2012). Test can be performed as bivariate analysis of BrdU incorporation and DNA content. 7-AAD is provided to stain the DNA. It allows the user to analyze BrdU incorporation and to assess amount of cellular DNA according to the 7-AAD staining intensity (cells in G0/1, S and G2/M phase of the cell cycle).

The method is compatible with the use of additional fluorescently labeled antibodies against other cell surface or intracellular molecules together with BrdU detection. The provided protocol was optimized for test tubes and 96-well plates.

5. Precautions

- For laboratory research use only, not for drug, diagnostic or other use.
- Wear protective gloves.
- Do not use reagents after their expiration date.
- Avoid contamination of reagents.
- The flow cytometer should be calibrated on a routine basis using fluorescent microbeads to ensure stable sensitivity of detectors.
- Any deviation from the recommended procedure can affect the results.

Warning:

BrdU component (ED7073-8) contains lyophilized 5-bromo-2-deoxyuridine.

DANGER

H phrases

H340 May cause genetic defects
H361 Suspected of damaging fertility or the unborn child.

P phrases

P201 Obtain special instructions before use.

P281 Use personal protective equipment as required.
P308 + P313 IF exposed or concerned: Get medical advice/ attention.
P501 Dispose of contents/container to authorized facility for dangerous wastes.

Fix and Lysing Solution component (ED7073-1) contains diethylene glycole; methanol and formaldehyde.

DANGER

H phrases

H302+312+332 Harmful if swallowed, in contact with skin or if inhaled.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction
H319 Causes serious eye irritation.
H335 May cause respiratory irritation.
H351 Suspected of causing cancer
H371 May cause damage to organs.
H373 May cause damage to organs (kidney) through prolonged or repeated exposure if swallowed.

P phrases

P270 Do not eat, drink or smoke when using this product.
P280 Wear protective gloves / protective clothing / eye protection / face protection.
P301+P312 IF SWALLOWED: Call a POISON Center or doctor/physician if you feel unwell.
P302+P352 IF ON SKIN: Wash with plenty of soap and water.
P305+P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

6. Reagents provided

- **ED7073-1 Fix and Lysing Solution**
1 vial containing 10x concentrated fixation and erythrocyte lysing solution, the vial is intended for 100 reactions.
- **ED7073-2 Permeabilizing Solution**, 1 vial containing 30 ml of ready-to-use cellular membrane permeabilizing solution, the vial is intended for 100 reactions.
- **ED7073-3 AS Solution 1** (the abbreviation AS stands for atomic scissors), 4 vials containing lyophilized sodium ascorbate, 1 vial is intended for 25 reactions.
- **ED7073-4 AS Solution 2** (the abbreviation AS stands for atomic scissors), 1 vial containing 10 ml of ready-to-use copper (II) solution, the vial is intended for 100 reactions.
- **ED7073-5 AS Solution 3** (the abbreviation AS stands for atomic scissors), 2 vials containing 30 ml of ready-to-use buffer, one vial is intended for 50 reactions.
- **ED7073-6 Anti-BrdU FITC**, 1 vial containing 1 ml of ready-to-use solution containing FITC labeled anti-BrdU antibody (clone MoBu-1, IgG1), the vial is intended for 100 reactions.
- **ED7073-7 7-AAD**, 1 vial containing 2 ml of ready-to-use 7-aminoactinomycin D solution, the vial is intended for 100 reactions.
- **ED7073-8 BrdU**, 4 vials containing lyophilized 5-bromo-2-deoxyuridine, 1 vial is intended to prepare 2 ml of 10 mM BrdU solution.

7. Necessary material not supplied

- Deionized water (dH₂O)
- PBS (phosphate buffered saline), see recipe on page 4.
- Test tubes (12 × 75 mm) or 96-well plates
- Rack for the test tubes
- Centrifuge
- Laboratory shaker
- Automatic pipettes with disposable tips
- Vortex
- Flow cytometer - blue laser excitation 488 nm, light emission at 525 nm (FITC channel) and 675 nm (PerCP channel).
- 50 ml polypropylene tubes or glass bottles for diluting and storing the 1x Fix and Lysing solution

8. Storage and expiration

Store the BrdU FlowEx® FITC Kit at 2-8 °C. Expiration date is stated on vial labels and on the box.

9. Preparation of reagents before assay

Fix and Lysing Solution

Dilute Fix and Lysing Solution ten times with deionized water, i.e. mix 1 part of the Fix and Lysing Solution with 9 parts of deionized water (e.g. 5 ml of Fix and Lysing Solution and 45 ml of dH₂O).

Unused diluted Fix and Lysing solution may be stored up to 4 weeks either at room temperature or at 2-8 °C.

AS Solution 1

Reconstitute the content of the vial in 4 ml of PBS. Unused reconstituted AS Solution 1 may be frozen and stored up to 3 months at -20 °C or lower. The solution can be thawed and refrozen repeatedly.

BrdU

Reconstitute the content of the vial in 2 ml of dH₂O. The reconstituted BrdU stock solution is in 10mM concentration (3.07 mg/ml).

Unused reconstituted BrdU solution may be frozen and stored up to 12 months at -20 °C or lower. The solution can be thawed and refrozen repeatedly.

10. BrdU labeling and detection procedure

- **BrdU labeling:**
Expose cells to BrdU by diluting the BrdU stock solution (10 mM) into the culture medium in final concentration ranging from 10-60 μM. The exact protocol depends on the experiment design (pulse, long time labeling).
- **Harvesting the cells**
 - a. **Cells grown in dishes or flasks:**
Bring the *cells growing attached* to plastic to suspension by trypsinization. If *cells are grown in suspension* proceed directly to the next step. Centrifuge the cells at 300 g for 10 min, discard supernatant.
Wash the cells by resuspending the pellet in 10 ml of PBS (room temperature), centrifuge at 300 g for 5-10 min and discard supernatant.

Resuspend the cells at $\sim 1 \times 10^7$ cells /ml of PBS.
Dispense 0.1 ml of the cell suspension to test tubes.

b. Cells grown in 96-well plates (cells in suspension):

Centrifuge the plate with cells at 300 g for 5 minutes at room temperature, discard supernatant.

Add 200 μ l PBS (room temperature), centrifuge the plate with cells at 300 g for 5 minutes, discard supernatant.

Resuspend the pellets in 20 μ l PBS.

c. Whole blood cultures

Add 20mM EDTA/PBS solution (not provided) to the whole blood cultures to final concentration 2 mM EDTA.

Incubate for 15 min at room temperature.

Mobilize the cells by vortexing/tapping the tube/flask and proceed to centrifugation step as with the cells grown in suspension (above).

• Staining the cellular surface antigens

Add fluorescent antibodies to the test tubes/wells containing the harvested cells and mix.

Incubate for 15-30 minutes according to the specific requirements.

Wash the cells by adding 2 ml PBS to tubes / 200 μ l PBS to wells in plates.

Centrifuge at 300 g for 5 minutes and discard the supernatant.

• BrdU detection

Two following protocols are optimized either for reactions performed in **test tubes** or **96-well plates**. The protocols differ in reagent volumes used in the procedure.

BrdU detection in test tubes

1. Add 1 ml of the diluted Fix and Lysing Solution to the tubes containing 0,1 ml of cell suspension and mix (i.e. use the diluted Fix and Lysing Solution in volume that is equal to ten volumes of the resuspended cells).
2. Incubate the mixture for 30 minutes at room temperature.
3. Wash the cells by adding 2 ml of PBS, centrifuge the cells at 500 g for 5 minutes at room temperature and discard supernatant.
4. Resuspend the cells in 0.3 ml of Permeabilizing Solution.
5. Incubate for 15 minutes at room temperature.
6. Wash the cells by adding 2 ml of PBS, centrifuge the tubes at 500 g for 5 minutes at room temperature and discard supernatant.
7. Add 0.1 ml of PBS.
8. Add 0.1 ml of AS Solution 1 and mix by gentle vortexing.
9. Add 0.1 ml of AS Solution 2.
10. Place the tubes on a shaker and incubate at 400 rpm for 10 minutes at room temperature.
11. Add 0.4 ml of AS Solution 3.
12. Incubate at 400 rpm for another 10 minutes at room temperature.

13. Wash the cells by adding 2 ml of PBS, centrifuge the tubes at 500 g for 5 minutes at room temperature and discard the supernatant.
14. Resuspend the pellet in 0.1 ml of AS Solution 3.
15. Add 10 μ l of Anti-BrdU FITC.
Add other fluorescently labeled antibodies against intracellular antigens and vortex gently.
16. Incubate for 30 minutes at room temperature protected from light.
17. Wash the cells by adding 2 ml of PBS, centrifuge the tubes at 500 g for 5 minutes at room temperature and discard the supernatant.
18. Resuspend the pellet in 0.2 ml of the diluted Fix and Lysing Solution.
19. Add 20 μ l of 7-AAD to stain cellular DNA (optional).
20. Incubate for 20 minutes at room temperature in the dark.
21. Measure with flow cytometer.

BrdU detection in 96-well plates

1. Add 200 μ l of the diluted Fix and Lysing Solution to well containing the resuspended cells (i.e. use the diluted Fix and Lysing Solution in volume that is equal to ten volumes of the resuspended cells).
2. Incubate the mixture for 30 minutes at room temperature.
3. Centrifuge the cells at 500 g for 5 minutes at room temperature. Discard supernatant.
Wash by adding 200 μ l of PBS. Centrifuge the cells at 500 g for 5 minutes at room temperature and discard supernatant.
4. Add 50 μ l of Permeabilizing Solution. Mix gently by tapping the plate several times.
5. Incubate for 15 minutes at room temperature.
6. Wash by adding 200 μ l of PBS. Centrifuge the cells at 500 g for 5 minutes at room temperature and discard supernatant.
7. Add 30 μ l of PBS.
8. Add 30 μ l of AS Solution 1.
9. Add 30 μ l of AS Solution 2.
10. Place the plate on a shaker and incubate at 400 rpm for 10 minutes at room temperature.
11. Add 150 μ l of AS Solution 3.
12. Incubate at 400 rpm for 10 minutes at room temperature.
13. Centrifuge the cells at 500 g for 5 minutes at room temperature. Discard the supernatant.
Wash by adding 200 μ l of PBS. Centrifuge the cells at 500 g for 5 minutes at room temperature and discard supernatant.
14. Add 50 μ l of AS Solution 3.
15. Add 5 μ l of Anti-BrdU FITC.
Add other fluorescently labeled antibodies to intracellular antigens and mix the content by tapping the plate.
16. Incubate at 300~400 rpm for 30 minutes at room temperature protected from light.
17. Wash by adding 200 μ l of PBS. Centrifuge the cells at 500 g for 5 minutes at room temperature and discard supernatant.
18. Resuspend the pellet in 200 μ l of the diluted Fix and Lysing Solution.

19. Add 20 μ l of 7-AAD to stain cellular DNA (optional).
20. Incubate for 20 minutes at room temperature in the dark.
21. Measure with flow cytometer.

11. Flow cytometric analysis

Bivariate analysis of DNA and BrdU

1. Adjust the PMT of SSC and FSC parameters so that the cell populations are on scale.
2. Set the PerCP channel to linear scale and adjust the PMT so that cell populations are on scale.
3. Adjust the PMT in FITC channel so that the negative populations fall into the lower third of the scale.
4. Acquire 10 000 – 20 000 events per sample.
5. Draw a gate around the target cell population (as shown in Figure 1).
6. Display the gated cells in FSC-Height versus FSC-Area and draw a diagonal gate that excludes doublets and aggregates (as shown in Figure 2).
7. Display singlets in FITC vs PerCP dot plot and draw gates around the G0/G1, S and G2/M populations, as shown in Figure 3. The subG1 population sometimes seen in cellular cultures represents the apoptotic cells.
8. Calculate the proportion of cells in each subpopulation.

12. Limitations of the assay

The method relies on unmasking of the incorporated BrdU by oxidative damage to DNA. Oxidation sensitive proteins such as R-Phycoerythrin and tandem conjugates of R-Phycoerythrin are affected as well. Therefore antibodies labeled with R-Phycoerythrin or with tandem conjugates of R-Phycoerythrin should only be used after the atomic scissors reaction, simultaneously with the Anti-BrdU FITC staining (step 15 of BrdU detection protocol).

13. Appendix

Recipe for PBS (other recipes may be used as well)

Dissolve:

- 8.0 g of Sodium Chloride (NaCl)
- 0.2 g of Potassium Chloride (KCl)
- 2.0 g of Potassium Phosphate monobasic (KH₂PO₄)
- 1.42 g of Sodium Phosphate dibasic (Na₂HPO₄·2H₂O)

in 800 ml of deionized H₂O.

Adjust the pH to 7.4 with HCl.

Add dH₂O to 1 liter.

Sterilize by autoclaving for 20 minutes at 15 psi or by filter sterilization.

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
Fig. 1 Side scatter (SSC) vs forward scatter (FSC) dot-plot depicting the region (R1) that delimits the target cell population (Jurkat cells).

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
Fig. 2 FSC-Height (Peak) vs FSC-Area (Integral) dot-plot of the gated target cell population depicting the diagonal region (R2) that delimits singlets.

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
Fig. 3 Jurkat cells incubated with 10 μ M BrdU for 1 hour. FITC (BrdU signal) vs PerCP (7-AAD signal) dot-plot of the gated singlets shows BrdU incorporation in relation to the 7-AAD signal. 7-AAD signal is proportional to the cellular DNA content (G1/G0, G2/M and S phase gates).

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Liboska R, Ligasová A, Strunin D, Rosenberg I, Koberna K. Most anti-BrdU antibodies react with 2'-deoxy-5-ethynyluridine -- the method for the effective suppression of this cross-reactivity. PLoS One. 2012;7(12)

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