

## LiveReceptor mGluR1 <Endogenous mGluR1 Labeling Reagent>


Catalog NO. FDV-0018C

Research use only, not for human or animal therapeutic or diagnostic use.

### Product Background


Neurotransmitter receptors including glutamate receptors and GABA receptors etc. located on post-synapse in neuronal cells play various roles in brain functions. To understand physiological roles of neurotransmitter receptors, live cell imaging is one of the powerful approaches. Conventional imaging methods on live cells rely on a genetically engineered proteins fused with fluorescent proteins such as GFP. However, one serious problem is that the functions and movement of over-expressed neurotransmitter receptors with non-physiological tags are not precisely correlated with endogenous native receptors. The labeling methods for endogenous receptors are desirable to observe physiological functions of receptors.

**LiveReceptor** is the world first reagent series for target-specific receptor labeling. The principle of **LiveReceptor** is based on **ligand-directed acyl imidazole (LDAI)** chemistry (ref.1,2). LDAI-based chemical labelling is driven by selective ligand-protein recognition, which facilitates an acyl substitution reaction of labeling reagents on nucleophilic amino acid residues including Lys, Ser and Tyr located near ligand-binding domain. After wash out, the labelled receptors which have free ligand-binding pockets are observed on live cells. Furthermore, based on pH-dependent fluorescent property of fluorescein, fluorescent signal of labeled receptors in endocytosis pathway are highly quenched and only cell surface receptors can be observed. LiveReceptors are powerful tools to monitor reduction of cell surface receptors by endocytosis upon extracellular stimulation.



**Figure 1. Principle of ligand-directed protein labeling**

“**LiveReceptor mGluR1**” is a specific labeling reagent for cell-surface metabotropic glutamate receptor 1 (mGluR1) which is the G-protein G<sub>q</sub> family coupled receptor, is strongly expressed in Purkinje cells in the cerebellum, and plays important roles in cerebellar functions. **LiveReceptor mGluR1** has three domains including FITM as an affinity ligand for mGluR1, fluorescein and acyl imidazole (a related compound was published in ref.3). Only when FITM binds to mGluR1, nucleophilic amino acid residues (Lys, Ser or Tyr) located near ligand-binding domain on mGluR1 are attacked acyl imidazole and fluorescein is transferred into mGluR1. After removing excess reagents and resultant ligand moiety, labeled mGluR1 can be observed in both live and fixed cells. The protocol is very simple, no genetic manipulation and additional treatment are required. Because **LiveReceptor mGluR1** shows no cell membrane permeability, only cell surface mGluR1 are labelled.



**Figure 2. Chemical structure of LiveReceptor mGluR1**

---

## Description

Catalog Number: FDV-0018C

Size : 10 µg

Formulation : C<sub>51</sub>H<sub>46</sub>FN<sub>9</sub>O<sub>11</sub>S

Molecular weight : 1012.04 g/mol

Visibility : Orange lyophilized powder

Solubility : Soluble in DMSO

\*This compound has water-solubility but it can be easily degraded in water and culture medium.

Please avoid store in the water.

---

## Spectrum

Excitation/ Emission: 495/515 nm

\*Compatible with FITC filter

---

## Application

- Live cell imaging
- Immunocytochemistry with specific antibodies
- Immunoprecipitation with anti-fluorescein antibody
- Immunoblotting with anti-fluorescein antibody

## Reconstitution and Storage

Reconstitution :

Reconstitute at 0.1 mM - 1 mM in 100% DMSO. Please optimize the final concentration of DMSO depended on your experiments. Before reconstitution, please spin down to collect the orange lyophilized powder on the bottom of a tube. Carefully add DMSO into the tube and vigorously mix to completely dissolve the powder.

Storage:

(powder) Store at -20°C. Protected from light.

(solution) DMSO stock solution is stable at least for 1 year at -80°C. Please make aliquots and avoid freeze and throw. Protected from light.

---

## How to use

General procedure for mGluR1 labeling


1. Prepare 10-100 nM of LiveReceptor mGluR1 in the appropriate medium
  - \* Note : Serum-free media are highly recommended. This compound is not stable in the medium. Please prepare assay solution at time of use.
2. Replace media of cultured cells to LiveReceptor mGluR1 containing medium.
3. Culture cells with LiveReceptor mGluR1 for 1-4 hours at 17-37°C
4. After labeling, wash cells several times or perfused continuously to remove excess reagents.
5. Labelled mGluR1 can be observed.

---

## Application data


### Live cell imaging of labelled mGluR1 in acute mouse cerebellum brain slice

Acute mouse cerebellum brain slices from 3 week-old mouse treated with 10 nM LiveReceptor mGluR1 in artificial cerebrospinal fluid (ACSF) for 4 hours at RT. After then, slices were washed three times by ACSF and fluorescent signal was observed by epi-fluorescent microscopy. Strong fluorescent signal was observed in molecular layer and Purkinje cells. When the slice was pretreated with FITM, an allosteric mGluR1 inhibitor, fluorescent signal was clearly suppressed.




### Validation of mGluR1 labeling in acute mouse cerebellum brain slice by western blotting

Acute mouse cerebellum brain slices from 3 week-old mouse treated with 10 nM LiveReceptor mGluR1 in artificial cerebrospinal fluid (ACSF) for 4 hours at RT in the absence or presence of 500 nM FITM. After then, slices were washed three times by ACSF and tissue lysate was prepared. The lysate was analyzed western blotting using anti-fluorescein (FL) or anti-mGluR1 antibodies. Two bands were observed by anti-FL antibody and these bands were well corresponded with anti-mGluR1. FITM-treatment clearly reduced the signal of anti-FL antibody.



## Live cell imaging of labelled mGluR1 in mGluR1-expressed HEK293

mGluR1-expressed HEK293 cells were treated with 100 nM of LiveReceptor mGluR1 in the absence or presence of 2  $\mu$ M FITM for 4 hour at 17°C and washed out two times with the basal medium. FITM-treatment clearly reduced the green fluorescent signal. (scale bars, 10  $\mu$ m).



## Reference

1. Fujishima *et al.*, *J. Am. Chem. Soc.*, **134**, 3961-3964 (2012). Ligand-directed acyl imidazole chemistry for labeling of membrane-bound proteins on live cells.
2. Miki *et al.*, *Chem. Biol.*, **21**, 1013-1022 (2014). LDAI-based Chemical Labeling of Intact Membrane Proteins and its Pulse-Chase Analysis under Live Cell Conditions.
3. Nonaka *et al.*, *bioRxiv.*, <https://doi.org/10.1101/2023.01.16.524180>. Bioorthogonal chemical labelling of endogenous neurotransmitter receptors in living mouse brain.

NOTE: Ref. 3 showed a similar mGluR1 labeling reagent called CmGlu1M which has CNITM, an analog of FITM, and Alexa Fluor 647 dye for *in vivo* imaging.

### Disclaimer/免責事項

This product has been commercialized by Funakoshi Co., Ltd. based on the results of academic research, and the advertisement text, figures and manuals (hereinafter “Product information”) have been prepared based on published research reports on June, 2023. The academic interpretation at the time of creation of the Product Information may change in accordance with future developments in the relevant research field and expansion of various scientific findings, and the latest version and certainty of the Product Information are not guaranteed. The specifications of this product and the Product Information are subject to change without notice. Please contact us for the latest information.

本製品は学術研究成果を基にフナコシ株式会社が製品化したもので、2023年6月時点における公開研究報告を基に広告文章およびマニュアル(以下、製品資料)を作成しています。今後の当該研究分野の発展および各種学術知見の拡大にともない、製品資料作成時の学術的解釈が変更になる可能性があり、最新性・確実性を保証するものではありません。また、本製品の仕様および製品資料を予告なく変更する場合がございます。最新の情報に関しましては、弊社までご確認いただけますようお願い申し上げます。

 E-mail Newsletter Sign Up

Japanese 

English 

## Related product

### LiveReceptor AMPAR <Endogenous AMPAR Labeling Reagent>

LiveReceptor AMPAR is a specific labeling reagent for AMPA-type glutamate receptor, AMPAR. Live imaging of cultured neuron and slice tissues were validated.


Catalog No. FDV-0018A

Size 10  $\mu$ g

#### Data examples

- Live cell imaging of labelled endogenous AMPARs in cultured neurons

Cultured hippocampal neurons were treated with 1  $\mu$ M of LiveReceptor AMPAR (in left) or Fluorescein-conjugated PFQX as negative control (in right) for 1 hour at 17°C and washed out three times with the basal medium. Dendritic spin-like punctual structures were observed on live cells by specifically LiveReceptor. (scal bars, 10  $\mu$ m)



### LiveReceptor GABA<sub>A</sub>R <GABA<sub>A</sub>R Labeling Reagent>

LiveReceptor GABA<sub>A</sub>R is a specific labeling reagent for ion channel-type GABA receptor, GABA<sub>A</sub>R.


Catalog No. FDV-0018B

Size 10  $\mu$ g

#### Data example:

Live cell imaging of labelled GABA<sub>A</sub>Rs in GABA<sub>A</sub>R-expressed HEK293

GABA<sub>A</sub>R ( $\alpha$ 1/ $\beta$ 3/ $\gamma$ 2)-expressed HEK293 cells were treated with 1  $\mu$ M of LiveReceptor GABA<sub>A</sub>R in the absence or presence of 100  $\mu$ M gabazine, a GABA<sub>A</sub>R selective inhibitor, for 3 hour and washed out three times with the basal medium. (scale bars, 20  $\mu$ m)



## Featured product

### LipiORDER <Membrane Lipid Order Imaging Dye>


LipiORDER is a solvatochromic dye for membrane lipid order imaging. LipiORDER exhibits green fluorescence with Lo phase and exhibits red fluorescence with Ld phase. The ratiometric analysis ( $F_{red}/F_{green}$ ) enables the quantitative visualization of membrane lipid order.

Catalog No. FDV-0041

Size 0.1 mg

#### Features

- Recommended Ex/Em: ~405 nm / 500-550 nm (Green channel) and 550-650 nm (Red channel)
- Enable to quantitatively monitor lipid order on plasma and inner membranes in live cells
- Highly photostable and cellularly stable compared with similar conventional dyes.



### LipiDye II <Lipid Droplet Live Imaging>


LipiDye II is a highly sensitive lipid droplet staining dye with extremely photostable property. This dye is the second generation of our previous reagent, LipiDye. This dye allows us to detect small lipid droplets (<1  $\mu\text{m}$ ) in non-adipocytes and to apply into long-term live cell imaging for dynamic lipid droplet movements.

Catalog No. FDV-0027

Size 0.1 mg

#### Features

- Recommended Ex/Em: 400-500 nm / 490-550 nm
- Enable to detect <1  $\mu\text{m}$  lipid droplets
- Suitable for long-term live cell imaging
- Extremely photostable compared with conventional dyes
- Compatible with both live and fixed cells



### FAOBlue <Fatty Acid Oxidation Detection Reagent>


FAOBlue is a cell-based fatty acid beta-oxidation (FAO) detection dye which emits blue fluorescence upon cellular FAO activity.

Catalog No. FDV-0033

Size 0.2 mg

#### Features

- Ex/Em: ~405 nm / 460 nm
- Enable to directly detect cellular FAO activity in live cells
- Apply quantitative comparison of FAO activity between different cell types
- Can monitor the drug-induced change of FAO activity



 **funakoshi**  
FRONTIERS IN LIFE SCIENCE

URL: <http://funakoshi.co.jp>  
9-7 Hongo 2-Chome, Bunkyo-ku, Tokyo 113-0033