

Data Sheet

IBA-lyse Bacterial Lysis Buffer

Cat. No.: 2-1017-050; 2-1017-250

Lot No.: 1017-

IBA Headquarters IBA GmbH

Rudolf-Wissell-Str. 28 37079 Goettingen Germany Tel. +49 (0) 551-5 06 72-0 Fax +49 (0) 551-5 06 72-181

IBA US Distribution Center

1328 Ashby Road Olivette, MO 63132 USA Tel. 1-877-IBA-GmbH (1-877-422-4624) Fax 1-888-531-6813

E-mail: info@iba-lifesciences.com http://www.iba-lifesciences.com

Last date of revision March 2012 Version 1017-8

| Description | IBA-lyse reagents are formulated for gentle disruption of <i>E. coli</i> and release of proteins. It consists of a reagent mix compatible with subsequent <i>Strep</i> -tag purification and provides a simple, rapid, inexpensive and most importantly more reproducible and milder alternative to mechanical methods such as sonication or French Press for preparation of cleared lysates to be submitted to affinity chromatography. |
|--------------------|--|
| Kit content | Buffer B which contains 1 mM EDTA, Tween 20 and avidin; store at 2-8 °C. (Cat. No. 2-1017-050: 50 ml; Cat. No. 2-1017-250: 2 x 125 ml) Lyophilized lysozyme; store at -20 °C. Lysozyme reconstitution buffer; store at 2-8 °C. (Cat. No. 2-1017-050: 0.1 ml; Cat. No. 2-1017-250: 0.5 ml) Lyophilized DNase I; store at -20 °C. DNase reconstitution buffer; store at 2-8 °C. (Cat. No. 2-1017-050: 0.1 ml; Cat. No. 2-1017-250: 0.5 ml) |
| Prior to first use | Dissolve provided lysozyme and DNase with the respective reconstitution buffer, store resulting stock solutions at -20°C. Determine how much cells are intended to be lysed in the next 3 months (2 ml Buffer B containing lysozyme will be needed for 100 ml bacterial culture). Add 2 µl lysozyme stock solution per ml Buffer B to the corresponding amount of Buffer B. Store Buffer B containing lysozyme at 2-8 °C for up to 3 months. DO NOT ADD DNase TO THE BUFFER! DNase is added after the lysis of the bacterial cells. |
| Stability | 12 months after shipping |
| Storage | Enzymes (reconstituted or lyophilized) at -20°C, buffers are stable at 2-8 °C. Store Buffer B containing lysozyme at 2-8 °C at which it will be stable for 3 months. |
| Shipping | Room temperature |
| Remarks | Toxicity not reported. Could provoke allergic skin reactions. |