

## Data Sheet

### pASK-IBA4

Cat. No.: 2-1303-000

Lot No.: 1303-

Last date of revision  
March 2012  
Version 1303-8

<b>Description</b>	Expression plasmid. The expression cassette is under transcriptional control of the tetracycline promoter/operator. The expressed recombinant protein will be secreted into the periplasm.
<b>Affinity tag</b>	<i>Strep</i> -Tactin <sup>®</sup> affinity tag ( <i>Strep</i> -tag II <sup>®</sup> ) for the purification of recombinant protein. The affinity tag is fused to the N-terminus of the recombinant protein.
<b>Secretion</b>	The ompA signal sequence directs the expressed protein into the periplasmic space and will be cleaved off during the translocation process
<b>Bacterial Expression</b>	Expression is induced upon addition of 200 µg anhydrotetracycline (order no.: 2-0401-001; 2-0401-002) per 1 liter <i>E. coli</i> shaking culture ( $A_{550} = 0.5$ ).
<b>Expression strain</b>	Any <i>E. coli</i> strain. The <i>tet</i> -promoter works independently from the genetic background of <i>E. coli</i> .
<b>Resistance</b>	Ampicillin
<b>Form</b>	5 µg, dissolved in 10 mM Tris/HCl pH 8.0, 1 mM EDTA; 20 µl
<b>Concentration</b>	250 ng/µl
<b>Storage</b>	4 °C for frequent usage, -20 °C for long-term storage

#### For research use only

##### Important licensing information

This product is based on *Strep*-tag and *tet* promoter technologies covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com) or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

##### Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

# Multiple Cloning Site of pASK-IBA4

1 CCATCGAATGGCCAGATGATTAATTCCTAATTTTGTGACACTCTATCATTGATAGAGTTATTTTACCACTCCCTATCA 80  
forward primer

M K K T A I A

81 GTGATAGAGAAAAGTGAATGAATAGTTTCGACAAAAATCTAGATAACGAGGGCAAAAAATGAAAAAGACAGCTATCGCGA 160  
XbaI

OmpA link Strep-tag

I A V A L A G F A T V A Q A A S W S H P Q F E K G A E

161 TTGCAGTGGCACTGGCTGGTTTCGCTACCGTAGCGCAGGCCGCTAGCTGGAGCCACCCGAGTTCGAAAAAGgcgCGAG 240  
NheI BbeI BsaI  
EheI PshAI  
KasI  
NarI

T A V P N S S S V P G D P S R S T C R G T M V S D I \*


241 ACCGCGGTCCCGAATTCGAGCTCGGTACCCGGGGATCCCTCGAGGTCGACCTGCAGGGGGACCATGGTCTCTgataTCTA 320  
SacII EcoRI KpnI BamHI SalI PstI BsmFI BsaI EcoRV  
BsmFI SstI SmaI XhoI PshAI  
NcoI

321 ACTAAGCTTGACCTGTGAAGTGAAAAATGGCGCACATTGTGCGACATTTTTTTGTCTGCCGTTTACCGCTACTGCGTCA 400  
HindIII reverse primer

**Please note:** Restriction enzymes in bold cut twice. The *BsaI* sites (isoschizomer of *Eco31I*) at each end of the multiple cloning site are useful for precise and oriented insertion of the recombinant gene by one cleavage reaction only. The “link” contains a restriction site which can be used e.g. for subcloning the recombinant gene into pEXPR-IBA vectors for mammalian expression. During secretion of the recombinant protein into the periplasmic space, the OmpA signal sequence will be cleaved off. The processed protein will start with the Ala-Ser-linker.

## Features of pASK-IBA4

	from bp	to bp
promoter	37	72
forward primer binding site	57	76
OmpA signal sequence	139	201
Strep-tag	202	231
multiple cloning site	231	313
reverse primer binding site	381	397
f1 origin	410	848
AmpR resistance gene	997	1857
tet-repressor	1867	2490
Col E1origin	2643	3231



### Cloning primers for the precise cloning using *BsaI* or *Eco31I*

Forward: 5'- NNNNNNGGTCTCNGC GCC <sup>(N<sub>20</sub>)</sup> NNN NNN...  
Reverse: 5'- NNNNNNGGTCTCNTA TCA NNN NNN...  
(N<sub>20</sub>)

### Sequencing primers:

Forward: 5'- GAGTTATTTTACCACTCCCT -3'  
Reverse: 5'- CGCAGTAGCGGTAAACG -3'