

Data Sheet

pASK-IBA44

Cat. No.: 2-1344-000

Lot No.: 1344 -

Last date of revision
March 2012
Version 1344-9

Description	Expression plasmid. The expression cassette is under transcriptional control of the tetracycline promoter/operator. The expressed recombinant protein will be secreted into the periplasm.
Affinity tag	The recombinant protein will contain two affinity tags: <ol style="list-style-type: none"> 1) <i>Strep-Tactin</i>[®] affinity tag (<i>Strep-tag II</i>[®]) for the purification of recombinant protein via <i>Strep-Tactin</i> resins. The <i>Strep-tag</i> is fused to the N-terminus of the recombinant protein. 2) 6xHistidine-tag for the purification of recombinant protein via Ni-NTA resins. The 6xHistidine-tag is fused to the C-terminus of the recombinant protein.
Secretion	The ompA signal sequence directs the expressed protein into the periplasmic space and will be cleaved off during the translocation process
Bacterial Expression	Expression is induced upon addition of 200 µg anhydrotetracycline (order no.: 2-0401-001; 2-0401-002) per 1 liter <i>E. coli</i> shaking culture ($A_{550} = 0.5$).
Expression strain	Any <i>E. coli</i> strain. The <i>tet</i> -promoter works independently from the genetic background of <i>E. coli</i> .
Resistance	Ampicillin
Form	5 µg, dissolved in 10 mM Tris/HCl pH 8.0, 1 mM EDTA; 20 µl
Concentration	250 ng/µl
Storage	4 °C for frequent usage, -20 °C for long-term storage

For research use only

Important licensing information

This product is based on *Strep-tag*, 6xHistidine-tag and tet promoter technologies covered by intellectual property (IP) rights and on completion of the sale IBA grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.


Multiple Cloning Site of pASK-IBA44

1	CCATCGAATGGCCAGATGATTAATTCCTAATTTTGTGACACTCTATCATTGATAGAGTTATTTTACCACTCCCTATCA	80
	forward primer	
	M K K T A I A	
81	GTGATAGAGAAAAGTGAAATGAATAGTTCGACAAAAATCTAGATAACGAGGGCAAAAATGAAAAAGACAGCTATCGCGA	160
	XbaI	
	OmpA link Strep-tag link	
	I A V A L A G F A T V A Q A A S W S H P Q F E K G A E	
161	TTGCAGTGGCACTGGCTGGTTTCGCTACCGTAGCGCAGGCCGCTAGCTGGAGCCACCCGCAGTTCGAAAAAGgcccCGAG	240
	NheI BbeI BsaI	
	EheI PshAI	
	KasI	
	NarI	
	T A V P N S S S V P G D P S R S T C R G T M V S G L R	
241	ACCGCGGTCCCGAATTCGAGCTCGGTACCGGGGATCCCTCGAGGTCGACCTGCAGGGGGACCATGGTCTCAGgcccTGAG	320
	SacII EcoRI KpnI BamHI SalI PstI BsmFI BsaI StuI	
	BsmFI SstI SmaI XhoI PshAI	
	NcoI	
	6xHistidine-tag	
	G S H H H H H H *	
321	AGGATCGCATCACCATCACCATCACTAATAAGCTTGACCTGTGAAGTGAAAAATGGCGCACATTGTGCGACATTTTTTTT	400
	HindIII	
401	GTCTGCCGTTTACCGCTACTGCGTCACGGATCTCCACGCGCCCTGTAGCGGCGCATTAAAGCGCGCGGGTGTGGTGGTTA	480
	reverse primer	

Please note: Restriction enzymes in bold cut twice. The *BsaI* sites (isoschizomer of *Eco31I*) at each end of the multiple cloning site are useful for precise and oriented insertion of the recombinant gene by one cleavage reaction only. The “link” contains a restriction site which can be used e.g. for subcloning the recombinant gene into pEXPR-IBA vectors for mammalian expression. During secretion of the recombinant protein into the periplasmic space, the OmpA signal sequence will be cleaved off. The processed protein will start with the Ala-Ser-linker.

Features of pASK-IBA44

	from bp	to bp
promoter	37	72
forward primer binding site	57	76
OmpA signal sequence	139	201
Strep-tag	202	231
multiple cloning site	232	318
6xHistidine-tag	319	348
reverse primer binding site	407	423
f1 origin	436	874
AmpR resistance gene	1023	1883
tet-repressor	1893	2516
Col E1origin	2669	3257



Cloning primers for the precise cloning using *BsaI* or *Eco31I*

Forward: 5'- NNNNNNGGTCTCNGC GCC ^(N₂₀) NNN NNN...
Reverse: 5'- NNNNNNGGTCTCNG GCC ^(N₂₀) NNN NNN...

Sequencing primers:

Forward: 5'- GAGTTATTTTACCACTCCCT -3'
Reverse: 5'- CGCAGTAGCGGTAAACG -3'