

CircAID-p-seq for Oxford Nanopore Technologies

| Product | Catalog no | Rxns |
|---------------|------------|------|
| CircAID-p-seq | #CA001 | 6 |

Shipping: Dry ice and 4°C

Storage Conditions: store components as indicated on data sheet

Shelf Life: 5 months

Description: CircAID-p-seq kit (**Circular Amplification and IDentification** of short RNA **sequences** bearing a 3 **Phosphate**) is designed for quick (1-day) high quality library preparation for short RNAs (20-50 nt) bearing a 3' -phosphate/2',3'-cyclic phosphate (3'-P/2',3'-cP) end. The protocol, suitable for the characterization of cP-forming endoribonucleases, is also applicable to ribosome profiling experiments and transcriptome analysis.

CircAID-p-seq is suitable for the Oxford Nanopore platform (Direct cDNA Sequencing Kit).

CircAID-p-seq for Oxford Nanopore Technologies



Reagents provided

| Product (label) | Cap Color | Cat. no. | Store condition | Quantity |
|--------------------------|-----------|---------------|----------------------------|----------------------|
| CircAID-p-seq kit | | #CA001 | according to manual | 1kit - 6 rxns |
| Buffer PK (BPK) | Red | #CA001-1 | -20°C | 40 µL |
| PK enzyme (PK) | Red | #CA001-2 | -20°C | 7 µL |
| ATP 10 mM | Red | #CA001-3 | -20°C | 45 µL |
| Buffer A (BA) | Blue | #CA001-4 | -20°C | 21 µL |
| Enzyme Mix A (mix A) | Blue | #CA001-5 | -20°C | 21 µL |
| MnCl ₂ | Blue | #CA001-6 | -20°C | 21 µL |
| GTP 1 mM | Blue | #CA001-7 | -20°C | 21 µL |
| Linker R™(R) | Blue | #CA001-8 | -80°C | 10.5 µL |
| Buffer B (BB) | yellow | #CA001-9 | -20°C | 14 µL |
| Enzyme Mix B (mix B) | yellow | #CA001-10 | -20°C | 7 µL |
| PEG 8000 | yellow | #CA001-11 | -20°C | 56 µL |
| Nuclease (N) | yellow | #CA001-12 | -20°C | 14 µL |
| Buffer N (BN) | yellow | #CA001-13 | -20°C | 18 µL |
| P1 oligo (P1) 20 µM | Green | #CA001-14 | -20°C | 19 µL |
| Buffer RT (BRT) | Green | #CA001-15 | -20°C | 14 µL |
| RT enzyme (RT) | Green | #CA001-16 | -20°C | 7 µL |
| dNTPs 10 mM | Green | #CA001-17 | -20°C | 7 µL |
| HI solution (HI) | Green | #CA001-18 | -20°C | 16 µL |
| Enhanced Buffer (EnB) | Green | #CA001-19 | RT | 50 µL |
| Buffer Taq (BT) | Clear | #CA001-20 | -20°C | 3.5 µL |
| Taq | Clear | #CA001-21 | -20°C | 70 µL |
| P2 oligo (P2) 20 µM | Clear | #CA001-22 | -20°C | 18 µL |

Shelf life: 5 months from the delivery date

CircAID-p-seq for Oxford Nanopore Technologies

kit code number.....

Number of samples (N).....



IMMAGINA[®]
• BioTECHNOLOGY •
The Relosome Company


Reagents and equipment to be supplied by user:

- RNA Clean & Concentrator™-5 (Zymo catalog. no. R1015 & R1016)
- Sodium Acetate 3M
- Nuclease-free water
- GlycoBlue (Ambion catalog no. AM9515)
- Isopropanol (Sigma catalog no. 278475)
- Microcentrifuge and nonstick RNase-free microfuge tubes (0.2 mL and 1.5 mL)
- Automatic wheel (rotator mixer)
- Magnetic stand for 1.5mL tube
- Qubit™ RNA HS Assay Kit (Thermo Fisher catalog. no. Q32852)
- Direct-cDNA Sequencing kit (SQK-DCS109)

Work always in an RNase-free environment!



Workflow



Sample Recommendations

Input RNA amount: it is recommended to start the protocol with at least 30 pmol of RNA

Input RNA quality: high RNA purity and integrity is recommended in order to ensure optimal downstream results

Step 1

5' phosphorylation

This step is essential when starting with short RNAs bearing 5'-OH ends. If RNA inputs already harbour 5'-P ends, the step can be omitted.

- **Step1:** Mix the following reagents in a 0.2 mL nuclease-free PCR tube:

| | |
|-------------------|-------------|
| Buffer PK | 5µL |
| 10mM ATP | 5 µL |
| PK | 1 µL |
| RNA 3'-P/2'3' -cP | 30-60 pmol |
| H ₂ O | Up to 50 µL |

- Incubate the reaction for 1h at 37 °C in a thermal cycler.
- Purify the reaction through the RNA Clean & Concentrator™-5 kit, following the protocol for small RNAs and performing the final elution in a volume of 6 µL of nuclease-free water.

- **Step2:** Mix the following reagents in a 0.2 mL nuclease-free PCR tube, scaling reagent volumes according to RNA input amounts as reported below:

| RNA (from step 1) | 30 pmol | 45 pmol | 60 pmol |
|-------------------|-------------|-------------|-------------|
| Buffer A | 1 µL | 1.5 µL | 2 µL |
| GTP 1mM | 1 µL | 1.5 µL | 2 µL |
| MnCl ₂ | 1 µL | 1.5 µL | 2 µL |
| Enzyme Mix A | 1 µL | 1.5 µL | 2 µL |
| Linker R™ | 0.5 µL | 1 µL | 1.5 µL |
| H ₂ O | Up to 10 µL | Up to 15 µL | Up to 20 µL |

- Incubate the reaction for 2h at 37 °C in a thermal cycler.
- Purify the reaction through the RNA Clean & Concentrator™-5 kit, following the protocol for small RNAs and performing the final elution in a volume of 8 µL of nuclease free water.

Step 3-4

Circularization

- **Step3:** Prepare 1 mM ATP by diluting the 10 mM stock in nuclease-free water. Assemble the following reaction in a 0.2 mL nuclease-free PCR tube:

| | |
|-------------------|------|
| Buffer A | 2 µL |
| ATP 1 mM | 1 µL |
| PEG8000 | 8 µL |
| Enzyme Mix B | 1 µL |
| RNA (from step 2) | 8 µL |

- Incubate the reaction for 2h at 25 °C in a thermal cycler.
- **Step 4:** at the end of the incubation, add the following reagents to the reaction mix:
 - 2.5 µL Buffer N
 - 2 µL Nuclease
 - H₂O up to 25 µL

- Incubate for 1 h at 37 °C in a thermal cycler.

Purify the reaction through RNA Clean & Concentrator™-5 kit, following the protocol for small RNAs and performing the final elution in a volume of 10 µL of nuclease free water. **OPTIONAL STOPPING POINT** (store at -80°C).

Step 5
Reverse Transcription

- **Recommended:** measure sample concentration using Qubit fluorimeter (Qubit™ RNA HS Assay Kit). Ideal input amounts for the following steps range between 10-100 ng.

- **Step 5:** For the generation of multimeric single-stranded cDNA, combine the following reagents:

| | |
|----------------------------|-------------|
| dNTPs 10 mM | 1 µL |
| Circular RNA (from step 4) | 10-100 ng |
| P1 | 2.5 µL |
| H ₂ O | Up to 13 µL |

- Heat the circular RNA-primer mix at 65°C for 5 minutes, then incubate on ice for at least 1 minute.
- Add the following reagents to the annealed RNA:

| | |
|-----------------|------|
| Buffer RT | 2 µL |
| Enhanced Buffer | 2 µL |
| RT enzyme | 1 µL |

- Incubate 4 h at 42 °C, then add 2.2 µL of HI and heat the mix for 20 min at 70 °C.
- Transfer the reaction to a new 1.5 mL tube.
- Add 156 µL nuclease-free water, 20 µL sodium acetate (3M), 300 µL isopropanol and 2 µL Glycoblue.
- Store at -80°C for at least 2 hours.
- Pellet the RNA by centrifugation (20000 g) for 30 min.
- Resuspend the pellet in 20 µL of nuclease-free water.

Step 6
Second strand synthesis

- **Step 6:** Set up the following PCR reaction in a 0.2 mL nuclease-free PCR tube:

| | |
|---------------------|-------------|
| Buffer T | 5 µL |
| dNTPs10mM | 1 µL |
| P2 | 2.5 µL |
| cDNA (from step 5) | 20 µL |
| Taq | 0.3 µL |
| Nuclease free water | Up to 50 µL |

Cycling conditions:

| Step | Temperature | Time |
|-----------------------------|-------------|---------|
| Initial denaturation | 94°C | 3 min |
| 1 Cycle | 94°C | 30 secs |
| | 51°C | 30secs |
| | 68°C | 1 min |
| Hold | 4°C | |

- Purify the reaction by adding 45µL of resuspended Agencourt AMPure XP beads and mix by flicking the tube.
- Incubate on a rotator mixer for 5 minutes at RT.
- Prepare 500 µL of fresh 70% ethanol in nuclease-free water.

- | ○ Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.
- | ○ Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.
- | ○ Repeat the previous step.
- | ○ Spin down and place the tube back on the magnet.
- | ○ Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.
- | ○ Remove the tube from the magnetic rack and resuspend pellet in 25 µl nuclease-free water. Incubate on a rotator mixer for 10 minutes at RT.
- | ○ Pellet beads on magnet until the eluate is clear and colourless.
- | ○ Remove and retain 25 µl of eluate into a clean nuclease-free 1.5 ml tube.

Step 7
ONT Library preparation

- | ○ **Step 7:** use the purified double-stranded cDNA for ONT library preparation, following the protocol Direct-cDNA Sequencing kit (SQK-DCS109), starting from **End Prep Step**.