Assay Performance Characteristics:

Standard range: 50-0.1ng/mL Limit of Detection: 0.2ng/mL Background: OD<0.05 at 450nm Coefficient of Determination: R-squared>0.98

Plate Template:

	1	2	3	4	5	6	7	8	9	10	11	12
А												
В												
С												
D												
Е												
F												
G												
Н												

References:

1. Shanti KN, Martin BM, Nagpal S, Metcalfe DD, Subba Rao PV. Identification of tropomyosin as the major shrimp allergen and characterization of its IgE-binding epitopes. J Immunol. 1993; 151 (10):5354-5363.

2. Daul CB, Slattery M, Reese G, Lehrer SB. Identification of the major brown shrimp (Penaeus aztecus) allergen as the muscle protein tropomyosin. Int Arch Allergy Immunol. 1994; 105(1):49-55.

y in



Indoor Biotechnologies, Inc. 700 Harris Street, Charlottesville, VA 22903 United States

Tel: (434) 984-2304 Fax:(434) 984-2709 Email: mail@inbio.com Indoor Biotechnologies Ltd. Vision Court Caxton Place Cardiff, Wales CF23 8HA United Kingdom

Tel: +44 (0) 29 2167 4640 Email: Info@indoorbiotech.co.uk

ail: Info@indoorbiotech.co.uk Email: info@inbioindia.com

RG

Indoor Biotechnologies,

Bangalore-560100, India

Tel: +91-9901722009

Bangalore Bioinnovation Centre.

BioTech Park, Electronic City Phase 1

India Private Limited

INDOOR and Variate Registered Trademarks of INDOOR Biotechnologies, Ltd., Manchester U.K.

© 2002 INDOOR biotechnologies LTD




Shrimp Tropomyosin ELISA 2.0

Pre-coated Plate Kit

Product Code: EPC-TPM-X Lot Number: xxxxx

Sample curve:



Contents:

Microtiter plate coated with anti-tropomyosin monoclonal antibody 1A6

Tropomyosin allergen standard (white cap) Concentration: 500ng/mL

Rabbit anti Shrimp Tropomyosin antiserum (brown cap)

Peroxidase-conjugated goat anti-rabbit IgG (blue cap)

Wash buffer (10x concentrate) Assay buffer (10x concentrate) TMB developing substrate Stop solution (0.5N sulfuric acid)

Store kit at 2-8°C Expiry: Xxx xx, xxxx

	Certificate of Analysis						
Pre-coated Plate:	96-well polystyrene microtiter plate coated with monoclonal antibody 1A6 and treated with stabilizing agent. Sealed in foil pouch with desiccant.						
Monoclonal Antibody: Immunogen: Isotype: Specificity:	 1A6 Mite (<i>D. pteronyssinus</i>) extract Mouse IgG1 Binds to specific epitope present on <i>D. pteronyssinus</i> tropomyosin allergen, Der p 10. Cross reactive with shellfish tropomyosin. Produced in ascites and purified by chromatography using Protein A. Single heavy and light chain bands on SDS-PAGE. 						
Purification:							
Lot Number:	XXXXX						
Detection Antibody:	Rabbit polyclonal antiserum						
Immunogen: Specificity: Activity:	Purified natural Shrimp Tropomyosin The pAb contains IgG antibodies to shellfish tropomyosin. Titrated for use in ELISA at 1/1000 dilution. Prepared in 1% BSA/50% glycerol/PBS, pH 7.4, 0.22 µm filtered,						
Lot Number:	preservative free. XXXX						
Allergen Standard:	Purified natural Shrimp Tropomyosin prepared in 1% BSA/50% glycerol/PBS, pH 7.4.						
Concentration: Lot Number:	500ng/mL (based on amino acid analysis) XXXXX						

Materials required, but not provided:

- Type I ultrapure water or 18.2MΩ de-ionized water
- Volumetric measuring equipment (e.g. serological pipette, graduated cylinder)
- Clean containers for buffer and reagent preparation
- Calibrated single and multi-channel micropipettes and tips
- Vortex mixer
- Plate reader capable of reading absorbance at 450nm
- Analysis software (recommended, but not required)

Protocol

Please read the entire protocol before starting the assay Bring all reagents to room temperature before use

- Prepare 1x working dilutions of the 10x wash and assay buffers in clean containers using 18.2MΩ de-ionized water or Type I ultrapure water. For one plate: Wash buffer: add 15mL concentrate to 135mL water Assay buffer: add 2.5mL concentrate to 22.5mL water Adjust volumes accordingly for multi-plate assays.
 *Diluted buffers may be stored at 4^oC for up to 1 week
- Remove the plate from the foil pouch and wash by adding 150µL wash buffer to each well. Empty the wells by inverting the plate and then tap on absorbent paper to remove residual buffer. Repeat the wash cycle two more times.
- 3. Add standards, samples, and blanks to the plate (final volume in all wells is 100µL).

Standards: add 180μ L assay buffer into wells A1 and B1, and 100μ L into remaining wells of rows A and B. Vortex the Can f 1 standard and add 20μ L to wells A1 and B1. Mix well by pipetting up and down 7-10 times and then transfer 100μ L into wells A2 and B2. Mix well and continue the serial doubling dilution scheme across the plate to column 10.

The assay buffer in wells A11, B11 and A12, B12 will serve as **Blanks**. **Samples:** dust extracts are routinely tested starting at 1/10 dilution and can be prepared directly on the pre-coated plate: add 20μ L sample to 180μ L assay buffer. Mix, then transfer 100μ L into 100μ L assay buffer in the next well. Continue across the plate for the desired number of dilutions. A minimum of three dilutions per sample should be tested; 6-12 dilutions are recommended.

Air filter extracts, allergen extracts, and other types of samples may require a different dilution scheme.

*Sample dilutions may also be prepared in tubes or on a 96-well dilution plate and transferred to the pre-coated plate.

- 4. Incubate the plate at room temperature (away from direct sunlight) for 1 hour.
- Wash the plate 3x with 150µL wash buffer per well. Vortex the polyclonal antibody and prepare a 1:1,000 detection antibody/conjugate mix by adding 10µL polyclonal antibody and 10µL peroxidase-conjugated goat anti-rabbit IgG to 10mL assay buffer.

Mix thoroughly and add 100µL to each well.

- 6. Incubate the plate at room temperature (away from direct sunlight) for 1 hour.
- 7. Pour the TMB substrate and stop solution into separate basins so they are ready to use in the next step. Wash the plate 3x with 150μL wash buffer per well.
- Use a <u>multi-channel</u> pipette to add 100µL TMB to each well and monitor the reaction as the blue color develops. Once OD450 reaches 0.08 for Standard 1, use a <u>multi-channel</u> pipette to add 50µL stop solution to each well (the color will change to yellow).
- 9. Read the plate at 450nm. The OD for Standard 1 should be between 1.2 and 3.5, with an ideal range of 2.0 2.5.

A list of frequently asked questions and troubleshooting guide can be found under the 'Support' tab on our web site: www.inbio.com.