

# Necrosis vs Apoptosis Assay

Catalog #9147, #9148

**FOR RESEARCH USE ONLY.**

Not for use in diagnostic procedures.


## 1. INTRODUCTION

Assessment of cellular cytotoxicity levels associated with cytolytic activity of T lymphocytes or natural killer cells is an important aspect of immunobiology related research<sup>1</sup>. Early assays developed to assess cytolytic activity utilized the release of <sup>51</sup>Cr from membrane compromised target cells, which were passively loaded with this radioactive indicator prior to exposure to the cytotoxic agent or cells<sup>2</sup>. A major draw-back to these <sup>51</sup>Cr release assays is their reliance on the use of a radioactive isotope with its associated disposal and safe handling issues. Additionally, chromium uptake methods suffer from variabilities in target cell preloading inconsistencies as well the tendency to spontaneously release the chromium isotope into the media in the absence of any cytotoxicity stimulus<sup>3</sup>.

Assessing potential cytotoxicity properties of chemical and biological agents is a mandatory requirement for the safe distribution of pharmaceuticals, vaccines, or additives associated with food product formulations. Early identification of unintended drug, vaccine, or chemical associated cytotoxicity properties is always an early priority of initial FDA approval testing protocols. With cellular cytotoxicity assessment playing a central role in countless research and environmental safety studies, there is an ever present need for simple, straightforward analysis methods like the Necrosis vs Apoptosis Assay developed by Immunochemistry Technologies, LLC (ICT).

ICT's Necrosis vs Apoptosis Assay simultaneously detects both apoptosis associated cytotoxicity events as well as cell death due to necrosis. Apoptotic cells are identified using ICT's Fluorescent Labeled Inhibitor of Caspases (FLICA) reagent probe<sup>4-7</sup>. The FAM-FLICA<sup>®</sup> probe covalently binds to active caspase enzymes, which are up-regulated during apoptosis, thus clearly labeling apoptotic cells for subsequent analysis<sup>8,9</sup>. Non-apoptotic cells will not contain the active caspase enzymes required for FAM-FLICA to remain covalently bound within the cell structure.

Loss of the integrity of the cell membrane, indicative of necrosis or late stage apoptosis, is detected using the vital staining dye, 7-aminoactinomycin D (7-AAD), a red fluorescing live/dead stain. This dye easily penetrates cell membrane-compromised cells, binding tightly to GC rich regions of the DNA<sup>10-14</sup>. 7-AAD will not label cells in early stages of apoptosis, as their cell membranes are still intact and are capable of excluding 7-AAD. As caspases are active during early apoptosis, combining FAM-FLICA with 7-AAD can provide a more detailed picture of the overall health of the cell population by revealing the percentage of cells that are 7-AAD-negative (membrane intact live cells) and yet FAM-FLICA positive (apoptotic) (Figures 3-6).



*Simultaneous  
detection of live,  
early apoptotic, late  
apoptotic, and  
necrotic cells*

FAM-FLICA probes optimally excite at 488-492 nm with maximal emission at 515-535 nm. The vital staining dye 7-AAD optimally excites at 546 nm and emits at 647 nm. This significant difference in fluorescence emission wavelength between the green FAM (carboxyfluorescein) label on the FAM-FLICA probe and the red 7-AAD vital dye simplifies flow cytometer gating and compensation. The FAM-FLICA probe (apoptosis) is monitored on the FL-1 channel, while 7-AAD (necrosis) is monitored on FL-3. Combining the use of ICT's FAM-FLICA apoptosis detection probe with a membrane integrity dye like 7-AAD makes it easy to distinguish between necrosis and apoptosis within a single sample.

## 2. KIT CONTENTS

### Catalog #9147, Trial Size, 50-100 Tests, Contains:

- 2 vials of FAM-FLICA poly caspase inhibitor reagent, #637
- 1 vial of 7-Aminoactinomycin D (vital dye), 0.26 mg, #6163
- 1 bottle of 10X Apoptosis Wash Buffer (60 mL), #634
- 1 bottle of Fixative (6 mL), #636

### Catalog #9148, Standard Size, 100-200 Tests, Contains:

- 4 vials of FAM-FLICA poly caspase inhibitor reagent, #637
- 2 vials of 7-Aminoactinomycin D (vital dye), 0.26 mg/vial, #6163
- 2 bottles of 10X Apoptosis Wash Buffer (60 mL/bottle), #634
- 1 bottle of Fixative (6 mL), #636

## 3. STORAGE

Store the unopened kit and each unopened component at 2-8°C until the expiration date. Once reconstituted with DMSO, use FAM-FLICA and 7-AAD immediately, or store at ≤-20°C for 6 months protected from light and thawed no more than twice during that time.

#### 4. SAFETY DATA SHEETS (SDS)

Safety data sheets are available at [www.immunochemistry.com](http://www.immunochemistry.com) or by calling 1-800-829-3194 or 952-888-8788.

#### 5. RECOMMENDED MATERIALS

- DMSO, up to 1,000  $\mu$ L (50  $\mu$ L per vial to reconstitute FAM-FLICA, 260  $\mu$ L per vial to reconstitute 7-AAD, and more to create controls)
- DiH<sub>2</sub>O, up to 1,080 mL (540 mL per bottle to dilute 10X Apoptosis Wash Buffer)
- Phosphate buffered saline (PBS) pH 7.4, up to 100 mL, to dilute FAM-FLICA and handle cells
- FBS and/or BSA to add to the buffer when handling cells
- Cultured cells treated with the experimental conditions ready to be labeled
- Reagents to induce caspase activity and create controls for FAM-FLICA staining (See Figure 1), such as staurosporine (catalog #6212) or camptothecin (catalog #6210)
- 90% ETOH or 3% formaldehyde to create dead cell population for 7-AAD staining control (See Figure 2)
- Hemocytometer
- Centrifuge at 200 x g
- 15 mL polypropylene centrifuge tubes (1 per sample)
- Ice bath

#### 6. DETECTION EQUIPMENT

The assay can be analyzed with a:

- Flow cytometer
- Fluorescence microscope

Use filter pairings that best approximate these settings:

- FAM-FLICA optimally excites at 488-492 nm and has a peak emission at 515-535 nm (use FL-1 channel).
- 7-AAD optimally excites at 546 nm. It has a peak emission at 647 nm (use FL-3 channel).

#### 7. EXPERIMENTAL PREPARATION

Staining apoptotic cells with FAM-FLICA and 7-AAD can be completed within a few hours. However, the experiment is performed on living cells, which require periodic maintenance and cultivation several days in advance. In addition, once the proper number of cells has been cultivated, time must be allotted for the experimental treatment or apoptosis induction, which typically requires a 2-6 hour incubation at 37°C based on the cell line and experimental conditions.

Create cell populations, such as:

- Cells that were exposed to the experimental treatment.
- A negative control population of cells that received a placebo treatment.

As FAM-FLICA detects the presence of catalytically active forms of caspase enzymes, plan the experiment so that FAM-FLICA will be diluted and administered at the time when caspases are expected to be activated in the cells.

Culture cells to a density optimal for the specific experimental conditions or apoptosis induction protocol. Cell density should not exceed 10<sup>6</sup> cells/mL. Cells cultivated in excess of this concentration may begin to naturally enter apoptosis. An initial experiment may be necessary to determine when and how much FAM-FLICA to use as the resulting positive signal is a direct measurement of caspase activity occurring during the incubation period.

#### 8. CONTROLS

Create experimental samples and control cell populations:


- Treated experimental population(s): cells exposed to the experimental condition(s).
- Negative control: non-treated cells grown in a normal cell culture environment.

#### FIGURE 1: APOPTOSIS CONTROLS STAINED WITH FAM-FLICA


ICT's green poly caspase inhibitor reagent, FAM-FLICA, was used to assess apoptosis in Jurkat cells. Using FAM-FLICA, non-apoptotic cells (unstained, left side of each histogram) can easily be distinguished from apoptotic cells exhibiting green fluorescence (right side of each histogram). The forward and side scatter graphs are also shown.

To create a positive control for FAM-FLICA, a population of apoptotic cells is needed (Section 9). If analyzing with a flow cytometer, it will be used to compensate the instrument.

In this example, Jurkat cells were grown to 5 x 10<sup>5</sup> cells/mL and split into two populations. One population (A, top) was treated with a placebo (DMSO, non-induced) while the other population (B, bottom) was treated with 1  $\mu$ M staurosporine for 4 hours to induce apoptosis. Cells were stained with FAM-FLICA for 1 hour at 37°C and analyzed using an Accuri C6 flow cytometer in FL-1. Only 7% of non-induced cells (A, right) are apoptotic compared with 96.5% of the induced cells (B, right). (ICT 226:17-19.)



A: Non-Induced Cells: 93% non-apoptotic (left)  
7% apoptotic (right)



B: Induced Cells: 3.5% non-apoptotic (left)  
96.5% apoptotic (right)

- c. Positive control for FAM-FLICA: cells induced to undergo apoptosis using a known caspase activation protocol (Section 9; Figure 1).
- d. Positive control for 7-AAD: cells treated with ethanol or formaldehyde to create a dead cell population (Section 10; Figure 2).

The induced positive cell population and negative control cell population tubes should come from a common pool of cells and contain similar quantities of cells. Create negative controls by culturing an equal volume of non-induced cells for every labeling condition. For example, if labeling with FAM-FLICA and 7-AAD stain, make 8 populations:

- 1&2. Unlabeled: induced and non-induced
- 3&4. FAM-FLICA-labeled: induced and non-induced
- 5&6. FAM-FLICA and 7-AAD-labeled: induced and non-induced
- 7&8. 7-AAD-labeled: induced and non-induced

## 9. APOPTOSIS INDUCTION (FAM-FLICA CONTROL)

Prior to commencing the experiment, determine a reproducible


### FIGURE 2: LIVE AND DEAD CONTROLS STAINED WITH 7-AAD

Detection of cell membrane integrity loss, indicative of necrosis or late stage apoptosis, can be done using the red fluorescent live/dead stain, 7-AAD. This vital dye works by penetrating cell membrane-compromised cells and tightly binding to GC rich regions of the DNA.

Using 7-AAD, live cells (unstained, left side of each histogram) can easily be distinguished from dead/membrane compromised cells exhibiting red fluorescence (right side of each histogram).

To create a positive control for 7-AAD, a population of dead or killed cell is needed (Section 10). If analyzing with a flow cytometer, it will be used to compensate the instrument.

In this example, Jurkat cells were grown to  $5 \times 10^5$  cells/mL and split into two populations. One population (A) was left untreated while the other population (B) was treated with 90% ethanol for 60 seconds. Cells that were exposed to ethanol (B) were treated with a 5-fold larger volume of PBS to stop the ethanol surface denaturation process. Cells were pelleted by centrifugation ( $200 \times g$  for 10 minutes) and resuspended in PBS. Cells were then stained with 7-AAD for 10 minutes on ice, and analyzed using an Accuri C6 flow cytometer in FL-3. Only 4% of untreated cells (A) are dead compared with 97.6% of the treated cells (B). Data courtesy of Dr. Kristi Strandberg, ICT 226:30-31.



A: Untreated Cells:  
96% live (left)  
4% dead (right)

B: Ethanol-Treated  
(Killed) Cells:  
2.4% live (left)  
97.6% dead (right)

method for obtaining a positive control by triggering caspase activity. This process can vary with each cell line. For example, apoptosis may be induced with 2-4  $\mu\text{g/mL}$  camptothecin, or 1-2  $\mu\text{M}$  staurosporine for >4 hours (Figure 1).

## 10. PREPARATION OF DEAD CELLS (7-AAD CONTROL)

Prior to commencing the experiment, determine a reproducible method for obtaining a population of dead/killed cells to use as a positive control for 7-AAD staining. This can easily be achieved using many different techniques. For example, briefly expose cells to 90% ethanol for 30-60 seconds at  $37^\circ\text{C}$ , or expose cells to 3% formaldehyde for 30 minutes on ice (Figure 2).

## 11. PREPARATION OF FAM-FLICA

FAM-FLICA is supplied as a lyophilized powder that may be slightly visible as an iridescent sheen inside the vial. Protect from light and use gloves when handling. Because the 30X FAM-FLICA solution must be used immediately, prepare it just before staining.

1. Reconstitute each vial of FAM-FLICA with 50  $\mu\text{L}$  DMSO to form the 150X stock. The stock solution should be colorless to light yellow. Once reconstituted, it may be stored at  $\leq 20^\circ\text{C}$  for 6 months protected from light and thawed no more than twice during that time.
2. Immediately prior to addition to the samples and controls, dilute FAM-FLICA 1:5 by adding 200  $\mu\text{L}$  PBS to each vial to form the 30X FAM-FLICA solution. Use 30X FAM-FLICA within 30 minutes of dilution into aqueous buffers.
  - The recommended volume of 30X FAM-FLICA for flow cytometry is 5 to 10  $\mu\text{L}$  per 300  $\mu\text{L}$  of cells at  $5 \times 10^5$  cells/mL. The recommended volume of 30X FAM-FLICA for microscopy is 10  $\mu\text{L}$  per 300  $\mu\text{L}$  of cells at  $5 \times 10^5$  cells/mL. These amounts are recommendations, however, the amount needed may vary based on the experimental conditions and the instrument used for analysis. Each investigator should adjust the amount of FAM-FLICA to accommodate the particular cell line and research conditions.

## 12. PREPARATION OF 1X APOPTOSIS WASH BUFFER

ICT's 10X Apoptosis Wash Buffer (catalog #634) is an isotonic solution used to wash cells following exposure to FAM-FLICA. It contains mammalian proteins to stabilize cells stained with FAM-FLICA, and sodium azide to retard bacterial growth (1X Apoptosis Wash Buffer contains 0.01% w/v sodium azide). Cell culture media containing FBS and other additives may be used to wash cells instead of Apoptosis Wash Buffer.

1. 10X Apoptosis Wash Buffer may form precipitates during cold storage. If this happens, gently warm it until all crystals have dissolved. Do not boil.
2. Dilute 10X Apoptosis Wash Buffer 1:10 in  $\text{diH}_2\text{O}$ . For example, add 60 mL 10X Apoptosis Wash Buffer to 540 mL  $\text{diH}_2\text{O}$  for a total of 600 mL.
  - 1X Apoptosis Wash Buffer may be stored at  $2-8^\circ\text{C}$  and used within 1 week or frozen and used within 6 months.



### 13. PREPARATION OF 7-AAD

Detection of cell membrane integrity loss, indicative of necrosis or late stage apoptosis, is detected using the red fluorescent live/dead stain, 7-AAD (catalog #6163). This vital dye works by penetrating cells with compromised membranes and tightly binding to GC rich regions of the DNA. 7-AAD will not label cells in the early stages of apoptosis, as cell membranes are still intact at this stage. Combining 7-AAD with the green fluorescent FAM-FLICA reagent allows for concurrent labeling of apoptotic cells. 7-AAD is supplied as a lyophilized powder that may be slightly visible as a red sheen inside the vial. Protect from light and use gloves when handling.

1. Reconstitute each vial of 7-AAD with 260  $\mu\text{L}$  DMSO to create a stock concentrate at 1 mg/mL. Mix by swirling or tilting the vial, allowing the DMSO to travel around the base of the amber vial until the reagent is completely dissolved. At room temperature, the reagent should be dissolved within a few minutes forming a red solution.
  2. If storing the stock concentrate for future use, prepare small aliquots (50  $\mu\text{L}$ ) to avoid freeze-thaw cycles. The stock concentrate will be stable for 6 months when protected from light and stored at or below  $-20^{\circ}\text{C}$ .
  3. When ready to stain cells, use 7-AAD at 1:200 dilution. For example, add 2  $\mu\text{L}$  per 400  $\mu\text{L}$  cells.
- **Danger!** 7-Aminoactinomycin D (7-AAD) is fatal if swallowed, may cause cancer if swallowed, and may damage the unborn child if swallowed. See SDS for further information.

### 14. FIXATIVE

ICT's Fixative (catalog #636) is a formaldehyde solution designed to cross-link and aggregate intracellular components. If the stained cell populations cannot be evaluated immediately after labeling with FAM-FLICA, add Fixative at a ratio of 1:5 or 1:10. For example, to use at a ratio of 1:10, add 100  $\mu\text{L}$  Fixative to 900  $\mu\text{L}$  cells. Never add Fixative until all the staining and final wash steps have been completed. Fixed cells may be stored on ice or at  $4^{\circ}\text{C}$  for up to 16 hours, protected from light.

ICT's Fixative will not interfere with the carboxyfluorescein (FAM) or 7-AAD labels. Do not use absolute ethanol- or methanol-based fixatives as they may inactivate the FAM-FLICA signal. Do not fix cells that will be stained later with 7-AAD (stain the cells before using Fixative).

- **Danger!** Fixative contains formaldehyde <10% and methanol <5% and is harmful. Avoid contact with skin, eyes, and clothing by wearing lab coat, gloves, and safety glasses. In case of exposure, immediately flush eyes or skin with water. See SDS for further information.

### 15. STAINING PROTOCOL FOR SUSPENSION CELLS

Prepare experimental and control cell populations. Ideally, the cell concentration should be  $3\text{-}5 \times 10^5$  cells/mL. The concentration should not exceed  $10^6$  cells/mL, as cells cultivated in excess of this concentration may begin to naturally enter apoptosis. Just prior to staining, cells may need to be concentrated to  $2\text{-}5 \times 10^6$  cells/mL as microscopy analysis methods (Section 18) require high cell concentrations. Start with a larger volume of cells at  $3\text{-}5 \times 10^5$  cells/mL (which is a typical density for cell culture) and then concentrate

cells and resuspend to 300  $\mu\text{L}$  per sample when ready for FAM-FLICA staining.


1. Expose cells to the experimental or control condition. If analyzing with a flow cytometer, set aside four populations to create instrument controls with 7-AAD-positive and 7-AAD-negative cells (Section 10), and FAM-FLICA-induced and FAM-FLICA-non-induced cells (See Figure 5).
2. If analyzing with a fluorescence microscope, concentrate cells to  $2\text{-}5 \times 10^6$  cells/mL just prior to FAM-FLICA staining. Fluorescence microscopy requires an excess of  $2 \times 10^6$  cells/mL to obtain 5-20 cells per image field. Flow cytometry can analyze samples at  $3\text{-}5 \times 10^5$  cells/mL. Sample fluorescent microscopy results are shown in Figure 3.
3. Transfer 290  $\mu\text{L}$  cells into fresh tubes.
4. Add 10  $\mu\text{L}$  30X FAM-FLICA solution, forming a final volume of 300  $\mu\text{L}$ . If different cell volumes were used, add 30X FAM-FLICA at a ratio of 1:30. Mix by gently flicking the tubes.
  - a. The amount of FAM-FLICA used for flow cytometry applications can be diluted 2-fold for some cell lines and experimental conditions, allowing the end user to test a greater number of samples per kit. If staining at a lower concentration, add 5  $\mu\text{L}$  of FAM-FLICA to 295  $\mu\text{L}$  of cells. The amount of FAM-FLICA should be optimized for each cell line and experimental condition.

#### FIGURE 3: MICROSCOPY ANALYSIS OF JURKAT CELLS STAINED WITH FAM-FLICA AND 7-AAD

Jurkat suspension cells were exposed to 1  $\mu\text{M}$  of staurosporine for 4 hours at  $37^{\circ}\text{C}$  to induce apoptosis. Cells were dually stained with the green fluorescent FAM-FLICA poly caspase probe to detect apoptosis via caspase activity, and the red fluorescent vital dye 7-AAD to detect necrosis.

The image shown below in panel A reveals cells in early stage apoptosis (which have active caspase enzymes and intact cell membranes and fluoresce green with FAM-FLICA), cells in mid-to-late stage apoptosis (which fluoresce green with FAM-FLICA and red with 7-AAD as these cells have active caspase enzymes and compromised cell membranes), and necrotic cells (which fluoresce red with 7-AAD). Panel B shows a corresponding differential interference contrast (DIC) image, which reveals cell morphology. In Figure 5, a similar experiment was performed and cells were analyzed via flow cytometry to quantitate each population.

Microscope images were obtained using an Olympus BH-2 photomicroscope equipped with bright field, DIC, and fluorescence optics. FAM-FLICA and 7-AAD were imaged using a 470-490 nm excitation filter plus  $>520$  nm long pass filter tandem. Data courtesy of Dr. Brian W. Lee (ICT 196:70).





5. Incubate cells at 37°C protected from light. The incubation period may range from 30 minutes to several hours and should be optimized for each cell line and experimental condition. As cells may settle on the bottom of the tubes, gently resuspend them by swirling cells every 10-20 minutes to ensure an even distribution of FAM-FLICA throughout the staining process.
6. Centrifuge at 200 x g for 10 minutes, and then discard the supernatant.
7. Add 2 mL 1X Apoptosis Wash Buffer and gently mix.
8. Centrifuge at 200 x g for 10 minutes at RT.
9. Carefully remove and discard supernatants. Gently vortex the pellets to disrupt clumping. If analyzing by fluorescence microscopy, repeat wash process a third time. If using a flow cytometer, two wash steps are generally sufficient. However, additional wash steps may be required if high FAM-FLICA background signal is observed.
10. Resuspend cells in 400 µL 1X Apoptosis Wash Buffer.
11. Stain with 7-AAD at a final concentration of 5 µg/mL. This can be accomplished by:
  - a. Add the stock solution directly to the cell culture at 1:200 (e.g. add 2 µL stock to 400 µL cell suspension).
  - b. Dilute the stock concentrate 1:10 to form the working solution, and then add the working solution to the cells at 1:20. For example:
    1. Add 50 µL 7-AAD stock concentrate to 450 µL PBS or sterile media to form the working solution.
    2. Mix by inverting or vortexing the vial at RT.
    3. Store on ice up to 2 hours.
    4. Add the working solution to the cell suspension at approximately 1:20; e.g. add 25 µL diluted 7-AAD working solution into 475 µL cell suspension.
12. Incubate for 10-30 minutes on ice while protected from light.
  - a. To analyze by flow cytometer, refer to Section 17.
  - b. To analyze using a microscope, refer to Section 18.

## 16. STAINING PROTOCOL FOR ADHERENT CELLS

Adherent cells need to be carefully washed to avoid the loss of any cells that round up and come off the plate surface. Loose cells may be harvested from the plate or slide surface and treated as suspension cells, while those remaining adherent to the surface should be washed as adherent cells. If the adherent cells are trypsinized, the loose cells can be recombined with the trypsinized pool; alternatively, the loose cells can be recombined with the adherent portion when the analysis is performed. If growing adherent cells in a plate, the entire plate may be gently spun as part of the wash process to

sediment any loose floating cells. If trypsin is necessary, it is preferable to trypsinize prior to FAM-FLICA staining. Avoid trypsinizing cells prior to labeling with a vital dye, like 7-AAD. Cell membranes exposed to trypsin could be transiently permeant to vital dyes for a variable time depending upon the cell line. Cells may be labeled with FAM-FLICA before or after trypsinization (Figure 6).

1. Culture cells in TC-flasks and expose to the experimental conditions.
2. Collect the cells. Note: Apoptotic cells may detach and begin to float into the media. Save and spin to pellet and include these cells in the analysis.
 

If trypsin is required to lift the adherent cells:

  - a. Rinse the adherent cells with 5-10 mL of PBS to remove any residual trypsin inhibitors present in FBS-containing medium.
  - b. Trypsinize by adding enough Trypsin/EDTA to cover the cell surface area.
  - c. Incubate cells with trypsin at 37°C for several minutes. Check frequently to determine when the cells have lifted.
  - d. Neutralize with trypsin inhibitor or cell culture media containing 10-20% FBS.
  - e. Pool cells with any pellets created in Step 2.
3. Centrifuge the collected cells at 200 x g for 5 minutes.
4. Remove all but ~100 µL supernatant. Resuspend cells in 300-500 µL cell culture media containing 10-20% FBS. If necessary, count cells and adjust the volume of cell suspension to fit the experiment. Transfer cells into a 15 mL tube.
5. Add 30X FAM-FLICA at 1:30.
  - a. The amount of FAM-FLICA used for flow cytometry applications can be diluted 2-fold for some cell lines and experimental conditions, allowing the end user to test a greater number of samples per kit. Use FAM-FLICA at a dilution between 1:30 and 1:60. The amount of FAM-FLICA should be optimized for each cell line and experimental condition.
6. Incubate 30 minutes to several hours at 37°C, mixing gently every 10-20 minutes.
7. Centrifuge at 200 x g for 10 minutes, and then discard the supernatant.
8. Add 2 mL 1X Apoptosis Wash Buffer, and then centrifuge at 200 x g for 10 minutes.
9. Aspirate supernatant and resuspend cells in 2 mL 1X Apoptosis Wash Buffer.
10. Incubate 10 minutes at 37°C to allow any unbound FAM-FLICA to diffuse out of the cells.
11. Centrifuge at 200 x g for 10 minutes.
12. Aspirate supernatant and resuspend cells in 400 µL 1X Apoptosis Wash Buffer.
13. Stain with 7-AAD at a final concentration of 5 µg/mL. This can be accomplished by:
  - a. Add the stock solution directly into the cell culture at 1:200 (e.g. add 2 µL 7-AAD stock solution to 400 µL cell suspension).
  - b. Dilute the stock concentrate 1:10 to form the working solu-

**FOR RESEARCH USE ONLY.**

Not for use in diagnostic procedures.

tion, and then add the working solution to the cells at 1:20. For example:

1. Add 50  $\mu\text{L}$  7-AAD stock concentrate into 450  $\mu\text{L}$  PBS or sterile media.
  2. Mix by inverting or vortexing the vial at RT.
  3. Store on ice up to 2 hours.
  4. Add the working solution to the cell suspension at approximately 1:20; e.g. add 25  $\mu\text{L}$  diluted 7-AAD working solution to 475  $\mu\text{L}$  cell suspension.
14. Incubate for 10-30 minutes on ice while protected from light.
- a. To analyze by flow cytometer, refer to Section 17.
  - b. To analyze using a microscope, refer to Section 18.

### 17. FLOW CYTOMETRY ANALYSIS

Follow Section 15 or Section 16.

To address compensation issues and set up the flow cytometer, prepare 4 instrument control populations:

Cells stained only with FAM-FLICA (refer to Section 9 and Figure 1):

1. Induced (apoptotic) cells
2. Non-induced (non-apoptotic) cells

Cells stained with only 7-AAD (refer to Section 10 and Figure 2)

3. Live cells (intact cell membranes)
4. Killed cells (cell membrane compromised)

These controls are needed to adjust the instrument PMT's to separate 7-AAD-positive from 7-AAD-negative samples and to compensate for bleed-over of the red 7-AAD signal from FL-3 into FL-1. They will also help to clearly differentiate the FAM-FLICA-positive population from the FAM-FLICA-negative population and compensate bleed-over of the green FAM-FLICA signal from FL-1 into FL-3.

While setting up the 7-AAD controls (Section 10), continue working with the test samples and experimental controls as prepared in Section 8.

1. Set up the instrument compensation.
  - a. Read the 7-AAD-positive and 7-AAD-negative controls to compensate bleed-over of the red 7-AAD signal from FL-3 into FL-1.
  - b. Read the FAM-FLICA-positive and FAM-FLICA-negative controls to compensate bleed-over of the green FAM-FLICA signal from FL-1 into FL-3.
2. To read the samples for bicolor analysis:
  - a. Measure green FAM-FLICA (FAM = carboxyfluorescein) on the FL-1 channel.
  - b. Measure red fluorescence (7-AAD) on the FL-3 channel.
  - c. Generate a log FL-1 versus log FL-3 dot plot (Figures 5 and 6).
  - d. This will reveal 4 populations of cells: (Figure 4)
    1. Live cells (unstained) only emit low level background fluoresce (lower left).
    2. Cells in early apoptosis fluoresce green with FAM-FLICA (lower right).

3. Cells in late apoptosis are dually stained with FAM-FLICA and 7-AAD; they fluoresce green (they have active caspases) and red (the cell membrane is compromised; upper right).
4. Necrotic membrane-compromised cells fluoresce red (upper left).


### 18. MICROSCOPY ANALYSIS

Follow Section 15 or Section 16.

1. If not viewing immediately, cells may be washed and then fixed for viewing up to 16 hours later.
  - a. Wash sample once with PBS prior to fixing to remove any unbound 7-AAD from the medium.
  - b. Add Fixative at a v/v ratio of 1:5-1:10.
  - c. Incubate 15 minutes at RT in the dark.
  - d. Place cells onto a microscope slide and allow to dry.
  - e. Briefly rinse cells with PBS.
  - f. Cover cells with mounting media and coverslip.
  - g. Store slides at 2-8°C for up to 16 hours.
2. To view cells immediately, place 1 drop of cell suspension onto a microscope slide and cover with a coverslip.
3. Observe cells under a fluorescence microscope using excitation at 470-490 plus a >520 nm long pass filter. Cells bearing active caspase enzymes covalently coupled to FAM-FLICA appear green. Necrotic or late-stage apoptotic cells containing nucleic acid-bound 7-AAD appear red (Figure 3).

**FIGURE 4:**  
**QUANTITATE 4 CELL POPULATIONS**

1. Live, unstained cells do not fluoresce (LL, lower left).
2. Early stage apoptotic cells fluoresce green with FAM-FLICA (LR, lower right).
3. Dually stained green and red fluorescing cells represent cells in mid to late stage apoptosis; these cells have active caspase enzymes and compromised cell membranes (UR, upper right).
4. Necrotic cells fluoresce red (UL, upper left).



Dot Plot Key




**FIGURE 5:  
FLOW CYTOMETRY ANALYSIS OF JURKAT CELLS TO  
QUANTIFY 4 POPULATIONS**

Jurkat cells were treated with a placebo (non-induced treatment with DMSO; A and B) or treated with 1 μM staurosporine for 4 hours to induce apoptosis via caspase activity (C and D). Cells were then dually stained with ICT's FAM-FLICA apoptosis reagent and 7-AAD, a vital dye. Cells were analyzed using an Accuri C6 flow cytometer. Compensation was set using cell populations stained individually with either FAM-FLICA or 7-AAD (data not shown). FAM-FLICA was analyzed on FL-1, and 7-AAD was analyzed on FL-3. Forward and side scatter density plots (A, non-induced; C, induced populations) and dual stained density plots (B, non-induced; D, induced populations) are shown. The density plot key is shown in Figure 4.

Flow cytometry can be used to quantitate 4 populations of cells. Live, unstained cells do not fluoresce (lower left, B and D). Early stage apoptotic cells fluoresce green with FAM-FLICA (lower right, B and D). Dually stained green and red fluorescing cells represent the population of Jurkat cells in mid to late stage apoptosis; these cells have active caspase enzymes and compromised cell membranes (upper right, B and D). Necrotic cells fluoresce red (upper left, B and D).

By including FAM-FLICA in the analysis, cells in early apoptosis (lower right, B and D) can be detected, which 7-AAD cannot detect alone. In the non-induced population (B), only 10.8% of cells were apoptotic (LR: 3.6% + UR: 2.3%) compared with 30.5% of the induced population (D; LR: 22.1% + UR: 8.4%). Cells in these quadrants fluoresce green (ICT 226:17-19). In Figure 3, a similar experiment was performed and images were taken using a fluorescence microscope.




**FIGURE 6:  
FLOW CYTOMETRY ANALYSIS OF U-2 OS  
ADHERENT CELLS TO QUANTIFY 4 POPULATIONS**

U-2 OS cells were treated with a placebo (non-induced, A and B) or treated with 1 μM staurosporine for 4 hours to induce apoptosis via caspase activity (induced, C and D). Cells were lifted off the tissue culture flask by trypsinization and stained with FAM-FLICA. Cells were then stained with 7-AAD and analyzed by flow cytometry using an Accuri C6 flow cytometer. Compensation was set using cell populations stained individually with either FAM-FLICA or 7-AAD (data not shown). FAM-FLICA was analyzed on FL-1, and 7-AAD was analyzed on FL-3. Forward and side scatter density plots (A, non-induced; C, induced populations) and dual stained density plots (B, non-induced; D, induced populations) are shown. The density plot key is shown in Figure 4.

Flow cytometry can be used to quantitate 4 populations of cells. Live, unstained cells do not fluoresce (lower left, B and D). Early stage apoptotic cells fluoresce green with FAM-FLICA (lower right, B and D). Dually stained green and red fluorescing cells represent the population of U-2 OS cells in mid to late stage apoptosis; these cells have active caspase enzymes and compromised cell membranes (upper right, B and D). Necrotic cells fluoresce red (upper left, B and D).

By including FAM-FLICA in the analysis, cells in early apoptosis (lower right, B and D) can be detected, which 7-AAD cannot detect alone. In the non-induced population (B), only 5.9% of cells were apoptotic (LR: 3.6% + UR: 2.3%) compared with 30.5% of the induced population (D; LR: 22.1% + UR: 8.4%). Cells in these quadrants fluoresce green (ICT 226:17-19).



**FOR RESEARCH USE ONLY.**  
Not for use in diagnostic procedures.



## 19. REFERENCES AND CITATIONS

1. Russell, J. H. & Ley, T. J. Lymphocyte-mediated cytotoxicity. *Annu Rev Immunol.* 20, 323-370, doi:10.1146/annurev.immunol.20.100201.131730 (2002).
2. Brunner, K. T., Mael, J., Cerottini, J. C. & Chapuis, B. Quantitative assay of the lytic action of immune lymphoid cells on 51-Cr-labelled allogeneic target cells in vitro; inhibition by isoantibody and by drugs. *Immunology.* 14, 181-196 (1968).
3. Jerome, K. R., Sloan, D. D. & Aubert, M. Measurement of CTL-induced cytotoxicity: the caspase 3 assay. *Apoptosis.* 8, 563-571, doi:10.1023/A:1026123223387 (2003).
4. Bedner, E., Smolewski, P., Amstad, P. & Darzynkiewicz, Z. Activation of caspases measured in situ by binding of fluorochrome-labeled inhibitors of caspases (FLICA): correlation with DNA fragmentation. *Exp Cell Res.* 259, 308-313, doi:10.1006/excr.2000.4955 (2000).
5. Ekert, P. G., Silke, J. & Vaux, D. L. Caspase inhibitors. *Cell Death Differ.* 6, 1081-1086, doi:10.1038/sj.cdd.4400594 (1999).
6. Darzynkiewicz, Z., Bedner, E., Smolewski, P., Lee, B. W. & Johnson, G. L. Detection of caspases activation in situ by fluorochrome-labeled inhibitors of caspases (FLICA). *Methods Mol Biol.* 203, 289-299, doi:10.1385/1-59259-179-5:289 (2002).
7. Smolewski, P., Grabarek, J., Halicka, H. D. & Darzynkiewicz, Z. Assay of caspase activation in situ combined with probing plasma membrane integrity to detect three distinct stages of apoptosis. *J Immunol Methods.* 265, 111-121 (2002).
8. Amstad, P. A. et al. Detection of caspase activation in situ by fluorochrome-labeled caspase inhibitors. *Biotechniques.* 31, 608-610, 612, 614, passim (2001).
9. Grabarek, J., Amstad, P. & Darzynkiewicz, Z. Use of fluorescently labeled caspase inhibitors as affinity labels to detect activated caspases. *Hum Cell.* 15, 1-12 (2002).
10. Murata, S., Herman, P. & Lakowicz, J. R. Texture analysis of fluorescence lifetime images of nuclear DNA with effect of fluorescence resonance energy transfer. *Cytometry.* 43, 94-100 (2001).
11. Gill, J. E., Jotz, M. M., Young, S. G., Modest, E. J. & Sengupta, S. K. 7-Amino-actinomycin D as a cytochemical probe. I. Spectral properties. *J Histochem Cytochem.* 23, 793-799 (1975).
12. Lecoeur, H., Ledru, E., Prevost, M. C. & Gougeon, M. L. Strategies for phenotyping apoptotic peripheral human lymphocytes comparing ISNT, annexin-V and 7-AAD cytofluorometric staining methods. *J Immunol Methods.* 209, 111-123 (1997).
13. Madhavarao, M. S., Chaykovsky, M. & Sengupta, S. K. N7-Substituted 7-aminoactinomycin D analogues. Synthesis and biological properties. *J Med Chem.* 21, 958-961 (1978).
14. Sengupta, S. K., Tinter, S. K., Lazarus, H., Brown, B. L. & Modest, E. J. 7-substituted actinomycin D analogs. Chemical and growth-inhibitory studies. *J Med Chem.* 18, 1175-1180 (1975).
15. Lee, B. W., Olin, M. R., Johnson, G. L. & Griffin, R. J. *In vitro* and *in vivo* apoptosis detection using membrane permeant fluorescent-labeled inhibitors of caspases. *Methods Mol Biol.* 414, 109-135 (2008).

**FOR RESEARCH USE ONLY.**

Not for use in diagnostic procedures.

**Thank you for using this kit! If you have any questions, or would like to share your data, please contact us at 1-800-829-3194 or 952-888-8788, or send an email to [help@immunochemistry.com](mailto:help@immunochemistry.com).**



**BRIGHT MINDS, BRIGHT SOLUTIONS.**

ImmunoChemistry Technologies, LLC gratefully acknowledges the significant contributions made by one of its founders, Brian W. Lee, Ph.D in the development of this product, including the creation and illustration of its strategy and protocol.

**ImmunoChemistry Technologies, LLC**

9401 James Ave. S. #155, Bloomington, MN 55431 USA 800-829-3194

Copyright © 2017 ImmunoChemistry Technologies, LLC. FLICA® is a registered trademark of ImmunoChemistry Technologies, LLC.

[immunochemistry.com](http://immunochemistry.com)