

Please note this data sheet has been changed effective October 24, 2014

Mayer's Hematoxylin

INTENDED USE

For In Vitro Diagnostic Use

<u>AVAILABILITY:</u>	<u>Catalog #</u>	<u>Volume</u>
	TA-060-MH	60 ml
	TA-125-MH	125 ml

<u>SPECIFICITY:</u>	N/A
<u>ENZYME:</u>	N/A
<u>CHROMOGEN/SUBSTRATE:</u>	N/A

DESCRIPTION

Mayer's Hematoxylin is a stain with histological applications. This formula is based on the original formulation developed by Mayer and requires slightly longer application periods.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store at room temperature. Each component is stable for 18 months.

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE


1. Rehydrate tissue if needed.
2. Stain nuclei 1-2 minutes with hematoxylin.
3. Place in Bluing Reagent until stain is blue.
4. Rinse in DI water.
5. Coverslip.

REFERENCES

N/A

TROUBLESHOOTING


Please contact Thermo Fisher Scientific Technical Support by phone (1-269-544-5600 or 1-800-522-7270) or by email (lab.reagents@thermofisher.com).



Please note this data sheet has been changed effective October 24, 2014



Lab Vision Corporation
46500 Kato Road.
Fremont, CA 94538-7310, USA
US Toll Free: 1 (800) 522-7270
Phone: +1 (269) 544-5600
Fax: +1 (269) 372-2674
www.thermoscientific.com/labvision



Thermo Fisher Scientific
Anatomical Pathology
Tudor Road, Manor Park
Runcorn, Cheshire WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049
sales.ap.uk@thermofisher.com