

Please note this data sheet has been changed effective October 24, 2014

Large Volume Tris Buffered Saline (25X)

INTENDED USE

For In Vitro Diagnostic Use

AVAILABILITY:

Catalog # Volume
TA-125-TB 125 ml

DESCRIPTION

Tris Buffered Saline has been formulated in a 25X concentration for use as a wash buffer for immunolabeling techniques.

WARNINGS & PRECAUTIONS

Refer to MSDS.

STORAGE & SHELF LIFE

Store Tris Buffered Saline at room temperature.

- Concentrate is stable for 24 months from date of manufacture.
- Diluted solution should be used within one week. If solution becomes cloudy, discard and make fresh solution.

NOTE: The salts in the 25X solution are highly concentrated.

- Temperatures below 44 degrees F (6 degrees C) may cause the salts to precipitate.
- If this occurs, stir the solution at room temperature for about 30 minutes to re-dissolve the salts (use an alcohol-cleaned stir-bar to prevent bacterial contamination).

MICROBIOLOGICAL STATE

Product(s) not sterile.

MATERIALS REQUIRED BUT NOT PROVIDED

N/A

SPECIMEN & REAGENT PREPARATION

Refer to Procedure.

PROCEDURE

To make a 1X TBS solution dilute concentrate 25X with distilled water.

- Use 40 ml TBS concentrate per 960 ml distilled water
- Measure and pour appropriate volume of 25X TBS concentrate into a mixing flask and add DI water to final volume.
- Stir briefly.

The 1X solution should be pH 7.6 ± 0.2 .

- The pH of the solution may require slight adjustment.
- Adjust pH with 1N HCl (if pH is high) or 1N NaOH (if pH is low).
- **NOTE:** Tris buffers require the use of a calomel pH probe or a double-junction pH probe to obtain the correct pH reading.

Diluted solution should be used within one week. If solution becomes cloudy, discard and make fresh solution.


REFERENCES

N/A

TROUBLESHOOTING



Lab Vision Corporation
46500 Kato Road
Fremont, CA 94538-7310, USA
US Toll Free: 1 (800) 522-7270
Phone: +1 (269) 544-5600
Fax: +1 (269) 372-2674
www.thermoscientific.com/labvision




Thermo Fisher Scientific
Anatomical Pathology
Tudor Road, Manor Park
Runcorn, Cheshire WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049
sales.ap.uk@thermofisher.com

Please note this data sheet has been changed effective October 24, 2014

Please contact Thermo Fisher Scientific Technical Support by phone (1-269-544-5600 or 1-800-522-7270) or by email (lab.reagents@thermofisher.com).



Lab Vision Corporation
46500 Kato Road
Fremont, CA 94538-7310, USA
US Toll Free: 1 (800) 522-7270
Phone: +1 (269) 544-5600
Fax: +1 (269) 372-2674
www.thermoscientific.com/labvision



Thermo Fisher Scientific
Anatomical Pathology
Tudor Road, Manor Park
Runcorn, Cheshire WA7 1TA, UK
Tel: +44 (0) 1928 534 050
Fax: +44 (0) 1928 534 049
sales.ap.uk@thermofisher.com