

Product information and content

AlumVax Hydroxide is wet gel (colloidal) of aluminum hydroxide 2%, provided as a ready-to-use suspension. It is sterilized and aseptically filled. AlumVax hydroxide is a crystalline aluminum oxyhydroxide that is positively charged at physiological pH (pI=11), suitable for adsorption of negatively charged, acidic proteins (such as albumin).

AlumVax Hydroxide is available in two quantities:
#AH0050: 50 mL and #AH0250: 250 mL.

Storage and stability

Shipping and storage: Room Temperature.
Product is stable for 3 years. DO NOT FREEZE (ice crystal may impair the product).

Chemical properties

CAS Number: 21645-51-2
Formulation: Al(OH)₃, Aluminium hydroxide gel
Appearance: White gelatinous precipitate
Aluminium content: 9.0 – 11.0 mg/ml; pH: ~6.5

Description

Aluminum hydroxide is the most common adjuvant used in approved prophylactic vaccines because of its excellent safety profile and ability to enhance protective humoral immune response. Since more than 80 years, it has been observed that aluminium compounds act by a depot effect and also by direct activation of the immune cells. Adsorption or entrapment of antigens in aggregates through hydrophobic and electrostatic interactions favors a high local antigen concentration and improved uptake by antigen presenting cells (APC).

AlumVax hydroxide stimulates Th2 response through the release of Th2-associated cytokines (IL4, IL-5, IL-13...) and Th2-associated antibodies (IgG1 & IgE). It increases Ag-specific CD4+ T Cell proliferation and promotes NALP3 inflammasome activation and caspase 1-mediated release of IL-1 and IL-18.

NOTE: Alum is frequently used as an alternative to Freund's adjuvants, as it is less hazardous and less likely to cause tissue necrosis at the injection site.

Results

Results presented below demonstrate the effect of Aluminum hydroxide adjuvant on immune system response:

AlumVax Hydroxide 2% Vaccine Adjuvant




Figure 1. Response of innate immune system cells to injection of Ag +/- Alum Adjuvant. The injection of ovalbumin (OVA) adsorbed onto Alum (OVA-Alum) in mice induces a marked increase in the levels of the monocyte chemotactic protein (MCP1; CCL2), the neutrophil chemotaxin KC (CXCL1) and the eosinophil chemotaxin eotaxin-1 (CCL11) compared to mice receiving OVA alone (adapted from Kool M. et al., J Exp Med. 2008; 205(4):869-82).




Figure 2. Innate inflammatory response induced by alum. Mice were injected i.p. with ovalbumin (OVA) or OVA-Aluminum hydroxide. 24H later, number of Neutrophils, Eosinophils and Monocytes and dendritic cells was determined in peritoneal lavage by flow cytometry (adapted from Kool M. et al., J Immunol. 2008; 181(6):3755-9).

Method/protocol

Recommendations before starting:

The inoculum should be free of extraneous microbial contamination; filtration of the antigen before mixing with the adjuvant is recommended.

Preparation of AlumVax hydroxide-immunogens mix is easier than other adjuvants as it does not require laborious emulsification.

1. Ensure complete re-suspension of AlumVax hydroxide adjuvant by vigorous shaking of the vial before use.

2. Dilute antigen mixture in saline buffer or phosphate buffer for a final immunogen concentration of 10-100 µg/100 µL.
3. Mix AlumVax hydroxide adjuvant with an equal volume of antigen solution for a 1:1 ratio:
 - a. Add AlumVax hydroxide dropwise with constant mixing to the immunogen solution.
 - b. Pipet up and down several times to ensure correct absorption of antigen by alum adjuvant and incubate 5 to 10 minutes.

NOTE: Ratio can be optimized from 1:1 (100µL adjuvant per 100µL antigen) to 1:9 (100µL adjuvant per 900µL antigen)

4. Inject into the animal according to the table below; the volume depends on the site of injection. Typical routes of administration include subcutaneous (SC), intramuscular (IM), intradermal (ID) or intraperitoneal (IP).

Species	Max vol/site	Primary injection	Subsequent Injection(s)
Mice, hamsters	0.1 mL	SC	SC
Mice, hamsters	0.05 mL	IM ^Δ	IM ^Δ
Mice	0.5 mL	IP [×]	SC, IM ^Δ
Guinea pigs, rats	0.2 mL	SC, IM ^Δ	SC, IM ^Δ
Rabbits	0.25 mL	SC, IM	SC, IM
Rabbits	0.025 mL	ID	SC, IM
Sheep, goats, donkeys, pigs, Chickens	0.5 mL	SC, IM	

^Δ Not recommended in general, in particular not for viscous adjuvants

[×] Not recommended for pAb production

Table 1: Maximum volumes for injection of immunogen/adjuvant mixtures per site of injection for different animal species (*Adapted from Leenars MPPA, Hendriksen CFM et al., 19*)

OZ Biosciences SAS

163 avenue de Luminy
Case 922, zone entreprise
13288 Marseille cedex 09 - FRANCE
Ph: +33 (0) 486 948 516
Fax: +33 (0) 486 948 515
contact@ozbiosciences.com
order@ozbiosciences.com

OZ Biosciences INC

4901 Morena Blvd,
Suite 901
San Diego CA 92117 - USA
Ph: + 1-858-246-7840
Fax: + 1-855-631-0626
contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com

Purchaser Notification

Limited License

The purchase of the AlumVax Hydroxide Vaccine Adjuvant grants the purchaser a non-transferable, non-exclusive license to use the included components. This reagent is intended for in-house research only by the buyer. Such use is limited to the transfection of nucleic acids as described in the product manual. In addition, research only use means that this formulation is excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences.

Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of the AlumVax Hydroxide Vaccine Adjuvant. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, contact the Director of Business Development at OZ Biosciences.

Buyers may end this License at any time by returning all AlumVax Hydroxide Vaccine Adjuvant material and documentation to OZ Biosciences, or by destroying AlumVax Hydroxide Vaccine Adjuvant components. Purchasers are advised to contact OZ Biosciences with the notification that a AlumVax Hydroxide Vaccine Adjuvant is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s).

This document covers entirely the terms of the AlumVax Hydroxide Vaccine Adjuvant research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

The AlumVax Hydroxide Vaccine Adjuvant is developed, designed, intended, and sold for research use only. It is not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the component by following proper research laboratory practices

