

CNE-CPO Vaccine Adjuvant

Cationic Nano Emulsion Vaccine Adjuvant

Product information

CNE-CPO vaccine adjuvant is a squalene oil-in-water emulsion containing cationic components.

CNE-CPO is available in two quantities. #CNE1000 1mL; #CNE5000 5mL.

Storage and stability

Shipping and storage: **CNE-CPO** is shipped at RT and stored at +4°C. **CNE-CPO** is stable for 1 year. DO NOT FREEZE.

Description

CNE-CPO is an oil-in-water Cationic Nano Emulsion made of squalene droplets and cationic polymers in a continuous aqueous phase. It is biodegradable, which is an important advantage over alternative oils that have been used in emulsion adjuvants, like Freund's adjuvant that contains mineral oil (paraffin oil) and has long term persistence in the organism. CNE induces local stimulation and recruitment of DCs and granulocytes, differentiation of monocytes into DCs and increased uptake of antigen by DCs. The emulsion acts more specifically on macrophages present at the site of injection. A local increase of chemokines released also influences the recruitment of immune cells from the blood to the site of vaccination, creating an amplification loop. This formulation enhances differentiation of monocytes towards a mature phenotype, thereby promoting migration of antigen-loaded cells to the draining lymph node. Compared to aluminum salts, a stronger immune response is elicited (e.g. higher antibody "humoral response, Th2" and T-cell response "cellular response Th1") with a mixed and more balanced Th1/Th2 cell phenotype.

The cationic components that compose this nano-emulsion make it a genetic adjuvant that allows the association with plasmid DNA to form an efficient nanoparticle delivery system (NPD). NPD are **non-viral gene delivery systems**, self-assembled from cationic entities and negatively charged immunogen that function as **vaccine carrier**.

CNE-CPO adjuvant is compatible with most immunization procedures: such as intramuscular, intraepidermal, intravenous, intraperitoneal or subcutaneous.

- **Plasmid DNA**

Cationic formulation-mediated antigen-coding plasmid DNA has been shown to greatly improve humoral and cell-mediated immunity. One of the possibilities is that these DNA vaccines could facilitate uptake of the plasmid by antigen-presenting cells (APC) and induce cytotoxic T lymphocyte response. Moreover, once entrapped into nanoparticles, DNA is protected from nucleases and depending on their size, some NPD may break down locally to release their vaccine content slowly; the accessibility of genetic material is thus prolonged.

Results

Results presented below demonstrate the effect of CNE adjuvant on immune system response:




Figure 1. Comparison of immune response in mice and rabbits injected with naked- or CNE-DNA at either 1, 10 or 500 µg. Sera were analyzed by ELISA. Adapted from G. Ott *et al.* *J. Control. Release*, 2002; 79(1-3), 1-5.

Method/protocol

Recommendations before starting:

The inoculum should be free of extraneous microbial contamination; use plasmid DNA as pure as possible. Adapt volumes according to the table 1 below.

1. Allow CNE adjuvant and immunogen solutions to reach room temperature before beginning.
2. Shake gently the CNE vial before opening.
3. Dilute immunogen/DNA/RNA in saline buffer or phosphate buffer for a final concentration of 10-100 µg/100 µL.

It is mandatory to not use buffer containing serum.

4. Mix CNE adjuvant with an equal volume of immunogen/DNA/RNA solution for a 1:1 ratio.
5. Pipette up and down several times to ensure correct mixing.
6. Incubate at room temperature for 20-30 min.
7. Inject into the animal according to the table below.

NOTE: do not store the complexes: discard solution after use. Prepare fresh LPD before each immunization

Volume (mL) for injection depends on the site of injection and the animal model. Typical routes of administration include intramuscular (IM), subcutaneous (SC), intradermal (ID) or intraperitoneal (IP).

Species	I.M.	S.C.	I.D.	I.P.
Mice, hamsters	0.05-0.1	0.1-0.2	0.025	0.5
Guinea pigs, rats	0.1-0.2	0.2-0.4	0.025	1.0
Rabbits	0.25	0.25	0.025	10
Pigs	0.25-0.5	0.5	0.5	50

Table 1: Recommended volumes (mL) for injection of immunogen/adjuvant mixtures per site of injection for different animal species (Adapted from Leenars MPPA, Hendriksen CFM *et al.*, 1999).

References and background reading

- S. W. Yi, T. Y. Yune *et al.* A cationic lipid emulsion/DNA complex as a physically stable and serum-resistant gene delivery system. *Pharmaceutical research*, **2000**, vol. 17, no 3, p. 314-320.
- T. W. Kim, H. Chung *et al.* In vivo gene transfer to the mouse nasal cavity mucosa using a stable cationic lipid emulsion. *Molecules and cells*, **2000**, vol. 10, no 2, p. 142-147.
- M. Singh, M. Briones *et al.* Cationic microparticles: a potent delivery system for DNA vaccines. *Proceedings of the National Academy of Sciences*, **2000**, vol. 97, no 2, p. 811-816.
- T. W. Kim, H. Chung *et al.* Optimization of lipid composition in cationic emulsion as in vitro and in vivo transfection agents. *Pharmaceutical research*, **2001**, vol. 18, no 1, p. 54-60.
- G. Ott, M. Singh *et al.* A cationic sub-micron emulsion (MF59/DOTAP) is an effective delivery system for DNA vaccines. *Journal of controlled release*, **2002**, vol. 79, no 1-3, p. 1-5.
- L. A. Brito, M. Chan *et al.* A cationic nanoemulsion for the delivery of next-generation RNA vaccines. *Molecular Therapy*, **2014**, vol. 22, no 12, p. 2118-2129.
- W. M. Bogers, H. Oostermeijer *et al.* Potent immune responses in rhesus macaques induced by nonviral delivery of a self-amplifying RNA vaccine expressing HIV type 1 envelope with a cationic nanoemulsion. *The Journal of infectious diseases*, **2015**, vol. 211, no 6, p. 947-955.
- K. Adlington, J. EL Harfi *et al.* Molecular design of squalene/squalane countertypes via the controlled oligomerization of isoprene and evaluation of vaccine adjuvant applications. *Biomacromolecules*, **2016**, vol. 17, no 1, p. 165-172.
- E. S. BesumbeS, C. Fornaguera *et al.* PLGA cationic nanoparticles, obtained from nano-emulsion templating, as potential DNA vaccines. *European Polymer Journal*, **2019**, vol. 120, p. 109229.
- A. S. A. Medeiros, M. Torres-rêgo *et al.* Self-assembled cationic-covered nanoemulsion as a novel biocompatible immunoadjuvant for antiserum production against Tityus serrulatus scorpion venom. *Pharmaceutics*, **2020**, vol. 12, no 10, p. 927.

Purchaser Notification

Limited License

The purchase of the CNE-CPO Vaccine Adjuvant grants the purchaser a non-transferable, non-exclusive license to use the included components. This reagent is intended for in-house research only by the buyer. Such use is limited to the transfection of nucleic acids as described in the product manual. In addition, research only use means that this formulation is excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of OZ Biosciences. Separate licenses are available from OZ Biosciences for the express purpose of non-research use or applications of the CNE-CPO Vaccine Adjuvant. To inquire about such licenses, or to obtain authorization to transfer or use the enclosed material, contact the Director of Business Development at OZ Biosciences.

Buyers may end this License at any time by returning all CNE-CPO Vaccine Adjuvant material and documentation to OZ Biosciences,

or by destroying CNE-CPO Vaccine Adjuvant components. Purchasers are advised to contact OZ Biosciences with the notification that a CNE-CPO Vaccine Adjuvant is being returned in order to be reimbursed and/or to definitely terminate a license for internal research use only granted through the purchase of the kit(s).

This document covers entirely the terms of the CNE-CPO Vaccine Adjuvant research only license, and does not grant any other express or implied license. The laws of the French Government shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

The CNE-CPO Vaccine Adjuvant is developed, designed, intended, and sold for research use only. It is not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the use of the component by following proper research laboratory practices.

OZ Biosciences SAS

163 avenue de Luminy
Case 922, zone entreprise
13288 Marseille cedex 09 - FRANCE
Ph: +33 (0) 486 948 516
Fax: +33 (0) 486 948 515
contact@ozbiosciences.com
order@ozbiosciences.com

OZ Biosciences INC

4901 Morena Blvd,
Suite 901
San Diego CA 92117 - USA
Ph: + 1-858-246-7840
Fax: + 1-855-631-0626
contactUSA@ozbiosciences.com
orderUSA@ozbiosciences.com

Rev 10/20AB&FP