## ANTICORPS SECONDAIRES

Fragments F(ab')<sub>2</sub> anti IgG (gamma) humaine  
conjugués à la fluorescéine (FITC)

CODE : BI 4118

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
## SECONDARY ANTIBODIES

Fluorescein labelled F(ab')<sub>2</sub>  
to human IgG (gamma)

CODE : BI 4118

for research use only

**Immunogen** : Highly purified immunoglobulin gamma heavy chain

**Host** : goat

**Quantity** : 1mg/1ml

**Format** : PBS, BSA 0.3%, NaN<sub>3</sub> 0.1%

For a better result in flow cytometry this product is glycerol free.

**Storage** : +4 C. For long storage this product should be distributed into aliquotes before freezing at -20 C.

**Applications** : Fab and F(ab')<sub>2</sub> are respectively mono and bivalent fragments of immunoglobulins, i.e. including one or two sites for the fixation of the antigen. They are prepared by enzymic digestion (papain or pepsin). Like the others P.A.R.I.S antibodies, they are purified by affinity chromatography and are thus particularly specific.

Fab and F(ab')<sub>2</sub> fragments have certain advantages with respect to the entire antibodies in immunohistochemistry, immunocytology and flow cytometry tests.

On the one hand their specificity is increased since they can no longer be bound in a non-specific way by their Fc part as entire Ig's are.

On the other hand their use also increase the sensitivity of the technique. In fact, their small size enhances their penetration into tissues and cells as well as their availability for the antigens (this is so particularly for Fab).

Finally, conjugate Fab's have a label/antigenic ratio that is twice as much compared to antibodies (except fluorescent conjugates).

**Conjugate properties** : Fluorescein labelling is carried out to obtain a good cue while at the same time minimising the background. The F/P (fluorochrome/protein) molar ratio is set between 3 and 7. FITC has the advantage of being easily detectable under a microscope. It can be easily used for carrying out multiple labelling in histochemistry or cytochemistry. FITC may rapidly lose its fluorescence.

Maximum Absorption : 495nm

Maximum Emission : 520nm

Color : green yellow

**Working dilutions** :

Dilutions must be determined according to user's experience. Only general bases are indicated here.

Indirect immunofluorescence (IIF) : 1/50 - 1/200

Immunocytology/Immunohistology : 1/10 - 1/100

Flow cytometry (FACS) : 1/10 - 1/100

Version janv 07