

## Instruction Manual

Product	Size	Catalog Number
NHEK.f single donor Pellet	1 million cells per pellet	C-14001
NHEK adult, single donor Pellet	1 million cells per pellet	C-14002
NHEK.f pooled Pellet	1 million cells per pellet	C-14003
NHEK adult pooled Pellet	1 million cells per pellet	C-14004
HFDPC Pellet	1 million cells per pellet	C-14005
HUVEC Growth Medium 2 single donor Pellet	1 million cells per pellet	C-14008
HUVEC Growth Medium 2 pooled Pellet	1 million cells per pellet	C-14009
HUVEC single donor Pellet	1 million cells per pellet	C-14010
HUVEC pooled Pellet	1 million cells per pellet	C-14011
HUVEC pre-screened Pellet	1 million cells per pellet	C-14012
HUAEC Pellet	1 million cells per pellet	C-14013
HDMEC Pellet	1 million cells per pellet	C-14015
HDMEC adult Pellet	1 million cells per pellet	C-14016
HDMEC pre-screened Pellet	1 million cells per pellet	C-14017
HDBEC Pellet	1 million cells per pellet	C-14018
HDBEC adult Pellet	1 million cells per pellet	C-14019
HDLEC Pellet	1 million cells per pellet	C-14020
HDLEC adult Pellet	1 million cells per pellet	C-14021
HCAEC Pellet	1 million cells per pellet	C-14022
HAoEC Pellet	1 million cells per pellet	C-14023
HPAEC Pellet	1 million cells per pellet	C-14024
HSaVEC Pellet	1 million cells per pellet	C-14025
HBdMEC Pellet	1 million cells per pellet	C-14026
HPMEC Pellet	1 million cells per pellet	C-14027
HUtMEC Pellet	1 million cells per pellet	C-14028
HCMEC Pellet	1 million cells per pellet	C-14029
NHDF Pellet	1 million cells per pellet	C-14030
NHDF adult Pellet	1 million cells per pellet	C-14031
HBMEC Pellet	1 million cells per pellet	C-14032
HVMF Pellet	1 million cells per pellet	C-14034
HPF Pellet	1 million cells per pellet	C-14035
HCF Pellet	1 million cells per pellet	C-14036
HAoAF Pellet	1 million cells per pellet	C-14037
HUF Pellet	1 million cells per pellet	C-14038
NHEM.f Pellet	1 million cells per pellet	C-14040
NHEM.f M2 Pellet	1 million cells per pellet	C-14042
NHEM M2 adult Pellet	1 million cells per pellet	C-14043

Product	Size	Catalog Number
HUASMC Pellet	1 million cells per pellet	C-14050
HPASMC Pellet	1 million cells per pellet	C-14051
HCASMC Pellet	1 million cells per pellet	C-14052
HAoSMC Pellet	1 million cells per pellet	C-14053
HBSMC Pellet	1 million cells per pellet	C-14055
HTSMC Pellet	1 million cells per pellet	C-14056
HBdSMC Pellet	1 million cells per pellet	C-14057
HPrSMC Pellet	1 million cells per pellet	C-14058
HUtSMC Pellet	1 million cells per pellet	C-14059
SkMC Pellet	1 million cells per pellet	C-14060
HNEpC Pellet	1 million cells per pellet	C-14062
HBepC Pellet	1 million cells per pellet	C-14063
HTEpC Pellet	1 million cells per pellet	C-14064
HSAepC Pellet	1 million cells per pellet	C-14065
HMEpC Pellet	1 million cells per pellet	C-14066
HREpC Pellet	1 million cells per pellet	C-14067
HRCEpC Pellet	1 million cells per pellet	C-14068
HPLEpC Pellet	1 million cells per pellet	C-14069
HCH Pellet	1 million cells per pellet	C-14070
HOB Pellet	1 million cells per pellet	C-14071
HWP subcutaneous Pellet	1 million cells per pellet	C-14072
HWP visceral Pellet	1 million cells per pellet	C-14073
HCM Pellet	1 million cells per pellet	C-14080
HHpC Pellet	1 million cells per pellet	C-14085
hMSC-BM Pellet	1 million cells per pellet	C-14090
hMSC-UC Pellet	1 million cells per pellet	C-14091
hMSC-AT Pellet	1 million cells per pellet	C-14092
hPC-PL Pellet	1 million cells per pellet	C-14095
hMNC-CB single donor Pellet	1 million cells per pellet	C-14096
hMNC-CB pooled Pellet	1 million cells per pellet	C-14097
hMNC-PB pooled Pellet	1 million cells per pellet	C-14099
hMoCD14 <sup>+</sup> -PB single donor Pellet	1 million cells per pellet	C-14110
hMo-PB single donor Pellet	1 million cells per pellet	C-14111
hMoCD14 <sup>+</sup> -PB single donor, pre-screened Pellet	1 million cells per pellet	C-14112
hMo-PB single donor, pre-screened Pellet	1 million cells per pellet	C-14113
hMNC-PB single donor Pellet	1 million cells per pellet	C-14098

## Product Description

---

PromoCell Cell Pellets are stored in RNAlater® - an aqueous protection solution that rapidly permeates the cells and stabilizes cellular DNA, RNA, and proteins and protects them from degradation. This makes the PromoCell Cell Pellets very useful for a wide range of research applications like gene expression analysis, cell characterization studies, and protein characterization experiments.\*

**Note:** RNAlater® denatures proteins. Therefore, proteins will be suitable for e.g. Western blotting but not for applications that require native proteins.

Because of the density of RNAlater Solution, greater centrifugal forces are required to pellet cells. We recommend centrifuging the cells at 5000 x g for 10 minutes at 4°C.

Each Cell Pellet is made from 1 million cells and is dissolved in 200 µl RNAlater®. The Cell Pellets are produced from released cell lots and have the same characteristics, i.e. expression profile or cell type specific antigens. On request, viable cells from the same lot are available for cultivation.

## Storage and Stability

---

The Cell Pellets can be stored indefinitely at -20°C, and for up to 1 month at 4°C.

**Note:** The Cell Pellets will not freeze at -20°C. Always centrifuge product briefly before opening vial.

## Quality Control

---

Rigid quality control tests are performed for each lot of the corresponding cells. They are tested for cell viability, cell mor-

phology, growth performance and for cell-type specific markers.

In addition, all cells have been tested for the absence of HIV-1, HIV-2, HBV, HCV, and microbial contaminants (fungi, bacteria, and mycoplasma).

## Intended Use

---

PromoCell Cell Pellets are for *in vitro* research use only and not for diagnostic or therapeutic procedures.

## Warning

---

Although tested negative for HIV-1, HIV-2, HBV, and HCV, the cells - like all products of human origin - should be handled as potentially infectious. No test procedure can completely guarantee the absence of infectious agents.

**Follow appropriate safety precautions!**

**PromoCell GmbH**

Sickingenstr. 63/65  
69126 Heidelberg  
Germany

Email: [info@promocell.com](mailto:info@promocell.com)  
[www.promocell.com](http://www.promocell.com)

**USA/Canada**

Phone: 1 – 866 – 251 – 2860 (toll free)  
Fax: 1 – 866 – 827 – 9219 (toll free)

**Deutschland**

Telefon: 0800 – 776 66 23 (gebührenfrei)  
Fax: 0800 – 100 83 06 (gebührenfrei)

**France**

Téléphone: 0800 90 93 32 (ligne verte)  
Téléfax: 0800 90 27 36 (ligne verte)

**United Kingdom**

Phone: 0800 – 96 03 33 (toll free)  
Fax: 0800 – 169 85 54 (toll free)

**Other Countries**

Phone: +49 6221 – 649 34 0  
Fax: +49 6221 – 649 34 40