

Sodium Dodecyl Sulfate (SDS)

Background:

Dodecyl sulfate, sodium salt (SDS), is a are potent protein denaturant and solubilizing agent which bind to proteins cancelling differences in intrinsic charge and converting their shape to rod-like particles whose lengths vary uniquely with molecular weights.

To meet the needs of researchers, Polysciences offers SDS, Electro Pure™ which contain less than 10 ppm of heavy metals. Heavy metals can precipitate or complex with heavy metals. These complexes can appear as false bands on electrophoresis gels. In addition, hexadecyl sulfate (C₁₆ alkyl sulfate), a common minor constituent of many grades of SDS, has been reported to inhibit the renaturation of enzyme activity after SDS polyacrylamide gel electrophoresis.¹ Polysciences SDS contains no C₁₆ alkyl sulfate.

Applications:

By binding with protein molecules, SDS reduces charge differences between proteins to make size differences of the complexes directly proportional to the molecular weight of the proteins. Thus molecular weights may then be determined by comparing the electrophoretic migration of unknown proteins with standards.²⁻⁴ Other applications of SDS include solubilization of proteins⁵ and enzyme substrates,⁶ and the fractionation of high and low molecular weight RNAs.⁷

Specifications:

	SDS, Electro/Pure™
C ₁₂ Assay	99% min.
Free fatty alcohol	0.2% max.
Phosphate	0.0001% max.
O.D., 220-350 nm (3% aqueous)	0.1% max.
Heavy Metals	10 ppm max.
pH of 1% solution	6.0-7.5

Caution:

These powders can cause irritation. Prolonged contact with skin may result in dermatitis.

First Aid:

In case of contact, flush with plenty of water for at least 15 minutes. Get medical attention for eyes. If ingested, wash out mouth with water and obtain medical attention.

Ordering Information:

Cat. #	Description	Size
03945	Sodium Dodecyl Sulfate, Electro Pure™	100g 1 kg
00019	Acrylamide, Electro Pure™, 99.9%, Recrystallized (Write for TDS #155)	100g 500g

03707	Coomassie Blue, G250	100g
00352	Coomassie Blue, R250	100g

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com**References:**

1. J. Biol. Chem., **244**, 4406 (1969).
2. Ibid., **245**, 1561 (1970).
3. Anal. Biochem., **44**, 246 (1971).
4. Meth. Enzymol., **48**, 3 (1978).
5. Anal. Biochem., **55**, 573 (1973).
6. Meth. Enzymol., **15**, 472 (1969).
7. Biochem., J., **102**, 251 (1967).
8. Noll, H., Methods in Enzymology, **XIIIB**, 130 (1968).
9. Delepelaire, P., et al., Proc. Nat. Acad. Science USA, **76**, 111 (1979).
10. Dam, V.T., et al., Anal. Biochem., **114**, 330 (1981).
11. Sinclair, J.F., et al., Anal. Biochem., **114**, 316 (1981).

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for his own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.