



Corporate Headquarters
400 Valley Road
Warrington, PA 18976
1-800-523-2575
FAX 1-800-343-3291
Email: info@polysciences.com
www.polysciences.com

Europe - Germany
Polysciences Europe GmbH
Handelsstr. 3
D-69214 Eppelheim, Germany
(49) 6221-765767
FAX (49) 6221-764620
Email: info@polysciences.de

TECHNICAL DATA SHEET 615

Page 1 of 2

Protein Conjugated Microspheres

Introduction:

Polysciences, Inc. offers antibodies, Protein A, and Protein G covalently coupled to dyed (fluorescent and non-fluorescent) and undyed microspheres. Antibody conjugated microspheres are used to detect trace amounts of antigens in solution via bead based ELISA & agglutination tests. Microspheres conjugated with Protein A and Protein G can be used to purify IgG by binding to the Fc portion of antibodies raised in most mammals.

Our protein conjugated microspheres have a diameter of 1.0 μ m and are offered as aqueous suspensions containing microspheres at a concentration of approximately 1.25%. They are packaged in a 0.02M Sodium Phosphate buffer, pH 7.4, containing 8mg/ml NaCl, 10mg/ml Bovine Serum Albumin, 0.1% Sodium Azide, and 5% glycerol. Each lot of the Protein Conjugated Microspheres is Quality Control tested for mean particle diameter, % solids, and the concentration of protein conjugated to the microspheres. This information is reported on the product label. The fluorescent yellow-green (YG) microspheres have an excitation max of 445nm and an emission max of 500nm, similar to FITC. For best results, these microspheres should be stored at 4°C, not frozen, and mixed before using.

These products are for research use only, not intended for use in humans or in vitro diagnostics use.

Antibody Conjugated Microspheres

Goat Anti-Mouse IgG (H&L) and Goat Anti-Rabbit IgG (H&L) are available on 1.0 μ m polystyrene microspheres. These antibodies are available on plain, visible blue dyed, or fluorescent yellow-green (YG) microspheres. The antibody concentration is 250-350 μ g/ml depending on the lot.

Protocol for Use:

(This protocol is offered as a guide. Specific situations may require one or more alterations of this protocol.)

1. In order to remove the sodium azide and transfer to the appropriate buffer, the microspheres must be washed before use. To do this, aliquot 0.5ml of 1.25% Antibody Conjugated Microspheres into an Eppendorf centrifuge tube (1.5-1.9ml capacity). Add sufficient amount of buffer or cell culture media to fill the tube. (Use buffer or culture media that is compatible with the current application.) Centrifuge in a micro-centrifuge at 10,000xG for 5-6 minutes. Carefully remove the supernatant using a Pasteur pipette. Discard supernatant and resuspend microspheres in fresh buffer. Repeat this procedure 3 times in order to sufficiently wash the microspheres.

2. Incubate the microspheres with an appropriate amount of the target antigen for a minimum of 20-30 minutes at 4°C. This can also be performed at room temperature if necessary. The amount of target antigen can vary widely. The researcher is strongly encouraged to optimize the antibody/microsphere ratio prior to using the microspheres in their applications. Antigen amounts may range from 10-300 μ g.
3. To remove any excess antigen that is not attached to the microspheres, wash the microspheres one time using the centrifugation procedure outlined in Step 1. The microsphere-antigen complex is now ready for analysis or experimentation.

Protein A and Protein G Conjugated Microspheres

Protein A and Protein G are available on 1.0 μ m polystyrene microspheres. These proteins are available on plain, visible blue dyed, or fluorescent yellow-green (YG) dyed microspheres. The Protein concentration is 150-250 μ g/ml depending on the lot.

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for his own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.

Protocol for Use:

(This protocol is offered as a guide. Specific situations may require one or more alterations of this protocol.)

Protein A/G Buffer: 0.1M Tris-HCl, 0.15M NaCl

1. Add 100µl of Protein A or Protein G microspheres to a 1.5ml microcentrifuge tube.
2. Wash microspheres three times with 750µl of Protein A/G Buffer by mixing with buffer, centrifuging in a micro-centrifuge for 5-6 minutes at 10,000xG and removing the supernatant.
3. To the washed microspheres, add 50µl of serum containing the target IgG and 50µl of Protein A/G Buffer. Incubate for 1 hour at 4°C, vortexing every five minutes.
4. Centrifuge microspheres 5-6 minutes at 10,000xG in a micro-centrifuge. Remove the supernatant and discard.
5. Wash microspheres with 750µl of Protein A/G Buffer three times.
6. Add 50µl of 0.1M glycine, pH 2.5, vortex and incubate for five minutes. Centrifuge to separate the microspheres and save supernatant, which contains the purified IgG.

Ordering Information:

Catalog #	Description	Size
Goat Anti-Mouse IgG (H&L)		
17697-1	Blue Dyed Polystyrene Microspheres	1ml
17694-1	Undyed Microspheres	1ml
17843-1	Fluoresbrite YG Polystyrene Microspheres	1ml
Goat Anti-Rabbit IgG (H&L)		
17696-1	Blue Dyed Polystyrene Microspheres	1ml
17693-1	Undyed Microspheres	1ml
17844-1	Fluoresbrite YG Polystyrene Microspheres	1ml
Protein A		
17699-1	Blue Dyed Polystyrene Microspheres	1ml
17698-1	Undyed Microspheres	1ml
17845-1	Fluoresbrite YG Polystyrene Microspheres	1ml

Protein G

21105-1	Blue Dyed Polystyrene Microspheres	1ml
21106-1	Undyed Microspheres	1ml
21107-1	Fluoresbrite YG Polystyrene Microspheres	1ml

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com

Should any of our materials fail to perform to our specifications, we will be pleased to provide replacements or return the purchase price. We solicit your inquiries concerning all needs for life sciences work. The information given in this bulletin is to the best of our knowledge accurate, but no warranty is expressed or implied. It is the user's responsibility to determine the suitability for his own use of the products described herein, and since conditions of use are beyond our control, we disclaim all liability with respect to the use of any material supplied by us. Nothing contained herein shall be construed as a recommendation to use any product or to practice any process in violation of any law or any government regulation.