

BioMag 15ml/50ml Tube Separator

Description

Designed to accommodate five 15ml tubes and three 50ml tubes for a variety of applications including cell sorting and couplings of ligands to BioMag. Each bar magnet, embedded in a plastic rack, has a magnetic strength of 35 megagauss Oersteads.

Use

Simply place the 15ml or the 50ml tubes in the appropriate holder. The rare earth magnets will draw the BioMag solution toward the wall of the separator to complete a separation in 3 to 10 minutes.

15ml tubes: Separations are optimized when volumes of BioMag solution are between 7 and 12mls. Once a separation is complete aspirate the supernatant without disturbing the separated BioMag.

50ml tubes: Separations are optimized when volumes of BioMag solution are between 25 and 40mls. Once a separation is complete aspirate the supernatant without disturbing the separated BioMag.

Care

The BioMag Tube Separator unit should be wiped clean with mild detergent.

DO NOT IMMERSERSE. Harsh chemicals and organic solvents may cause damage to the surface of the unit.

Warning

Strong Magnetic Field. This separation unit should be used and stored at a safe distance from magnetic media such as computers, videotapes and credit cards etc.

Ordering Information:

Cat. #	Description	Size
84102S	BioMag 15ml/50ml Tube Separator	1 each

To Order:

In The U.S. Call: 1-800-523-2575 • 215-343-6484

In The U.S. FAX: 1-800-343-3291 • 215-343-0214

In Germany Call: (49) 6221-765767

In Germany FAX: (49) 6221-764620

Order online anytime at www.polysciences.com

