

# q•Blue Colorant

for use in PCR and qPCR

*for in vitro use only*

## q•Blue Colorant

Cat. No.	Description	Volume	No. Rxns
QBC-18	q•Blue Colorant	18 µL	100 x 20 µL
QBC-45	q•Blue Colorant	45 µL	250 x 20 µL
QBC-100	q•Blue Colorant	100 µL	555 x 20 µL

### Protocol:

Master Mix  
concentration      Protocol

2X	Mix 18 µL of q•Blue with 1 mL of 2X PCR or qPCR master mix
5X	Mix 45 µL of q•Blue with 1 mL of 5X PCR or qPCR master mix

*For use with ABI instruments, additional instrument settings may be required. Please contact [technical.support@qartabio.com](mailto:technical.support@qartabio.com) or your local sales representative for more information.*

### Description

q•Blue Colorant is an inert colorant for coloring one or more of the components of a PCR or qPCR reaction, such as the DNA template or master mix, for visible detection. It does NOT interfere with PCR amplification or the detection of PCR product using fluorescent dyes.

### Applications

- Coloring DNA templates for visible detection
- Coloring PCR reagents for visible detection
- Coloring qPCR reagents for visible detection

### Eco-friendly Shipping & Storage

- Shipped at ambient temperature without ice, foam, or other wasteful packaging
- Store at 4°C to -20°C
- Avoid repetitive freeze/thawing
- Once q•Blue is added to reagent, the mix can be stored at 4°C for one month

*Some applications of this product are covered by patents issued to parties other than qARTA Bio, and may require a license which is not provided by the purchase of this product. User should obtain a patent license if appropriate.*

T\_+1-310-326-1542    F\_+1-310-606-2220    E\_ [sales@qartabio.com](mailto:sales@qartabio.com)    W\_ [www.qartabio.com](http://www.qartabio.com)



Revision date: 2011.10.01  
© 2011 qARTA Bio, Inc.