


## DATA EXAMPLES

Differentiation of naïve CD4<sup>+</sup> T cells into Th1 cells is confirmed by flow cytometry (Figure 1) and ELISA (Figure 2). These data demonstrate that Th1-polarized cells express and secrete IFN- $\gamma$ . The corresponding tests for IL-4 (Th2 cell marker) and IL-17 (Th17 cell marker) are low/negative.




**Figure 1: Intracellular Cytokine Staining of Differentiated Mouse Th1 Cells.** Flow cytometry data of mouse naïve CD4<sup>+</sup> T cells (A, C) and differentiated Th1 cells (B, D) generated using reagents included in this kit. After 6 days of differentiation, naïve CD4<sup>+</sup> cells or differentiated Th1 cells were stimulated with Cell Activation Cocktail (Tocris®, Catalog # 5476) and stained with conjugated Anti-Mouse IFN- $\gamma$ , Anti-Mouse IL-4 (Clone 11B11), and Anti-Mouse IL-17A Monoclonal Antibodies. Quadrants were set based on isotype-stained controls. All R&D Systems® and Novus Biologicals® antibodies and corresponding catalog numbers used in this figure are shown in the table below.

## SUGGESTED REAGENTS FOR FLOW CYTOMETRY

CATALOG #	DESCRIPTION
IC485P, and IC006P	Mouse IFN- $\gamma$ PE-conjugated Antibody (R&D Systems®), and Rat IgG <sub>2a</sub> PE-conjugated Antibody (R&D Systems®)
NBP1-72027, and IC006P	Mouse IL-17/IL-17A PE-conjugated Antibody (Novus Biologicals®), and Rat IgG <sub>2a</sub> PE-conjugated Antibody (R&D Systems®)
FAB554P, and IC005P	Mouse CD4 PE-conjugated Antibody (R&D Systems®), and Rat IgG <sub>2a</sub> PE-conjugated Antibody (R&D Systems®)
FC009	Flow Cytometry Fixation & Permeabilization Buffer (R&D Systems®)
5476	Cell Activation Cocktail 500X (Tocris®)

## DATA EXAMPLES CONTINUED



**Figure 2: Th1-differentiated Mouse CD4<sup>+</sup> Cells Secrete High Levels of IFN- $\gamma$ .** Mouse naïve CD4<sup>+</sup> T cells were differentiated for 6 days using the reagents included in this kit. On day 6 of differentiation cells were harvested and re-stimulated with Anti-Mouse CD3 and Anti-Mouse CD28 overnight. The cell culture supernatant was collected and cytokine secretion was determined using the Mouse IFN- $\gamma$  Quantikine® ELISA Kit, the Mouse IL-4 Quantikine® ELISA Kit, and the Mouse IL-17 Quantikine® ELISA Kit. All relevant Quantikine® and DuoSet® ELISA kits and corresponding R&D Systems® catalog numbers are listed in the table below.

## SUGGESTED REAGENTS FOR ELISA

CATALOG #	DESCRIPTION
M1700, or DY421	Mouse IL-17 Quantikine® ELISA Kit, or Mouse IL-17 DuoSet® ELISA
MIF00, or DY485	Mouse IFN- $\gamma$ Quantikine® ELISA Kit, or Mouse IFN- $\gamma$ DuoSet® ELISA
M4000B, or DY404	Mouse IL-4 Quantikine® ELISA Kit, or Mouse IL-4 DuoSet® ELISA

## REFERENCES

- Zhu, J. and W.E. Paul (2010) Immunol. Rev. **238**:247.
- Dardalhon, V. *et al.* (2008) J. Autoimmun. **31**:252.
- Seder, R.A. and W.E. Paul (1994) Annu. Rev. Immunol. **12**:635.

CellXVivo™

## Mouse Th1 Cell Differentiation Kit

Catalog Number: CDK018

## BACKGROUND

T helper type 1 (Th1) cells are a lineage of CD4<sup>+</sup> effector T cells that promote cell-mediated immune responses against intracellular viral and bacterial pathogens (1). Enhanced Th1 responses are associated with autoimmune diseases including rheumatoid arthritis, multiple sclerosis, inflammatory bowel disease, and type-1 diabetes (2). Differentiation of CD4<sup>+</sup> effector T cells into the Th1 lineage is promoted by cytokines such as IL-12 and IFN- $\gamma$  (3). Th1 cells secrete IFN- $\gamma$ , IL-10, and TNF- $\alpha$ . The CellXVivo™ Mouse Th1 Cell Differentiation Kit contains optimized reagents for Th1 differentiation from naïve CD4<sup>+</sup> cells. The quantity of components in this kit is sufficient to differentiate approximately  $8 \times 10^6$  naïve CD4<sup>+</sup> T cells, and generate  $8 \times 10^7$  CD4<sup>+</sup> cells of which  $\geq 70\%$  are IFN- $\gamma$ <sup>+</sup> Th1 polarized cells.

This package insert must be read in its entirety before using this product. For research use only. Not for use in diagnostic procedures.

### MANUFACTURED AND DISTRIBUTED BY:

**USA & Canada | R&D Systems, Inc.**

614 McKinley Place NE, Minneapolis, MN 55413, USA  
TEL: (800) 343-7475 (612) 379-2956 FAX: (612) 656-4400  
E-MAIL: info@RnDSystems.com

### DISTRIBUTED BY:

**UK & Europe | R&D Systems Europe, Ltd.**

19 Barton Lane, Abingdon Science Park, Abingdon OX14 3NB, UK  
TEL: +44 (0)1235 529449 FAX: +44 (0)1235 533420  
E-MAIL: info@RnDSystems.co.uk

**China | R&D Systems China Co., Ltd.**

24A1 Hua Min Empire Plaza, 726 West Yan An Road, Shanghai PRC 200050  
TEL: +86 (21) 52380373 FAX: +86 (21) 52371001  
E-MAIL: info@RnDSystemsChina.com.cn

## MATERIALS PROVIDED & STORAGE CONDITIONS

Store the unopened kit at  $\leq -20^{\circ}\text{C}$ . Do not use past kit expiration date.

COMPONENTS	PART #	# VIALS	STORAGE OF OPENED/ RECONSTITUTED MATERIAL
Hamster Anti-Mouse CD3, Th1	968107	1 vial	May be stored at $2-8^{\circ}\text{C}$ under sterile conditions for up to 1 month or at $-20^{\circ}\text{C}$ to $-70^{\circ}\text{C}$ in a manual defrost freezer for up to 3 months.*
Mouse Th1 Reagent 1	968108	1 vial	
Mouse Th1 Reagent 2	968109	1 vial	
Mouse Th1 Reagent 3	968110	1 vial	
Mouse Th1 Reagent 4	968111	1 vial	
Reconstitution Buffer 1	967552	2 vials	May be stored under sterile conditions for up to 3 months at $2-8^{\circ}\text{C}$ .*
Reconstitution Buffer 2	967553	1 vial	
20X Wash Buffer	967557	1 vial	

\* Provided this is within the expiration date of the kit.

## OTHER MATERIALS & SUPPLIES REQUIRED

- Laboratory mice
- MagCollect™ Mouse Naïve CD4<sup>+</sup> T Cell Isolation Kit (R&D Systems®, Catalog # MAGM205, or equivalent).
- X-VIVO™ 15 Chemically Defined, Serum-free Hematopoietic Cell Medium (Lonza, or equivalent)
- Penicillin/Streptomycin (optional)
- Sterile deionized water
- Cell Activation Cocktail 500X (Tocris, Catalog #5476)
- Tissue culture flasks and/or plates
- Pipettes and pipette tips
- Inverted microscope
- Hemocytometer
- $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$  incubator
- Centrifuge

## REAGENT PREPARATION

**Mouse Th1 Reagent 1 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 1 to Mouse Th1 Reagent 1 to produce Mouse Th1 Reagent 1 (400X).

**Mouse Th1 Reagent 2 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 1 to Mouse Th1 Reagent 2 to produce Mouse Th1 Reagent 2 (400X).

**Mouse Th1 Reagent 3 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 1 to Mouse Th1 Reagent 3 to produce Mouse Th1 Reagent 3 (400X).

**Mouse Th1 Reagent 4 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 2 to Mouse Th1 Reagent 4 to produce Mouse Th1 Reagent 4 (400X).

**Mouse Th1 Differentiation Media** - Add Mouse Th1 Reagent 1 (400X), Mouse Th1 Reagent 2 (400X), Mouse Th1 Reagent 3 (400X), and Mouse Th1 Reagent 4 (400X) to a final concentration of 1X in the desired amount of X-VIVO™ 15 Medium. Addition of Penicillin (100 units/mL) and Streptomycin (100  $\mu\text{g}/\text{mL}$ ) is recommended. Make fresh as needed.

**Hamster Anti-Mouse CD3 (50X)** - Add 100  $\mu\text{L}$  of Reconstitution Buffer 2 to Hamster Anti-Mouse CD3, Th1 Antibody to produce Hamster Anti-Mouse CD3 (50X).

**Wash Buffer (1X)** - Dilute 20X Wash Buffer 1:20 in sterile deionized water to produce Wash Buffer (1X).

**Hamster Anti-Mouse CD3 (1X)** - Dilute Hamster Anti-Mouse CD3 (50X) stock 1:50 with Wash Buffer (1X) to produce Hamster Anti-Mouse CD3 (1X).

## PROTOCOL FOR Th1 DIFFERENTIATION

The quantity of components in this kit is sufficient to differentiate approximately  $8 \times 10^6$  naïve CD4<sup>+</sup> T cells, and generate  $8 \times 10^7$  CD4<sup>+</sup> cells of which 70-90% are IFN- $\gamma$ <sup>+</sup> Th1 polarized cells.

**Note:** *If starting with fewer cells, adjust starting volumes/number of wells accordingly.*

- Coat the desired tissue culture plate with Hamster Anti-Mouse CD3 (1X).
  - Add Hamster Anti-Mouse CD3 (1X) to plate using the suggested coating volumes below.
  - Incubate overnight at  $2-8^{\circ}\text{C}$  or 2-3 hours at  $37^{\circ}\text{C}$ .
  - Wash plate or flask twice with Wash Buffer (1X) just prior to adding cells.

PLATE	SUGGESTED COATING VOLUME	# OF WELLS PER KIT
24-well plate	0.4 mL/well	8
96-well plate	0.1 mL/well	40

- Prepare a single cell suspension of mouse splenocytes and isolate mouse naïve CD4<sup>+</sup> T cells according to the product insert for the MagCollect™ Mouse Naïve CD4<sup>+</sup> T Cell Isolation Kit. Perform a cell count.
 

**Note:** *1 mouse spleen will provide roughly enough naïve CD4<sup>+</sup> T cells for 1 well of a 24-well plate. The quantity of spleens needed may vary based on mouse strain, age, and/or health.*
- Suspend mouse naïve CD4<sup>+</sup> T cells at  $1 \times 10^6$  cells/mL in Mouse Th1 Differentiation Media.
- Add the cells to a Hamster Anti-Mouse CD3 antibody-coated plate using the suggested volumes below.

PLATE	SUGGESTED VOLUME
24-well plate	1.0 mL/well
96-well plate	0.2 mL/well

- Centrifuge the plate at 300 x g for 1 minute and incubate the cells in a  $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$  humidified incubator for 3 days.
- On day 3 of differentiation, harvest cells and dilute them 1:10 by adding fresh Mouse Th1 Differentiation Media in an appropriate sized conical tube.
- Transfer diluted cells to a new, uncoated plate or flask using the volumes indicated in the table below. Incubate the cells in a  $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$  humidified incubator for an additional 3 days.

PLATE OR FLASK	SUGGESTED VOLUME
24-well plate	1 mL/well
6-well plate	4 mL/well
T-25	8-10 mL/well
T-75	20 mL/well

## PROTOCOL FOR Th1 DIFFERENTIATION CONTINUED

- On day 6 of differentiation, the differentiated mouse Th1 cells are ready to be used for downstream applications.
- To verify Th1 cell differentiation via flow cytometry, collect the cells and wash with X-VIVO™ 15 Medium once, resuspend the cells in 1 mL X-VIVO™ 15 Medium and Cell Activation Cocktail (1X). Addition of Penicillin (100 units/mL) and Streptomycin (100  $\mu\text{g}/\text{mL}$ ) is recommended. Incubate the cells in a  $37^{\circ}\text{C}$ , 5%  $\text{CO}_2$  humidified incubator for 4-5 hours. Analyze cytokine expression via flow cytometry.

## PROTOCOL OUTLINE


**Coat** the desired tissue culture plate with Hamster Anti-Mouse CD3, Th1 antibody.




**Isolate** mouse splenocytes.




**Isolate** mouse naïve CD4<sup>+</sup> T cells from splenocytes (e.g., using magnetic cell selection).




**Perform** a cell count.



**Suspend**  $1 \times 10^6$  naïve CD4<sup>+</sup> T cells/mL in Mouse Th1 Differentiation Media. Culture the cells on plates pre-coated with Hamster Anti-Mouse CD3 antibody.



**Harvest** cells on day 3. **Dilute** cells 1:10 with fresh Mouse Th1 Differentiation Media. **Culture** cells in a new plate or flask for 3 days.



**Verify** Th1 cell differentiation by analyzing cytokine expression via flow cytometry or ELISA (optional).

