


## DATA EXAMPLES

Differentiation of naïve CD4<sup>+</sup> T cells into Th2 cells is confirmed by flow cytometry (Figure 1) and ELISA (Figure 2). These data demonstrate that Th2-polarized cells express and secrete IL-4. The corresponding tests for IFN- $\gamma$  (Th1 cell marker) and IL-17 (Th17 cell marker) are low/negative.




**Figure 1: Intracellular Cytokine Staining of Differentiated Mouse Th2 Cells.** Flow cytometry data of mouse naïve CD4<sup>+</sup> T cells (A, C) and differentiated Th2 cells (B, D) generated using reagents included in this kit. After 6 days of differentiation, naïve CD4<sup>+</sup> cells or differentiated Th2 cells were stimulated with Cell Activation Cocktail (Tocris®, Catalog # 5476) and stained with conjugated Anti-Mouse IL-4 (Clone 11B11), Anti-Mouse IFN- $\gamma$ , and Anti-Mouse IL-17A Monoclonal Antibodies. Quadrants were set based on isotype-stained controls. All R&D Systems® and Novus Biologicals® antibodies and corresponding catalog numbers used in this figure are shown in the table below.

## SUGGESTED REAGENTS FOR FLOW CYTOMETRY

CATALOG #	DESCRIPTION
IC485P, and IC006P	Mouse IFN- $\gamma$ PE-conjugated Antibody (R&D Systems®), and Rat IgG <sub>2a</sub> PE-conjugated Antibody (R&D Systems®)
NBP1-72027, and IC006P	Mouse IL-17/IL-17A PE-conjugated Antibody (Novus Biologicals®), and Rat IgG <sub>2a</sub> PE-conjugated Antibody (R&D Systems®)
FAB554P, and IC005P	Mouse CD4 PE-conjugated Antibody (R&D Systems®), and Rat IgG <sub>1</sub> PE-conjugated Antibody (R&D Systems®)
FC009	Flow Cytometry Fixation & Permeabilization Buffer (R&D Systems®)
5476	Cell Activation Cocktail 500X (Tocris®)

## DATA EXAMPLES CONTINUED



**Figure 2: Th2-differentiated Mouse CD4<sup>+</sup> Cells Secrete High Levels of IL-4.** Mouse naïve CD4<sup>+</sup> T cells were differentiated for 6 days using the reagents included in this kit. On day 6 of differentiation cells were harvested and restimulated with Anti-Mouse CD3 and Anti-Mouse CD28 overnight. Cell culture supernatant was collected and cytokine secretion determined using the Mouse IL-4 Quantikine® ELISA Kit, the Mouse IFN- $\gamma$  Quantikine® ELISA Kit, and the Mouse IL-17 Quantikine® ELISA Kit. All relevant Quantikine® and DuoSet® ELISA kits and corresponding R&D Systems catalog numbers are listed in the table below.

## SUGGESTED REAGENTS FOR ELISA

CATALOG #	DESCRIPTION
M1700, or DY421	Mouse IL-17 Quantikine® ELISA Kit, or Mouse IL-17 DuoSet® ELISA
MIF00, or DY485	Mouse IFN- $\gamma$ Quantikine® ELISA Kit, or Mouse IFN- $\gamma$ DuoSet® ELISA
M4000B, or DY404	Mouse IL-4 Quantikine® ELISA Kit, or Mouse IL-4 DuoSet® ELISA

## REFERENCES

- Li, Z. *et al.* (2013) Prot. Cell **2**:604.
- Luckheeram, R.V. *et al.* (2012) Clin. Dev. Immunol. **2012**:925135.
- Hirahara, K. *et al.* (2011) Immunology **134**:235.3.

CellXVivo™

## Mouse Th2 Cell Differentiation Kit

Catalog Number: CDK019

## BACKGROUND

CD4<sup>+</sup> T cells differentiate into separate T helper cells under the influence of various cytokines and cellular interactions that induce expression of specific transcription factors. T helper type 2 (Th2) cells are a lineage of CD4<sup>+</sup> effector T cells that provide host protection against intestinal parasites and extracellular bacteria. In addition, they provide support for B cell-dependent humoral responses. Pathological Th2 cell activity is also a hallmark of allergic inflammation and asthma (1). Differentiation of naïve CD4<sup>+</sup> T cells into the Th2 lineage is promoted by cytokines such as IL-4 in combination with either IL-2, IL-7, or TSLP (2, 3). Th2 cells secrete IL-4, IL-5, IL-9, IL-13, and IL17E/IL-25. The CellXVivo™ Mouse Th2 Cell Differentiation Kit contains the necessary components to differentiate mouse naïve CD4<sup>+</sup> T cells into Th2-polarized cells that are IL-4<sup>+</sup> IFN- $\gamma$  IL-17. The quantity of the components in the kit is sufficient to differentiate approximately 4x10<sup>6</sup> naïve CD4<sup>+</sup> T cells, and generate 80 x 10<sup>6</sup> CD4<sup>+</sup> T cells of which > 40% are IL-4<sup>+</sup> Th2 polarized cells.

This package insert must be read in its entirety before using this product. For research use only. Not for use in diagnostic procedures.

### MANUFACTURED AND DISTRIBUTED BY:

**USA & Canada | R&D Systems, Inc.**

614 McKinley Place NE, Minneapolis, MN 55413, USA  
TEL: (800) 343-7475 (612) 379-2956 FAX: (612) 656-4400  
E-MAIL: info@RnDSystems.com

### DISTRIBUTED BY:

**UK & Europe | R&D Systems Europe, Ltd.**

19 Barton Lane, Abingdon Science Park, Abingdon OX14 3NB, UK  
TEL: +44 (0)1235 529449 FAX: +44 (0)1235 533420  
E-MAIL: info@RnDSystems.co.uk

**China | R&D Systems China Co., Ltd.**

24A1 Hua Min Empire Plaza, 726 West Yan An Road, Shanghai PRC 200050  
TEL: +86 (21) 52380373 FAX: +86 (21) 52371001  
E-MAIL: info@RnDSystemsChina.com.cn

## MATERIALS PROVIDED & STORAGE CONDITIONS

Store the unopened kit at  $\leq -20^{\circ}\text{C}$ . Do not use past kit expiration date.

COMPONENTS	PART #	# VIALS	STORAGE OF OPENED/ RECONSTITUTED MATERIAL
Hamster Anti-Mouse CD3, Th2	968098	1 vial	May be stored at 2-8 °C under sterile conditions for up to 1 month or at -20 °C to -70 °C in a manual defrost freezer for up to 3 months.*
Rat Anti-Mouse CD28, Th2	968099	1 vial	
Mouse Th2 Reagent 1	968100	1 vial	
Mouse Th2 Reagent 2	968101	1 vial	
Mouse Th2 Reagent 3	968102	1 vial	
Reconstitution Buffer 1	967552	1 vial	May be stored under sterile conditions for up to 3 months at 2-8 °C.*
Reconstitution Buffer 2	967553	2 vials	
20X Wash Buffer	967557	3 vials	

\* Provided this is within the expiration date of the kit.

## OTHER MATERIALS & SUPPLIES REQUIRED

- Laboratory mice
- MagCollect™ Mouse Naive CD4<sup>+</sup> T Cell Isolation Kit (R&D Systems®, Catalog # MAGM205, or equivalent).
- X-VIVO™ 15 Chemically Defined, Serum-free Hematopoietic Cell Medium (Lonza, or equivalent)
- Penicillin/Streptomycin (optional)
- Sterile deionized water
- Cell Activation Cocktail 500X (Tocris, Catalog #5476)
- Tissue culture flasks and/or plates
- Pipettes and pipette tips
- Inverted microscope
- Hemocytometer
- 37 °C, 5% CO<sub>2</sub> incubator
- Centrifuge

## REAGENT PREPARATION

**Mouse Th2 Reagent 1 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 1 to Mouse Th2 Reagent 1 to produce Mouse Th2 Reagent 1 (400X).

**Mouse Th2 Reagent 2 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 1 to Mouse Th2 Reagent 2 to produce Mouse Th2 Reagent 2 (400X).

**Mouse Th2 Reagent 3 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 2 to Mouse Th2 Reagent 3 to produce Mouse Th2 Reagent 3 (400X).

**Rat Anti-Mouse CD28 (400X)** - Add 250  $\mu\text{L}$  of Reconstitution Buffer 2 to Rat Anti-Mouse CD28, Th2, to produce Rat Anti-Mouse CD28 (400X).

**Mouse Th2 Differentiation Media** - Add Rat Anti-Mouse CD28 (400X), Mouse Th2 Reagent 1 (400X), Mouse Th2 Reagent 2 (400X), and Mouse Th2 Reagent 3 (400X) to a final concentration of 1X in the desired amount of X-VIVO™ 15 Medium. Addition of Penicillin (100 units/mL) and Streptomycin (100  $\mu\text{g}/\text{mL}$ ) is recommended. Make fresh as needed.

**Hamster Anti-Mouse CD3 (100X)** - Add 150  $\mu\text{L}$  of Reconstitution Buffer 2 to Hamster Anti-Mouse CD3, Th2, to produce Hamster Anti-Mouse CD3 (100X).

**Wash Buffer (1X)** - Dilute 20X Wash Buffer 1:20 in sterile deionized water to produce Wash Buffer (1X).

**Hamster Anti-Mouse CD3 (1X)** - Dilute Hamster Anti-Mouse CD3 (100X) stock 1:100 with Wash Buffer (1X) to produce Hamster Anti-Mouse CD3 (1X). Make fresh as needed.

## PROTOCOL FOR Th2 DIFFERENTIATION

The quantity of the components in the kit is sufficient to differentiate approximately  $4 \times 10^6$  naive CD4<sup>+</sup> T cells, and generate  $80 \times 10^6$  CD4<sup>+</sup> T cells of which > 40% are IL-4<sup>+</sup> Th2 polarized cells.

**Note:** If starting with fewer cells, adjust starting volumes/number of wells accordingly.

- Coat the desired tissue culture plate with Hamster Anti-Mouse CD3 (1X).
  - Add Hamster Anti-Mouse CD3 (1X) to plate using the suggested coating volumes below:

TISSUE CULTURE PLATE	SUGGESTED COATING VOLUME	# OF WELLS PER KIT
24-well plate	0.4 mL/well	Up to 3
96-well plate	0.1 mL/well	Up to 20

- Incubate overnight at 2-8 °C or 2-3 hours at 37 °C.
  - Wash plate twice with Wash Buffer (1X) just prior to adding cells.
- Prepare a single cell suspension of mouse splenocytes and isolate the mouse naive CD4<sup>+</sup> T cells according to the product insert for the MagCollect™ Mouse Naive CD4<sup>+</sup> T Cell Isolation Kit. Perform a cell count.
 

**Note:** 1 mouse spleen will provide roughly enough naive CD4<sup>+</sup> T cells for 1 well of a 24-well plate. The quantity of spleens needed may vary based on mouse strain, age, and/or health.
  - Suspend mouse naive CD4<sup>+</sup> T cells at  $1 \times 10^6$  cells/mL in Mouse Th2 Differentiation Media.
  - Add the cells to a Hamster Anti-Mouse CD3 antibody-coated plate using the suggested volumes below.

TISSUE CULTURE PLATE	SUGGESTED VOLUME
24-well plate	1.3 mL/well
96-well plate	0.2 mL/well


- Centrifuge the plate at 300 x g for 1 minute and incubate the cells in a 37 °C, 5% CO<sub>2</sub> humidified incubator for 3 days.
- On day 3 of differentiation, harvest cells from each well and dilute around 1:20 in fresh Mouse Th2 Differentiation Media in an appropriate sized conical tube.
- Transfer diluted cells to a new, uncoated flask using the volumes indicated in the table below (ie., 1 well from a 24-well plate can be expanded into 1 T-75 flask or 3 T-25 flasks). Incubate the cells in a 37 °C, 5% CO<sub>2</sub> humidified incubator for 3 days.

FLASK	SUGGESTED VOLUME
T-25	8 mL/flask
T-75	25 mL/flask


- On day 6 of differentiation, the differentiated mouse Th2 cells are ready to be used for downstream applications.
- To verify Th2 cell differentiation via flow cytometry, collect the cells and wash with X-VIVO™ 15 Medium once, resuspend the cells in 1 mL X-VIVO™ 15 Medium and Cell Activation Cocktail (1X). Addition of Penicillin (100 units/mL) and Streptomycin (100  $\mu\text{g}/\text{mL}$ ) is recommended. Incubate the cells in a 37 °C, 5% CO<sub>2</sub> humidified incubator for 4-5 hours. Analyze cytokine expression via flow cytometry (See Figure 1).

## PROTOCOL OUTLINE


**Coat** the desired tissue culture plate with Hamster Anti-Mouse CD3 antibody.




**Isolate** mouse splenocytes.




**Isolate** mouse naive CD4<sup>+</sup> T cells from splenocytes (e.g., using magnetic cell selection).




**Perform** a cell count.



**Suspend**  $1 \times 10^6$  naive CD4<sup>+</sup> T cells/mL in Mouse Th2 Differentiation Media. Culture the cells on plates pre-coated with Hamster Anti-Mouse CD3 antibody.



**Harvest** cells on day 3. **Dilute** cells 1:20 with fresh Mouse Th2 Differentiation media. **Culture** cells in a new flask for an additional 3 days.



**Verify** Th2 cell differentiation by analyzing cytokine expression via flow cytometry or ELISA (optional).

