

## PRODUCT DESCRIPTION

This kit contains three conjugated antibodies (and a rabbit IgG-PE isotype control) that can be used for the single-step staining of mouse regulatory T cells (1-9).

## MATERIALS PROVIDED & STORAGE

Store the unopened kit at 2-8 °C **in the dark**. Use within 6 months of receipt.

	PART	CATALOG #/PART #	DESCRIPTION
Positive Markers	Human/Mouse FoxP3-PE Rabbit IgG; Clone 1054C	IC8214P	250 µL of a phycoerythrin-conjugated antibody specific for human/mouse FoxP3.
	Mouse CD4-FITC Rat IgG <sub>2b</sub> ; Clone GK1.5	FAB554F	250 µL of a fluorescein-conjugated antibody specific for mouse CD4.
	Mouse IL-2 Ra/CD25-APC Rat IgG <sub>2a</sub> ; Clone 280406	FAB2438A	250 µL of an allophycocyanin-conjugated antibody specific for mouse IL-2 Ra/CD25.
Isotype Control	Rabbit IgG Control-PE	IC105P	250 µL of a phycoerythrin-conjugated control specific for rabbit IgG.
Staining Buffers	FoxP3 Fixation Concentrate (4X)	894065	8 mL of a formaldehyde solution.
	FoxP3 Fixation Diluent	894068	25 mL of a buffered detergent.
	FoxP3 Permeabilization and Wash Buffer (10X)	894356	25 mL of a buffered protein base with detergent and preservatives. <b>Note:</b> <i>May contain a precipitate but will not affect product performance.</i>
	Staining Buffer (1X)	895068	50 mL of a 1X Staining Buffer


## INTENDED USE

This product is designed for the flow cytometric analysis of mouse regulatory T cells using three fluorochrome-conjugated antibodies.

## STAINING PROTOCOL

1. Wash mouse splenocytes ( $1 \times 10^6$  cells per sample) with 2 mL of Flow Cytometry Staining Buffer (R&D Systems, Catalog # FC001) or other BSA-containing buffer, by spinning at 300 x g for 5 minutes, using 5 mL flow cytometry tubes.
2. Fc-block cells with blocking IgG ( $1 \mu\text{g}$  IgG/ $10^6$  cells) for 10 minutes at 2-8 °C.
3. Add 10 µL of CD4-FITC and 10 µL of CD25-APC antibodies or isotype controls (R&D Systems, Catalog # IC013F and IC006A).
4. Incubate the mixture for 30-45 minutes at 2-8 °C **in the dark**.
5. Wash the cells two times with **cold** 1X PBS. During the washes, make up fresh 1X FoxP3 Fixation Buffer by diluting FoxP3 Fixation Concentrate (4X) with FoxP3 Fixation Diluent (*i.e.* 100 µL FoxP3 Fixation Concentrate (4X) + 300 µL FoxP3 Fixation Diluent).
6. Resuspend the cells in fresh 1X FoxP3 Fixation Buffer using 0.5 mL/tube. Incubate at 2-8 °C for 30 minutes. During this incubation, make up 1X FoxP3 Permeabilization and Wash Buffer by diluting FoxP3 Permeabilization and Wash Buffer (10X) with distilled water (*i.e.* 100 µL FoxP3 Permeabilization and Wash Buffer (10X) + 900 µL diH<sub>2</sub>O) and keep at 2-8 °C.
7. Wash two times with fresh, cold, 1X FoxP3 Permeabilization and Wash Buffer.
8. Add 10 µL of FoxP3 antibody or the rabbit IgG-PE isotype control included in this kit to the cells and incubate for 30 minutes at 2-8 °C.
9. Wash the cells one time with cold 1X FoxP3 Permeabilization and Wash Buffer.
10. Resuspend the cells in Flow Cytometry Staining Buffer and run on a flow cytometer.

## DATA EXAMPLES



**Figure 1: Detection of FoxP3 in CD57/B6 Mouse Splenocytes Natural Tregs by Flow Cytometry.** CD57/B6 mouse splenocyte natural regulatory T cells (Tregs) were surface stained with **(A)** Rat Anti-Mouse CD4 Fluorescein-conjugated Monoclonal Antibody (R&D Systems, Catalog # FAB554F) and **(B)** Rat Anti-Mouse IL-2 R $\alpha$ /CD25 APC-conjugated Monoclonal Antibody (R&D Systems, Catalog # FAB2438A), followed by intracellular staining using Rabbit Anti-Human/Mouse FoxP3 PE-conjugated Antigen Affinity-purified Monoclonal Antibody (R&D Systems, Catalog # IC8214P). To facilitate intracellular staining, cells were fixed and permeabilized with FlowX FoxP3 Fixation & Permeabilization Buffer Kit (R&D Systems, Catalog # FC012). Cells were gated on lymphocytes.

## PRECAUTIONS

FoxP3 Staining Buffer contains formaldehyde which is a suspected carcinogen. Avoid contact with skin, eyes, and mucous membranes, and avoid inhaling fumes. In case of contact, wash immediately with water and seek medical advice.

The Staining Buffers in this kit contain sodium azide which may react with lead and copper plumbing to form explosive metallic azides. Flush with large volumes of water during disposal.

## REFERENCES

1. Tang, Q. and J.A. Bluestone (2008) *Nature Immunol.* **9**:239.
2. Zheng, Y. and A.Y. Rudensky (2007) *Nature Immunol.* **8**:457.
3. Hori, S. *et al.* (2003) *Science* **299**:1057.
4. Kim, J.M. and A.Y. Rudensky (2006) *Immunol. Rev.* **212**:86.
5. Fontenot, J.D. *et al.* (2003) *Nat. Immunol.* **4**:330.
6. Khatri, R. *et al.* (2003) *Nat. Immunol.* **4**:337.
7. Bluestone, J.A. and A.K. Abbas (2003) *Nat. Rev. Immunol.* **3**:253.
8. Sakaguchi, S. (2004) *Ann. Rev. Immunol.* **22**:531.
9. Josefowicz, S.Z. and A.Y. Rudensky (2009) *Immunity* **30**:616.