



# Magnetic Luminex® Performance Assay Human CCL4/MIP-1 $\beta$ Kit

**Catalog Number:** LUHM271

**Pack Size:** 100 Tests

## SPECIFICATIONS AND USE

### Recommended Sample Types

### Microparticle Region

### Components

### Other Supplies Required


### Storage

### Instructions for Use

- Cell culture supernates, serum, EDTA plasma, and heparin plasma.
- Region-35
- Microparticle Concentrate (Part 894446) is supplied as a 100X concentrated stock (0.075 mL) with preservatives.
- Biotin-Antibody Concentrate (Part 892633) is supplied as a 100X concentrated stock solution (0.075 mL) with preservatives.
- Magnetic Luminex Performance Assay Human Base Kit A (Catalog Number LUHM000).
- Store the unopened kit at 2-8 °C. Do not use past the expiration date on the label.
- **Avoid freezing microparticles.**
- **Protect microparticles from light.**
- Refer to the Base Kit insert for the Luminex Performance Assay procedure.

## TYPICAL DATA

This human MIP-1 $\beta$  standard curve is provided only for demonstration. A standard curve must be generated each time an assay is run, utilizing values from the Standard Value Card included in the Base Kit.



Standard	pg/mL	MFI	Average	Corrected
Blank	0	18 18	18	—
1	6200	17,067 17,148	17,108	17,090
2	2067	10,821 10,832	10,827	10,809
3	689	4228 4291	4260	4242
4	230	1463 1476	1469	1451
5	77	403 406	405	387
6	26	114 115	115	97
7	9	52 52	52	34

## PERFORMANCE CHARACTERISTICS

**All data were collected with assays run as a multiplex.**

**Data obtained with polystyrene and magnetic beads were equivalent.**

**Sensitivity** - The Minimum Detectable Dose (MDD) was determined by adding two standard deviations to the MFI of twenty zero standard replicates and calculating the corresponding concentration.

Fourteen assays were evaluated, and the MDD of human MIP-1 $\beta$  ranged from 0.23-0.72 pg/mL. The mean MDD was 0.44 pg/mL.

**Intra-assay Precision (precision within an assay)** - Three samples of known concentration were tested twenty times on one plate to assess precision within an assay.

**Inter-assay Precision (precision between assays)** - Three samples of known concentration were tested in twenty-five separate assays to assess precision between assays.

	Intra-assay Precision				Inter-assay Precision		
Sample	1	2	3		1	2	3
n	20	20	20		20	20	20
Mean (pg/mL)	442	929	1672		440	921	1687
Standard Deviation	19	40	100		22	72	269
% CV	4.3	4.4	6.0		5.0	7.8	15.9

**Recovery and Linearity** – Samples containing and/or spiked with high concentrations of MIP-1 $\beta$  were evaluated for recovery and were serially diluted to evaluate assay linearity.

Recovery			Linearity					
Sample Type	Average % Recovery	Range (%)			Cell culture supernates	Serum	EDTA Plasma	Heparin Plasma
Cell culture supernates	101	87-119	1:2	Average % of Expected	101	107	98	83
				Range (%)	83-123	94-127	70-121	66-98
Serum	108	95-120	1:4	Average % of Expected	101	104	99	82
				Range (%)	79-120	92-137	87-122	71-96
EDTA plasma	106	86-131	1:8	Average % of Expected	106	104	94	85
				Range (%)	90-128	91-133	71-134	72-97
Heparin plasma	96	76-121						

**Specificity** - This assay recognizes natural and recombinant human MIP-1 $\beta$ . The assay was tested for cross-reactivity and interference with the following factors. Less than 0.5% cross-reactivity and interference was observed.

Recombinant human:			Recombinant mouse:		Recombinant rat:	Recombinant porcine:	Recombinant human multiplex partners:	
6Ckine	IL-1 RII	IL-17	G-CSF	IL-8	GM-CSF	GM-CSF	ENA-78	IL-6
CNTF	IL-2 R $\alpha$	IL-18	GM-CSF	IL-10	IFN- $\gamma$	IL-1 $\alpha$	FGF basic	IL-8
$\beta$ -ECGF	IL-2 R $\beta$	LIF	IFN- $\gamma$	IL-17	IL-1 $\alpha$	IL-1 $\beta$	G-CSF	IL-10
FGF acidic	IL-2 R $\gamma$	LIF R	IL-1 $\alpha$	MIP-1 $\alpha$	IL-1 $\beta$	IL-2	GM-CSF	IL-17
FGF-4	IL-3 R $\alpha$	MIP-1 $\alpha$	IL-1ra	MIP-1 $\beta$	IL-2	IL-4	IFN- $\gamma$	MCP-1
FGF-5	IL-4 R	MIP-3 $\alpha$	IL-1	RANTES	IL-4	IL-6	IL-1 $\alpha$	MIP-1 $\alpha$
FGF-6	IL-5 R $\alpha$	MIP-3 $\beta$	IL-2	Tpo	IL-6	IL-8	IL-1 $\beta$	RANTES
FGF-9	IL-6 R	MCP-2	IL-4	TNF- $\alpha$	IL-10	IL-10	IL-1ra	Tpo
FGF-10	IL-10 R	MCP-3	IL-5	VEGF	TNF- $\alpha$	Leptin	IL-2	TNF- $\alpha$
FGF-18	IL-3	MCP-4	IL-6			TNF- $\alpha$	IL-4	VEGF
GCP-2	IL-7	M-CSF					IL-5	
GRO $\alpha$	IL-9	TNF RI						
GRO $\beta$	IL-11	TNF- $\alpha$						
GRO $\gamma$	IL-12 p40	VEGF <sub>121</sub>						
I-309	IL-12 p70	VEGF <sub>165</sub>						
IGF-I	IL-13	VEGF-D						
IGF-II	IL-15							
IL-1 RI	IL-16							

## TECHNICAL HINTS

- Protect the microparticles and streptavidin-PE from light at all times.
- Refer to the Base Kit Standard Value Card for reconstitution volume and values of the reconstituted standard.
- Diluted microparticles cannot be stored. Make a fresh dilution of microparticles each time the assay is run.
- The use of a magnetic device made to accommodate a microplate is necessary for washing.
- Discrepancies may exist in values obtained for the same analyte utilizing different technologies.

Luminex Performance Assays afford the user the benefit of multianalyte analysis of biomarkers in a complex sample. For each sample type, a single, multipurpose diluent is used to optimize recovery, linearity, and reproducibility. Such a multipurpose diluent may not optimize any single analyte to the same degree that a unique diluent selected for analysis of that analyte can optimize conditions. Therefore, some performance characteristics may be more variable than those for assays designed specifically for single analyte analysis.