



MagCollect™ Human NK Cell Isolation Kit

Catalog Number: MAGH109

INTENDED USE

The MagCollect Human NK Cell Isolation Kit is designed to isolate Natural Killer cells via a negative selection principle. The resulting cell preparation is highly enriched with NK cells. Typical purity of recovered NK cells ranges from 80-90%.

BACKGROUND

R&D Systems MagCollect products are designed for the isolation of cells in a "liquid phase". MagCollect technology is based on the use of ferrofluids or magnetic nanoparticles that have no magnetic memory (superparamagnetic) and behave like colloidal particles. This feature allows the ferrofluids to remain in solution without the need for mixing and additionally allows for efficient diffusion kinetics during the binding reaction. The proprietary manufacturing technology of MagCollect Ferrofluids generates particles with higher ligand binding capacity per mass compared to many other larger diameter magnetic particles.

PRINCIPLE OF SELECTION

A mononuclear cell suspension is first incubated with the MagCollect Antibody Cocktail which targets unwanted cells. MagCollect Streptavidin Ferrofluid is added to the reaction and the streptavidin-coated nanoparticles interact with the biotinylated antibody tagged cells. The tube containing the cell suspension is then placed within a magnetic field. Magnetically tagged cells will migrate toward the magnet (unwanted cell fraction), leaving the untagged cells or desired cell population in suspension. This population of cells can then be harvested by aspiration while the tube remains in the magnetic field. The enriched cell preparation for a variety of applications including tissue culture, immune status monitoring and flow cytometry.

MATERIALS PROVIDED & STORAGE CONDITIONS

Store the unopened kit at 2-8 °C. **DO NOT FREEZE.**

This kit contains sufficient reagents to process 1×10^9 total cells; up to 25 isolations.

| PART | PART # | DESCRIPTION | STORAGE OF OPENED/DILUTED MATERIAL |
|--|--------|--|---|
| Human NK Cell Biotinylated Antibody Cocktail | 860124 | 1.0 mL of a phosphate buffered solution containing BSA. | May be stored at 2-8 °C when handled aseptically.* |
| Streptavidin Ferrofluid | 860127 | 2 vials (1.25 mL/vial) of a solution containing BSA and preservatives. | |
| 10X Buffer | 860040 | 10 mL of a 10-fold concentrated buffer. | May be stored for up to 24 hours at 2-8 °C after dilution.* |

* Provided this is within the expiration date of the kit.

OTHER SUPPLIES REQUIRED

- MagCollect Magnet (R&D Systems, Catalog # MAG997) or equivalent
- Human Erythrocyte Lysing Kit (R&D Systems, Catalog # WL1000)
- 12 x 75 mm (5 mL) polystyrene round bottom tubes
- Sterile Pasteur pipettes or transfer pipettes (ThermoFisher, Catalog # 13-711-9B) or equivalent
- Sterile deionized or distilled water

PRECAUTION

Some components of this kit contain sodium azide, which may react with lead and copper plumbing to form explosive metallic azides. Flush with large volumes of water during disposal.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

REAGENT PREPARATION

Prepare 5 mL of 1X MagCollect Buffer for each 5×10^7 cells to be processed by mixing 0.5 mL of 10X Buffer with 4.5 mL sterile deionized or distilled water. **The buffer must be kept cold (2-8 °C) for the following procedure.**

CELL PREPARATION

1. Process cells on a density gradient, like Ficoll Hypaque to enrich for mononuclear cells.
2. Recover the "buffy coat" containing the mononuclear cells and wash the cells two times with excess PBS to remove any residual separation media. This can be done by spinning the cells for 10 minutes at 200 x g.
3. After the second wash step, disrupt the cell pellet by "racking" the tube, resuspend the cells in H-Lyse Buffer from R&D Systems' Human Erythrocyte Lysing Kit (Catalog # WL1000) that has been diluted to 1X strength with sterile distilled water. Quickly vortex the tube (10 mL of 1X H-Lyse solution per 250 million cells is recommended).
4. Incubate the cells for 10 minutes at room temperature and then fill the tube with 1X Wash Buffer from the Lysing kit.
Note: *The Wash Buffer must also be diluted with sterile water to 1X strength prior to use.*
5. Spin the cells for 10 minutes at 200 x g and then resuspend the cells in a small volume of 1X MagCollect Buffer.
6. Perform a cell count and then adjust the cell concentration to 1×10^8 cells per mL with cold 1X MagCollect Buffer.
7. Continue the cell selection by referring to step 1 of the Cell Selection Procedure.

CELL SELECTION PROCEDURE

This procedure is for processing 5×10^7 total cells using 5 mL tubes and the MagCollect Magnet. For processing other cell numbers, please refer to the Technical Hints section of the insert. Cells and reagents should be kept cold using an ice bath or a refrigerator. **Reaction incubations must be carried out at 2-8 °C in a refrigerator and not in an ice bath to avoid excessively low temperatures that can slow the kinetics of the optimized reactions.**


1. Prepare a single cell suspension of human leukocytes by traditional methods or by following the instructions outlined in the Cell Preparation section on this insert. Cells must be suspended in cold 1X MagCollect Buffer prior to beginning the procedure and be at a cell density of 1×10^8 cells/mL.
2. Transfer 5×10^7 cells (0.5 mL volume) into a 5 mL polystyrene tube. Add 50 μ L of Human NK Cell Biotinylated Antibody Cocktail. Gently mix the cell-antibody suspension, avoiding bubble formation, and incubate at 2-8 °C in a refrigerator for 15 minutes.
3. Add 125 μ L of Streptavidin Ferrofluid to the cell suspension, mix gently and incubate at 2-8 °C in a refrigerator for 15 minutes.
4. At the end of the incubation period bring the volume of the reaction in the tube to 2 mL by adding 1.35 mL of 1X MagCollect Buffer. Mix gently to ensure that all reactants in the tube are in suspension.
5. Place the reaction tube in the MagCollect Magnet that has been positioned horizontally to accommodate 5 mL tubes and incubate for 8 minutes at room temperature. Magnetically tagged cells will migrate toward the magnet (these are the unwanted cells), leaving the untouched desired cells in suspension in the supernatant.
6. Recovery of desired cells is achieved as follows: While the tube is firmly held in the magnet, using a sterile Pasteur pipette or transfer pipette, carefully aspirate all of the reaction supernatant and place it in a new 5 mL tube. Remove the tube containing the magnetically trapped cells from the magnet, and discard.
7. To ensure that all of the magnetic nanoparticles have been removed, repeat the magnetic depletion (steps 5 and 6) with the new tube containing the recovered cells. The supernatant obtained at the end of these steps is the final depleted cell fraction containing the desired enriched NK cells. The cells are now ready for counting and further downstream applications.

TECHNICAL HINTS

- If sterile cells are required following cell selection, the entire procedure should be carried out in a laminar flow hood to maintain sterile conditions. Use sterile equipment when pipetting reagents that will be reused at a later date.
- Avoid antibody capping on cell surfaces and non-specific cell tagging by working quickly, keeping cells and solutions cold through the use of pre-cooled solutions and by adhering to the incubation times and temperatures specified in the protocol. Increased temperature and prolonged incubation times may lead to non-specific cell labeling thus lowering cell purity and yield.
- When processing different numbers of cells, follow the recommendations in the following table or observe the following guidelines: keep antibody cocktail and ferrofluid incubation times and temperatures the same; keep the cell density at 1×10^8 cells/mL; add 10 μ L of the antibody cocktail per 1×10^7 cells being processed; add 20 μ L of Streptavidin Ferrofluid per 1×10^7 cells being processed. A minimum of 75 μ L and a maximum of 300 μ L of Streptavidin Ferrofluid is required per isolation.
- Use the following table for recommended quantities to be used in steps 2-4 of the Cell Selection Procedure:

| Number of Cells in Starting Preparation | 5×10^7 | 1×10^8 |
|--|-----------------|-----------------|
| Reaction Volume | 0.5 mL | 0.5 mL |
| Human NK Cell Biotinylated Antibody Cocktail | 50 μ L | 100 μ L |
| Streptavidin Ferrofluid | 125 μ L | 200 μ L |

DATA EXAMPLES



Enrichment of NK cells from PBMCs using this MagCelect Human NK Cell Isolation Kit. Cells, before and after enrichment, were double-stained with APC-conjugated Anti-Human CD3 Antibody (R&D Systems, Catalog # FAB100A) and PE-conjugated Anti-Human CD56 Antibody (R&D Systems, Catalog # FAB2408P).