



Human TGF- β 1 ELISA

650 010 096 1 x 96 tests
650 010 192 2 x 96 tests

INTENDED USE

The TGF- β 1 ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of human transforming growth factor beta-1 levels in cell culture supernatants, human serum, plasma, or other body fluids. **The TGF- β 1 ELISA is for research use only. Not for use in diagnostic or therapeutic procedures.**

SUMMARY

Transforming growth factor- β (TGF- β) belongs to a family of dimeric 25 kDa polypeptides that are ubiquitously distributed in tissues and synthesized by many different cells (13). Three isoforms of transforming Growth Factor- β (TGF-beta1, beta-2 and beta-3) exist in mammals. They play critical roles in growth regulation and development. Each isoform is encoded by a unique gene on different chromosomes. All three of these growth factors are secreted by most cell types, generally in a latent form, requiring activation before they can exert biological activity. The TGF-betas possess three major activities: they inhibit proliferation of most cells, but can stimulate the growth of some mesenchymal cells; they exert immunosuppressive effects and they enhance the formation of extracellular matrix. Two types of membrane receptors possessing kinase activity are involved in signal transduction. The TGF-betas are involved in wound repair processes and in starting inflammatory reaction and then in the resolution through chemotactic attraction of inflammatory cells and fibroblast (15).

TGF- β 1 is the first recognized transforming growth factor (5), its subunits of each 12.5 kDa are bound via disulphide bridges. TGF- β 1 is inhibitive to T- and B cell proliferation as well as to maturation and activation of macrophages. It furthermore inhibits activity of natural killer cells and lymphokine activated killer cells and blocks production of cytokines.

Measurement of TGF- β 1 in blood has been advocated for diagnosis of various diseases. TGF- β 1 has been shown to be an organizer of responses to neurodegeneration (10).


In this context, it turned out to be interesting in monitoring Alzheimer's disease (18), Down's syndrome, AIDS and Parkinson's disease (11). Serum and cerebrospinal fluid levels of Multiple Sclerosis patients were shown to be of great value to monitor remission and acute phases (4, 21). TGF- β 1 is thought to play an important role in bone metabolism (22), it is considered a putative regulator of osteoclastic-osteoblastic interaction, thus it can be regarded as a marker for osteoporosis (14). TGF- β 1 is involved in the pathogenesis of glomerular diseases (3, 23) such as diabetic nephropathy and glomerulosclerosis (28). TGF- β 1 has been described to be functionally connected to major immune system abnormalities as in autoimmunity (SLE) (8). Serum levels have been shown to correlate with disease activity in autoimmune hepatitis (2). Elevated serum levels of TGF- β 1 are determined in Chronic fatigue syndrome patients (6) and in Guillain-Baire syndrome patients (24). An inverse correlation with disease activity was described for TGF- β 1 levels in Kawasaki disease (17) and patients with IgA deficiency (19).

TGF- β 1 has been confirmed to promote fibrotic processes, thus it is implicated in the myelofibrosis with myeloid metaplasia (16). Increased serum levels of TGF- β 1 in patients affected by thrombotic thrombocytopenic purpura implicate its function on bone marrow haematopoiesis (29, 25). Determination of circulating TGF- β 1 turned out to reflect the various stages in solid tumors as has been shown for cervical cancer (7), elevations were furthermore found in prostatic cancer (27), bladder cancer (9), and liver cancer (20).


Decreased levels of TGF- β 1 in the serum of sepsis and acute stroke patients (1, 12) may reflect the changing immunological-inflammatory status of these patients. Decreased TGF- β 1 serum levels were described for patients with acute Plasmodium falciparum malaria (26).

PRINCIPLES OF THE TEST


An anti-TGF- β 1 coating antibody is adsorbed onto microwells.




TGF- β 1 present in the sample or standard binds to antibodies adsorbed to the microwells; a HRP-conjugated monoclonal anti-TGF- β 1 antibody is added and binds to TGF- β 1 captured by the first antibody.



Following incubation unbound enzyme conjugated anti-TGF- β 1 is removed during a wash step and substrate solution reactive with HRP is added to the wells.



A coloured product is formed in proportion to the amount of TGF- β 1 present in the sample. The reaction is terminated by addition of acid and absorbance is measured at 450nm. A standard curve is prepared from seven TGF- β 1 standard dilutions and TGF- β 1 sample concentration determined.



REAGENTS PROVIDED

REAGENTS (store at 2-8°C)	QTY 1 plate	QTY 1 plates	RECONSTITUTION
96-wells precoated microtiter plate	1	2	Ready-to-use
Plate covers	2	4	
TGF- β 1 Standard:	2 vials	4 vials	See label on the vial
HRP conjugate monoclonal antibody to TGF- β 1	1 vial	2 vials	(0.1 ml)
Assay Buffer Concentrate*	1 vial	2 vials	(5 ml) 20X concentrate. Dilute in distilled water
Wash Buffer Concentrate*	1 bottle	2 bottles	(50 ml) 20X concentrate. Dilute in distilled water
Substrate Solution	1 vial	2 vials	(15 ml) Ready-to-use
Stop Solution (1 M Phosphoric acid)	1 vial	2 vials	(12 ml) Ready-to-use
1N HCl (pretreatment of samples)	1 vial	2 vials	(3 ml) Ready to use
1N NaOH (pretreatment of samples)	1 vial	2 vials	(3 ml) Ready to use
Blue Dye	1 vial	2 vials	(0.4 ml) Make a 1/250 dilution in the appropriate diluent
Green Dye	1 vial	2 vials	(0.4 ml) Make a 1/100 dilution in the appropriate diluent

MATERIALS REQUIRED BUT NOT PROVIDED

- 5 ml and 10 ml graduated pipettes
- 5 µl to 1,000 µl adjustable single channel micropipettes with disposable tips
- 50 µl to 300 µl adjustable multichannel micropipette with disposable tips
- Multichannel micropipette reservoir
- Beakers, flasks, cylinders necessary for preparation of reagents
- Device for delivery of wash solution (multichannel wash bottle or automatic wash system)
- Microwell strip reader capable of reading at 450 nm (620 nm as optional reference wave length)
- Glass-distilled or deionized water
- Statistical calculator with program to perform linear regression analysis

SAFETY

- Reagents are intended for research use only and are not for use in diagnostic or therapeutic procedures.
- Do not mix or substitute reagents with those from other lots or other sources.
- Do not use kit reagents beyond expiration date on label.
- Do not expose kit reagents to strong light during storage or incubation.
- Do not pipette by mouth.
- Do not eat or smoke in areas where kit reagents or samples are handled.
- Avoid contact of skin or mucous membranes with kit reagents or specimens.
- Rubber or disposable latex gloves should be worn while handling kit reagents or specimens.
- Reagents containing thimerosal as preservative may be toxic if ingested.
- Avoid contact of substrate solutions with oxidizing agents and metal.
- Avoid splashing or generation of aerosols.
- In order to avoid microbial contamination or cross-contamination of reagents or specimens which may invalidate the test use disposable pipette tips and/or pipettes.
- Use clean, dedicated reagent trays for dispensing the conjugate and substrate reagents.
- Exposure to acids will inactivate the conjugate.
- Glass-distilled water or deionized water must be used for reagent preparation.
- Substrate solutions must be at room temperature prior to use.
- Decontaminate and dispose specimens and all potentially contaminated materials as if they could contain infectious agents. The preferred method of decontamination is autoclaving for a minimum of 1 hour at 121.5°C.
- Liquid wastes not containing acid and neutralized waste may be mixed with sodium hypochlorite in volumes such that the final mixture contains 1.0 % sodium hypochlorite. Allow 30 minutes for effective decontamination. Liquid waste containing acid must be neutralized prior to the addition of sodium hypochlorite.

STORAGE INSTRUCTIONS

Store kit reagents between 2° and 8°C. Immediately after use remaining reagents should be returned to cold storage (2° to 8°C). Expiry of the kit and reagents is stated on box front labels.

The expiry of the kit components can only be guaranteed if the components are stored properly, and if, in case of repeated use of one component, the reagent is not contaminated by the first handling.

SPECIMEN COLLECTION, PROCESSING AND STORAGE

Cell culture supernatants, human serum, plasma or other biological samples will be suitable for use in the assay. Remove serum from the clot or red cells, respectively, as soon as possible after clotting and separation.

Pay attention to a possible “**hook effect**” due to high sample concentrations (see chapter 11).

Samples containing a visible precipitate must be clarified prior to use in the assay. Do not use grossly hemolyzed or lipemic specimens.

Clinical samples should be kept at 2° to 8°C and separated rapidly before storing at -20°C to avoid loss of bioactivity. Avoid repeated freeze-thaw cycles.

PREPARATION OF REAGENTS

1. Wash Buffer

If crystals have formed in the Wash Buffer Concentrate, warm it gently until they have completely dissolved.

Pour entire contents (50 ml) of the **Wash Buffer Concentrate** into a clean 1,000 ml graduated cylinder. Bring final volume to 1,000 ml with glass-distilled or deionized water. Mix gently to avoid foaming. The pH of the final solution should adjust to 7.4.

Transfer to a clean wash bottle and store at 2° to 25°C. Please note that the Wash Buffer is stable for 30 days. Wash Buffer may be prepared as needed according to the following table.

Number of Strips	Wash Buffer Concentrate (ml)	Distilled Water (ml)
1 - 6	25	475
1 - 12	50	950

2. Assay Buffer

Mix the contents of the bottle well. Add contents of **Assay Buffer Concentrate** (5.0 ml) to 95 ml distilled or deionized water and mix gently to avoid foaming. Store at 2° to 8°C. Please note that the Assay Buffer is stable for 30 days. Assay Buffer may be prepared as needed according to the following table:

Number of Strips	Assay Buffer Concentrate (ml)	Distilled Water (ml)
1 - 6	2.5	47.5
1 - 12	5.0	95.0

3. Preparation of HRP-Conjugate

The **HRP-Conjugate** must be diluted 1:100 with Assay Buffer just prior to use in a clean plastic test tube.

Please note that the HRP-Conjugate should be used within 30 minutes after dilution. HRP-Conjugate may be prepared as needed according to the following table:

Number of Strips	HRP-Conjugate (ml)	Assay Buffer (ml)
1 - 6	0.03	2.97
1 - 12	0.06	5.94

4. Preparation of Standard

Dilute TGF-β1 **Standard** by addition of **Assay Buffer**. Dilution volume is stated on the label of the standard vial. Mix gently to ensure complete homogeneity.

5. Addition of Colour-giving Dyes

In order to help our customers to avoid any mistakes in pipetting, DIACLONE now offers a new tool that helps to monitor the addition of even very small volumes of a solution to the reaction well by giving distinctive colours to each step of the ELISA procedure.

This procedure is optional, does not in any way interfere with the test results, and is designed to help the customer with the performance of the test, but can also be omitted, just following the instruction booklet. Alternatively, the dye solutions from the stocks provided (**Blue-Dye, Green-Dye**) can be added to the reagents according to the following guidelines:

A. Diluent:

Before sample dilution add the **Blue-Dye** at a dilution of 1:250 (see table below) to the appropriate diluent (1x) according to the test protocol. After addition of **Blue-Dye**, proceed according to the instruction booklet.

5 ml Assay Buffer	20 µl Blue-Dye
-------------------	-----------------------

	1	2	3	4	5	6	7	8	9	10	11	12
D	3.75	3.75										
E	1.88	1.88										
F	0.94	0.94										
G	0.5	0.5										
H	Blank	Blank										

- f. Add 100 µl of **Assay Buffer**, in duplicate, to the blank wells.
- g. Add 100 µl of each pretreated **Sample**, in duplicate, to the designated wells.
- h. Prepare **HRP-Conjugate**. (Refer to preparation of reagents)
- i. Add 50 µl of diluted **HRP-Conjugate** to all wells.
- j. Cover with a **Plate Cover** and incubate at room temperature (18°C to 25°C) for 4 hours on a **rotator set at 100 rpm**.
- k. Remove **Plate Cover** and empty wells. Wash microwell strips 3 times according to point d. of the test protocol. Proceed immediately to the next step.
- l. Pipette 100 µl of **TMB Substrate Solution** to all wells, including the blank wells.
- m. Incubate the microwell strips at room temperature (18° to 25°C) for about 10 minutes, if available on a rotator set at 100 rpm. Avoid direct exposure to intense light. The point at which the substrate reaction is stopped is often determined by the ELISA reader being used. Many ELISA readers record absorbance only up to 2.0 O.D. **The O.D. values at the plate should be monitored and the substrate reaction stopped before positive wells are no longer properly recordable.**
- n. Stop the enzyme reaction by pipetting 100 µl of **Stop Solution** into each well, including the blank wells. It is important that the Stop Solution is spread quickly and uniformly throughout the microwells to completely inactivate the enzyme. Results must be read immediately after the Stop Solution is added or within one hour if the microwell strips are stored at 2 - 8°C in the dark.
- o. Read absorbance of each microwell on a spectro-photometer using 450 nm as the primary wave length (optionally 620 nm as the reference wave length; 610 nm to 650 nm is acceptable). Blank the plate reader according to the manufacturer's instructions by using the blank wells. Determine the absorbance of both, the samples and the TGF-β1 standards.

CALCULATION OF RESULTS

- Calculate the average absorbance values for each set of duplicate standards and samples. Duplicates should be within 20 per cent of the mean.
- Create a standard curve by plotting the mean absorbance for each standard concentration on the ordinate against the TGF-β1 concentration on the abscissa. Draw a best fit curve through the points of the graph.
- To determine the concentration of circulating TGF-β1 for each sample, first find the mean absorbance value on the ordinate and extend a horizontal line to the standard curve. At the point of intersection, extend a vertical line to the abscissa and read the corresponding TGF-β1 concentration.


For samples which have been diluted according to the instructions given in this manual (e.g. 1:12) the concentration read from the standard curve must be multiplied by the dilution factor (e.g. x12).

Note: Calculation of samples with an O.D. exceeding 2.0 may result in incorrect, low TGF-β1 levels. Such samples require further dilution with Standard / Sample Diluent in order to precisely quantitate the actual TGF-β1 level.

It is suggested that each testing facility establishes a control sample of known TGF-β1 concentration and runs this additional control with each assay. If the values obtained are not within the expected range of this control, the assay results may be invalid.

- A representative standard curve is shown in Figure 3. This curve cannot be used to derive test results. Every laboratory must prepare a standard curve for each group of microwell strips assayed.

Figure 3. Representative standard curve for TGF-β1 ELISA. Recombinant soluble TGF-β1 was diluted in Assay Buffer; symbols represent the mean of three parallel titrations. Do not use this standard curve to derive test results. A standard curve must be run for each group of microwell strips assayed.



Typical data using the TGF-β1ELISA

Measuring wavelength: 450 nm
 Reference wavelength: 620 nm

Standard	TGF-β1 Concentration (ng/ml)	O.D. Mean
1	30 30	1.867
2	15 15	1.037
3	7.5 7.5	0.5665
4	3.75 3.75	0.287
5	1.88 1.88	0.152
6	0.94 0.94	0.0835
7	0.47 0.47	0.0555
Blank	0 0	0.008

LIMITATIONS OF THE PROCEDURE

- Since exact conditions may vary from assay to assay, a standard curve must be established for every run.
- Bacterial or fungal contamination of either samples or reagents or cross-contamination between reagents may cause erroneous results.
- Disposable pipette tips, flasks or glassware are preferred, reusable glassware must be washed and thoroughly rinsed of all detergents before use.
- Improper or insufficient washing at any stage of the procedure will result in either false positive or false negative results. Completely empty wells before dispensing fresh Wash Buffer, fill with Wash Buffer as indicated for each wash cycle and do not allow wells to sit uncovered or dry for extended periods.

PERFORMANCE CHARACTERISTICS

Sensitivity

The limit of detection for recombinant TGF- β 1, spiked into normal human serum, defined as the analyte concentration resulting in an absorption significantly higher than the absorption of the dilution medium (mean plus two standard deviations) was determined to be 23.76pg/ml (mean of three independent assays).

Reproducibility

a. Intra-assay

Intra-assay variability was determined by 6 replicates of 8 serum samples. The average coefficient of variation was 6.7%.

b. Inter-assay

Inter-assay variability was determined by 18 replicates of 8 serum samples. The average coefficient of variation was 8.5%.

Specificity

The ELISA was shown to specifically detect human TGF- β 1. No cross reactivity was apparent with human serum factors spiked into human serum, notably there was no cross reactivity with human TGF- β 2 and TGF- β 3.

Spike Recovery

Spike recovery was determined to be 102% in average of recombinant TGF- β 1 spiked into various samples.

Dilution Linearity

Linearity of dilution was measured in various samples. Recovery was 110% in average.

Expected values

A panel of 40 serum samples from randomly selected apparently healthy donors (males and females) was tested for human TGF- β 1.

The detected human TGF- β 1 levels ranged between not detectable and 3.46 ng/ml with a mean level of 1.96 ng/ml and a standard deviation of 0.79 ng/ml.

BIBLIOGRAPHY

1. Astiz M, Saha D, Lustbader D, Lin R, Rackow E. Monocyte response to bacterial toxins, expression of cell surface receptors, and release of anti-inflammatory cytokines during sepsis. J Lab Clin Med 1996 Dec;128(6):594-600.
2. Bayer EM, Herr W, Kanzler S, Waldmann C, Meyer Zum Buschenfelde KH, Diense HP; Lohse AW. Transforming growth factor-beta1 in autoimmune hepatitis: correlation of liver tissue expression and serum levels with disease activity. J Hepatol 1998 May;28(5):803-811.
3. Border WA. Transforming growth factor-beta and the pathogenesis of glomerular diseases. Curr Opin Nephrol Hypertens 1994 Jan;3(1):54-58.
4. Carrieri PB, Provitera V, Bruno R, Perrella M, Tartaglia G, Busto A, Perrella O. Possible role of transforming growth factor-beta in relapsing-remitting multiple sclerosis. Neurol Res 1997 Dec;19(6):599-600.
5. Chambaz EM, Souchevnik S, Pellerin S, Defaye G, Cochet C, Feige JJ. Transforming growth factors-beta s: a multifunctional cytokine family. Implication in the regulation of adrenocortical cell endocrine functions. Horm Res 1996;45(3-5):222-226.
6. Chao CC, Janoff EN, Hu SX, Thomas K, Gallagher M, Tsang M, Peterson PK. Altered cytokine release in peripheral blood mononuclear cell cultures from patients with the chronic fatigue syndrome.

Cytokine 1991 Jul;3(4):292-298.

7. Chopra V, Dinh TV, Hannigan EV. Circulating serum levels of cytokines and angiogenic factors in patients with cervical cancer. *Cancer Invest* 1998;16(3):152-159.
8. Del Giudice G, Crow MK. Role of transforming growth factor beta (TGF beta) in systemic autoimmunity. *Lupus* 1993 Aug;2(4):213-220.
9. Eder IE, Stenzl A, Hobisch A, Cronauer MV, Bartsch G, Klocker H. Expression of transforming growth factors beta-1, beta 2 and beta 3 in human bladder carcinomas. *Br J Cancer* 1997;75(12):1753-1760.
10. Finch CE, Laping NJ, Morgan TE, Nichols NR, Pasinetti GM. TGF-beta 1 is an organizer of responses to neurodegeneration. *J Cell Biochem* 1993 Dec;53(4):314-322.
11. Flanders KC, Ren RF, Lippa CF. Transforming growth factor-betas in neurodegenerative disease. *Prog Neurobiol* 1998 Jan;54(1):71-85.
12. Kim JS, Yoon SS, Kim YH, Ryu JS. Serial measurement of interleukin-6, transforming growth factor-beta, and S-100 protein in patients with acute stroke. *Stroke* 1996 Sep;27(9):1553-1557.
13. Kropf J, JO Schurek, A Wollner, and AM Gressner. Immunological measurement of transforming growth factor-beta I (TGF- β 1) in blood; assay development and comparison. *Clinical Chemistry* 1997;43(10):1965-1974.
14. Langdahl BL, Knudsen JY, Jensen HK, Gregersen N, Eriksen EF. A sequence variation: 713-8delC in the transforming growth factor-beta 1 gene has higher prevalence in osteoporotic women than in normal women and is associated with very low bone mass in osteoporotic women and increased bone turnover in both osteoporotic and normal women. *Bone* 1997 Mar;20(3):289-294.
15. Lawrence DA. Transforming growth factor-beta: a general review. *Eur Cytokine Netw* 1996 Sep;7(3):363-374.
16. Martyre MC, Magdelenat H, Bryckaert MC, Laine-Bidron C, Calvo F. Increased intraplatelet of platelet-derived growth factor and transforming growth factor-beta in patients with myelofibrosis with myeloid metaplasia. *Br J Haematol* 1991 Jan;77(1):80-86.
17. Matsubara T, Umezawa Y, Tsuru S, Motohashi T, Yabuta K, Furukawa S. Decrease in the concentrations of transforming growth factor-beta 1 in the sera of patients with Kawasaki disease. *Scand J Rheumatol* 1997;26(4):314-317.
18. Mattson MP, Barger SW, Furukawa K, Bruce AJ, Wyss-Coray T, Mark RJ, Mucke L. Cellular signaling roles of TGF beta, TNF alpha and beta APP in brain injury responses and Alzheimer's disease. *Brain Res Rev* 1997 Feb;23(1-2):47-61.
19. Muller F, Aukrust P, Nilssen DE, Froland SS. Reduced serum level of transforming growth factor-beta in patients with IgA deficiency. *Clin Immunol Immunopathol* 1995 Aug;76(2):203-208.
20. Murawaki Y, Ikuta Y, Nishimura Y, Koda M, Kawasaki H. Serum markers for fibrosis and plasma transforming growth factor-beta 1 in patients with hepatocellular carcinoma in comparison with patients with liver cirrhosis. *J Gastroenterol Hepatol* 1996 May;11(5):443-450.
21. Ossege LM, Sindern E, Voss B, Malin JP. Corticosteroids induce expression of transforming-growth-factor-beta1 mRNA in peripheral blood mononuclear cells of patients with multiple sclerosis. *J Neuroimmunol* 1998 Apr 1;84(1):1-6.
22. Pfeilschifter J, Diel I, Scheppach B, Bretz A, Krempien R, Erdmann J, Schmid G, Reske N, Bismar H, Seck T, Krempien B, Ziegler R. Concentration of transforming growth factor beta in human bone tissue: relationship to age, menopause, bone turnover, and bone volume.

23. Shankland SJ, Johnson RJ. TGF-beta in glomerular disease. *Miner Electrolyte Metab* 1998;24(2-3):168-173.
24. Sindern E, Schweppe K, Ossege LM, Malin JP. Potential role of transforming growth factor-beta 1 in terminating the immune response in patients with Guillain-Barre syndrome. *J Neurol* 1996 Mar;243(3):264-268.
25. Tornquist SJ, Oaks JL, Crawford TB. Elevation of cytokines associated with the thrombocytopenia of equine infectious anaemia. *J Gen Virol* 1997 Oct;78(Pt 10):2541-2548.
26. Wenisch C, Parschalk B, Burgmann H, Looareesuwan S, Graninger W. Decreased serum levels of TGF-beta in patients with acute *Plasmodium falciparum* malaria. *J Clin Immunol* 1995 Mar;15(2):69-73.
27. Wolff JM, Fandel T, Borchers H, Brehmer B Jr, Jakse G. Transforming growth factor-beta 1 serum concentration in patients with prostatic cancer and benign prostatic hyperplasia. *Br J Urol* 1998 Mar;81(3):403-405.
28. Yokoyama H, Deckert T. Central role of TGF-beta in the pathogenesis of diabetic nephropathy and macrovascular complications: a hypothesis. *Diabet Med* 1996 Apr;13(4):313-320.
29. Zauli G, Gugliotta L, Catani L, Vianelli N, Borgatti P, Belmonte MM, Tura S. Increased serum levels of transforming growth factor beta-1 in patients affected by thrombotic thrombocytopenic purpura (TTP): its implications on bone marrow haematopoiesis. *Br J Haematol* 1993 Jul;84(3):381-386.

REAGENT PREPARATION SUMMARY

- A. Wash Buffer** Add **Wash Buffer Concentrate** 20 X (50 ml) to 950 ml distilled water
- B. Assay Buffer**
- | Number of Strips | Assay Buffer Concentrate (ml) | Distilled Water (ml) |
|------------------|-------------------------------|----------------------|
| 1 - 6 | 2.5 | 47.5 |
| 1 - 12 | 5.0 | 95.0 |
- C. HRP-Conjugate**
- | Number of Strips | HRP-Conjugate (ml) | Assay Buffer (ml) |
|------------------|--------------------|-------------------|
| 1 - 6 | 0.03 | 2.97 |
| 1 - 12 | 0.06 | 5.94 |
- D. Standard** Reconstitute TGF- β 1 Standard by addition of distilled water as stated on vial label.

TEST PROTOCOL SUMMARY

- Pretreat Samples
- Wash Microwell Strips twice with Wash Buffer
- Add 100 μ l Assay Buffer, in duplicate, to all standard wells
- Pipette 100 μ l TGF- β 1 Standard into the first standard wells and create standard dilutions ranging from 30 to 0.5 ng/ml by transferring 100 μ l from well to well; discard 100 μ l from the last wells.
- Add 100 μ l Assay Buffer, in duplicate, to the blank wells
- Add 100 μ l pretreated Sample to designated wells
- Prepare HRP-Conjugate
- Add 50 μ l diluted HRP-Conjugate to all wells
- Cover microwell strips and incubate 4 hours at room temperature (18° to 25°C) on a rotator set at 100 rpm
- Empty and wash microwell strips 3 times with Wash Buffer
- Add 100 μ l of TMB Substrate Solution to all wells including blank wells
- Incubate the microwell strips for about 15 minutes at room temperature (18° to 25°C)
- Add 100 μ l Stop Solution to all wells including blank wells
- Blank microwell reader and measure colour intensity at 450 nm
- **Note: Calculation of samples with an O.D. exceeding the range of the standard curve may result in incorrect, low TGF- β 1 levels. Such samples require further dilution with Standard/Sample Diluent in order to precisely quantitate the actual TGF- β 1 level.**