Instruction

Wieslab[®] hCOMP quantitative kit

Enzyme immunoassay for detection of human COMP
(Cartilage Oligomeric Matrix Protein)

Microtitration strips (12x8) 96 wells
Store the kit at +2-8° C

For Research Use Only. Not for use in diagnostic procedures.


Document No. E-23-0179-04 RUO
October, 2013

Wieslab[®] hCOMP quantitative kit

English: page 2
Français: page 10
Español: página 16
Deutsch: Seite 22
Italiano: pagina 28
Svenska: sida..... 34

REF

COMP 200 RUO

PURPOSE OF RESEARCH PRODUCT

The Wieslab[®] hCOMP is an enzyme immunoassay for quantitation of intact or fragmented human Cartilage Oligomeric Matrix Protein (COMP) in serum or synovial fluid. The assay serves as an indicator of COMP turnover. The analysis should be performed by trained laboratory professionals. The result shall not be used for clinical diagnosis or patient management.

Summary and explanation

COMP is a high molecular weight, multi-subunit protein originally isolated from cartilage (Hedbom et al., 1992; Mörgelin et al., 1992). The protein is abundant in cartilage but is also found in tendon and other tissues (DiCesare et al., 1994). The protein belongs to the thrombospondin family (Oldberg et al., 1992).

Measurement of intact COMP and fragments thereof in synovial fluid or serum have been shown to correlate to cartilage destruction in rheumatoid arthritis and osteoarthritis patient studies (Månsson et al., 1995; Saxne and Heinegård, 1992; Sharif et al., 1995). In synovial fluid COMP is present to some extent as an apparently intact protein but the majority is found as several different fragments (Neidhart et al., 1997). In serum the fragments are 50-90 kDa. The exact nature of the fragments is not known.

Principle of the Wieslab[®] COMP assay

The assay utilises native human articular cartilage COMP coated to 96-well microtiter plates and a rabbit polyclonal antiserum directed to human COMP. It is in standard ELISA inhibition format with an overnight preincubation step with sample and primary antiserum. After the overnight preincubation the solution is transferred to the COMP coated plate. Bound antibody is detected using an alkaline phosphatase labelled anti rabbit IgG conjugate.

Warnings and precautions

- **For Research Use Only.**
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of COMP in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.
- Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

Warning

| | | |
|------|------|-----|
| BUF | WASH | 30X |
| DIL | | |
| CONJ | | |
| CAL | | |

| | |
|-----------|------|
| CONTROL | C1 |
| CONTROL | C2 |
| SUBS | pNPP |
| Anti-COMP | |

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

- H317: May cause an allergic skin reaction.
P264: Wash hands thoroughly after handling.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P302+352: IF ON SKIN: Wash with plenty of soap and water.
P333+313: If skin irritation or rash occurs: Get medical advice/attention.

Kit components and storage of reagents

- One frame with strips (12x8) coated with COMP and one lid sealed in a foil pack with a dry pack.
- Six calibrators 600 µL/each, (10, 20, 30, 40, 60, 80 ng/mL), containing COMP in diluent, sufficient for three calibrator curves.
- Two prediluted control sera (C1 and C2), 600 µL each.
- 2 x 32 mL diluent i.e. sample dilution buffer (red colour).
- 8 mL anti-COMP reagent for preincubation (red colour).
- 13 mL conjugate buffer with alkaline phosphatase labelled conjugate directed against rabbit IgG (blue colour).
- 30 mL wash buffer 30x concentrated.
- 13 mL pNPP-substrate solution

All reagents, except the wash solution, are ready for use. For storage keep all components refrigerated at 4° C. Remove only the number of COMP coated strips needed for testing, resealing the foil package carefully.

Materials or equipment required but not provided

- 96-well round bottom preincubation plate.
- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- NaCl, Tween 20 and distilled water for preparing more washing solution.
- Washer for strips, absorbent tissue, tubes and a timer.

Specimen collection

Serum samples are quite stable with normal handling and refrigerated storage for shorter periods. However, for extended periods of storage samples should be stored frozen at –20° C. Freezing and thawing of serum samples have little effect on COMP levels but should always be avoided since precipitations may occur. Synovial fluid samples appear to be more sensitive to degradation. Avoid using sera which are icteric, lipemic and hemolyzed.

The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

PROCEDURE

All calibrators, controls and samples should be tested in duplicate.
All solutions should be allowed to reach room temperature before use.

- FIRST DAY (AFTERNOON) -

Dilution of serum and preincubation with first antibody

1. Dilute serum samples 1/50 with the dilution buffer. Dilute synovial fluid 1/500.
(Serum samples should be diluted 1/30 or more. Recommended dilution of serum samples is 1/50. Synovial fluids should be diluted 1/200 or more. Recommended starting dilution is 1/500.)

2. Pipet 75 μ L of calibrators C10, C20, C30, C40, C60, C80 and controls C1, C2 and the diluted samples (S) into a preincubation plate (round bottom plate) according to the diagram below.

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|---|---|---|---|---|---|----|----|----|
| A | C80 | C20 | S1 | | | | | | | | | |
| B | C80 | C20 | S1 | | | | | | | | | |
| C | C60 | C10 | S2 | | | | | | | | | |
| D | C60 | C10 | S2 | | | | | | | | | |
| E | C40 | C1 | etc | | | | | | | | | |
| F | C40 | C1 | | | | | | | | | | |
| G | C30 | C2 | | | | | | | | | | |
| H | C30 | C2 | | | | | | | | | | |

3. Add 75 μ L of anti-COMP reagent to all wells, mix on a shaker for 5 min. Incubate overnight (12-18 hrs) at about 4° C

- SECOND DAY -

All solutions should be allowed to reach room temperature before use.

Preparation of washing solution

Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

If necessary, more wash solution can be prepared by dissolving 9.0 g NaCl in 1.0 L distilled water, add 0.5 mL Tween 20, and mix.

Incubation in the antigen coated plate

1. Transfer 100 μ L from each well in the preincubation plate to the antigen coated plate.
2. Incubate at room temperature for 60 minutes.

After sample incubation

Wash 3 times with washing solution (approximately 300 μ L/well), by filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding conjugate

Add 100 μ L conjugate to each well. Incubate for 60 minutes.

After conjugate incubation

Wash as before.

Adding substrate solution

Add 100 μ L pNPP-substrate solution to each well, incubate for 60 minutes \pm 10 minutes. Read the absorbance at 405 nm on a microplate reader.

Analysis of results

Calculate amount of COMP in sample by using appropriate software. Best curve fitting results are obtained using a curve fitting algorithm like spline. To calculate levels in serum or other samples remember to adjust for dilutions.

Quality Control

The control values of C1 and C2 and the limits are shown in the QC certificate. Actual control ranges should be established in your laboratory.

Expected results

Serum levels of COMP in 80 unselected blood donors age 22-59 were found to be 1.35 ± 0.40 μ g/mL (range 0.99 – 2.54). No normal values are available for synovial fluids but values are generally about ten to twenty times higher than serum levels. In a clinical material of knee synovial fluids COMP ranged from 14 to 48 μ g/mL.

Performance characteristics

Intra-assay precision was determined by testing two samples in duplicate. Each sample was assayed 20 times on the same plate.

| Sample | Mean ($\mu\text{g/mL}$) | S.D. | C.V.% |
|--------|---------------------------|------------|------------|
| 1 | 1.03 | ± 0.11 | ± 10.7 |
| 2 | 1.71 | ± 0.06 | ± 3.5 |

Inter-assay precision was determined by testing two samples in duplicate. Results were obtained from fifteen different runs.

| Sample | Mean ($\mu\text{g/mL}$) | S.D. | C.V.% |
|--------|---------------------------|------------|------------|
| 1 | 0.97 | ± 0.10 | ± 10.3 |
| 2 | 1.82 | ± 0.12 | ± 6.6 |

Troubleshooting

| Problem: | Possible causes: | Solution: |
|------------------------------|--|--|
| Control values out of range. | <ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. | <ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread. |
| All test results negative. | <ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. | <ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagents. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread. |
| All test results yellow. | <ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. | <ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water solution used to prepare. 3. Repeat test. |
| Poor precision. | <ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. | <ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips. |

References

1. **DiCesare, P., Hauser, N., Lehman, D., Pasumarti, S., and Paulsson, M. (1994).** Cartilage oligomeric matrix protein (COMP) is an abundant component of tendon. *FEBS Lett.* 354, 237-240.
2. **Hedbom, E., Antonsson, P., Hjerpe, A., Aeschlimann, D., Paulsson, M., Rosa Pimentel, E., Sommarin, Y., Wendel, M., Oldberg, A., and Heinegård, D. (1992).** Cartilage matrix proteins. An acidic oligomeric protein (COMP) detected only in cartilage. *J.Biol.Chem.* 267, 6132-6136.
3. **Mörgelin, M., Heinegård, D., Engel, J., and Paulsson, M. (1992).** Electron microscopy of native cartilage oligomeric matrix protein purified from the Swarm rat chondrosarcoma reveals a five-armed structure. *J.Biol.Chem.* 267, 6137-6141.
4. **Månsson, B., Carey, D., Alini, M., Ionescu, M., Rosenberg, L.C., Poole, A.R., Heinegård, D., and Saxne, T. (1995).** Cartilage and bone metabolism in rheumatoid arthritis. Differences between rapid and slow progression of disease identified by serum markers of cartilage metabolism. *J.Clin.Invest.* 95, 1071-1077.
5. **Neidhart, M., Hauser, N., Paulsson, M., DiCesare, P., Michel, B.A., and Häuselmann, H.J. (1997).** Small fragments of cartilage oligomeric matrix protein in synovial fluid and serum as markers for cartilage degradation. *Br.J.Rheum.* 36, 1151-1160.
6. **Oldberg, Å., Antonsson, P., Lindblom, K., and Heinegård, D. (1992).** COMP (cartilage oligomeric matrix protein) is structurally related to the thrombospondins. *J.Biol.Chem.* 267, 22346-22350.
7. **Saxne, T. and Heinegård, D. (1992).** Cartilage oligomeric matrix protein: a novel marker of cartilage turnover detectable in synovial fluid and blood [published erratum appears in *Br J Rheumatol* 1993 Mar;32(3):247]. *Br.J.Rheumatol.* 31, 583-591.
8. **Sharif, M., Saxne, T., Shepstone, L., Kirwan, J.R., Elson, C.J., Heinegård, D., and Dieppe, P.A. (1995).** Relationship between serum cartilage oligomeric matrix protein levels and disease progression in osteoarthritis of the knee joint. *Br.J.Rheumatol.* 34, 306-310.
9. **Volck, B., Johansen, J.S., Stoltenberg, M., Garbarsch, C., Price, P.A., Østergaard, M., Østergaard, K., Løvgreen-Nielsen, P., Sonne-Holm, S. and Lorenzen, I. (2001).** Studies on YKL-40 in knee joints of patients with rheumatoid arthritis and osteoarthritis. Involvement of YKL-40 in the joint pathology.
10. **Garnero, P., Piperno, M., Gineyts, E., Christgau, S., Delmas P.D. and Vignon, E. (2001).** Cross sectional evaluation of biochemical markers of bone, cartilage, and synovial tissue metabolism in patients with knee osteoarthritis: relations with disease activity and joint damage. *Ann. Rheum. Dis.* 60, 619-626.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Förklaringar till symboler.

| | |
|---|---|
|  | Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. Usare entro. Använd före |
|  | Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Biologisk risk. |
|  | Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Förvaringstemperatur. |
|  | Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Tillverkare. |
|  | Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Satsnummer. |
|  | Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Katalognummer. |
|  | Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Läs instruktionsmanualen. |
|  | Warning. Attention. Atención. Achtung. Attenzione. Varning. |
|  | Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Innehåller tillräckligt för 96 test. |

| | |
|---------------------------------------|---|
| Ag | Antigen. Antigène, Antigeno. Antigene. L'antigene. |
| DIL | Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Spädningsbuffert. |
| CONJ | Conjugate. Conjugué. Conjugado. Conjugato. Konjugat. |
| BUF WASH 30X | Washsolution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x conc. Soluzione di lavaggio 30x conc. Vaskebuffert/Tvättbuffert 30x conc. |
| SUBS pNPP | Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP. |
| CAL X | Calibrator. Etalon. Calibrador. Calibratore. Kalibrator. |
| CONTROL X | Control. Kontrolle. Contrôle. Controllo. Kontroll. |
| Anti-COMP | Anti-COMP reagent, Anti-COMP Reagenz, Anti-COMP Reactivo, Reagente anti-COMP, Anti-COMP reagens |

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EURO DIAGNOSTICA AB